BD BBL™ Medium Supplement for the Selection of Pathogenic *Neisseria*

8810281(03)
2019-02

English: pages 1 – 3

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Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyks lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviseletétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL™ V-C-N Inhibitor is an antibiotic mixture of vancomycin, colistin and nystatin which is incorporated into culture media to permit selective isolation of *Neisseria gonorrhoeae* and *N. meningitidis*.

BD BBL™ V-C-N-T Inhibitor contains the same antibiotic mixture as BD BBL V-C-N Inhibitor with the addition of trimethoprim to improve recovery of pathogenic *Neisseria* by increasing selectivity of isolation media.

SUMMARY AND EXPLANATION

Media containing BD BBL V-C-N Inhibitor, e.g., Thayer-Martin Selective Agar,¹ are used for isolation of pathogenic *Neisseria* from the throat, vagina, rectum or urethra.¹⁻³ Not only is the overgrowth of gonococci and meningococci by contaminants minimized, but the “saphrophytic” *Neisseria* species are almost completely suppressed.

Martin et al. modified Thayer-Martin Selective Agar by adding trimethoprim to produce Transgrow Medium in bottles with a carbon dioxide-enriched atmosphere and Modified Thayer-Martin Agar in plates.^{4,5} A significantly greater number of positive gonococcal isolates from clinical specimens was reported as compared with Thayer-Martin Selective Agar due to the inhibition of swarming *Proteus* species.⁵⁻⁷ Because of its improved performance, it is recommended over earlier formulations for the isolation of *N. gonorrhoeae*.⁸⁻¹⁰

PRINCIPLES OF THE PROCEDURE

BD BBL V-C-N Inhibitor contains vancomycin to inhibit gram-positive contaminants, colistin to inhibit gram-negative bacteria, including *Pseudomonas* species, and nystatin to suppress the growth of yeasts.

BD BBL V-C-N-T Inhibitor contains the same antibiotic mixture as BD BBL V-C-N Inhibitor with the addition of trimethoprim, which inhibits *Proteus* species.

REAGENTS

Approximate Formula* per 1 mL Reconstituted Solution:

BD BBL V-C-N Inhibitor: Vancomycin300 µg
Colistin750 µg
Nystatin1250 units

BD BBL V-C-N-T Inhibitor: Vancomycin300 µg
Colistin750 µg
Nystatin1250 units
Trimethoprim Lactate500 µg

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions

For Laboratory Use.

BD BBL V-C-N Inhibitor and BD BBL V-C-N-T Inhibitor are for use in culture media and not for use in human or animal therapy. Observe aseptic techniques in the reconstitution and addition of these media supplements. Read “Instructions.”

This product contains dry natural rubber.

Danger

H301 Toxic if swallowed.

P264 Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P310** IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. **P321** Specific treatment (see on this label). **P330** Rinse mouth.

P405 Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions and Reconstitution: On receipt, store at -20 to +8 °C. After reconstitution, use immediately or store below -20 °C and use within two weeks. Avoid repeated freezing and thawing.

Using a sterile syringe and needle, reconstitute each lyophilized vial by aseptically adding the following volumes of sterile purified water: 212227–2 mL; 212228 and 212408–10 mL.

The expiration date applies to the lyophilized product, stored as directed.

Product Deterioration: Examine reconstituted reagents at the time of use for evidence of contamination, evaporation, or other signs of deterioration.

PROCEDURE

Material Provided: Depending upon which product is ordered, one of the medium supplements listed above is provided.

Materials Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Instructions: Preparation of Thayer-Martin Agar

1. Prepare a double-strength base by suspending 7.2 g of BBL™ GC Agar Base or GC II Agar Base in 100 mL of purified water, using a 500 mL flask. Mix thoroughly, heat with frequent agitation and boil for about 1 min to assure complete solution of ingredients.
2. Suspend 2 g of BBL Hemoglobin Powder in 100 mL purified water to make a 2% solution. (Mix 2 g of Hemoglobin Powder with 2 to 3 mL purified water until a smooth paste is achieved. Gradually add the balance of the water until the solution is homogeneous. If larger volumes are required, use the same method, maintaining the same ratio of Hemoglobin to purified water.) Alternatively, use BBL Hemoglobin Solution 2% warmed to approximately 50 °C.
3. Sterilize the GC Agar Base or GC II Agar Base and BBL Hemoglobin Solution, if prepared from the powder, by autoclaving at 121 °C for 15 min.
4. Cool the sterile solutions to approximately 50 °C.
5. Reconstitute BBL IsoVitaleX™ Enrichment, 2 mL (see insert for directions).
6. To reconstitute BD BBL V-C-N Inhibitor, see “Storage Instructions and Reconstitution.”
7. Aseptically add the 100 mL of Hemoglobin, 2 mL of BBL IsoVitaleX Enrichment, 2 mL of BD BBL V-C-N Inhibitor to the 100 mL of sterile cooled GC Agar Base or GC II Agar Base.
8. Mix gently but thoroughly and distribute into sterile Petri dishes or other sterile containers.

The BD BBL V-C-N Inhibitor, 10 mL, is used similarly, by adding the reconstituted contents of one vial to 500 mL of sterile cooled GC Agar Base or GC II Agar Base (36.0 g of the base in 500 mL purified water to make a double-strength base), 500 mL of sterile 2% BBL Hemoglobin Solution (approximately 50 °C), and 10 mL of BBL IsoVitaleX Enrichment.

Preparation of Modified Thayer-Martin Agar and Transgrow Medium with Trimethoprim

1. For preparation of Transgrow Medium, add extra agar to GC Agar Base or GC II Agar Base to bring the final agar concentration to 2% and sterilize by autoclaving at 121 °C for 15 min.
2. After mixing the sterile, cooled (approximately 50 °C) GC Agar Base or GC II Agar Base, sterile and cooled BBL Hemoglobin Solution, BBL IsoVitaleX Enrichment and BD BBL V-C-N-T Inhibitor, add 6 mL of sterile 25% dextrose solution to each liter of medium.*
3. Aseptically dispense into Petri dishes (Modified Thayer-Martin Agar) or (for preparation of Transgrow Medium with Trimethoprim) into horizontally positioned sterile 1 oz. prescription bottles (8 to 10 mL volumes) and loosely apply rubber-lined screw caps.
4. After the medium has cooled and solidified in the bottles, introduce a carbon dioxide atmosphere into the bottles by placing a group of bottles in a vacuum chamber, exhaust the air with a vacuum pump (15 lb negative pressure) and refill the chamber with a filtered mixture of 10% CO₂ – 90% air until the chamber is at 5 lb positive pressure; repeat this step three times. After the third gassing, leave the chamber at positive pressure for 2 to 3 h before returning it to atmospheric pressure.
5. Upon opening the chamber, tighten screw caps to make an airtight seal.

* In the BBL formulation for Modified Thayer-Martin Agar the extra dextrose has been eliminated for improved growth of *N. gonorrhoeae*.

User Quality Control

Examine lyophilized and reconstituted supplement for signs of deterioration as noted under “Product Deterioration.” Check performance of the finished medium by inoculation with pure cultures of stable control organisms, producing known, desired reactions. The following cultures are recommended:

Thayer-Martin Agar

Neisseria gonorrhoeae, ATCC® 43069

Growth

Staphylococcus epidermidis, ATCC 12228

Inhibition (partial)

Modified Thayer-Martin Agar and Transgrow Medium with Trimethoprim

In addition to the cultures listed above:

Proteus mirabilis, ATCC 43071

Inhibition (partial)

RESULTS

Typical colonial morphology on these media is as follows:

Neisseria gonorrhoeae

Small, grayish-white to colorless, mucoid

Neisseria meningitidis

Medium to large, blue-gray, mucoid

LIMITATIONS OF THE PROCEDURE

Media into which BD BBL V-C-N Inhibitor and BD BBL V-C-N-T Inhibitor are incorporated are selective media that may inhibit other pathogenic bacteria, e.g., *Haemophilus*. Also, the existence of strains of *N. gonorrhoeae* inhibited by vancomycin and trimethoprim lactate have been reported.^{11,12} It is suggested that Chocolate Agar and blood agar plates be used in conjunction with selective media to aid in the isolation of other pathogens that may be present in the specimen.

Nystatin has been reported to be relatively ineffective in inhibiting the growth of yeast contaminants at the recommended concentration.¹³ Gonococci may be inhibited in the presence of *Candida albicans*.^{14,15}

While "saprophytic" *Neisseria* are generally suppressed by selective media, the occasional recovery of *N. lactamica* on Thayer-Martin Selective Agar has been reported.¹⁶

Some strains of *Capnocytophaga* species may grow on these selective media when inoculated with oropharyngeal specimens.¹⁷

Since there is no such entity as a perfect medium, some strains of microorganisms are encountered that grow poorly on a particular medium; the nature of the specimens or samples themselves and the physiologic state of the organisms on isolation can influence recovery of desired species, as well as modify the effects of inhibitory characteristics of a selective medium for undesired species.

It should be noted that situations are relatively rare when a single medium will suffice both for detection and for enumeration of specific microorganisms. Each selective medium represents a compromise in that selective agents, while being inhibitory to many undesired species, may also be somewhat inhibitory to specific strains of the desired species for which the medium was designed.

Appropriate texts should be consulted for further information.^{17,18}

AVAILABILITY

Cat. No.	Description
212227	BD BBL™ V-C-N Inhibitor, Ten vials, lyophilized; each reconstitutes to 2 mL
212228	BD BBL™ V-C-N Inhibitor, Ten vials, lyophilized; each reconstitutes to 10 mL
212408	BD BBL™ V-C-N-T Inhibitor, Ten vials, lyophilized; each reconstitutes to 10 mL

REFERENCES

1. Thayer, J.D., and J.E. Martin, Jr. 1966. Improved medium selective for cultivation of *N. gonorrhoeae* and *N. meningitidis*. Public Health Rep. 81:559–562.
2. Mitchell, M.S., D.L. Rhoden, and B.B. Marcus. 1966. Immunofluorescence techniques for demonstrating bacterial pathogens associated with cerebrospinal meningitis. III. Identification of meningococci from the nasopharynx of asymptomatic carriers. Am. J. Epidemiol. 83:74–85.
3. Martin, J.E., T.E. Billings, J.F. Hackney, and J.D. Thayer. 1967. Primary isolation of *N. gonorrhoeae* with a new commercial medium. Public Health Rep. 82:361–363.
4. Martin, J.E., Jr., and A. Lester. 1971. Transgrow, a medium for transport and growth of *Neisseria gonorrhoeae* and *Neisseria meningitidis*. HSMHA Health Rep. 86:30–33.
5. Martin, J.E., J.H. Armstrong, and P.B. Smith. 1974. New system for cultivation of *Neisseria gonorrhoeae*. Appl. Microb. 27:802–805.
6. Center for Disease Control. January 2, 1975. Memorandum: recommendation to use the same medium, Modified Thayer-Martin (MTM), in both plates and bottles for the GC culture screening program. U.S. Public Health Service, Atlanta.
7. Seth, A. 1970. Use of trimethoprim to prevent overgrowth by *Proteus* in the cultivation of *N. gonorrhoeae*. Br. J. Vener. Dis. 46:201–202.
8. Center for Disease Control. 1975. Criteria and techniques for the diagnosis of gonorrhoeae. U.S. Public Health Service, Atlanta.
9. Evangelista, A.T., and H.R. Beilstein. 1993. Cumitech 4A, Laboratory diagnosis of gonorrhoeae, Coord. ed., C. Abramson. American Society for Microbiology, Washington, D.C.
10. Morello, J.A., W.M. Janda, and G.V. Doern. 1991. *Neisseria* and *Branhamella*, p. 258–276. In A. Balows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), Manual of clinical microbiology. 5th ed. American Society for Microbiology, Washington, D.C.
11. Cross, R.C., M.B. Hoger, R. Neibaur, B. Pasternack, and F.J. Brady. 1971. VCN-inhibited strains of *Neisseria gonorrhoeae*. HSMHA Health Rep. 86:990–992.
12. Phillips, I., D. Humphrey, A. Middleton, and C.S. Nicol. 1972. Diagnosis of gonorrhoeae by culture on a selective medium containing vancomycin, colistin, nystatin, and trimethoprim (VCNT). A comparison with Gram-staining and immunofluorescence. Br. J. Vener. Dis. 48:287–292.
13. Faur, Y.C., M.H. Weisburd, and M.E. Wilson. 1973. A new medium for the isolation of pathogenic *Neisseria* (NYC Medium). II. Effect of amphotericin B and trimethoprim lactate on selectivity. Health Lab. Sci. 10:55–60.
14. Hipp, S.S., W.D. Lawton, N.C. Chen, and H.A. Gaafar. 1974. Inhibition of *Neisseria gonorrhoeae* by a factor produced by *Candida albicans*. Appl. Microbiol. 27:192–196.
15. Hipp, S.S., W.D. Lawton, M. Savage, and H.A. Gaafar. 1975. Selective interaction of *Neisseria gonorrhoeae* and *Candida albicans* and its possible role in clinical specimens. J. Clin. Microbiol. 1:476–477.
16. Edberg, S.C. 1974. The growth of *Neisseria lactamica* on media selective for pathogenic *Neisseriaceae*. Am. J. Clin. Pathol. 62:445.
17. Baron, E.J., and S.M. Finegold. 1990. Bailey and Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.
18. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Change History

Revision/Date	Section	Change Summary
(03) 2019-02	All p14 Warnings and Precautions	Updated BD branding. Updated trademark line. Updated the "Warnings and Precautions" section with new GHS requirements for P codes.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Becton Dickinson France S.A.S.
11 rue Aristide Bergès
38800 Le Point de Claix, France

ATCC is a trademark of the American Type Culture Collection.

© 2019 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.