BD CDC Anaerobe Agar + 5% Sheep Blood

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LITTERATUR

1. Dowell, V.R., Jr., G.L. Lombard, F.S. Thompson, and A.Y. Armfield. 1977. Media for isolation, characterization, and identification of obligately anaerobic bacteria. CDC laboratory manual. Center for Disease Control, Atlanta.
2. Dowell, V.R., Jr., and T.M. Hawkins. 1987. Laboratory methods in anaerobic bacteriology. CDC laboratory manual. HHS Publication No. (CDC) 87-8272. Centers for Disease Control, Atlanta.
3. Rodloff, A.C., P.C. Appelbaum, and R.J. Zabransky. 1991. Cumitech 5A, Practical anaerobic bacteriology. Coordinating ed., A.C. Rodloff. American Society for Microbiology, Washington, D.C.
4. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Starr, S.E., G.E. Killgore, and V.R. Dowell, Jr. 1971. Comparison of Schaedler agar and Trypticase soy- yeast extract agar for the cultivation of anaerobic bacteria. Appl. Microbiol. 22:655-658.
6. Gibbons, R.J., and J.B. MacDonald. 1960. Hemin and vitamin K compounds as required factors for the cultivation of certain strains of *Bacteroides melaninogenicus*. J. Bacteriol. 80:164-170.

7. Wilkins, T.D., S.L. Chalgren, F. Jimenez-Ulate, C.R. Drake, Jr., and J.L. Johnson. 1976. Inhibition of *Bacteroides fragilis* on blood agar plates and reversal of inhibition by added hemin. *J. Clin. Microbiol.* 3:359-363.
8. Holdeman, L.V., E.P. Cato, and W.E.C. Moore (ed.). 1977. *Anaerobe laboratory manual*, 4th ed. Virginia Polytechnic Institute and State University, Blacksburg.
9. Engelkirk, P.G., J. Duben-Engelkirk, and V.R. Dowell, Jr. 1992. *Principles and practice of clinical anaerobic bacteriology*. Star Publishing Co., Belmont, Calif.
10. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. *Wadsworth anaerobic bacteriology manual*, 5th ed. Star Publishing Co., Belmont, Calif.
11. Miller, J.M., and H.T. Holmes. 1995. Specimen Collection, transport, and storage. *In: Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
12. Murray, P.R., and D.M. Citron. 1991. General processing of specimens for anaerobic bacteria. *In: A. Balows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), Manual of clinical microbiology*, 5th ed. American Society for Microbiology, Washington, D.C.
13. Forbes, B.A. and P.A. Granato. 1995. Processing specimens for bacteria. *In: Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
14. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
15. Thomson Jr., R.B., and J.M. Miller. 2003. Specimen Collection, transport, and processing: bacteriology. *In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
16. Chapin, K.C., and T.-L. Lauderdale. 2003. Reagents, stains, and media: bacteriology. *In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.

EMBALLERING/BESTILLING

BD CDC Anaerobe Agar with 5% Sheep Blood

Katalognummer 256506

Plademedier klar til brug, cpu 20

YDERLIGERE OPLYSNINGER

For further information please contact your local BD representative.

Becton Dickinson GmbH

Tullastrasse 8 – 12

D-69126 Heidelberg/Germany

Phone: +49-62 21-30 50 Fax: +49-62 21-30 52 16

Reception_Germany@europe.bd.com

<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection

BD, BD logo, Trypticase, Stacker, Port-A-Cul and GasPak are trademarks of Becton, Dickinson and Company

© 2013 BD