BD CultureSwab Max V

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
Referanslar

1. Stuart R.D. The diagnosis and control of gonorrhea by bacteriological cultures. Glasgow M. J. 27:131–142, May 1946.
2. Moffett M. Young J.L. and Stuart R.D. Centralized gonococcus culture for dispersed clinics. Brit. M. J.. 2: 421–424. Aug. 28, 1948.
3. Stuart R.D., Toshach S.R. and Patsula T.M. The problem of transport of specimens for culture of gonococci. Cand. J. Health, 45: 73–83, February, 1954.
4. Stuart R. D. Transport Medium for specimens in Public Health Bacteriology. Public Health Reports 74: No. 5, 431–438, May, 1959.
5. Amies C.R. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296–300.
6. Isenberg H. D., Schoenkencht F.D. and Von Graeventiz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
7. Balows A., Hausler, Jr. W. J., Herrmann K.L, Isenberg H. D., Shadomy H.J. Manual of Clinical Microbiology. Fifth Edition. American Society for Microbiology, Washington DC, 1991.
8. Isenberg H. D. (Editor in Chief). Clinical Microbiology Procedures Handbook. American Society for Microbiology, Washington DC, 1992.
9. Van Horn K., Toth C. and Wegienek J. Viability of aerobic microorganisms in four swab systems. Poster Session 249/C Abstract C-436.98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
10. Warshauer D. M. and J.T. Paloucek. Comparison of Bacterial Recovery Rates of Four Commercial Swab Transport Systems. Poster Session 249/C. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
11. Perry J. L. Inhibitory properties of a commercially available swab transport device. 97th General Meeting of American Society for Microbiology, Miami Beach, Florida, May 1997.
12. Perry J. L. and Ballou D. R. Inhibitory properties of a swab transport device. J. Clin. Micro, Dec 1997, 3367–3368.
13. Krepel J. The effect of different sponge (pledget) material on the survival of Group A Streptococcus (GAS) during swab transport. Abstracts from Conjoint Meeting on Infectious Diseases, Canadian Association for Clinical Microbiology and Infectious Diseases, October 1997.
14. Perry J. L. Evaluation of fastidious organism survival in swab transport systems. Abstracts of 39th Interscience Conference on Antimicrobial Agents and Chemotherapy (ICAAC), American Society for Microbiology, San Francisco, CA. September 1999.
15. Campos J. M., Ruthman L, and Tshimanga M. Survival of fastidious bacteria on specimen collection swabs stored at room temperature. Abstract C-150.100th General Meeting of the American Society for Microbiology, Los Angeles, CA. May 2000.
16. Tóth C, and Van Horn K. Comparison of four swab transport systems for the recovery of aerobic microorganisms. Abstract C-153.100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.

© 2017 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.

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R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

CE 0123

Identification number of notified body / Identification de l'organisme notifié / Identifizierung der benannten Stelle / Identificación del organismo notificado / Identificazione dell'organismo notificato / Identificação do organismo notificado / Identifikační číslo notifikovaného subjektu / Identifikasjonsnummer ved påvist organisme / Identifieringsnummer av certifierad myndighet

STERILE R

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismeetod: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация әдісі – сәуле түсіру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestråling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metóda sterilizácie: ožiarovanie / Metoda sterilizacije: ozračavanje / Steriliseringssmetod: strålning / Sterilizasyon yöntemi: iradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射

Do not reuse / Не используйте повторно / Ne používejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Ne používať opakovaně / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullanna / Не використовувати повторно / 请勿重复使用

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperatuurlimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Peel / Обелете / Otevřete zde / Åbn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ыстіңі қабатын алып таста / 벗기다 / Pléștiția / Atīmēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipște / Отклеить / Odtrhnite / Oljuštiti / Dra isär / Ayırma / Відклеїти / 撕下

Manufactured by:
Copan Italia SpA
Via F. Perotti, 10
25125 Brescia Italy

Distributed by:
Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2017 BD, BD, and the BD Logo are trademarks of Becton, Dickinson and Company.