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REFERENCE

1. Hall, C.B, Weinberg, G.A., Iwane, M.K., et al. 2009. The Burden of Respiratory Syncytial Virus Infection in Young Children. N Engl J Med 2009; 360:588–598.
2. Nair, H., Nokes, D.J., Gessner, B.D., et al., 2010. Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. Lancet 375 (9725), 1545–1555.
3. Falsey, A.R. and E. E. Walsh. 2000. Respiratory syncytial virus infection in adults. Clin. Microbiol. Rev. 13: 371–384.
4. Murata Y. and A. R. Falsey. 2006. RSV Infection in Elderly Adults, In: Patricia Cane, Editor(s), Perspectives in Medical Virology, Elsevier, Volume 14, Pages 163–182.
5. Crowcroft, N.S., F. Cutts and M.C. Zambon. 1999. Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. Commun Dis Public Health. 2: 234–241.
6. Thompson WW, Shay DK, Weintraub E, et al. 2003. Mortality associated with influenza and respiratory syncytial virus in the United States. JAMA. 2003. 289(2):179–186.
7. Henrickson KJ and C.B.Hall, Diagnostic assays for respiratory syncytial virus disease. 2007. Pediatr Infect Dis J. Nov; 26(11 Suppl):S36–40.
8. Therese Popow-Kraupp, T and J. H. Aberle. 2011. Diagnosis of Respiratory Syncytial Virus Infection. Open Microbiol J.; 5:128–134. Published online 2011 December 30.
9. Barenfanger, J., C. Drake, N. Leon, T. Mueller and T. Troutt. 2000. Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. J. Clin. Microbiol. 38:2824–2828.
10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, PA.
11. Garner, J. S. 1996. Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53–80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021–0045.

Tehnička služba: obratite se lokalnom predstavniku kompanije BD ili bd.com.

Istorija promena

Revizija	Datum	Change Summary
(11)	2020-04	Ime BD Synapsys™ Microbiology Informatics Solution promenjeno je u BD Synapsys™ Informatics Solution. Moguće su i dodatne opcije za dokumentaciju rezultata implementacijom rešenja BD Synapsys Microbiology Informatics Solution, kao i dodavanjem modula BD Veritor InfoScan i BD Veritor Plus Connect (koji zamenjuju BD Veritor System Reader (Čitač sistema BD Veritor) i modul BD InfoSync). U odeljcima C i D, revidirani su koraci 10C i 10D koji se odnose na analizator Veritor Plus Analyzer sa modulom BD Veritor InfoScan. Dodate su informacije o pristupu kako bi se pribavio dokument bd.com/e-labeling.
(12)	2020-05	Ispravljene su greške u tipografiji, bez uticaja na sadržaj.

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Атқарушы / 제조업체 / Gamintojas / Ražotājs / Tilviker / Producent / Producător / Производител / Výrobca / Proizvođač / Tilverkare / Üretici / Виробник / 生产厂商

Use by / Използвайте до / Spottebujte do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейин пайдаланура / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použít do / Upotrebiti do / Använd före / Son kullanna tarhi / Використати до / 使用截止日期
 YYYY-MM-DD / YYYY-MM (MM = end of month)
 TTTT-MM-DD / TTTT-MM (MM = края на месеца)
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce)
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
 EEEE-HH-NN / EEEE-HH (HH = hónap utolsó napja)
 AAAA-MM-GG / AAAA-MM (MM = fine mese)
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)
 YYYY-MM-DD/YYYY-MM (MM = 월말)
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
 AAAA-MM-DD / AAAA-MM (MM = sluten av måneden)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
 TTTT-MM-DD / TTTT-MM (MM = конец месяца)
 RRRR-MM-DD / RRRR-MM (MM = koniec mesiacu)
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
 YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номер / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号

Authorized Representative in the European Community / Оторизиран представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitadud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizuirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közösségben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유럽 공동체의 위임 대표 / Įgaliotasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autorisert representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Автура Топлулуğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 欧洲共同体授权代表

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisinskhes in-vitro-Diagnostikum / In vitro diagnostic medical device / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiiniaparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnostikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жүзisetін медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietaisai / Medicinas ierices, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Уполномоченный прибор для диагностики ин витро / Medicinska romčka na diagnostiku in vitro / Medicinski uređaji za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Dijagnostik Tibbi Cihaz / Медицинский пристрій для діагностики in vitro / 体外诊断医疗设备

Temperature limitation / Температури ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hörmersékleti határ / Limiti di temperatura / Температураны шекте / 온도 제한 / Laikumo temperatūra / Temperatūras ierobežojumi / Temperaturlimit / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ochránenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(로트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-code (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (lot) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партии / 批号 (並批)

Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečně množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenido suficiente per <n> test / <n> тесттері үшін жеткілікті / <n> 테스트가 충분한 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor <n> testen / Innholder tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достаточное для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli miktarda içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostdutujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

Do not reuse / Не используйте повторно / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolosiți / Не использовать повторно / Nepoužívajte opakovane / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullannayin / Не використовувати повторно / 请勿重复使用

Serial number / Серийн номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / N° de serie / Seerianumber / Numéro de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Sérijas numurs / Serie nummer / Numer serijny / Número de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号

For IVD Performance evaluation only / Само за оценка качеството на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Μόνο για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réserve à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietaisų veikimo charakteristikoms tikrinti / Vienīgi IVD darbības novērtēšanai / Utlslutend voor doeltreffendheidsonderzoek / Kun for evaluering av IVD-ytelse / Tylko do oceny wydajności IVD / Uso esclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики ин витро / Určené iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirmesi için / Тільки для оцінювання якості діагностики in vitro / 仅限 IVD 性能评估

For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najnižia dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураны төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurgrense / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controlo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / Контроль / 对照

	Positive control / Положительный контроль / Pozitívny kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitívna kontrola / Pozitiv kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Negativa kontrolé / Negativá kontrola / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂
	Negative control / Отрицательный контроль / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negatiivne kontroll / Contrôle négatif / Negativa kontrola / Negativ kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativa kontrolé / Negativá kontrola / Negatieve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂
	Method of sterilization: ethylene oxide / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Sterilisierungsmethode: Ethenoxid / Steriliseringmethode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Sterilisierungsmethode: Ethenoxid / Steriliseringmethode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Стерилизация эдиси – этилен тотыгы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringmethode: etylenoksid / Metoda sterylizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metoda sterilizácie: etýlenoxid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷
	Method of sterilization: irradiation / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisierungsmethode: bestrahlung / Steriliseringmethode: bestraling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisierungsmethode: bestrahlung / Steriliseringmethode: bestraling / Sterilisationsmethode: Bestrahlung / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация эдиси – сәулә тәсірі / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmethode: bestraling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metoda de sterilizare: iradiere / Метод стерилизации: облучение / Metoda sterilizácie: ožiarenie / Metoda sterilizacije: ozračevanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: radyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐照
	Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Bioloogilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологический тәуекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Biologisk risiko / Zagrozenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologiskie riziko / Biološki rizici / Biologisk risk / Biyolojik risk / Биологічна небезпека / 生物学风险
	Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatus! Lugea kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristite prateću dokumentaciju / Fysgeleem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dămesio, ziurăteci pridedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsigtig, se vedlagt dokumentation / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vysťah, pozri spríevodné dokumenty / Paźnia! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увага: див. супутню документацию / 小心, 请参阅附带文档。
	Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrense / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Górná granica temperatura / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranice teploty / Gornja granica temperatura / Øvre temperaturgrens / Sicaklık üst sınırı / Максимальна температура / 温度上限
	Keep dry / Пазете сухо / Skladujte v suchém prostředí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausai / Uzglabāt sausu / Droog houden / Hoides tørt / Przechowywać w suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Берегти від вологості / 请保持干燥
	Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жынау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de coheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamanı / Час забору / 采集时间
	Peel / Obeneri / Otevěte zde / Abn / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Устірті қабатын алып таста / 벗기기 / Plešti čia / Atliimt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipse / Отклеить / Odtrhnite / Oljuštiti / Dra isär / Ayırma / Відклеїти / 撕下
	Perforation / Перфорация / Perforace / Perforering / Διάτρηση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Tecik tetsy / 절취선 / Perforacija / Perforácia / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔
	Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιό. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Ерер пакет бузылган болса, пайдаланба / पैकिज کا 손상된 경우 사용 금지 / Jei pakutė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包裝破損, 請勿使用
	Keep away from heat / Пазете от топлина / Nevystavujte přilísnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қаңаңыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo do calor / A se feri de căldură / Не нагревать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тепла / 请远离热源
	Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lögata / Découper / Reži / Vágja ki / Tagliare / Keciş / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciać / Cortar / Descupați / Отрезать / Odstihnite / Iseći / Klipp / Kesme / Pozpazati / 剪下
	Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жынаған тізбекүні / 수집 날짜 / Paėmimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de coheita / Data colectării / Дата сбора / Datum odboru / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期
	µL/test / µL/тест / µL/Test / µL/εξέταση / µL/prueba / µL/teszt / µL/테스트 / мкл/тест / µL/tyrimas / µL/pārbaude / µL/teste / мкл/анализ / µL/检测
	Keep away from light / Пазете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қаңаңыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Feriti de lumină / Хранити в темноті / Uchovávať mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Išiktan uzak tutun / Берегти від дії світла / 请远离光线
	Hydrogen gas generated / Образован е водород газ / Možnost úniku plynného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadrži hydrogen vodik / Hidrogén gázt fejleszt / Produzione di gas idrogeno / Газетекс сутери пайда болды / 수소 가스 생성됨 / Išskiria vandenilio dujas / Rodas idėnradis / Waterstofgas gegenereerd / Hydrogengas generert / Powoduje powstawanie wodoru / Produção de gás de hidrogênio / Generare gaz de hydrogen / Выделение водорода / Vyrobené puzitím vodíka / Oslobada se vodonik / Genererad vätegas / Aciqa çikan hidrojen gazı / Реакция з виділенням водню / 会产生氢气
	Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikációs broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттің идентификациялық нөмірі / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificatienummer van de patiënt / Patientens ID-nummer / Numer ID paciența / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienťa / ID broj pacijenta / Patientnummer / Hasta kimlik numarasi / Идентификатор пациента / 患者标识号
	Fragile, Handle with Care / Чупливо, Роботете с необходимото внимание. / Khehé. Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Ευθραusto. Χειρτοτείτε το με προσοχή. / Frágil. Manipular con cuidado. / Őrn, kásitsege ettevaatlikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сығыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkités atargiai. / Trausis; rtkolies uzmanigi / Breekbaar, voorzichtig behandelen. / Őmtálig, händter försiktigt. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Khehé, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşıyın. / Тендітна, звертатися з обережністю / 易碎, 小心轻放

 bd.com/e-labeling

Europe, CH, GB, NO:		+800 135 79 135
International:		+31 20 794 7071
AR	+800 135 79 135	LT 8800 30728
AU	+800 135 79 135	MT +31 20 796 5693
BR	0800 591 1055	NZ +800 135 79 135
CA	+1 855 805 8539	RO 0800 895 084
CO	+800 135 79 135	RU +800 135 79 135
EE	0800 0100567	SG 800 101 3366
GR	00800 161 22015 7799	SK 0800 606 287
HR	0800 804 804	TR 00800 142 064 866
IL	+800 135 79 135	US +1 855 236 0910
IS	800 8996	UY +800 135 79 135
LI	+31 20 796 5692	VN 122 80297

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

BD, the BD logo, Synapsys, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.