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REFERENSER

- James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Kandry, David W. Warnock, eds. Manual of Clinical Microbiology. 10th edition. Washington, DC: ASM; 2011.
- Lynne S. Garcia. (Editor in Chief). Clinical Microbiology Procedures Handbook. 3rd edition. American Society for Microbiology, Washington, DC, 2007.
- Isenberg HD. (Editor in Chief). Essential Procedures for Clinical Microbiology. American Society for Microbiology, Washington DC, 1998.
- C.H. Collins, Patricia M. Lyne, J.M. Grange. Microbiological Methods. Sixth Edition. Butterworths, 1989.
- Miller JM. A Guide to Specimen Management in Clinical Microbiology. Second Edition. American Society for Microbiology. Washington, DC. 1999.

6. Quality Control of Microbiology Transport Systems; Approved Standard-Second Edition. M40-A2 Vol 34 No. 9. June 2014. Clinical and Laboratory Standards Institute.
 7. Kevin Gough, MD, Michelle Alfa, PhD, and Godfrey Harding, MD. Evaluation of routine enteric pathogens in hospitalized patients: A Canadian perspective. The Canadian Journal of Infectious Diseases Pulsus Group.
 8. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr., et al. A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013.
 9. Public Health England. UK Standard for Microbiology Investigations: Investigation of Faecal Specimens for Enteric Pathogens. Standard Unit, Microbiology Services, PHE. Bacteriology B30, Issue no: 8.1, 24 April 2014.
 10. M. Wasfy, B. Oyoyo, A. Elgindy, A. Churilla. Comparison of Preservation Media for Storage of Stool Samples. Journal of Clinical Microbiology, Aug. 1995.
 11. W. Nugent, M.J. Russell, S. Beck and G.P. Leonardi. Nassau University Medical Center (NY). Evaluation of the physical and diagnostic nature of swabs. 23rd Annual Clinical Virology Symposium, Clearwater Beach (FL), 2007.
 12. Fleming D. Biological Safety: Principles and Practices. January 2000. ASM, Washington DC.
 13. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA.
<http://www.cdc.gov/od/ohs/symp5/jrtext.htm>.
 14. Richardson JH. Biosafety in Microbiological and Biomedical Laboratories. December 1994. Diane Publishing Company.
 15. Hansen DJ. Healthcare, Laboratories and Biosafety. Vol 2., 1992. CRC Press.
 16. Peter H. Gilligan, J. Michael Janda, Mohamed A. Karmali, J. Michael Miller, Laboratory Diagnosis of Bacterial Diarrhea CUMITECH 12°, ASM April 1992.
 17. Collection and Transport of Fecal specimens, Chapters 1–3 V.2. p65 Centers for Disease Control and Prevention.
 18. JOB AIDS: How To Collect a Fecal Specimen And Transfer To Transport Medium Centers for Disease Control and Prevention.
 19. Bailey and Scott, Diagnostic Microbiology Chapter 11 o1st edition of Bailey and Scott's Diagnostic Microbiology (1962).
 20. Patzi-Vargas S1, Zaidi M, Bernal-Reynaga R, León-Cen M, Michel A, Estrada-Garcia TJ. Persistent bloody diarrhoea without fever associated with diffusely adherent Escherichia coli in a young child. Med Microbiol. 2013 Dec;62 (Pt 12):1907–10.
 21. Marchou B. Traveller's diarrhea: epidemiology, clinical practice guideline for the prevention and treatment Presse Med. 2013 Jan;42(1):76–81.
 22. Diagnostic Procedure for Stool Collection, Centers for Disease Control and Prevention Richard K. Root, Clinical Infectious Disease 1999 p. 18.
- Teknisk service: Kontakta närmaste BD-representant eller besök bd.com

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante /
 Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirkir / Producent / Producător / Производители / Výrobca / Proizvođač / Tillverkare / Üretici / Виробник /
 生产厂商

Identification number of notified body/ Numero di identificazione dell'organismo notificato/ Identification de l'organisme notifié / Identifizierung der benannten Stelle/ Identificación del organismo notificado / Identificação do organismo notificado / Identifikasjon for godkjenningsorganet / Identifikasjonsnummer ved påvist organisme/Identifikationsnummer av certifierad myndighet

Use by / Используйте до / Spotbefujete do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 /
 Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izljetiot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Использовать до / Použít do / Upotrebiti do / Använd före / Son kullanna tarihi / Використати до/line / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce)
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)
 AAAA-MM-GG / AAAA-MM (MM = fine mese)
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)
 YYYY-MM-DD/YYYY-MM(MM = 월말)
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
 AAAA-MM-DD / AAAA-MM (MM = slutten av månaden)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
 YYYY-AA-GG / YYYY-AA (AA = айн сону)
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
 YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue /
 Kataloški broj / Katalogusszám / Numero di catalogo / Каталог номери / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus number / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numerasi / Номер за каталогом / 目录号

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsinaparatuur / Dispositif médical de diagnostic in vitro / Medicinska romagala za In Vitro Dijagnostiku / In vitro diagnostikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жұризетін медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietaisais / Medicinas ierices, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispozitiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro / Medicinska romôčka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Diyagnostik Tibbi Cihaz / Медицинский пристрій для діагностики in vitro / 体外诊断医疗设备

Temperature limitation / Температури ограничєния / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Temperaturuuri shekteu / 온도 제한 / Laikumi temperatūra / Temperaturas ierobežojumi / Temperaturilimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ograničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-code (lot) / Batch-Code (Charge) / Κυβικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-code (parti) / Kod partii (serie) / Código do lote / Cod de serie (lot) / Код партии (порт) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)

Contains sufficient for <n> tests / Содержанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> testhez elegendő / Contenido suficiente per <n> test / <n> testleri için yeterli / <n> 테스트가 충분하다 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor <n> testen / Innholder tilstrækkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достатчно для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli malzeme içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeada kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i brugsanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultati instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se brugsanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

Do not reuse / He използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 제사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolosiți / Не использовать повторно / Nepoužívať opakovane / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullannayin / Не використовувати повторно / 请勿重复使用

STERILE R

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisieremismetode: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация адісі – сәуле тусып / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestrålning / Metoda sterilizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metodă sterilizácie: ožiarenie / Metoda sterilizacije: ozračivanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: ırradyasyon / Метод стерилизації: опромінення / 灭菌方法: 辐射

Peel / Обелете / Otevřete zde / Åbn / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстіңі қабатын алып таста / 벗기기 / Plești ăia / Altimët / Schillen / Trekk av / Oderwać / Destacar / Se dezlipște / Отклеить / Odrhniť / Oljušiti / Dra isär / Ayırma / Відклеїти / 撕下

Do not use if package damaged / He използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhl beschädigter Packung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Eper paket бузылган болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayin / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用

R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Manufactured by:
Copan Italia SpA
Via F. Perotti, 10
25125 Brescia Italy

Distributed by:
Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

FecalSwab™ is a trademark of Copan Group.

ATCC is a trademark of the American Type Culture Collection.

© 2019 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.