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LITTERATUR

1. Ruoff, K.L., R.A. Whiley, and D. Beighton. *Streptococcus*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
2. Schuchat, A.. 1998. Epidemiology of group B streptococcal disease in the United States. Clin. Microbiol. Rev. 11: 497-513.
3. Juncosa, T., et al. 1998. Infección neonatal por *Streptococcus agalactiae*. Estudio multicéntrico en el área de Barcelona. Enferm. Infecc. Microbiol. Clin. 16: 312-315.
4. Centers for Disease Control and Prevention. 1996. Prevention of perinatal group B streptococcal disease: a public health perspective. Morbid. Mortal. Weekly Rep. 45: 1-24.
5. Islam, A.K.M. 1977. Rapid recognition of group B streptococci. Lancet I: 256-257.
6. de la Rosa, M., et al. 1983. Granada medium for detection and identification of group B streptococci. J. Clin. Microbiol. 18: 779-785.
7. de la Rosa, M., et al. 1992. New Granada medium for detection and identification of group B streptococci. J. Clin. Microbiol. 30: 1019-1021.
8. Reardon, E.P., et al. 1984. Evaluation of a rapid method for the detection of vaginal group B streptococci in women in labor. Am. J. Obstet. Gynecol. 184: 575-578.
9. Overman, S.B., et al. 2002. Evaluation of methods to increase the sensitivity and timeliness of detection of *Streptococcus agalactiae* in pregnant women. J. Clin. Microbiol. 40: 4329-4331.
10. de la Rosa-Fraile, M. 2003. Granada agar sensitivity and detection of group B Streptococcus. Letter to the Editor. J. Clin. Microbiol. 41: 4007.
11. Thomson, R.B., and J.M. Miller. 2003. Specimen collection, transport, and processing: bacteriology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
12. Spellerberg, B., et al. 1999. Identification of genetic determinants for the hemolytic activity of *Streptococcus agalactiae* by ISS1 transposition. J. Bacteriol. 181: 3212-3219.
13. Pritzlaff, C.A., et al. 2001. Genetic basis for the β -haemolytic/cytolytic activity of group B *Streptococcus*. Molec. Microbiol. 39: 236-247.
14. Edwards, M. and J. Baker. 2000. *Streptococcus agalactiae*. In: Mandell, G.I., et al. (eds.) Principles and practice of infectious diseases. New York Medical Publications, New York, p. 2156-2167.
15. Data i arkiv. 2003; 2010. BD Diagnostic Systems Europe. Heidelberg, Tyskland.
16. Rosa-Fraile, M., et al. 2006. Granadaene: proposed structure of the Group B Streptococcus polyenic pigment. Appl. Environm. Microbiol. 72: 6367-6370.

EMBALLERING/BESTILLING

BD Group B Streptococcus Differential Agar (Granada Medium)

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2013 BD

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Becton Dickinson GmbH