BD GeneOhm™ MRSA Assay
Test BD GeneOhm™ MRSA
BD GeneOhm™ MRSA Test
Prueba BD GeneOhm™ MRSA
Saggio BD GeneOhm™ MRSA

REF 441242

IVD

200 Tests / Pruebas / Test

REF 441244

IVD

48 Tests / Pruebas / Test

Table of Contents / Table des Matières / Inhaltstabelle / Índice / Indice

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English

Intended use

BD GeneOhm™ MRSA assay is a qualitative *in vitro* diagnostic test for the direct detection of nasal colonization by methicillin-resistant *Staphylococcus aureus* (MRSA) to aid in the prevention and control of MRSA infections in healthcare settings. The test performed on the SmartCycler® instrument with a nasal swab specimen from patients at risk for colonization, utilizes polymerase chain reaction (PCR) for the amplification of MRSA DNA and fluorogenic target-specific hybridization probes for the detection of the amplified DNA. BD GeneOhm™ MRSA assay is not intended to diagnose MRSA infections nor to guide or monitor treatment for MRSA infections. Concomitant cultures are necessary only to recover organisms for epidemiological typing or for further susceptibility testing.

Summary and explanation of the test

A nasal specimen is collected and transported to the laboratory using the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided). For testing, the swab is placed in sample buffer. The specimen is concentrated and lysed. An aliquot of the lysate is added to PCR reagents which contain the MRSA-specific primers used to amplify the genetic target, if present. The assay also includes an internal control (IC) to detect PCR inhibitory specimens and to confirm the integrity of assay reagents. Amplified targets are detected with hybridization probes labeled with quenched fluorophores (molecular beacons). The amplification, detection and interpretation of the signals are done automatically by the Cepheid SmartCycler® software. The whole procedure takes about 60 to 75 minutes, depending on the number of specimens processed. For the recovery of MRSA for epidemiological typing or for further antibiotic susceptibility testing, appropriate culture media can be inoculated during specimen preparation or up to 24 hours after its preparation.

S. aureus is a major cause of nosocomial infections. Most transmissions occur through the contaminated hands of a person carrying *S. aureus*. While *S. aureus* causes infections with clinical manifestations ranging from pustules to sepsis and death¹, it is commonly found in the nose or on the skin of healthy individuals (asymptomatic carriers). Treatment of *S. aureus* infections has become a real challenge with the emergence of strains resistant to previously effective antimicrobial agents. Methicillin-resistant strains of *S. aureus* are frequently encountered in health-care settings, and represent over 50% of isolates from hospital acquired *S. aureus* in some North American hospitals². Risk factors for infection with MRSA in health-care settings include prolonged hospital stay, proximity to patients infected with MRSA, exposure to multiple and prolonged broad-spectrum antibiotic treatments, and MRSA nasal carriage. Early screening of patients for MRSA nasal carriage to identify those patients that require isolation can be part of an effective infection control program for MRSA. Current techniques used to detect MRSA all require a culture step and the isolation of pure colonies followed by either oxacillin susceptibility testing, detection of the *mecA* gene or detection of the penicillin binding protein (PBP 2a) encoded by the *mecA* gene. This increases the time of resolution of MRSA carrier status to a minimum of 16 hours; with a median time of more than 48 hours. With the rapidity at which *S. aureus* infections can spread, especially in health-care settings where carriers are common, the capability of providing results of MRSA nasal carriage on the day of admission represents a definite advantage for infection control programs.

Principle of the procedure

Following lysis, amplification of the target [sequence near the insertion site of Staphylococcal Cassette Chromosome *mec* (SCC*mec*)] will occur. Amplification of the IC, a DNA fragment of 335-bp including a 277-bp sequence not found in MRSA, will also take place unless there are PCR inhibitory substances.

The amplified DNA targets are detected with molecular beacons, a hairpin-forming single-stranded oligonucleotides labelled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the absence of target, the fluorescence is quenched. In the presence of target, the hairpin structure opens upon beacon/target hybridization, resulting in emission of fluorescence. For the detection of MRSA amplicons, the molecular beacon contains the fluorophore FAM at the 5' end and the non-fluorescent quencher moiety DABCYL at the opposite end of the oligonucleotide. For the detection of the IC amplicons, the molecular beacon contains the fluorophore TET at the 5' end and the quencher DABCYL at the 3' end. Each beacon-target hybrid fluoresces at a wavelength characteristic of the fluorophore used in the particular molecular beacon. The amount of fluorescence at any given cycle, or following cycling, depends on the amount of specific amplicons present at that time. The SmartCycler® simultaneously monitors the fluorescence emitted by each beacon, interprets all data and at the end of the cycling program provides a final result (see Interpretation of Results).

Reagents

BD GeneOhm™ MRSA Assay	48 Tests	200 Tests
Sample Buffer	60 X 1mL	240 X 1mL
Tris-EDTA buffer		
Lysis tube	50 tubes	200 tubes
Glass beads		
Master Mix (8 reactions each)	8 tubes	34 tubes
< 0.0005 % DNA polymerase complex		
< 0.001% Internal Control - non-infectious DNA containing MRSA-primer binding sequences and a unique sequence for probe hybridization		
< 0.002 % primers		
< 0.002 % molecular probes		
< 0.05 % dATP, dCTP, dGTP, dTTP		
Bovine serum albumin		
Carbohydrate		
MgCl ₂		
< 0.001 % non-infectious <i>Staphylococcus epidermidis</i> genomic DNA (ATCC 14990)		

Control DNA	8 tubes	34 tubes
Tris-EDTA buffer		
Carbohydrate		
< 0.001 % non-infectious genomic MRSA DNA (ATCC 43300)		
Diluent	8 X 700 µL	34 X 700 µL
Tris-HCl buffer		
MgCl ₂		
(NH ₄) ₂ SO ₄		
KCl		

Precautions

This test is for *in vitro* diagnostic use only.

- Do not use the kit if the outer carton safety seal is broken.
- Do not use reagents if the protective pouches are open or torn upon arrival.
- Close protective pouches of Master Mix and Control DNA quickly with the zip seal after each use.
- Do not use reagents if desiccant is not present inside Master Mix and Control DNA pouches.
- Do not remove desiccant from Master Mix and Control DNA pouches.
- Reagents are not interchangeable between lots.
- Never pool reagents from different tubes even if they are from the same lot.
- Do not use the reagents after their expiration date.
- Do not interchange caps among reagents as contamination may occur and compromise test results.
- Avoid microbial and deoxyribonuclease (DNase) contamination of reagents when removing aliquots from tubes. The use of sterile DNase-free disposable filter-blocked or positive displacement pipettor tips is recommended.
- To avoid contamination of the environment with MRSA amplicons, do not open the reaction tubes post-amplification.
- Use a new tip for each specimen or reagent.
- Performing the assay outside of the time ranges recommended can produce invalid results. Assays not done within specified time ranges should be repeated.
- Additional controls may be tested according to guidelines or requirements of local, state, provincial and/or federal regulations or accrediting organizations.
- In cases where open-tube PCR tests are also conducted by the laboratory, separated and segregated working areas should be used for specimen preparation and amplification/detection activities. Supplies and equipment should be dedicated to each area and should not be moved from one area to another. Gloves must always be worn and must be changed before going from one area to another. Gloves must be changed before manipulating lyophilized reagents.
- Do not freeze the collection device. Store at room temperature. Unopened, the collection device is stable until the expiration date indicated.
- Always handle specimens as if they are infectious and in accordance with safe laboratory procedures such as those described in *Biosafety in Microbiological and Biomedical Laboratories*³ and in the CLSI Document M29⁴.
- Wear protective clothing and disposable gloves while handling kit reagents. Wash hands thoroughly after performing the test.
- Do not pipet by mouth.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of unused reagents and waste in accordance with country, federal, provincial, state and local regulations.

Materials provided

- Sample Buffer
- Lysis tube
- Master Mix
- Control DNA
- Diluent
- SmartCycler® reaction tubes, 25 µL
- Specimen identification labels

Storage, handling and stability

Collected specimens

Specimens should be kept between 2°C and 30°C during transport. Protect against freezing or exposure to excessive heat.

Specimens with evidence of potentially interfering substances (e.g. blood, excessive nasal secretions, etc.) should be tested within 24-36 hours. Otherwise, specimens can be stored up to 5 days at 2-8°C before testing. Specimens that can be tested within 36 hours can be kept at room temperature (15-30°C).

Reagents

Note: Storage conditions must follow the specifications written on each pouch. Tubes outside of their protective bag and unused within specified time limit should be discarded.

Kit Component		Master Mix and Control DNA (white and red strip labels)	Lysis tube (yellow cap)	Sample Buffer, and Diluent (blue cap, and black strip label respectively)
Sealed pouch	Temperature	2-25 °C	2-25 °C	2-25 °C
	Stability	Expiration date	Expiration date	Expiration date
Opened pouch	Temperature	2-8 °C ¹	2-25 °C	2-25 °C ²
	Stability	1 month ³	Expiration date	2 months ³

¹ Once the original seal on the pouch is broken, carefully close the pouch with the zip seal after each use and store at appropriate temperature.

² Although these reagents can be stored at room temperature, they should be kept with their accompanying reagents of the same lot at 2-8°C.

³ Provided the bag is properly closed with the zip seal after each use.

Kit Component outside of their protective pouch			Master Mix and Control DNA (white and red strip labels)
Not reconstituted tubes		Temperature	15-25 °C
		Stability	2 hours
Reconstituted tubes ¹	Original container	Temperature	2-8 °C
		Stability	3 hours
	SmartCycler® tube	Temperature	2-8 °C
		Stability	1 hour

¹ Discard unused tubes after expiration of the delay.

Materials required but not provided

- **BBL™ CultureSwab™ Liquid Stuart** (Becton Dickinson catalog no. 220099), **Copan Transystem™ Liquid Stuart** (Copan Italia International catalog no. 141C.USE), **Copan Venturi Transystem™ Liquid Stuart** (Copan Diagnostics Inc. catalog no. 141C.US), **HealthLink TransPorter™ single Liquid Stuart** (HealthLink catalog no. 4432)
- **BBL™ CHROMagar™ Staph aureus** catalog no. 214982, **Mannitol Salt Agar (MSA)** catalog no. 221773 or 221271 or equivalent media (optional),
- **Vortex** Genie 2 (Fisher) with 1.5 mL microtube holder or equivalent; for processing multiple samples, adapter with multiple holding sites can be used
- **Micropipettors** (accurate range between 1-10 µL, 10-100 µL and 100-1000 µL)
- **Sterile DNase-free** filter-blocked or positive displacement **micropipettor tips**
- **Sterile fine tip transfer pipets** (e.g. VWR catalog nos. 14670-329, 14670-332) or extended-length sterile DNase-free filter-blocked tips of 1250 µL (e.g. MATRIX catalog no. 8255)
- **Scissors**
- **Gauzes**
- Disposable **gloves**, powderless
- **Microcentrifuge** for high (must reach 14 000 x g) and low speed centrifugation
- **Dry heating block** for 1.5 mL tubes or water bath
- Ice or **cooling block** for 1.5 mL tubes
- Cap removal tool (e.g. MATRIX catalog no. 4469) (optional)
- Stopwatch or **timer**
- **SmartCycler® starter system** with **Dx Software** (processing block, user manual, accessory kit, and desktop computer) (Cepheid, Sunnyvale, CA, USA)

Instructions for use

Specimen collection

In order to obtain an adequate specimen, the procedure for specimen collection must be followed closely.

Using the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided), nasal specimens are collected according to the following procedure:

1. Moisten the swab with two drops (about 50 µL) of sterile physiological serum (saline) or use it dry.
2. Carefully insert the swab into the patient's nostril (**the swab tip must be inserted up to 2.5 cm (1 inch) from the edge of the nares**).
3. Roll the swab 5 times.
4. Insert the same swab into the second nostril and repeat sampling as 2 and 3.
5. Place the swab in its container.
6. Label the container.
7. Ship the swab to the laboratory according to hospital standard operating procedures.
8. Refer to the section entitled Storage, handling and stability – Collected specimens for storage and handling.

Specimen preparation

Note: One Sample Buffer tube (**blue cap**) and one lysis tube (**yellow cap**) are required for **each specimen** to be tested. An additional Sample Buffer tube (**blue cap**) is also required for each group of 20 specimens to be tested. Remove the required number of tubes from their protective pouch, **remove the excess air, and close the pouch quickly with the zip seal.**

For culturing clinical specimens (swab) prior to performing the BD GeneOhm™ MRSA assay, refer to section Culturing of clinical specimens - Streak-plate method for details.

1. **Place the collection device (swab) in a Sample Buffer tube (blue cap).**
Identify the Sample Buffer tube on the cap and/or the tube label.
2. **Break the swab stem and close the tube tightly.**
Hold the swab by the stem near the rim of the tube (use gauze to minimize risks of contamination). Lift the swab a few millimeters from the bottom and bend the stem against the edge of the tube to break it. Alternative method: use clean scissors to cut the stem. Make sure the cap will close tightly.
3. **Vortex at high speed for one minute.**
For processing multiple samples, an adapter with multiple holding sites can be used.
4. **Transfer the entire cell suspension to a lysis tube (yellow cap).**
Use a sterile fine tip transfer pipet or a P-1000 micropipettor with an extended-length tip.
For culturing clinical specimens with the remaining swab, refer to section Culturing of clinical specimens - Enrichment broth for details.
5. **Centrifuge at high speed (between 14 000 x g and 21 000 x g) for 5 minutes at room temperature.**
6. **Remove the supernatant and discard it.**
Use a sterile fine tip transfer pipet; take care not to touch the pellet. Use a new transfer pipet for each specimen.
7. **Add 50 µL of sample buffer to the lysis tube; close tightly.**
Use a new pipettor tip for each specimen. Up to 20 specimens can be prepared with one Sample Buffer tube. For more than 20 specimens, use as much Sample Buffer tubes as required (**keep the unused buffer for later steps of the BD GeneOhm™ MRSA assay procedure**).
8. **Vortex for 5 minutes at high speed.**
For processing multiple samples, an adapter with multiple holding sites can be used.
9. **Centrifuge the lysis tube briefly (quick spin).**
At low speed for 2 to 5 seconds; to bring the content at the bottom of the tube.
10. **Heat at 95 ± 2 °C for two (2) minutes.**
Use a dry heating block for 1.5 mL tubes or a water bath.
11. **Keep the lysis tube on ice or on a cooling block.**
Lysates are stable for up to 4 hours at 2-8 °C. If the lysates have not been used at the end of this time frame, store them at -20 ± 5 °C for later use.

BD GeneOhm™ MRSA assay procedure

Note: One reconstituted Master Mix tube (**white label**) will yield enough reagents to run 8 reactions. Allow one SmartCycler® tube per specimen to be tested and 2 additional SmartCycler® tubes for the positive and the negative control. **One (1) positive and one (1) negative control must be included in each BD GeneOhm™ MRSA assay run.** One control DNA (**red strip label**) is required per assay run. One diluent tube (**black strip label**) is required for the reconstitution of up to 3 Master Mix tubes. Remove the required number of tubes from their protective pouch, **remove the excess air, and close the pouch quickly with the zip seal.**

Prepare only enough SmartCycler® tubes to fill available I-CORE® modules on the SmartCycler® instrument.

1. **Place the required number of Master Mix tube(s) on ice or on a cooling block for 1.5 mL tubes.**
2. **Add 225 µL of diluent (black strip label) to each Master Mix tube.**
Insert the micropipettor tip through the septum of the cap of the Master Mix tube. Do not insert the tip too deeply into the cap. Deliver the diluent. Discard the unused diluent afterward.
3. **Vortex the tube(s) for 5-10 seconds.**
4. **Place the required number of SmartCycler® tubes on the SmartCycler® cooling block.**
Allow one SmartCycler® tube per specimen and two more SmartCycler® tubes for the controls. **Avoid touching the optical detection windows at the bottom edges of the tube and the lower diamond-shaped area.**

The following steps MUST be performed within one hour time frame:

5. **Add 25 µL of reconstituted Master Mix to the SmartCycler® tubes.**
Remove the septum cap before pipetting the reagent. Deliver the liquid into the reservoir (upper part) of the SmartCycler® tubes. Identify the SmartCycler® tubes on the cap. Specimen identification labels can be used (provided with the kit). Unused Master Mix **immediately stored at 4-8 °C for a maximum of 3 hours away from light** can be used for immediate retesting of unresolved results if necessary. Discard the unused Master Mix after this period.
6. **Add 2.8 µL of each lysed specimen to a different SmartCycler® tube previously filled; close the tubes.**
Take care not to aspirate beads when pipetting into the lysis tube. After addition of the specimen, pipette up and down 2-3 times in the reservoir to ensure transfer of the complete volume. Use a new micropipettor tip for each specimen.
7. **Place a control DNA tube (red strip label) on ice or on a cooling block for 1.5 mL tubes.**
8. **Add 225 µL of Sample Buffer (blue cap) to the control DNA tube.**
Use the Sample Buffer tube from step 7 of the Specimen Preparation protocol. Insert the micropipettor tip through the septum of the cap of the control DNA tube. Do not insert the tip too deeply into the cap. Deliver the sample buffer.
9. **Vortex the tube for 5-10 seconds.**
10. **Add 2.8 µL of the reconstituted control DNA to the next to last SmartCycler® tube (Positive Control); close the tube.**
After addition of the DNA, pipette up and down 2-3 times in the reservoir to ensure transfer of the complete volume. Identify as the positive control. Unused control DNA **immediately stored at 4-8 °C for a maximum of 3 hours away from light** can be used for the preparation of a new positive control in case of immediate retesting of unresolved results. Discard the unused control DNA after this period.
11. **Add 2.8 µL of Sample Buffer (blue cap) to the last SmartCycler® tube (Negative Control); close the tube.**
Use the Sample Buffer tube from step 7 of the Specimen Preparation protocol. This will monitor for PCR contamination that might occur during the manipulation of the specimens. Identify as the negative control. Discard the unused sample buffer afterward.
12. **Centrifuge all reaction tubes for 5-10 seconds.**
Use the specially adapted microcentrifuge provided with the SmartCycler® instrument.
13. **Keep the tubes at 2-8 °C on the SmartCycler® cooling block before loading on the instrument.**
The remaining lysates should be frozen at -20 ± 5 °C for later use, if necessary.
14. **Create a run with the BD GeneOhm™ MRSA assay protocol.**
Refer to the SmartCycler® Dx Software Operator Manual if needed. You should enter the identification parameters for the specimens before starting the run.
15. **Insert each reaction tube into an I-CORE® module of the SmartCycler® and close the I-CORE® lid.**
Place the positive and negative controls at their appropriate position (see the section entitled "Quality control"). Press down all the tubes firmly into place.
16. **Start the run.**

Quality control

Positive and negative controls

Quality control procedures are designed to monitor assay performance. The positive control is intended to monitor for substantial reagent failure. The negative control is used to detect reagent or environmental contamination (or carry-over) by either MRSA DNA or MRSA amplicons. Positive and negative controls are assay controls (run controls). An invalid control invalidates the run. Finally, an internal control incorporated into each reaction mixture is intended to monitor PCR inhibition in each specimen and reagent integrity.

One positive control and one negative control must be included in each assay run on the SmartCycler®. The software automatically assigns the position of the controls on the instrument (refer to the SmartCycler® Dx Software Operator Manual).

Specimen processing controls

Control strains may be tested according to guidelines or requirements of local, state and/or federal regulations or accreditation organizations. A reference MRSA strain (e.g. American Type Culture Collection, ATCC 43300) or a well characterized MRSA clinical isolate may be used as a specimen processing control while a culture of methicillin sensitive *Staphylococcus aureus* (e.g. ATCC 29213) or of any other non-*Staphylococcus aureus* (e.g. *Staphylococcus epidermidis* ATCC 14990) may be used as an external negative control.

Resuspend isolated colonies from an 18- to 24-h 5% sheep blood agar plate (eg. BBL™ Trypticase Soy Agar (TSA II) with 5% Sheep Blood, BD catalog no. 221239 or 221261) in saline to a turbidity of 0.5 McFarland ($\sim 1.5 \times 10^8$ CFU/mL). Dilute with saline to obtain a suspension of $\sim 10^6$ CFU/mL. Dip the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided) into the bacterial suspension, press out the excess fluid, place the swab in its container (to allow contact with the transport medium), and let stand at room temperature for 5 minutes. Process and test as a clinical specimen (refer to the sections entitled "Specimen preparation" and "BD GeneOhm™ MRSA" assay procedure"), including controls. All specimens and controls should yield valid results (no invalid positive or negative control and no failed internal control).

For general QC guidance, the user may wish to refer to CLSI MM3⁵ and C24⁶.

Culturing of clinical specimens

To perform antimicrobial susceptibility testing or epidemiological typing, clinical specimens can be cultured at two different steps of the BD GeneOhm™ MRSA assay procedure. It can be done from the collection device (swab) before conducting the analysis with BD GeneOhm™ MRSA using the streak-plate method or it can be done with the swab following the sample preparation procedure using an enrichment step.

Streak-plate method

This culture method can be performed prior to the specimen preparation procedure for BD GeneOhm™ MRSA.

1. Remove the collection device (swab) from its container.
2. Inoculate a mannitol salt agar (MSA) plate by streaking onto the first quadrant of the plate.
3. Return the swab in its container or break it in a Sample Buffer tube (blue cap) of BD GeneOhm™ MRSA and continue according to instructions in section "Specimen preparation".
4. With a loop needle, streak the inoculum into the remaining quadrants of the MSA plate.
5. Incubate the MSA plate for 24-48 hours at 35 °C.
6. Identify and confirm *S. aureus* colonies and test for methicillin resistance according to standard methods.

Enrichment broth

This culture method can be performed with swabs left behind in the Sample Buffer tube following the transfer of the cell suspension into the lysis tube. Swabs can be stored at 2-8°C in closed Sample Buffer tube for up to 24 hours before culturing; recovery may not be successful with clinical specimens with low amounts of MRSA.

1. Add 1.0 mL of enrichment broth to Sample Buffer tubes containing the swab.
TSB (tryptic soy broth) supplemented with 6.5% NaCl is suggested.
2. Incubate for 24-48 hours at 35 °C.
3. Subculture on an appropriate solid media for 24-48 hours at 35 °C.
5% sheep blood agar plate is suggested.
4. Identify and confirm *S. aureus* colonies and test for methicillin resistance according to standard methods.

Interpretation of results

The decisional algorithm for BD GeneOhm™ MRSA assay is embedded in the SmartCycler® software. The interpretation of assay results is done according to the following criteria:

Assay result reported	IC result reported	Interpretation of result ¹
NEG	PASS	No MRSA DNA detected, MRSA nasal colonization unlikely
POS	NA	MRSA DNA detected, MRSA nasal colonization
Unresolved	FAIL	Unresolved—inhibitory specimen or reagent failure
ND	ND	Not determined due to I-CORE® Module failure (with Warning or Error Codes ²)

IC - Internal Control; NA – not applicable; ND – not determined

¹BD GeneOhm™ MRSA results may be used to guide isolation and level of precautions in accordance with institutional programs and practices.

² Refer to the SmartCycler® Dx Software Operator Manual for interpretation of warning and error codes.

An invalid positive or negative control invalidates the assay run. In such cases, assay results obtained for that run are invalid and must not be reported. Invalid assay run or instrument error codes or warnings are flagged on-screen and on reports. Before reporting MRSA results, always verify that the assay run is valid.

Refer to the SmartCycler® Dx Software Operator Manual for printing of results.

Invalid Assay Run

Using frozen lysate(s), prepare new reaction tubes for all clinical specimens within that assay run along with new control tubes.

Unresolved specimen

Repeat testing with the corresponding specimen frozen lysate. The effect of the freeze-thaw cycle has been shown to reduce PCR inhibitory substances.

Specimen not determined due to I-CORE® module failure

Repeat testing with the corresponding specimen frozen lysate. For the interpretation of warning or error code messages, refer to the SmartCycler® Dx Software Operator Manual.

Limitations of the procedure

- Performance of this test has been established with the SmartCycler® instrument, with nasal specimens obtained from hospitalized patients or from patients at admission and with specimens collected using the Copan Venturi Transystem® with Liquid Stuart Medium. Therefore, this product can only be used with the SmartCycler®; moreover, the use of specimen collection and transport system other than those listed in the Material required but not provided section is not recommended. Other clinical sources have not been assessed and performance characteristics of this test are unknown on other specimen types or patient populations.
- Negative test results may also occur from improper specimen collection, handling or storage, presence of inhibitor, technical error, sample mix-up or because the number of organisms in the specimen is below the analytical sensitivity of the test. Careful compliance to the instructions given in this insert and in the SmartCycler® Dx Software Operator Manual is necessary to avoid erroneous results. Use of this test should be limited to personnel trained on the procedure and on the use of the SmartCycler®.
- Although there is no need for reagent preparation and the main technical operation is pipetting, good laboratory technique is essential to the proper performance of this assay. Due to the high analytical sensitivity of this test, extreme care should be taken to preserve the purity of all reagents, especially in cases where multiple aliquots are taken from a tube.
- Screening determines the colonization status at a given time and could vary depending on patient treatment (e.g. decolonization regime), patient status (e.g. not actively shedding MRSA) or exposure to high risk environments (e.g. contact with MRSA carrier, prolonged hospitalization). Monitoring colonization status should be done according to hospital policies.
- Results from BD GeneOhm™ MRSA assay should be used as an adjunct to nosocomial infection control efforts to identify patients needing enhanced precautions. The test is not intended to identify patients with staphylococcal infection. Results should not be used to guide or monitor treatment for MRSA infections.
- A BD GeneOhm™ MRSA positive result does not necessarily indicate treatment eradication failure since DNA may persist. A negative result following a previous positive test result may indicate treatment eradication success or may occur due to intermittent shedding.
- BD GeneOhm™ MRSA assay results may sometimes be unresolved or invalidated due to an invalid control, and require retesting that can result in a delay for obtaining the results.
- A positive test result does not necessarily indicate the presence of viable organisms. It is however presumptive for the presence of MRSA as BD GeneOhm™ MRSA detects simultaneously the *SCCmec* cassette (carrying the *mecA* gene) and a *S. aureus* specific sequence located within the *orfX* gene. The BD GeneOhm™ MRSA does not detect the *mecA* gene directly nor the penicillin binding protein (PBP 2a) encoded by this gene.
- BD GeneOhm™ MRSA was shown to detect various types of MRSA that have been associated with hospital acquired and community acquired infections; evaluations did not assess phenotypic or genotypic traits of MRSA isolates to identify strains that may be more virulent in the hospital environment.
- Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown MRSA variants resulting in a false negative result with the BD GeneOhm™ MRSA assay.

Interfering substances

Potentially interfering substances include, but are not limited to the following: blood, excessive nasal secretions/mucus, decongestant and substances occasionally used to relieve nasal dryness and/or irritation. The presence of excessive nasal secretions and blood may inhibit PCR and give unresolved results.

In an investigational study conducted on 786 nasal specimens, potentially interfering substances were reported for 43% of collected specimens. In total, there were 35 (4.5%) unresolved specimens. Of these 35 specimens, there was no potentially interfering substance reported for 24 of them, nasal secretions for 10 specimens and a combination of substances for one specimen. Twenty-seven (27) were resolved following a freeze-thaw cycle. For the 8 specimens that could not be resolved, no potentially interfering substance was observed for 4 specimens, nasal secretions was observed on three specimens and a combination of substances on one specimen.

Expected values

Humans are a natural reservoir for *S. aureus*. Colonization may be transient or persistent and can last for years. *S. aureus* nasal carriage rates of 25 to 30% have been reported for the general population and rates of 10 to 40% have been reported for outpatient population or on admission⁷. Methicillin-resistant *S. aureus* now represent over 50% of isolates from hospital acquired *S. aureus* in some North American hospitals².

In the investigational study for BD GeneOhm™ MRSA assay, the overall *S. aureus* nasal carriage rate determined by culture (*S. aureus* isolated in specimens by either culture technique) was 36.1% with a range of 33-40%. Fifty one percent (51%) of *S. aureus* isolates were methicillin resistant by the oxacillin screen agar method [Mueller-Hinton agar plate supplemented with oxacillin (6 µg/mL) and NaCl (4% w/v)] for an overall MRSA nasal carriage rate of 18.6%. With BD GeneOhm™ MRSA assay, the overall MRSA nasal carriage rate was 20.3%.

Performance characteristics

Clinical performance

Performance characteristics of BD GeneOhm™ MRSA assay were determined in a multi-site prospective investigational study: 4 medical centers, two in Canada and two in the U.S with MRSA culture-based screening programs in place. To be enrolled in the study, patients had to provide written consent and be eligible for MRSA screening according to hospital policies and not have received antibiotic therapy towards MRSA. Screening criteria included but were not necessarily limited to: systematic screening of all patients at admission, prior MRSA infection or carriage, transfer from another institution, prolonged hospital stay or history of prolonged hospitalization and contact with an MRSA carrier.

The reference method consisted of an initial analysis with the oxacillin screen agar test after selective growth on mannitol salt agar. Specimens negative for MRSA were subjected to an additional analysis consisting of an enrichment step in trypticase soy broth (TSB) containing 6.5% NaCl followed by the oxacillin screen agar test. An MRSA culture-positive specimen was defined as a specimen positive for MRSA by either culture technique. An MRSA culture-negative specimen was defined as a specimen negative for MRSA by both culture techniques.

For the screening method with selective growth on mannitol salt agar (MSA), plates were directly inoculated with nasal swab specimens followed by incubation at 35°C for 24 to 48 hours. When necessary, presumptive colonies of *S. aureus* were subcultured to 5% sheep blood agar plates and incubated for 24-48 hours at 35°C. Otherwise isolated colonies were tested directly from MSA plates. Presumptive staphylococcal colonies were confirmed with a latex agglutination assay or by a tube coagulase test. Confirmed *S. aureus* colonies were tested on Mueller-Hinton agar plate supplemented with oxacillin (6 µg/mL) and NaCl (4% w/v) incubated at 33 to 35°C (not exceeding 35°C) for a full 24 hours. MRSA screening with an enrichment step in tryptic soy broth (TSB) containing 6.5% NaCl was also performed in cases where a negative result for MRSA was obtained with the MSA screening method. For this purpose, TSB broths were inoculated after direct plating onto MSA, incubated overnight at 35°C, subcultured to 5% sheep blood agar plates for 24-48 hours at 35°C and the plates stored at 2-8°C for testing if needed. Confirmation of presumptive colonies and determination of methicillin resistance were done as described above. In one investigational site, every specimen was tested using both culture screening methods.

In total, 786 nasal specimens collected with the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided) were screened for MRSA with the culture method of reference described above and with BD GeneOhm™ MRSA assay. Compared to the culture method of reference, the BD GeneOhm™ MRSA identified 92.5% of the specimens positive for MRSA by either culture techniques and 96.4% of the specimens negative by both culture techniques (Tables 1 and 2). For the population tested, this results in a negative predictive value of 98.2% and a positive predictive value of 85.4%.

Table 1. Results obtained with BD GeneOhm™ MRSA assay in comparison to the reference method¹.

		BD GeneOhm™ MRSA		Total
		Positive	Negative	
Culture techniques	Positive	135	11	146
	Negative	23	609	632
	Total	158	620	778

¹ Eight (8) specimens that gave initially unresolved results remained unresolved upon retesting with BD GeneOhm™ MRSA assay and are not included in the table. All 8 were culture negative.

Fourteen (14) of the 23 culture-negative specimens but BD GeneOhm™ MRSA positive were found to be MRSA culture-positive upon further investigation, resulting in a total of 149 culture-positive and BD GeneOhm™ MRSA positive specimens out of a total of 160 culture positive specimens. For 2 of the culture-positive specimens but BD GeneOhm™ MRSA negative, none of the isolates that exhibited methicillin resistance on oxacillin agar plates could be shown to carry the *mecA* gene when tested with the *mecA*-specific PCR assay described by Martineau *et al.*⁸.

Table 2. Performance of BD GeneOhm™ MRSA assay obtained by the investigational sites when compared to the reference method

	Sensitivity (95% CI) ¹	Specificity (95% CI) ¹	No. of unresolved specimens ²	Invalid/ total no. of runs
Site 1	86.8% (n=38) (71.9-95.6%)	99.6% (n=261) (97.9-100%)	6	0/26
Site 2	100% (n=30) (88.4-100%)	98.0% (n=102) (93.1-99.8%)	16	0/21
Site 3	89.3% (n=28) (71.8-97.7%)	90.8% (n=119) (84.1-95.3%)	11	1/15
Site 4	94.0% (n=50) (83.5-98.7%)	94.0% (n=150) (88.9-97.2%)	2	2/27
Total	92.5% (n=146) (86.9-96.2%)	96.4% (n=632) (94.6-97.27%)	35	3/89

¹ Binomial 95% confidence intervals.

² All specimens were unresolved due to failed internal controls indicative of inhibition or reagent failure. Twenty-seven (27) of the 35 were resolved upon re-testing.

SCCmec typing of MRSA culture-positive specimens according to Oliveira and de Lencastre⁹ revealed specimens of types I, II and IV. There were no SCCmec type III specimens isolated from the study. However, in a separate study, reference strains and other clinical isolates of all SCCmec types were tested and detected with BD GeneOhm™ MRSA assay.

Performances obtained by the investigational sites for BD GeneOhm™ MRSA, and the individual culture techniques as compared to the culture reference method (both culture techniques) are presented in Tables 3 and 4.

Table 3. Results obtained for BD GeneOhm™ MRSA assay and each culture screening technique with specimens positive for MRSA by the culture reference method.

Site	MRSA prevalence ¹	Sensitivity (95% CI) ²		
		BD GeneOhm™ MRSA	OxaMSA ³	OxaTSB ⁴
Site 1	12.7% (38/300)	86.8% (71.9-95.6%)	81.6% (65.7-92.3%)	ND
Site 2	21.7% (30/138)	100% (88.4-100%)	93.3% (77.9-99.2%)	ND
Site 3	18.9% (28/148)	89.3% (71.8-97.7%)	64.3% (44.1-81.4%)	ND
Site 4	25.0% (50/200)	94.0% (83.5-98.7%)	80.0% (66.3-90.0%)	78.0% (64.0-88.5%)
Total	18.6% (146/786)	92.5% (86.9-96.2%)	80.1% (72.7-86.3%)	ND

¹ Determined from results obtained with the culture reference method (combination of OxaMSA and OxaTSB)² Binomial 95% confidence intervals.³ Oxacillin screen agar test after selective growth on MSA.⁴ Oxacillin screen agar test after enrichment in TSB with 6.5% NaCl. Site 4 tested all specimens with both culture techniques while the remaining sites tested only specimens negative with the OxaMSA culture method.**Table 4.** Results obtained for BD GeneOhm™ MRSA assay and each culture screening technique with specimens negative for MRSA by the culture reference method.

Site	MSSA ¹	Specificity (95% CI) ²		
		BD GeneOhm™ MRSA	OxaMSA ³	OxaTSB ⁴
Site 1	73/262	99.6% (97.9-100%)	100% (98.6-100%)	ND
Site 2	25/108	98.0% (93.2-99.8%)	100% (96.4-100%)	ND
Site 3	24/120	90.8% (84.1-95.3%)	100% (96.9-100%)	ND
Site 4	16/150	94.0% (88.9-97.2%)	100% (97.6-100%)	100% (97.6-100%)
Total	138/640	96.4% (94.6-97.7%)	100% (99.4-100%)	ND

¹ Number of MSSA in total MRSA culture-negative population² Binomial 95% confidence intervals.³ Oxacillin screen agar test after selective growth on MSA.⁴ Oxacillin screen agar test after enrichment in TSB with 6.5% NaCl. Site 4 tested all specimens with both culture techniques while the remaining sites tested only specimens negative with the OxaMSA culture method.**Analytical specificity**

Genomic DNA from 37 ATCC strains representing species phylogenetically related to *S. aureus* and members of the nasal commensal flora, 27 strains (15 reference strains and 12 clinical isolates) of methicillin-sensitive coagulase negative staphylococci and 44 strains (two reference strains and 42 clinical isolates) of methicillin-resistant coagulase negative staphylococci were tested. The specificity was 100%.

As indicated in Table 4, 138 specimens analyzed with culture techniques were deemed by investigational sites to contain methicillin-sensitive *S. aureus*. Seven (7) yielded positive results with BD GeneOhm™ MRSA assay. Upon further investigation with more enhanced culture techniques, 4 were shown to actually contain MRSA.

Analytical sensitivity

The analytical sensitivity (limit of detection or LOD) of BD GeneOhm™ MRSA assay was determined with 7 strains of MRSA. Quantitated culture and purified genomic DNA diluted in the Sample Buffer of BD GeneOhm™ MRSA assay were tested in 5 replicates. The LOD is defined as the smallest concentration at which at least 92.5% of all replicates test positive.

The LOD of BD GeneOhm™ MRSA assay is 15 genome copies per reaction. The LOD in CFU is 5 CFU/reaction. Taking into account the dilution factor due to specimen processing, this translates into approximately 325 CFU/swab.

Reproducibility

A panel of 10 simulated specimens with varying concentrations of MRSA and the two controls (positive and negative) of the BD GeneOhm™ MRSA assay were tested in triplicate on three different days at each of 3 sites (10 specimens plus two controls tested X 3 X 3 days X 3 sites). This was repeated with three lots of reagents.

Table 5. Cumulative data of reproducibility study.

Specimen ID	Lot 1	Lot 2	Lot 3	Total agreement	Total % agreement
Negative	27/27	27/27	27/27	81/81	100%
Negative	27/27	27/27	27/27	81/81	100%
Weak positive ¹	19/27	18/27	16/27	53/81	65%
Strong positive	27/27	27/27	27/27	81/81	100%
Strong positive	27/27	27/27	27/27	81/81	100%
Weak positive ¹	14/27	13/27	18/27	45/81	56%
Positive	27/27	27/27	27/27	81/81	100%
Positive	27/27	27/27	27/27	81/81	100%
Weak positive	27/27	27/27	27/27	81/81	100%
Positive	27/27	27/27	27/27	81/81	100%
Positive control	26/27	27/27	27/27	80/81	99%
Negative control	26/27	27/27	27/27	80/81	99%
Total agreement	301/324	301/324	304/324	906/972	93%
Total % agreement	93%	93%	94%	93%	

¹ Specimens with CFUs below the limit of detection of the assay.

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References / Références / Referenzen / Referencias / Riferimenti

- ¹ Centers for Disease Control and Prevention. Methicillin-resistant *Staphylococcus aureus* skin or soft tissue infection in a state prison. Mississippi, 2000. MMWR 2001; 50:919-922.
- ² National Nosocomial Infections Surveillance (NNIS) System Report, data summary from January 1992 to June 2002, issued August 2002. Am J Infect Control. 2002; 30:458-475.
- ³ Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Richmond JY and McKinney RW (eds) (1993). HHS Publication number (CDC) 93-8395.
- ⁴ Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline – Second edition. Document M29 (Refer to the latest edition).
- ⁵ Clinical and Laboratory Standards Institute. Molecular Diagnostic Methods for Infectious Diseases; Approved Guideline. CLSI document MM3-A2 (Refer to the latest edition).
- ⁶ Clinical and Laboratory Standards Institute. Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline – Second Edition, Document C24 (Refer to the latest edition).
- ⁷ Jernigan, J.A. *et al.* Prevalence of and risk factors for colonization with methicillin-resistant *Staphylococcus aureus* at the time of hospital admission. Infect Control Hosp Epidemiol. 2003; 24:409-414.
- ⁸ Martineau, F. *et al.* Correlation between the resistance genotype determined by multiplex PCR assays and the antibiotic susceptibility patterns of *Staphylococcus aureus* and *Staphylococcus epidermidis*. Antimicrob. Agents Chemother. 2000; **44**(2) : 231-238
- ⁹ Oliveira, D.C. and H. de Lencastre. Multiplex PCR strategy for rapid identification of structural types and variants of the *mec* element in methicillin-resistant *Staphylococcus aureus*. Antimicrob Agents Chemother. 2002; 46:2155-2161.

This kit is sold under license from the Public Health Research Institute of the City of New York, Inc. and may be used under the PHRI patent rights only for human *in vitro* clinical diagnostics.

The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences for providing human *in vitro* diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

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
BENEX Limited, Bay K 1a/d,
Shannon Industrial Estate,
Shannon, County Clare,
Ireland

BD Diagnostics,
2555 boul. Parc-Technologique
Québec, Qc, Canada, G1P 4S5