BD MacConkey Agar with Sorbitol

USO PREVISTO

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BIBLIOGRAFIA

1. Riley, L.W. et al. 1983: Hemorrhagic colitis associated with a rare *Escherichia coli* serotype. *New Engl. J. Med.* 308: 681-685.
2. Kaper, J.B., O'Brien, A.D. (eds.). 1998: *Escherichia coli* O157:H7 and Other Shiga Toxin-Producing *E. coli* strains. American Society for Microbiology, Washington, DC, USA.
3. Dorn, C.R., and E.J. Angrick. 1991. Serotype O157:H7 *Escherichia coli* from bovine and meat sources. *J. Clin. Microbiol.* 29: 1225-1231.
4. Wells, J.G. et al. 1983. Laboratory investigation of hemorrhagic colitis outbreaks associated with a rare *Escherichia coli* serotype. *J. Clin. Microbiol.* 18: 512-520.
5. Willshaw, G.A., et al. 1994: Vero cytotoxin-producing *Escherichia coli* O157 in beefburgers linked to an outbreak of diarrhea, haemorrhagic colitis and haemolytic uraemic syndrome in Britain. *Letters Appl. Microbiol.* 19: 304-307.
6. Ewing, W. H., and P. R. Edwards. 1954. Isolation and preliminary identification of *Escherichia coli* serotypes associated with cases of diarrhea of the newborn. *Public Health Lab.* 12:75-81.
7. March, S.B., and S. Ratman. 1986. Sorbitol-MacConkey medium for detection of *Escherichia coli* O157:H7 associated with hemorrhagic colitis. *J. Clin. Microbiol.* 23: 869-872.
8. Rappaport, F., and E. Henig. 1952. Media for the isolation and differentiation of pathogenic *Escherichia coli* (serotypes 0111 and 055). *J. Clin. Pathology.* 5:361-362.
9. Bopp, C. A., F. W. Brenner, P. I. Fields, J. G. Wells, and N. A. Stockbrine. 2003. *Escherichia, Shigella, and Salmonella*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
10. Adams, S. 1991. Screening for verotoxin-producing *Escherichia coli*. *Clinical Lab Science* 4(1):19-20.
11. March, S. B., and S. Ratnam. 1986. Sorbitol-MacConkey medium for detection of *Escherichia coli* O157:H7 associated with hemorrhagic colitis. *J. Clin. Microbiol.* 23:869-872.
12. Meng, J., Feng, P., and M.P. Doyle. 2001. Pathogenic *Escherichia coli*. In: Downes, F.P., and K. Ito (ed.). *Compendium of methods for the microbiological examination of foods.*, 4th edition. American Public Health Association, Washington. D.C.
13. Hitchins, A. D., P. Feng, W. D. Watkins, S. R. Rippey, and L. A. Chandler. 1995. *Escherichia coli* and the coliform bacteria. p. 4.01-4.29. In *Bacteriological analytical manual*, 8th ed. AOAC International, Gaithersburg, MD.
14. Mortlock, S. 1994. Recovery of *Escherichia coli* O157:H7 from mixed suspensions: evaluation and comparison of pre-coated immunomagnetic beads and direct plating. *Brit. J. Biomed. Sci.* 51: 207-214.

15. Ogden, I.D., Hepburn, N.F., and M. MacRae. 2001. The optimization of isolation media used in immunomagnetic separation methods for the detection of *Escherichia coli* O157 from foods. J. Appl. Microbiol. 91: 373-379.
16. Kist, M., et al. 2000. Infektionen des Darmes. In: Mauch, H., Lüttiken, R., and S. Gatermann (eds.): MiQ - Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik, vol. 9. Urban & Fischer, Munich, Germany.

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38800 Le Pont de Claix/France

Tel: +33-476 68 3636 Fax: +33-476 68 3292 <http://www.bd.com>

BD, BD logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company

Dynabeads and Dynal are trademarks of Dynal

CHROMagar is a trademark of Dr. A. Rambach

ATCC is a trademark of the American Type Culture Collection

© 2003 Becton, Dickinson and Company