MICROTROL DISCS

Becton Dickinson France S.A.
 BD Diagnostic Systems Eur.
 11 rue Aristide Bergès
 38800 Le Pont de Claix
 France
 Tel +33 04 76 68 36 36
 Fax +33 04 76 68 32 92

ISS. NO. 8/10/2003

2157/8/02 FRONT

1

MICROTROL – ENGLISH

Microtrol discs are 1st generation derivatives from the original freeze dried vial of the ATCC, NCTC, NCYC or CIP culture providing a quick and reliable source of control micro-organisms of known reaction. They are freeze dried, water-soluble discs that are offered in pack sizes of 10 or 25 per vial depending on the strain. The viability of the organisms is stabilised by the incorporation of activated charcoal in the disc and silica gel within the vials.

Microtrol discs are intended for use in microbiology laboratories for the control of test methods. Being 1st generation derivatives traceable to vials of recognised national type culture strains, Microtrol discs are acceptable in accredited laboratories for the production of working stock cultures (see limitations).

Removing the disc from the vial - A disc is easily removed from the vial by using sterile forceps or a sterile 10µl inoculating loop.

For use on solid media - Place a disc on appropriate solid medium. Allow the disc to soften for 10–15 minutes. The plate may be placed in an incubator to accelerate the process. Spread the softened disc around the plate and incubate under optimum conditions for the strain used.

For use in liquid media - Place the disc in 1-10mls of the appropriate broth. Avoiding the production of aerosols, shake the broth gently to dissolve the disc and incubate under optimum conditions to achieve a stock culture.

Alternatively, for rapid use, dissolve the disc, allow the culture to recover by incubating at 35–37°C for 1 hour and then use immediately. By experimentation you will be able to find the best dilution for your own methods.

Limitations - Repeated sub-culture can cause the characteristics of a strain to change. We recommend that a fresh disc be used for each set of tests performed.

In accredited laboratories, the discs may only be used for the production of working stocks. Further subculture will invalidate their use.

Precautions & storage

1. Microtrol discs contain viable pathogenic organisms and must only be used in appropriately equipped laboratories by microbiologists or persons under the supervision of microbiologists qualified by training or experience to work with micro-organisms.

2. Store at the temperature indicated on the label. Loss of viability may occur if the discs are left above this temperature range for longer than necessary. **Product to be used within 9 months of opening.**

3. Allow the vial to reach room temperature before opening and replace the cap on each vial immediately after a disc has been removed. Moisture in air entering the vial will cause a reduction in the number of viable organisms, eventually leading to complete loss of viability.

4. Do not use discs that are past their expiration date or if silica gel colour changes. These may show a loss of viability and changes in antibiotic and chemical reactions.

Breakages, spillage & disposal

- All exposed discs, contaminated packaging and broken glassware should be placed in a suitable container and either incinerated or autoclaved at 121°C for 30 minutes.
- Unbroken vials may be removed with forceps, washed in a suitable bactericidal solution, rinsed, dried and retained for further use.
- All contaminated surfaces must be disinfected with a suitable bactericidal solution.

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