
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	8408-16

NOTES:

- BD Catalog Number: 220093, 220097, 220099, 220105, 220109, 220115, 220116, 220129, 220130, 220132, 220133, 220134, 220135, 220146, 220147, 220148, 220149
- Blank (Sheet) Size: Length: 23 5/8" Width: 3 15/32"
- Number of Pages: 2 Number of Sheets: 1
- Page Size: Length: N/A Width: N/A Final Folded Size: N/A
- Ink Colors: No. of Colors: 1 PMS#: 2755 Blue
- Printed two sides: ☒ Yes ☐ No *Peel-off insert / Label. Insert perforated 3/8" from bottom with 3/16" notch for easy peel-off. Top insert size 4 3/16" x 3 15/32" includes 3/8" bottom portion of glued area. Insert is folded under with an horizontal fold and followed by two parallel folds to form the size of 3 1/4" x 3 15/32". Insert to be glued on top 5/16" and bottom 3/8" to roll label, size 3 15/53" (web) x 4 3/8".
- Style (see illustrations below): N/A

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- Material specification controled by Copan.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 9:49 am, Aug 24, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Nichole Graham at 1:43 pm, Aug 24, 2016		
Checked By	THIRD EYE BY By Mary Schmidt at 1:40 pm, Aug 25, 2016		
Part Number:	L000021	Category and Description Peel-Off Insert / Label, CultureSwab	Sheet: 1 of 3 Scale: N/A A

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкциите. / Pokyny vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жерінікі BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL™ CultureSwab™ devices are sterile ready-to-use systems intended for the collection, transport and preservation of clinical specimens for bacteriological examination.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This can be accomplished using the BD BBL CultureSwab collection and transport device. Each BD BBL CultureSwab unit is comprised of a sterile peel pouch containing a rayon- or polyester-tipped swab applicator used to collect the sample and a tube containing transport medium into which the swab applicator is placed after sampling.

The BD BBL CultureSwab transport media (Amies Liquid Medium, Liquid Stuart Medium and Cary-Blair Transport Medium) are non-nutritious, buffered with phosphate and provide a reduced environment due to their formulation with sodium thioglycollate.¹ Organisms in the sample material are protected from drying by moisture in the transport medium. The medium is designed to maintain the viability of organisms during transit to the laboratory. BD BBL CultureSwab pouches are made of a plastic film which retards the penetration of atmospheric air into the product.

BD BBL CultureSwab media are contained in a constricted (venturi) tube. Nitrogen gas is flushed into the transport tube during the media filling and capping process. During final packaging of the swab and tube, air is removed from the pouch by vacuum and nitrogen gas is flushed inside.

REAGENTS

Approximate Formula Per Liter Distilled Water

Amies Liquid Medium		Stuart Liquid Medium		Cary-Blair Agar Medium	
Calcium Chloride	0.10 g	Calcium Chloride	0.10 g	Bacteriological Agar	5.60 g
Disodium Phosphate	1.15 g	Mercaptoacetic Acid	1.0 mL	Calcium Chloride	0.09 g
Magnesium Chloride	0.10 g	Sodium		Disodium Hydrogen Phosphate	1.00 g
Monopotassium Phosphate	0.20 g	Glycerophosphate	10.00 g	Sodium Chloride	5.10 g
Potassium Chloride	0.20 g			Sodium Thioglycollate	1.15 g
Sodium Chloride	3.00 g				
Sodium Thioglycollate	1.00 g				

Precautions: For *in vitro* Diagnostic Use.

It must be assumed that all specimens contain infectious microorganisms; therefore, all specimens should be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

Ⓢ BD BBL CultureSwab is for single use only; reuse may cause a risk of infection and/or inaccurate results.

Storage: Store BD BBL CultureSwab devices at 5–25 °C.

Product Deterioration: Contents sterile if unopened or not damaged. Do not use if they show evidence of damage, dehydration or contamination. Do not use if past expiration date.

SPECIMEN COLLECTION AND HANDLING

The BD BBL CultureSwab devices are available with different applicator shafts which facilitate the collection of specimens from various sites of the patient. For specific recommendations about collection of specimens for microbiological analysis and primary isolation techniques, consult appropriate references.²⁻⁴

Once a swab sample is collected, it should be placed in the tube of medium, transported to the laboratory as soon as possible and cultured onto appropriate primary isolation media.

PROCEDURE

Materials Provided: Fifty (50) units of sterile BD BBL CultureSwab devices, which contain media, are contained in each Vi-Pak Pouch.

One hundred (100) units of sterile BD BBL CultureSwab swab and transport tubes are contained in each Vi-Pak Pouch.

Materials Required But Not Provided: Appropriate materials for isolating, differentiating and culturing aerobic and anaerobic bacteria. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations.

Directions For Use:

The directions for use are printed on each BD BBL CultureSwab unit, along with descriptive diagrams, and are summarized as follows:

1. Peel open the BD BBL CultureSwab pouch.
2. Remove cap from transport tube.
3. Remove applicator swab and collect specimen.
During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.
4. Place applicator swab in transport tube.
5. Record patient's name and information on tube label.
6. Send specimen to the laboratory for immediate analysis.

EXPECTED RESULTS

The survival of bacteria in a transport medium depends on many factors. These include the types of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium. BD BBL CultureSwab devices maintain viability of many microorganisms for 24–48 h. For fastidious bacteria, such as *Neisseria gonorrhoeae* and *Streptococcus pneumoniae*, swab specimens should be plated directly onto culture medium or transported immediately to the laboratory and cultured within 24 h.

LIMITATIONS OF THE PROCEDURE

BD BBL CultureSwab Amies Liquid Medium and Stuart Liquid Medium are intended for the collection and transport of bacteriological samples only. Preferred samples for anaerobic investigations are: tissue samples obtained during surgical procedures, biopsies from tissue or bone, fluid, pus, or aspirates collected using a syringe. For detailed information and recommendations for transporting fluid and tissue specimens for anaerobic culture, refer to specific publications.³⁻⁷ Samples containing viruses or chlamydia should be collected and transported using alternative specific transport systems.

PERFORMANCE CHARACTERISTICS

Recovery studies were performed using BD BBL CultureSwab products with a variety of aerobic organisms. Swabs were dosed with inoculum and inserted into the transport tube containing transport media. The tubes were stored at room temperature prior to subculturing onto appropriate culture media. Listed are the organisms tested with the various transport media.

Organism	Amies Liquid form	Transport Media Stuart Liquid form	Cary-Blair Agar
<i>Neisseria gonorrhoeae</i> ATCC® 43069	*	*	
<i>Streptococcus pyogenes</i> ATCC 19615	*	*	
<i>Haemophilus influenzae</i> ATCC 19418	*	*	
<i>Campylobacter jejuni</i> ATCC 33291			*
<i>Shigella flexneri</i> ATCC 12022			*
<i>Yersinia enterocolitica</i> ATCC 9610			*

REFERENCES

1. Amies CR. 1967. A modified formula for the preparation of Stuart's transport medium. Can. J. Public Health. 58: 296-300.
2. Isenberg H.D., F.D. Schoenkencht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
3. Balows, A., W.J. Hausler, Jr, K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
4. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Zavala, M.K., D.M. Citron, E.J.C. Goldstein. 1998. Evaluation of a novel specimen transport system for anaerobic bacteria. Clin. Infect. Dis. 25 (supplement 2): S132-133.
6. Perry, J.L. 1997. Assessment of swab transport systems for aerobic and anaerobic organism recovery. J. Clin. Microbiol. 35: 1269-1271.
7. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.

TECHNICAL INFORMATION: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

ATCC is a trademark of the American Type Culture Collection.
© 2016 BD. BD, the BD Logo, and all other trademarks are property of Becton, Dickinson and Company.

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