
Rev from	Rev to	ECO #
01	02	6954-14

NOTES:

1. BD Catalog Number: 220144, 220145
2. Blank (Sheet) Size: Length: 11" Width: 25.5"
3. Number of Pages: 6 Number of Sheets: 1
4. Page Size: Length: 11" Width: 8.5" Final Folded Size: 5.5" x 8.5"
5. Ink Colors: No. of Colors: 1 PMS#: Standard Black
6. Printed two sides: Yes ☒ No ☐
7. Style (see illustrations below): # 4

8. See specification control no. n/a for material information.
9. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design *[Signature]* Date *8/25/14*
 Proofer *[Signature]* Date *8/28/14*
 Checked By *[Signature]* Date *10-10-14*
 Part Number: L007980

COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.

Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA

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INTENDED USE

The BBL™ CultureSwab™ EZ/BBL™ CultureSwab™ EZ II Single and Dual Swab Specimen Collection and Transport Systems are sterile devices for collecting and transporting aerobic microbiological specimens.

SUMMARY AND EXPLANATION

The BBL CultureSwab EZ/BBL CultureSwab EZ II Systems are self-contained, ready-to-use specimen collection systems that provide a simplified method for the transport and maintenance of viable microorganisms. The devices incorporate one or two polyurethane-tipped swabs on plastic shafts secured to a cap. The swab(s) contained within the cap is inserted into a tube. Studies have shown that BBL CultureSwab EZ/BBL CultureSwab EZ II Systems yield equivalent results as compared to traditional fiber-tipped collection swabs placed in transport media.¹

For isolation and identification of aerobic microorganisms, the specimen must be removed from the transport container and cultured on appropriate media or processed as appropriate for other microbiological procedures.^{2,4}

PRINCIPLES OF THE PROCEDURE

The non-toxic nature and unique characteristics of the polyurethane-tipped swab enable the transport of specimens without using a liquid transport medium. Thus, the polyurethane-tipped swab does not dilute the specimen as occurs with media-containing devices, and overgrowth by organisms which may utilize transport medium components for growth cannot occur.

PRODUCT DESCRIPTION

Each individually wrapped device consists of one or two sterile polyurethane-tipped swabs on plastic shafts attached to a cap within a tube.

Precautions: For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁵⁻⁸ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. When collecting swab samples from patients, care should be taken not to use excessive force or pressure which may result in breakage of the swab shaft.

⊗ BBL CultureSwab EZ/BBL CultureSwab EZ II is for single use only; reuse may cause a risk of infection and/or inaccurate results.

Storage Instructions: Store at room temperature (15 to 30°C). Do not freeze. Do not use after expiration date.

Production Deterioration: Product sterility guaranteed if primary container is intact. Do not use if there is evidence of damage or contamination.

SPECIMEN COLLECTION AND TRANSPORT

The procedure is used directly with clinical specimens or other sources containing microorganisms.

Whenever possible, the laboratory should provide collection devices and instructions for their proper use. Specimens should be obtained before antibiotics or other antimicrobial agents have been administered. If therapy was initiated prior to collection of the specimen, this should be noted on the specimen container or on forms accompanying the specimen. Material should be collected from the site where the suspected organism is most likely to be found, with as little external contamination as possible, and at the clinical stage of disease most likely to yield a positive culture. Specimens should be of sufficient quantity to permit complete examination. To help ensure that pathogens are viable upon receipt, provision must be made for prompt delivery to the laboratory (see "Expected Results" section as follows).

PROCEDURE

Material Provided: One hundred (100) units of sterile BBL CultureSwab EZ/ BBL CultureSwab EZ II devices are contained in each box.

Materials Required But Not Provided: Materials required for isolation, identification and for other microbiological procedures.

DIRECTIONS FOR USE:

1. Peel open the sterile pouch at the point marked "Peel Here," and remove device from package.
2. Remove the swab and take the sample. To minimize potential contamination, the applicator tip should only touch the area where the infection is suspected.
3. After taking the sample, return the swab to tube.
4. Complete patient identification label.
5. Send to laboratory for immediate processing.

EXPECTED RESULTS

The BBL CultureSwab EZ/BBL CultureSwab EZ II Systems will transport and maintain viability of aerobic microorganisms for 24 to 48 h. Transport time for fastidious microorganisms, such as *Neisseria gonorrhoeae* and *Streptococcus pneumoniae* should be 24-h or less. All organisms may exhibit some loss of viability over time, therefore, transport time should be minimized as much as possible.

LIMITATIONS OF THE PROCEDURE

Optimum recovery is achieved by direct specimen plating and/or smear preparation at the time of collection. Chlamydiae, mycoplasma, viruses and anaerobes require special transport systems.

Recovery of microorganisms is influenced by numerous factors; e.g., specimen composition, organism type and concentration, transport temperature and time in transit. Temperature extremes during transport should be avoided. Excessive transport time may result in loss of viability. Transport time for fastidious organisms, such as *N. gonorrhoeae* and *S. pneumoniae*, should not exceed 24 h.

PERFORMANCE CHARACTERISTICS

Recovery studies were performed using BBL CultureSwab EZ /BBL CultureSwab EZ II Systems with a variety of aerobic organisms. Swabs were dosed with a standard inoculum and inserted into the transport tubes. The tubes were stored at room temperature prior to subculturing onto appropriate media. Aerobic organisms evaluated were:

Organisms	ATCC™	Organisms	ATCC™
<i>Candida albicans</i>	60193	<i>Pseudomonas aeruginosa</i>	27853
<i>Enterococcus faecalis</i>	29212	<i>Salmonella choleraesuis</i>	14028
<i>Escherichia coli</i>	25922	<i>Shigella sonnei</i>	9290
<i>Haemophilus influenzae</i>	10211	<i>Staphylococcus aureus</i>	25923
<i>Haemophilus influenzae</i>	49247	<i>Streptococcus agalactiae</i>	12386
<i>Neisseria gonorrhoeae</i>	31426	<i>Streptococcus pyogenes</i>	12379
<i>Neisseria gonorrhoeae</i>	43069	<i>Streptococcus pyogenes</i>	19615
<i>Neisseria gonorrhoeae</i>	49226	<i>Streptococcus pneumoniae</i>	6305
<i>Neisseria meningitidis</i>	13090	<i>Streptococcus pneumoniae</i>	27336
<i>Neisseria meningitidis</i>	13102	<i>Streptococcus pneumoniae</i>	49619
<i>Proteus mirabilis</i>	43071		

All organisms tested remained viable for 24 h when maintained at room temperature.

AVAILABILITY

Cat. No.	Description
220144	BBL™ CultureSwab™ EZ Collection and Transport System, Single Swab, Carton of 100.
220145	BBL™ CultureSwab™ EZ II Collection and Transport System, Double Swab, Carton of 100.

REFERENCES

1. Gosnell, C.M., J.F. Monthony, D.D. Hardy, L.G. Wood, E. Bergogne-Berezin and C. Muller. 1991. New aerobic transport system, abstr. C-79, p. 355. Abstr. 91st Gen. Meet. Am. Soc. Microbiol. 1991.
2. Forbes, B.A., D.F. Sahn, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
3. Isenberg, H.D., F.D. Schoenkecht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Reuben. American Society for Microbiology, Washington, D.C.
4. Miller, J.M., and H.T. Holmes. 1999. Specimen collection, transport and storage, p.33-63. In Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
5. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, Pa.
6. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
7. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton Dickinson France S.A.S.
Becton Dickinson Pty Ltd.