Part Number:	L008319	BALTSO0191 Version 13.0 Template 4 Inserts
Category and Description: Package Insert, BD ProbeTec CT/GC Urethral Dry Transport Kit	Rev from: 02	Rev to: 03
	Job Number:	1220-19

Catalog Number: 440461

Blank (Sheet) Size: Length: 8.5" Width: 27.5"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 8.5" X 5.5"

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. N/A for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Nancy Carlsen at 11:25 am, Jul 17, 2019

Proofing
Approved By: **PROOFING APPROVED BY**
By Paul Meyers at 9:49 am, Jul 18, 2019

Third Eye By: **THIRD EYE BY**
By MZ Addy at 5:37 pm, Jul 23, 2019

 BD ProbeTec™ ET *Chlamydia trachomatis* / *Neisseria gonorrhoeae* (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit

English:	pages	1 – 2	Italiano:	pagine	5 – 6
Dansk:	side	2 – 3	Español:	páginas	6 – 7
Français :	pages	3 – 4	Svenska:	sidan	7 – 8
Deutsch:	seiten	4 – 5			

L008319(03)
2019-07

Europe: Medical Device
93/42/EEC

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Die Packungsbeilage erhalten Sie bei Ihrer örtlichen BD-Vertretung. / Póngase en contacto con su representante local de BD para instrucciones. / Contacter le représentant local de BD pour les instructions. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Rivolgersi al rappresentante BD di zona per istruzioni. / Нұсқаулар үшін жерінікі BD өкілімен хабарласыңыз. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit offers a convenient way to collect and transport specimens from the patient to the laboratory. It is designed for the collection of male urethral specimens. This transport system is for use with the **BD ProbeTec ET *Chlamydia trachomatis*** and **Neisseria gonorrhoeae** amplified DNA assays on the **BD ProbeTec ET System**.

SUMMARY AND EXPLANATION

The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit is a sterile, self-contained, disposable system. The rayon-tipped collection swab consists of a flexible aluminum wire in a white plastic shaft that is scored. The rayon-tipped swab facilitates obtaining the specimen and maintaining the organisms for nucleic acid amplification on the **BD ProbeTec ET System**. The plastic portion of the shaft remains attached to the cap of the transport tube.

REAGENTS

Materials Provided: The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit contains 100 units. Each unit contains one rayon-tipped swab on a combination aluminum/plastic shaft. The plastic portion of the shaft is attached to the cap of the transport tube.

Materials Required But Not Provided: Gloves.

Storage Instructions: Store swabs at 15–30 °C. Do not use past the expiration date. Product sterility guaranteed if primary container is intact.

Precautions: For *in vitro* Diagnostic Use

Warnings:

1. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
2. Optimal performance of the **BD ProbeTec ET** amplified DNA assays requires proper specimen collection, handling and transport.
3. Use only collection and transport systems labeled for use with the **BD ProbeTec ET *Chlamydia trachomatis*** and **Neisseria gonorrhoeae** Amplified DNA Assays.
4. Proper labeling should accompany each specimen to the laboratory.
5.  **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit is for single use only; reuse may cause a risk of infection and/or inaccurate results.

SPECIMEN COLLECTION AND TRANSPORT

Male Urethral Swab Specimen Collection

1. Insert the male urethral specimen collection swab 2–4 cm into the urethra and rotate for 3–5 s.
2. Withdraw the swab and place the cap/swab into the transport tube. Make sure the cap is tightly secured to the tube.
3. Label the tube with patient information and date/time collected.

Swab Storage and Transport

The male urethral specimen collection swab must be stored and transported to the laboratory and/or test site at 2–27 °C within 4–6 days of collection. Storage up to 4 days has been validated with clinical specimens; storage up to 6 days has been demonstrated with seeded specimens. In addition, storage up to 30 days at 2–8 °C has been demonstrated with seeded specimens.

NOTE: If specimens cannot be transported directly to the testing laboratory under ambient temperatures (15–27 °C) and must be shipped, an insulated container with ice should be used with either an overnight or 2-day delivery vendor.

For domestic and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53–80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021–0045.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton Dickinson Pty Ltd.
Becton Dickinson France S.A.S.