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REFERENCES

1. Kent, HL. Am J. Obstet Gynecol; 165 (4 Pt 2):p1168-76;Oct 1991.
2. Vulvovaginitis. 1989. ACOG Technical Bulletin. American College of Obstetricians and Gynecologists
3. Sobel, JD. Med Clin North Am; 74 (6) p1573-602; Nov 1990.
4. Gravett MG, HP Nelson, T DeRouen, et al. 1986. J Am Med Assoc 256:1899-1903.
5. Faro S. 1991. J Reprod Med 34:602-604.
6. Gravett MG, D Hummel, DA Eschenbach, et al. 1986. Obstet Gynecol 67:229-237.
7. McGregor JA, JI French, K Seo. 1991. Am J Obstet Gynecol 165:867-875.
8. McGregor JA, JI French, R Richter, et al. 1990. Am J Obst Gynecol 163:1465-1473
9. Hillier SL, NB Kiviat, C. Critchlow, et al. Abstract presented at IDSOG, Aug. 6-8, 1992. San Diego, pg. 12.
10. Spiegel, CA, R Amse, and KK Holmes. 1983. J Clin Micro 18:170-177.
11. Krieger JN. 1981. Invest Urog 18:411-417.
12. Rein MF. 1990. Trichomonas vaginalis. In: Principles and Practice of Infectious Diseases, 3rd ed. Mandell GL, RG Douglas, Jr, and JE Bennett (eds). New York: Churchill Livingstone Inc., p.2115-2118.
13. Martius J, MA Krohn, SL Hillier, et al. 1988 Obstet Gynecol 71:89-95
14. Nugent, RP, MA Krohn and SL Hillier. 1991. J Clin Microbiol. 29:297-301.
15. Morton K, L Regan, S Sprige, and E Houang. 1990. Eur J Obstet Gynecol Reprod Biol 37:231-236
16. Soper DE, RC Bump, and WG Hurt. 1990. Am J Obstet Gynecol 163:1016-23
17. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
18. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
19. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
20. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
21. Data is on file at Becton, Dickinson and Company.
22. Isenberg, Henry D, and Richard F D'Amato. 1991. Manual of Clinical Microbiology, 5th ed, Ballows, Hausler, Herrman, Isenberg, Shadomy (ed), ASM, Washington, DC, p.8.
23. Bergey's Manual of Systematic Bacteriology, Vol. 2. PHA Sneath, NS Mair, ME Sharpe, JG Holt (eds), Baltimore, Williams and Wilkens. 1989. p. 1427.

BD Diagnostics teknisk service: kontakta närmaste BD-representant eller besök www.bd.com/ds.

Tabell 1

JUSTERADE

Symtomatiska patienter

ORGANISM	Ref.metod	Affirm VPIII			Utstryk med hängande droppe			Gramfärgning		
		Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet
<i>Candida</i> sp.	Odling	82,3 %	98,4 %	95,2 %	71,9 %	92,2 %	88,1 %	45,6 %	99,6 %	90,3 %
		79/96	377/383	456/479	69/96	353/383	422/479	26/57	271/272	297/329
<i>G. vaginalis</i>	Gramfärgning	95,2 %	100 %	95,3 %	—	—	—	—	—	—
		118/124	5/5	123/129	—	—	—	—	—	—
	Odling	98,3 %	100 %	98,4 %	—	—	—	—	—	—
		118/120	9/9	127/129	—	—	—	—	—	—
Alla patienter										
ORGANISM	Ref.metod	Affirm VPIII			Utstryk med hängande droppe			Gramfärgning		
		Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet
<i>Candida</i> sp.	Odling	80,6 %	98,2 %	95,3 %	60 %	94,0 %	88,2 %	44,4 %	99,3 %	92,1 %
		100/124	605/616	705/740	74/124	579/616	653/740	28/63	418/421	446/484
<i>G. vaginalis</i>	Gramfärgning	83,8 %	100 %	89,0 %	—	—	—	—	—	—
		171/204	95/95	266/299	—	—	—	—	—	—
	Odling	89,0 %	99,1 %	92,6 %	—	—	—	—	—	—
		170/191	107/108	277/299	—	—	—	—	—	—
<i>T. vaginalis</i>	Utstryk med hängande droppe	92,8 %	99,9 %	98,9 %	—	—	—	—	—	—
		103/111	740/741	843/852	—	—	—	—	—	—
	Odling 5-7 dagar	89,6 %	99,9 %	98,5 %	82,6 %	99,6 %	97,3 %	—	—	—
		103/115	736/737	839/852	95/115	734/737	829/852	—	—	—

8

EJ JUSTERADE

Symtomatiska patienter

ORGANISM	Ref.metod	Affirm VPIII			Utstryk med hängande droppe			Gramfärgning		
		Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet
<i>Candida</i> sp.	Odling	79,3 %	95,0 %	92,3 %	67,1 %	88,9 %	85,2 %	42,6 %	98,5 %	89,4 %
		65/82	377/397	442/479	55/82	353/397	408/479	23/54	271/275	294/329
<i>G. vaginalis</i>	Gramfärgning	95,1 %	83,3 %	94,6 %	—	—	—	—	—	—
		117/123	5/6	122/129	—	—	—	—	—	—
	Odling	98,1 %	40,9 %	88,4 %	—	—	—	—	—	—
		105/107	9/22	114/129	—	—	—	—	—	—
Alla patienter										
ORGANISM	Ref.metod	Affirm VPIII			Utstryk med hängande droppe			Gramfärgning		
		Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet	Sensitivitet	Specificitet	Noggrannhet
<i>Candida</i> sp.	Odling	78,0 %	95,9 %	93,2 %	54,1 %	91,8 %	86,2 %	41,7 %	98,6 %	91,5 %
		85/109	605/631	690/740	59/109	579/631	638/740	25/60	418/424	443/484
<i>G. vaginalis</i>	Gramfärgning	83,5 %	96,0 %	87,6 %	—	—	—	—	—	—
		167/200	95/99	262/299	—	—	—	—	—	—
	Odling	87,5 %	81,7 %	84,9 %	—	—	—	—	—	—
		147/168	107/131	254/299	—	—	—	—	—	—
<i>T. vaginalis</i>	Utstryk med hängande droppe	91,8 %	98,1 %	97,4 %	—	—	—	—	—	—
		90/98	740/754	830/852	—	—	—	—	—	—
	Odling 5-7 dagar	89,2 %	99,3 %	98,0 %	82,0 %	99,1 %	96,8 %	—	—	—
		99/111	736/741	835/852	91/111	734/741	825/852	—	—	—

*Sensitivitet = TP/(TP+FN); Specificitet = TN/(TN+FP); Noggrannhet (TP+TN)/(TP+FP+FN+TN)

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Атқарушы / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производитель / Výrobca / Proizvodač / Tillverkare / Üretici / Виробник

Use by / Исполняйте до / Spotřebujte do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейин пайдаланура / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použite do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до line

YYYY-MM-DD / YYYY-MM (MM = end of month)
 TTGG-MM-DD / TTGG-MM (MM = края на месеца)
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce)
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)
 AAAA-MM-GG / AAAA-MM (MM = fine mese)
 ЖОЖОК-АА-КК / ЖОЖОК-АА / (АА = айдың соңы)
 MMMM-MM-DD / MMMM-MM (MM = mėnesio pabaiga)
 GGGG-MM-DD / GGGG-MM (MM = mėneša beigās)
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
 AAAA-MM-DD / AAAA-MM (MM = slutten av måneden)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
 TTGG-MM-DD / TTGG-MM (MM = конец месяца)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог нөмөрі / Katalogo numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер на каталог / Katalógové číslo / Kataloški broj / Katalog numarası / Номер за каталог

Authorized Representative in the European Community / Оторизирани представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentante autorisée pour la Communauté européenne / Autorizuirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Községszénben / Représentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / Įgaliotasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Reprezentante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Avrupa Topluluğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС

In vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku *in vitro* / *In vitro* diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / *In vitro* διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico *in vitro* / *In vitro* diagnostisk meditsiiniparatuur / Dispositif médical de diagnostic *in vitro* / Medicinska pomagala za *In vitro* Diagnostiku / *In vitro* diagnosztikai orvosi eszköz / Dispositivo medicale per diagnostica *in vitro* / Жасанды жағдайда жүргізілетін медициналық диагностика аспабы / *In vitro* diagnostikos prietais / Medicīnas ierīces, ko lieto *in vitro* diagnostikā / Medisch hulpmiddel voor *in-vitro* diagnostiek / *In vitro* diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki *in vitro* / Dispositivo médico para diagnóstico *in vitro* / Dispozitiv medical pentru diagnostic *in vitro* / Медицинский прибор для диагностики *in vitro* / Medicīnska romēška na diagnostiku *in vitro* / Medicinski uređaj za *in vitro* dijagnostiku / Medicinteknisk produkt för *in vitro*-diagnostik / In Vitro Diagnotisk Tibbi Cihaz / Медицинский прибор для диагностики *in vitro*

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμός θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектегі / Laikymo temperatūra / Temperatūras ierobežojumi / Temperaturulimiet / Temperaturbegrensnig / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії

Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenuto sufficiente per <n> test / <n> тесттегі үшін жеткілікті / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor <n> testen / Innholder tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достаточно для <n> тестов(а) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli miktarda içerir / Вистачить для аналізів: <n>

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanim Talimatları'na başvurun / Див. інструкції з використання

Do not reuse / Не используйте повторно / Ne používejte znovu / Neponožte znovu / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolosiți / Не использовать повторно / Neponožujte opakovane / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kulanmayın / Не використовувати повторно

Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA

Benex Limited
 Pottery Road, Dun Laoghaire
 Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

Proclin is a trademark of Rohm and Haas Co.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD