

BD BACTEC™ MGIT™ 960 SIRE Kits

For the Antimycobacterial Susceptibility Testing of *Mycobacterium tuberculosis*

R_x Only  

8008200(05)
2019-09

INTENDED USE

The BD BACTEC™ MGIT™ 960 SIRE Kit is a rapid qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*, from culture, to streptomycin (STR), isoniazid (INH), rifampin (RIF) and ethambutol (EMB). The BD BACTEC MGIT 960 STR 4.0 Kit and the BD BACTEC MGIT 960 INH 0.4 Kit are for testing at higher drug concentrations.

The BD BACTEC MGIT 960 susceptibility test kits are used with the BD BACTEC MGIT 960 and BD BACTEC MGIT 320 Systems.

SUMMARY AND EXPLANATION

Antimycobacterial susceptibility testing is valuable in the proper treatment of patients with tuberculosis. The treatment of tuberculosis is commonly through a multiple drug regimen that includes the antimycobacterial drugs streptomycin, isoniazid, rifampin and ethambutol. It is important that the antimycobacterial drugs prescribed show appropriate activity against *Mycobacterium tuberculosis*, i.e., susceptibility of the isolate to the drug.

Multi-drug resistant *Mycobacterium tuberculosis* (MDR-TB) has recently become a serious public health problem.¹ Resistance to any of the primary drugs, streptomycin (STR), isoniazid (INH), rifampin (RIF) and ethambutol (EMB), makes the disease more difficult and expensive to treat. The rapid detection of these resistant isolates is critical to effective patient management.

A widely used method for antimycobacterial susceptibility testing, known as the Method of Proportion,² uses Middlebrook and Cohn 7H10 Agar. It compares colony counts on drug-containing and drug-free media. Resistance to a drug is detected when 1% or more of the bacterial population is resistant to the drug concentration under test. Results are generally available after 21 days of incubation.

Historically, the Method of Proportion (MOP) procedure has included susceptibility testing of *M. tuberculosis* using two concentrations of antimicrobials. The Clinical and Laboratory Standards Institute (CLSI) continues to recommend that the MOP test procedure include two concentrations of the primary drugs for testing except rifampin. The recommended low concentrations for the MOP procedure are the critical concentrations for these drugs. The critical concentration is defined as the drug concentration that allows the interpretation of a result as either resistant or susceptible. An isolate is determined resistant if 1% or more of the test population grows in the presence of the critical concentration of the drug. The high drug concentration is used to profile the degree of resistance within the population. This result provides information to the physician to assist in determining whether a modification to the therapy regimen is necessary.

The BD BACTEC MGIT 960 SIRE test provides the susceptibility results to be reported earlier, in most cases, than the MOP procedure.

The BD BACTEC MGIT 960 SIRE test was developed with critical concentrations for streptomycin, isoniazid, rifampin and ethambutol that are slightly lower than the critical concentrations used in the MOP in order to avoid false susceptibility. This is most apparent for streptomycin where many isolates are close to the recommended critical concentration as performed by the MOP. For this reason, a second, higher drug concentration was developed for streptomycin and isoniazid. A susceptible result at the critical concentration can be reported and no other testing is necessary. Isolates that are resistant at the critical concentration of streptomycin, isoniazid and/or ethambutol, should be tested at a higher drug concentration either in the BD BACTEC MGIT system or using an alternate method. In this case, a final result of resistant at the critical concentration may be reported, with notification that an additional test at a higher concentration is being performed.

Testing of resistant isolates at a higher concentration is important to identify those that exhibit low-level resistance, i.e., resistant at the critical concentration and susceptible at the high concentration. The high concentrations in the BD BACTEC MGIT system were designed to be lower than the concentrations used in the MOP. This design of the BD BACTEC MGIT system is such that a resistant result, especially for streptomycin, may not always correlate to a resistant result at the high concentration in MOP. In the event that a streptomycin result is obtained that is resistant at the high concentration, an alternate method of testing at this concentration should be performed.

PRINCIPLES OF THE PROCEDURE

The BD BBL™ MGIT 7 mL Mycobacteria Growth Indicator Tube is a tube containing a modified Middlebrook 7H9 Broth which supports the growth and detection of mycobacteria (see BD BBL MGIT 7 mL package insert). The MGIT tube contains a fluorescent compound embedded in silicone on the bottom of a 16 x 100 mm round-bottom tube. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. The initial concentration of dissolved oxygen quenches the emission from the compound, and little fluorescence can be detected. Later, actively growing and respiring microorganisms consume the oxygen that allows the compound to fluoresce.

The BD BACTEC MGIT 960 SIRE Kit is a 4–13 day qualitative test. The test is based on growth of the *M. tuberculosis* isolate in a drug-containing tube compared to a drug-free tube (Growth Control). The BD BACTEC MGIT instrument monitors tubes for increased fluorescence. Analysis of fluorescence in the drug-containing tube compared to the fluorescence of the Growth Control tube is used by the instrument to determine susceptibility results.

The BD BACTEC MGIT instrument automatically interprets these results, and reports a susceptible or resistant result.

REAGENTS

BD BACTEC MGIT 960 SIRE Kit contains one each lyophilized vials of streptomycin, isoniazid, rifampin and ethambutol and eight vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Streptomycin 332 µg
Approximate Formula* Per Vial Lyophilized drug: Isoniazid 33.2µg
Approximate Formula* Per Vial Lyophilized drug: Rifampin 332 µg
Approximate Formula* Per Vial Lyophilized drug: Ethambutol 1,660 µg

BD BACTEC MGIT 960 IR Kit contains one each lyophilized vials of isoniazid and rifampin and four vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Isoniazid 33.2µg
Approximate Formula* Per Vial Lyophilized drug: Rifampin 332 µg

BD BACTEC MGIT 960 STR 4.0 Kit contains one vial lyophilized streptomycin and two vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Streptomycin 664 µg

BD BACTEC MGIT 960 INH 0.4 Kit contains one vial lyophilized isoniazid and two vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Isoniazid 66.4µg

BD BACTEC MGIT 960 SIRE Supplement contains 20 mL Middlebrook OADC enrichment

Approximate Formula* Per L Purified Water

Bovine albumin 50.0 g Catalase 0.03 g
Dextrose 20.0 g Oleic Acid 0.6 g

**Adjusted and/or supplemented as required to meet performance criteria.*

Storage and reconstitution of reagents: BD BACTEC MGIT 960 SIRE Drug vials – On receipt, store the lyophilized drug vials at 2–8 °C. Once reconstituted, the antibiotic solutions may be frozen and stored at -20 °C or colder up to six months, not to exceed the original expiration date. Once thawed, use immediately. Discard unused portions.

Product Deterioration:

Some variation in the appearance of the lyophilized SIRE drugs may occur. This results from the lyophilization process and does not affect performance of the products.

BD BACTEC MGIT SIRE Supplement – On receipt, store in dark at 2–8 °C. Avoid freezing or overheating. Open and use prior to the expiration date. Minimize exposure to light.

Directions For Use:

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Streptomycin lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Isoniazid lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 8.3 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Rifampin lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Ethambutol lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 415 µg/mL.

NOTE: The following are reconstituted with a different volume. Failure to use the appropriate volume of sterile distilled water for reconstitution of the higher drug concentrations will invalidate those test results.

Reconstitute each BD BACTEC MGIT 960 STR 4.0 Kit Streptomycin lyophilized drug vial with **2 mL** of sterile distilled/deionized water to make a stock solution of 332 µg/mL.

Reconstitute each BD BACTEC MGIT 960 INH 0.4 Kit Isoniazid lyophilized drug vial with **2 mL** of sterile distilled/deionized water to make a stock solution of 33.2 µg/mL.

Warnings and Precautions

For *in vitro* Diagnostic Use.

POTENTIALLY INFECTIOUS TEST SPECIMEN: Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. “Standard Precautions”³⁻⁶ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

BD BACTEC MGIT 960 SIRE – Catalog number 245123

BD BACTEC MGIT 960 Rifampin

Danger



H301 Toxic if swallowed. **H332** Harmful if inhaled.

P261 Avoid breathing dust/fume/gas/mist/vapors/spray. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P271** Use only outdoors or in a well-ventilated area. **P301+P310** IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. **P330** Rinse mouth. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P321** Specific treatment (see on this label). **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

BD BACTEC MGIT 960 Ethambutol

Danger



H360 May damage fertility or the unborn child.

P280 Wear protective gloves/protective clothing/eye protection/face protection. **P201** Obtain special instructions before use. **P202** Do not handle until all safety precautions have been read and understood. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

BD BACTEC MGIT 960 IR – Catalog number 245157

BD BACTEC MGIT 960 IR Rifampin

Danger



H301 Toxic if swallowed. **H332** Harmful if inhaled.

P261 Avoid breathing dust/fume/gas/mist/vapors/spray. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P271** Use only outdoors or in a well-ventilated area. **P301+P310** IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. **P330** Rinse mouth. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Working with *M. tuberculosis* growth in culture requires Biosafety Level (BSL) 3 practices, containment equipment and facilities. Read and follow directions contained in all appropriate package inserts including the BD BBL MGIT 7 mL Mycobacteria Growth Indicator Tube.

Prior to use, the user should examine the tubes and vials for evidence of contamination or damage. Discard any tubes or vials if they appear unsuitable. Dropped tubes should be examined carefully. If damage is seen, the tube should be discarded.

In the event of tube breakage: 1) Close the instrument drawers; 2) Turn off the instrument; 3) Vacate the area immediately; 4) Consult your facility/CDC guidelines. An inoculated leaking or broken tube may produce an aerosol of mycobacteria; appropriate handling should be observed.

Autoclave all inoculated BD MGIT tubes prior to disposal.

INOCULUM PREPARATION

All preparations detailed below must be from pure cultures of *M. tuberculosis*. The laboratory should confirm, by appropriate identification techniques, that the isolate to be tested is a pure culture of *M. tuberculosis*.

Inoculum can be prepared from solid media or from a positive BD BACTEC MGIT 7 mL tube. In addition, cultures grown in liquid and on solid media can be used to prepare a seed BD MGIT tube, which can then be used to prepare the inoculum. Each of these options is described below.

Preparation of the Inoculum from Solid Media:

NOTE: It is important to prepare the inoculum according to the following instructions to obtain the appropriate organism concentration for the susceptibility test.

1. Add 4 mL of BD BBL Middlebrook 7H9 Broth (or BD BBL MGIT broth) to a 16.5 x 128 mm sterile tube with cap containing 8–10 glass beads.
2. Scrape with a sterile loop as many colonies as possible from growth no more than 14 days old, trying not to remove any solid medium. Suspend the colonies in the Middlebrook 7H9 Broth.
3. Vortex the suspension for 2–3 min to break up the larger clumps. The suspension should exceed a 1.0 McFarland standard in turbidity.
4. Let the suspension sit for 20 min without disturbing.
5. Transfer the supernatant fluid to another 16.5 x 128 mm sterile tube with cap (avoid transferring any of the sediment) and let sit for another 15 min.
6. Transfer the supernatant fluid (it should be smooth, free of any clumps) to a third 16.5 x 128 mm sterile tube. **NOTE:** The organism suspension should be greater than a 0.5 McFarland Standard at this step.
7. Adjust the suspension to a 0.5 McFarland standard by visual comparison to a 0.5 McFarland turbidity standard. Do not adjust below a 0.5 McFarland Standard.
8. Dilute 1 mL of the adjusted suspension in 4 mL of sterile saline (1:5 dilution). Proceed to “Inoculation Procedure for Susceptibility Test.”

Preparation of the Inoculum from a Positive BD BACTEC MGIT 7 mL Tube:

NOTE: It is important to prepare the inoculum using the following time references to obtain the appropriate organism concentration for the susceptibility test.

1. The first day of an instrument positive BD MGIT tube is considered Day 0.
2. For the preparation of the test inoculum, a positive 7 mL BD MGIT tube should be used the day **after** it first becomes positive on the BD BACTEC MGIT instrument (Day 1), up to and including the fifth day (Day 5) after instrument positivity. A tube which has been positive longer than five days should be subcultured to a fresh 7 mL BD MGIT tube containing BD BACTEC MGIT Growth Supplement and tested on the BD BACTEC MGIT instrument until positive, and used from one to five days following positivity. See “Preparation of a Seed BD MGIT Tube from Liquid Media.”
3. If the tube is a Day 1 or Day 2 positive, use the BD MGIT broth suspension for the inoculation procedures. Mix well. Proceed to “Inoculation Procedure for Susceptibility Test.”
4. If the tube is a Day 3, Day 4, or Day 5 positive, mix well then dilute 1 mL of the positive broth in 4 mL of sterile saline (1:5 dilution). Mix tube thoroughly. Use the diluted suspension for the inoculation procedures. Proceed to “Inoculation Procedure for Susceptibility Test.”

Preparation of a Seed BD MGIT Tube from Liquid Media

1. Mix the tube by inversion or vortexing.
2. Make a 1:100 dilution by adding 0.1 mL of the culture into 10 mL of BD BBL Middlebrook 7H9 Broth or BD BBL MGIT Broth. Mix well.
3. Add 0.5 mL of this suspension into a 7 mL BD MGIT tube supplemented with 0.8 mL of BD BACTEC MGIT Growth Supplement.
4. Cap tightly and gently mix by inverting 2–3 times.
5. Enter the tube into a BD BACTEC MGIT instrument and test until positive.
NOTE: Time to positivity **must** be ≥ 4 days for use as AST inoculum. If tube becomes positive in < 4 days, return to step 1 and prepare a new seed tube.
6. This tube may now be used from one to five days following positivity. Proceed to “Preparation of the Inoculum from a Positive BD BACTEC MGIT 7 mL Tube” above.

Preparation of a Seed BD MGIT Tube from Solid Media

1. Using a sterile loop, scrape growth from a slant and add to a 7 mL BD MGIT tube supplemented with 0.8 mL of BD BACTEC MGIT Growth Supplement.
2. Cap tightly and gently mix by inverting 2–3 times.
3. Enter the tube into a BD BACTEC MGIT instrument and test until positive.
NOTE: Time to positivity **must** be ≥ 4 days for use as AST inoculum. If tube becomes positive in < 4 days, return to step 1 and prepare a new seed tube.
4. This tube may now be used from one to five days following positivity. Proceed to “Preparation of the Inoculum from a Positive BD BACTEC MGIT 7 mL Tube” above.

PROCEDURE

Materials Provided: BD BACTEC MGIT 960 SIRE Kit containing one vial each lyophilized drug and eight vials of SIRE Supplement (approximately 40 tests per drug per kit). BD BACTEC MGIT 960 STR 4.0 Kit containing one vial each lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit) and BD BACTEC MGIT 960 INH 0.4 Kit containing one vial each lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit).

Materials Required But Not Provided: BD BBL MGIT 7 mL Mycobacteria Growth Indicator Tubes, ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Inoculation Procedure for BD BACTEC MGIT 960 SIRE Kit Susceptibility Test:

1. Label five 7 mL BD MGIT tubes for each test isolate. Label one as GC (Growth Control), one as STR, one as INH, one as RIF and one as EMB. Place the tubes in the correct sequence in the appropriate size AST Set Carrier (see BD BACTEC MGIT Instrument User's Manual).
2. Aseptically add 0.8 mL of BD BACTEC MGIT SIRE Supplement to each tube. NOTE: It is important to use the supplement supplied with the kit.
3. Aseptically pipet, using a micropipet, 100 µL of 83 µg/mL MGIT STR solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 8.3 µg/mL MGIT INH solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 83 µg/mL MGIT RIF solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 415 µg/mL MGIT EMB solution to the appropriately labeled BD MGIT tube. It is important to add the correct drug to the corresponding tube. No antibiotics should be added to the BD MGIT GC tube.

Drug	Concentration of Drug after Reconstitution*	Volume Added to BD MGIT Tubes for Test	Final Concentration in BD MGIT Tubes
MGIT STR	83 µg/mL	100 µL	1.0 µg/mL
MGIT INH	8.3 µg/mL	100 µL	0.1 µg/mL
MGIT RIF	83 µg/mL	100 µL	1.0 µg/mL
MGIT EMB	415 µg/mL	100 µL	5.0 µg/mL

* These drugs must be reconstituted using 4 mL sterile/deionized water to achieve concentrations indicated.

4. **Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see "INOCULUM PREPARATION") into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate 0.5 mL of the 1:100 Growth Control suspension into the BD MGIT tube labeled "GC."
5. **Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see "INOCULUM PREPARATION") into each of the FOUR remaining drug tubes (STR, INH, RIF, EMB).
6. Tightly recap the tubes. Mix tubes thoroughly by gentle inversion three to four times.
7. Enter the AST set into the BD BACTEC MGIT instrument using the AST set entry feature (refer to the BD BACTEC MGIT Instrument User's Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the set carrier definitions selected when performing the AST set entry feature.
8. Streak 0.1 mL of the organism suspension to a BD Trypticase™ Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35–37 °C.
9. Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the BD BACTEC MGIT Instrument User's Manual) and repeat testing with pure culture.

Inoculation Procedure for BD BACTEC MGIT 960 STR 4.0 and INH 0.4 Kits Susceptibility Test:

It is recommended if resistance occurs at the critical concentration, a susceptibility profile test be performed which, at a minimum, tests the high concentration of the drug to which the isolate was originally resistant.

Isolate Source: The isolate used for this testing must have been prepared as described in "INOCULUM PREPARATION." A seed tube may be prepared from the drug-free Growth Control tube from the previously tested AST set of the isolate, by inoculating 0.5 mL to a fresh BD MGIT 7 mL tube containing BD BACTEC MGIT Growth Supplement. Once the seed tube is instrument positive, proceed as described in "INOCULUM PREPARATION: Preparation of the Inoculum from a Positive BD MGIT tube."

1. Label enough BD MGIT 7 mL tubes for the test isolate to have a MGIT GC (Growth Control) and a MGIT drug tube for each antimicrobial tested. Place the tubes in the correct sequence in the appropriate size AST Set Carrier (see BD BACTEC MGIT Instrument User's Manual).
2. Aseptically add 0.8 mL of BD BACTEC MGIT SIRE Supplement to each tube. NOTE: It is important to use the supplement supplied with the kit.
3. Aseptically pipet, using a micropipet, 100 µL of the drug solution to the appropriately labeled BD MGIT tube. It is important to add the correct drug to the corresponding tube. No antibiotics should be added to the BD MGIT GC tube.

Drug	Concentration of Drug after Reconstitution*	Volume Added to BD MGIT Tubes for Test	Final Concentration in BD MGIT Tubes
MGIT STR 4.0	332 µg/mL	100 µL	4.0 µg/mL
MGIT INH 0.4	33.2 µg/mL	100 µL	0.4 µg/mL

* These drugs must be reconstituted using 2 mL sterile/deionized water to achieve concentrations indicated.

4. **Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see "INOCULUM PREPARATION") into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate 0.5 mL of the 1:100 Growth Control suspension into the BD MGIT tube labeled "GC."
5. **Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see "INOCULUM PREPARATION") into each of the drug tubes.
6. Tightly recap the tubes. Mix tubes thoroughly by gentle inversion three to four times.
7. Enter the AST set into the BD BACTEC MGIT instrument using the AST set entry feature (refer to the BD BACTEC MGIT Instrument User's Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the set carrier definitions selected when performing the AST set entry feature.

8. Streak 0.1 mL of the organism suspension to a BD Trypticase Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35–37 °C.
9. Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the BD BACTEC MGIT Instrument User's Manual) and repeat testing with pure culture.

NOTE: The susceptibility test may be configured in a variety of formats. For example, a five tube carrier set containing only the critical concentrations can be configured into the system. A variety of other tube carrier sets can be configured depending on the optional profile tests being run (refer to BD BACTEC MGIT Instrument User's Manual).

User Quality Control: Upon receipt of a new shipment or lot number of BD BACTEC MGIT 960 SIRE Kit vials, it is recommended that the control organism shown below be tested. The control organism should be a pure culture and the culture should be prepared according to the "INOCULUM PREPARATION" instructions.

The quality control (QC) AST Set should be prepared according to the "Inoculation Procedure for Susceptibility Test" instructions for the drug kits being tested. Important considerations when preparing the QC AST Set are the proper reconstitution of the lyophilized drugs and the proper dilution of the QC organism for the Growth Control and drug tubes.

It is important to add the appropriate drug to the corresponding labeled tube. The use of the pan-sensitive QC organism will not detect incorrect drug pipetted into AST Set tubes.

Observation of the proper results, as shown below, within 4–13 days indicates that the BD BACTEC MGIT 960 SIRE Kits are ready for use in testing patient isolates. If the proper results are not observed, repeat the test. If, after repeating the test, the proper results are still not observed, do not use the product until you have contacted BD Technical Service and Support at 1.800.638.8663 (United States only).

Strain	GC	MGIT STR	MGIT INH	MGIT RIF	MGIT EMB
<i>M. tuberculosis</i> ATCC® 27294	Positive	Susceptible	Susceptible	Susceptible	Susceptible

Strain	GC	MGIT STR 4.0	MGIT INH 0.4
<i>M. tuberculosis</i> ATCC 27294	Positive	Susceptible	Susceptible

The same control organism should be run as batch QC once each week when susceptibility testing is performed. If the batch QC fails, do not report patient results for the drug(s) that failed for that testing period. Repeat the QC for the drug(s) and patient isolates affected by the initial QC failure. If the repeat QC does not perform as expected, do not report patient results. Do not use the product until you have contacted BD Technical Service and Support at 1.800.638.8663 (United States only).

During the external evaluation of the BD BACTEC MGIT 960 SIRE Kits, the most common causes of QC failure were contaminated QC cultures, over/under inoculated AST Sets, drug not added to appropriate tubes and instrument error conditions.

RESULTS

The BD BACTEC MGIT instrument will monitor AST Sets until a susceptible or resistant determination is made. Once the set testing is completed, the results are reported by the BD BACTEC MGIT instrument (refer to the BD BACTEC MGIT Instrument User's Manual). The BD BACTEC MGIT instrument will report an AST Set result as an Error(X), no susceptibility interpretation, when certain conditions occur that may affect the test results. Conditions that may result in an Error(X) result are described in Section 7-Troubleshooting of the BD BACTEC MGIT Instrument User's Manual.

When reporting results, it is important to include the test method, drug name and concentration, whether the result is obtained with the BD BACTEC MGIT System or an alternate method. The Pulmonary and/or Infectious Disease specialist in TB control should be consulted concerning the appropriate therapeutic regimen and dosages.

In the event of unexpected resistant results, verify identification of isolate tested as *M. tuberculosis*. Ensure that only a pure culture was used (rule-out presence of mixed mycobacteria, etc.) Mono-resistance to ethambutol is uncommon and should be verified.^{2,7}

BD BACTEC MGIT 960 SIRE critical concentration result reporting

Drug (concentration)	MGIT System Result	Recommended Report	Action
STR (1.0 µg/mL)	Susceptible (SIRE)	Isolate tested with BD BACTEC MGIT System [drug/concentration] and result is susceptible.	No report.
INH (0.1 µg/mL)	Resistant (SIR)	Isolate tested with BD BACTEC MGIT System [drug/concentration] and result is resistant. Results of testing [drug] at a higher concentration to follow (if tested).	Recommend testing at higher concentration (STR and/or INH).
RIF (1.0 µg/mL)	Resistant (E)	If resistant to more than ethambutol (EMB) Isolate tested with BD BACTEC MGIT System [Ethambutol 5.0 µg/mL] and result is resistant. Consult laboratory for testing EMB at a higher concentration. If mono-resistant to ethambutol (EMB) Isolate tested with BD BACTEC MGIT System [Ethambutol 5.0 µg/mL] and result is resistant. Mono-resistance to ethambutol is uncommon. Consult laboratory for confirmation.	Recommend testing EMB at higher concentration using alternate method. Recommend testing EMB at both critical concentration and higher concentration with an alternate method.
EMB (5.0 µg/mL)	Error (X)	No report.	Repeat test.

BD BACTEC MGIT 960 STR 4.0 and INH 0.4 result reporting

Drug (concentration)	MGIT System Result	Recommended Report	Action
STR (4.0 µg/mL)	Susceptible	Isolate tested with BD BACTEC MGIT System streptomycin 4.0 µg/mL and result is susceptible. This isolate, with a resistant result at 1.0 µg/mL and a susceptible result at 4.0 µg/mL, indicates low-level resistance to streptomycin.	No action.
	Resistant	Isolate tested with BD BACTEC MGIT System streptomycin 4.0 µg/mL and result is resistant. Consult laboratory for confirmation.	Isolate should be tested by an alternate method to verify result.
	Error (X)	No report.	Repeat test.
INH (0.4 µg/mL)	Susceptible	Isolate tested with BD BACTEC MGIT System isoniazid 0.4 µg/mL and result is susceptible. This isolate, with a resistant result at 0.1 µg/mL and a susceptible result at 0.4 µg/mL, indicates low-level resistance to isoniazid.	No action.
	Resistant	Isolate tested with BD BACTEC MGIT System isoniazid 0.4 µg/mL and result is resistant.	No action.
	Error (X)	No report.	Repeat test.

LIMITATIONS OF THE PROCEDURE

The BD BACTEC MGIT System susceptibility test does not interpret the degree of susceptibility of the isolate being tested. Results are reported as either S, susceptible, or R, resistant, for the drug and concentration tested.

The BD BACTEC MGIT 960 SIRE test was developed with critical concentrations for streptomycin, isoniazid, rifampin and ethambutol that are slightly lower than the critical concentrations used in the MOP in order to avoid false susceptibility. Testing of the higher concentrations, as recommended, will enhance the ability to detect isolates with low-level resistance.

The BD BACTEC MGIT System susceptibility tests can only be performed using a BD BACTEC MGIT instrument. The AST Sets cannot be read manually.

Use only pure cultures of *M. tuberculosis*. Cultures that are contaminated or that may contain multiple species of mycobacteria may give erroneous results and should not be tested. Direct testing from clinical specimens is not recommended.

Suspensions made from solid media must be allowed to settle for the prescribed times prior to standardization. Inoculum preparations made from solid media should be visually compared to a 0.5 McFarland turbidity standard; failure to do so may give inaccurate results or cause an AST Set error.

Failure to use the 1:5 dilution of the organism suspension, when indicated, to inoculate the drug containing tubes may give inaccurate results.

Failure to use a 1:100 dilution of the organism suspension for the inoculation of the Growth Control tube may give inaccurate results or cause an AST Set error.

Failure to reconstitute the drugs with the appropriate volume of sterile distilled/deionized water may give inaccurate results.

Thorough mixing of inoculated tubes is important. Failure to mix the tubes adequately can lead to false resistant results.

Failure to load the tubes of the AST Set into the AST Set Carrier in the proper sequence may give inaccurate results. Failure to select the appropriate set carrier drug definition may result in invalid or inaccurate results.

Failure to load the AST Set into the instrument correctly will result in an anonymous condition that must be resolved within eight hours. If condition is not resolved within eight hours, the AST Set must be discarded and set up again.

Failure to use the SIRE Supplement in the AST Set may give inaccurate results. DO NOT add BD BACTEC MGIT Growth Supplement to the AST Set.

EXPECTED VALUES

A total of 106 clinical isolates of *M. tuberculosis* were tested with the BD BACTEC MGIT 960 SIRE susceptibility test at four geographically diverse sites. The testing included both fresh clinical and stock isolates from both liquid and solid culture sources. A total of 200 susceptibility tests (liquid and solid) were performed at the critical concentrations of streptomycin (STR), isoniazid (INH) and rifampin (RIF) and a total of 223 susceptibility tests (liquid and solid) were performed at the critical concentration of ethambutol (EMB) during separate testing. The overall average time-to-result for the BD BACTEC MGIT 960 SIRE susceptibility test is seven to eight days with a range from four to fourteen days. The data are shown in Figure 1 at end of insert.

PERFORMANCE CHARACTERISTICS

ANALYTICAL STUDIES

Liquid and Solid Media AST Inoculum Ranges:

Liquid media – The recommended procedure for preparing an AST Set from a positive BD MGIT 7 mL tube uses a direct inoculum on Day 1 and Day 2 post-positivity and a dilute (1:5) inoculum on Day 3 to Day 5 post-positivity. Internal studies show that inocula prepared from a Day 1 to Day 5 positive BD MGIT 7 mL tube range between 0.8×10^5 to 3.2×10^5 CFU/mL.

Solid media – The recommended procedure for preparing an AST Set from growth on solid media (up to 14 days) uses a 1:5 dilution of an organism suspension equivalent to a 0.5 McFarland Standard. Internal studies show that inocula prepared from solid medium culture range between 1.4×10^5 to 2.4×10^6 CFU/mL.

Lot Reproducibility:

Lot reproducibility was evaluated using twenty-five *M. tuberculosis* isolates (to include five ATCC strains). Each BD BACTEC MGIT 960 SIRE test at the critical drug concentration was performed in triplicate for a total of seventy-five results per drug. Each replicate represented a separate test condition differentiated by lot of SIRE drug and SIRE Supplement used (three lots each).

Those isolates that were determined resistant to streptomycin, isoniazid or ethambutol in the initial test were then tested with the high drug concentration, except the ATCC strains. In addition to the resistant isolates tested, two sensitive isolates to STR (critical concentration), two sensitive isolates to INH (critical concentration) and two sensitive isolates to EMB (critical concentration) were included in the susceptibility profile test. Observed results were compared to the expected results.

The overall reproducibility for each drug at the critical concentration is 96% for STR, 100% for INH, 100% for RIF and 100% for EMB. The overall reproducibility for each drug at the high concentration is 96% for STR 4.0 and 100% for INH 0.4.

CDC Challenge Panel Testing:

The performance of the BD BACTEC MGIT 960 SIRE susceptibility test was evaluated using a panel of challenge isolates obtained from the Centers for Disease Control and Prevention (CDC), GA, USA. The panel consisted of thirty isolates of *M. tuberculosis* with known susceptibility patterns (using MOP). The panel was tested twice with the BD BACTEC MGIT 960 SIRE susceptibility test and both results were in agreement. The BD BACTEC MGIT 960 SIRE results were compared to the CDC expected results.

The overall agreement with CDC expected results for each drug at the critical concentration is 93% for STR, 100% for INH, 100% for RIF and 100% for EMB. The overall agreement with CDC expected results for each drug at the high concentration is 100% for STR 4.0 and 100% for INH 0.4.

CLINICAL EVALUATION

The BD BACTEC MGIT 960 SIRE susceptibility test was evaluated at four geographically diverse clinical sites, composed of regional reference centers and university hospital-based laboratories, including one foreign site. The BD BACTEC MGIT 960 SIRE susceptibility test was compared to the Method of Proportion (MOP)² susceptibility test method. The initial evaluation included the drugs streptomycin, isoniazid and rifampin. A separate evaluation was performed for the drug ethambutol.

Reproducibility Testing:

The reproducibility of the BD BACTEC MGIT 960 SIRE test was evaluated at the clinical sites using a panel of ten qualified isolates, including several isolates resistant to each of the drugs. The BD BACTEC MGIT 960 SIRE test results were compared to the expected results. The overall reproducibility for each drug at the critical concentration is 98.9% for STR, 99.7% for INH, 99.2% for RIF and 97.5% for EMB. Individual site reproducibility ranged from 89.9% to 100% for the combined critical concentration drug results. The overall reproducibility for each drug at the high concentration is 99.7% for STR 4.0 and 95.6% for INH 0.4. Individual site reproducibility ranged from 92.2% to 100% for the combined high concentration drug results.

CDC Challenge Panel Testing:

The performance of the BD BACTEC MGIT 960 SIRE susceptibility test was evaluated using a panel of challenge isolates obtained from the Centers for Disease Control and Prevention (CDC), GA, USA. The panel consisted of thirty isolates of *M. tuberculosis* with known susceptibility patterns (using MOP) tested by each clinical site.

Table 1 shows the agreement of the BD BACTEC MGIT 960 SIRE susceptibility test for each drug compared to the CDC expected results.

Table 1: CDC Challenge Panel – BD BACTEC MGIT 960 Clinical Site Testing

MGIT 960	Number tested	Number Correct	% Correct
STR 1.0	120	111	92.5
INH 0.1	120	119	99.2
RIF 1.0	120	120	100
EMB 5.0	119	111	93.3
STR 4.0	29*	29	100
INH 0.4	87*	82	94.3

* Only isolates resistant at critical concentrations tested at STR 4.0 and INH 0.4.

Clinical Isolate Testing:

A total of 106 clinical isolates of *M. tuberculosis* were tested with the BD BACTEC MGIT 960 SIRE susceptibility test and the MOP susceptibility test. This included testing of both fresh clinical and stock isolates from both liquid and solid culture sources. This generated a total of 195 test results for the initial susceptibility test performed for streptomycin, isoniazid and rifampin (critical concentration). A separate evaluation of ethambutol was performed from frozen aliquots of the original clinical and stock isolates as well as prospective clinical isolates from both liquid and solid culture sources. This generated a total of 223 test results for the ethambutol test at the critical concentration.

Table 2 presents the results from clinical isolate testing for each drug (critical concentration) from liquid source cultures. Table 3 presents the results from clinical isolate testing for each drug (critical concentration) from solid source cultures.

Table 2: Clinical Isolate Results – BD BACTEC MGIT 960 AST Compared to Method of Proportion from Liquid Source Cultures

Method of Proportion				MGIT 960 AST System				
DRUG	Concentration	S	R	Concentration	Susceptible Results		Resistant Result	
					# Agreement	% Agreement (95% CI)	# Agreement	% Agreement (95% CI)
STR	2.0 µg/mL	69	27	1.0 µg/mL	62	90 (80–96)	26	96 (81–100)
INH	0.2 µg/mL	59	37	0.1 µg/mL	57	97 (88–100)	36	97 (86–100)
RIF	1.0 µg/mL	72	24	1.0 µg/mL	71	99 (93–100)	24	100 (95–100)
EMB	5.0 µg/mL	91	20	5.0 µg/mL	88	97 (91–99)	17	85 (62–97)

All isolates with discordant MGIT results were tested by MOP at two independent sites. Of the seven discordant STR resistant (R-960, S-MOP) isolates, three had resistant results from both sites and one had susceptible results from both sites. The remaining three had resistant results from one site and susceptible results from the other site. The discordant STR susceptible (S-960, R-MOP) isolate had susceptible results from both sites. The two discordant INH resistant (R-960, S-MOP) isolates had susceptible results from both sites. The discordant INH susceptible (S-960, R-MOP) isolate had susceptible results from both sites. The discordant RIF resistant (R-960, S-MOP) isolate had resistant results from both sites. The three discordant EMB resistant (R-960, S-MOP) isolates had susceptible results from both sites. Of the three discordant EMB susceptible (S-960, R-MOP) isolates, two had susceptible results from both sites and one had a resistant result from one site and a susceptible result from the other site.

Table 3: Clinical Isolate Results – BD BACTEC MGIT 960 AST Compared to Method of Proportion from Solid Source Cultures

Method of Proportion				MGIT 960 AST System				
DRUG	Concentration	S	R	Concentration	Susceptible Results		Resistant Result	
					# Agreement	% Agreement (95% CI)	# Agreement	% Agreement (95% CI)
STR	2.0 µg/mL	70	29	1.0 µg/mL	65	93 (84–98)	28	97 (82–100)
INH	0.2 µg/mL	63	36	0.1 µg/mL	62	98 (92–100)	35	97 (86–100)
RIF	1.0 µg/mL	70	29	1.0 µg/mL	70	100 (95–100)	26	90 (73–98)
EMB	5.0 µg/mL	87	25	5.0 µg/mL	86	99 (94–100)	20	80 (59–93)

All isolates with discordant MGIT results were tested by MOP at two independent sites. Of the five discordant STR resistant (R-960, S-MOP) isolates, two had resistant results from both sites and one had susceptible results from both sites. The remaining two had resistant results from one site and susceptible results from the other site. The discordant STR susceptible (S-960, R-MOP) isolate had resistant results from both sites. The discordant INH resistant (R-960, S-MOP) isolate had susceptible results from both sites. The discordant INH susceptible (S-960, R-MOP) isolate had resistant results from both sites. The three discordant RIF susceptible (S-960, R-MOP) isolates had susceptible results from both sites. The one discordant EMB resistant (R-960, S-MOP) isolate had resistant results from both sites. Of the five discordant EMB susceptible (S-960, R-MOP) isolates, four had susceptible results from both sites. The remaining isolate had a resistant result from one site and a susceptible result from the other site.

Table 4 presents the results from clinical isolate testing for streptomycin and isoniazid (high concentration) from liquid source cultures. Table 5 presents the results from clinical isolate testing for streptomycin and isoniazid (high concentration) from solid source cultures.

Table 4: Clinical Isolate Results – BD BACTEC MGIT 960 AST Compared to Method of Proportion from Liquid Source Cultures

Method of Proportion				MGIT 960 AST System				
DRUG	Concentration	S	R	Concentration	Susceptible Results		Resistant Result	
					# Agreement	% Agreement (95% CI)	# Agreement	% Agreement (95% CI)
STR	10.0 µg/mL	77	19	4.0 µg/mL	73*	95 (87–99)	17	90 (67–99)
INH	1.0 µg/mL	65	31	0.4 µg/mL	65*	100 (95–100)	29	94 (74–99)

* Assumes MGIT high drug S result for all isolates with MGIT low drug S result.

All isolates with discordant MGIT results were tested by MOP at two independent sites. The four discordant STR resistant (R-960, S-MOP) isolates had susceptible results from both sites. Of the two discordant STR susceptible (S-960, R-MOP) isolates, one had susceptible results from both sites and one had resistant results from both sites. Of the two discordant INH susceptible (S-960, R-MOP) isolates, one had susceptible results from both sites and one had resistant results from both sites.

Table 5: Clinical Isolate Results – BD BACTEC MGIT 960 AST Compared to Method of Proportion from Solid Source Cultures

Method of Proportion		MGIT 960 AST System						
DRUG	Concentration	S	R	Concentration	Susceptible Results		Resistant Result	
					# Agreement	% Agreement (95% CI)	# Agreement	% Agreement (95% CI)
STR	10.0 µg/mL	78	21	4.0 µg/mL	73*	94 (86–98)	17	81 (58–95)
INH	1.0 µg/mL	68	31	0.4 µg/mL	68*	100 (95–100)	30	87 (83–100)

* Assumes MGIT high drug S result for all isolates with MGIT low drug S result.

All isolates with discordant MGIT results were tested by MOP at two independent sites. The five discordant STR resistant (R-960, S-MOP) isolates had susceptible results from both sites. Of the four discordant STR susceptible (S-960, R-MOP) isolates, three had susceptible results from both sites and one had resistant results from both sites. The discordant INH susceptible (S-960, R-MOP) isolate had resistant results from both sites.

AVAILABILITY

Cat. No. Description


- 245123 BD BACTEC™ MGIT™ 960 SIRE Kit, carton of 4 lyophilized drug vials, and 8 SIRE Supplements.
- 245125 BD BACTEC™ MGIT™ 960 STR 4.0 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.
- 245157 BD BACTEC™ MGIT™ 960 IR Kit, carton of 2 lyophilized drug vials, and 4 SIRE Supplements.
- 245126 BD BACTEC™ MGIT™ 960 INH 0.4 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.

REFERENCES

- Barenfanger, J. 1993. Making your lab safe against multi-drug resistant *Mycobacterium tuberculosis*. Clin. Microbiol. Newsl. 15: 76–80.
- Clinical and Laboratory Standards Institute. 2003. Approved Standard M24-A. Susceptibility testing of mycobacteria, nocardiae, and other aerobic actinomycetes. CLSI, Wayne, Pa.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53–80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021–0045.
- Ridderhof, J. 2001. Multicenter evaluation of 3.75 µg/mL ethambutol (EMB) in BD BACTEC vials for susceptibility testing of *Mycobacterium tuberculosis*. Abstract C-244, American Society for Microbiology Abstracts 2001.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

Figure 1: Distribution of BD BACTEC MGIT 960 AST Time in Protocol



Change History

Revision	Date	Change Summary
(05)	2019-09	Converted printed instructions for use to electronic format and added access information to obtain the document from BD.com/e-labeling. Per Safety Data Sheet for catalog number 245123 existing precautionary codes and statements updated for BD BACTEC MGIT 960 Rifampin; added health hazard pictogram, signal word "Danger", all hazard & precautionary codes and statements for BD BACTEC MGIT 960 Ethambutol. Per Safety Data Sheet for catalog number 245157 existing hazard & precautionary codes and statements updated for BD BACTEC MGIT 960 Rifampin.

US Customers only: For symbol glossary, refer to www.bd.com/symbols-glossary



Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvođač / Tillverkare / Üretici / Виробник / 生产厂商



Use by / Използвайте до / Spotføjbejdt do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Използovať до / Použít do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
AAAA-MM-DD / AAAA-MM (MM = fin del mes)
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
EEEE-HH-NN / EEEE-HH (HH = hónap utolsó napja)
AAAA-MM-GG / AAAA-MM (MM = fine mese)
ЖОЖОЖ-АА-КК / ЖОЖОЖ-АА / (АА = айдың соңы)
YYYY-MM-DD/YYYY-MM (MM = 월말)
MMMM-MM-DD / MMMM-MM (MM = mėnesio pabaiga)
GGGG-MM-DD/GGGG-MM (MM = mėneša beigas)
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
AAAA-MM-DD / AAAA-MM (MM = slutten av månaden)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = fim do mês)
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
YYYY-MM-DD / YYYY-MM (MM = 月末)



Catalog number / Каталоген номер / Katalógové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номер / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalógové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号



Authorized Representative in the European Community / Оторизирани представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret representant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizuirani predstavnik u Evropskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유럽 공동체의 위임 대표 / Įgaliotasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Аврпа Топлuluğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 欧洲共同体授权代表



In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiiniparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnosztikai orvosi eszköz / Dispositivo medical per diagnostica in vitro / Жасанды жағдайда жүргізілетін медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietais / Medicinas ierīces, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispozitiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro / Medicinska pomůcka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Diagnostik Tibbi Cihaz / Медицинский прибор для диагностики in vitro / 体外诊断医疗设备



Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperaturlimit / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制



Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)



Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenido suficiente per <n> test / <n> тесттері үшін жеткінікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Continut suficient pentru <n> teste / Достаточное для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli malmaze içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测



Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaityti naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明



Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Nepoužívejte opakovaně / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kulanmayin / Не використовувати повторно / 请勿重复使用



Serial number / Серийн номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / N° de serie / Seerianumber / Numéro de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Sérijas numurs / Serie nummer / Numer seryjny / Número de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号

	For IVD Performance evaluation only / Само за оценка качества работы на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Μόνο για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réserve à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietaisų veikimo charakteristikoms tikrinti / Vienīgi IVD darbības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluering av IVD-ytelse / Tylko do oceny wydajności IVD / Uso exclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики in vitro / Určenie iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirmesi için / Тільки для оцінювання якості діагностики in vitro / 仅限 IVD 性能评估
	For US: "For Investigational Use Only"
	Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūros žemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrense / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Мінімальна температура / 温度下限
CONTROL	Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / 对照
CONTROL +	Positive control / Положителен контрол / Pozitivni kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitivna kontrola / Pozitív kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívá kontrolé / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивний контроль / 阳性对照试剂
CONTROL -	Negative control / Отрицателен контрол / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negatiivne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo / Негативті бақылау / 음성 컨트롤 / Neigiama kontrolė / Negatívá kontrolé / Negatieve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативний контроль / 阴性对照试剂
STERILE EO	Method of sterilization: ethylene oxide / Метод на стерилизация: етиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringensmetode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseringensmetode: etüleenoksiid / Méthode de stérilisation : oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metoda di sterilizzazione: ossido di etilene / Sterilizacija edici – etilen totyгы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksisds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringensmetode: etylenoksid / Metoda sterylizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metodá sterilizácie: etylenoxid / Metoda sterilizacije: etilen oksid / Steriliseringensmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷
STERILE R	Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringensmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseringensmetode: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metoda di sterilizzazione: irradiazione / Sterilizacija edici – сәуле тысы / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringensmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metodá sterilizácie: ožiarenie / Metoda sterilizacije: ozračevanje / Steriliseringensmetod: strålning / Sterilizasyon yöntemi: ırradyasyon / Метод стерилизації: опроміннення / 灭菌方法: 辐射
	Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologilised riskid / Risques biologiques / Biološki rizik / Biológiai veszélyes / Rischio biologico / Биологический тауекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Biologische опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Біологічна небезпека / 生物学风险
	Caution, consult accompanying documents / Внимание, направете справка в придружаващите документи / Pozor! Pristudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatust! Lugeka kaasnevat dokumentatsiooni! / Attention, consulter les documents joints / Urozorenje, koristi prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dămesio, žiūrėkite priedamus dokumentus / Piesardziba, skatit pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vystraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увага: див. супутню документацію / 小心, 请参阅附带文档。
	Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Ausgėjá temperatūros robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Góna granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sıcaklık üst sınırı / Міксимальна температура / 温度上限
	Keep dry / Пазете сухо / Skladujte v suchém prostředi / Opbevarer tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausiai / Uzglabāt sausu / Droog houden / Holdes tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Бергетти від вологості / 请保持干燥
	Collection time / Време на събиране / Čas odběru / Opsamlingsstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingsstid / Toplama zamanı / Час забору / 采集时间
	Peel / Обелете / Otevfete zde / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстиғи қабатын алып таста / 벗기 / Plęsti čia / Atfimt / Schillen / Trek av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odrhните / Oljuštiti / Dra isär / Ayırma / Відкрити / 撕下
	Perforation / Перфорация / Perforace / Perforering / Διάτρηση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 찔뽀침 / Perforacija / Perforácia / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔
	Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά / No usar si el paquete está dañado / Mit kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Ереп пакет бұзылған болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用
	Keep away from heat / Пазете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valguusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Сапқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не награвать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Бергетти від дії тепла / 请远离热源
	Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Reži / Vágja ki / Tagliare / Keciңiz / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciąć / Cortar / Decupaji / Отрезать / Odstřihněte / Iseći / Klipp / Kesme / Pozpizati / 剪下



Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuurpäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paemimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Datum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期



μL/test / μL/тест / μL/Test / μL/εξέταση / μL/prueba / μL/teszt / μL/테스트 / мкл/тест / μL/tyrimas / μL/pårbaude / μL/teste / мкл/анализ / μL/检测



Keep away from light / Παзете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródła światła / Manter ao abrigo da luz / Feriți de lumină / Хранить в темноте / Uchovávajte mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Беретти від дії світла / 请远离光线



Hydrogen gas generated / Образуван е водород газ / Možnost úniku plynného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadrží hydrogen vodik / Hidrogén gázt fejleszt / Produzione di gas idrogeno / Газтөктес сутері пайда болды / 수소 가스 생성됨 / Išskiria vandenilio dujas / Rodas ūdenradis / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Osloбаda se vodonik / Genererad vätgas / Açığa çıkan hidrojen gazı / Реакция з виділенням водню / 会产生氢气



Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттің идентификациялық нөмірі / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificationnummer van de patiënt / Pasientens ID-nummer / Numer ID pacjenta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентифікатор пацієнта / 患者标识号



Fragile, Handle with Care / Чупливо, Работете с необходимото внимание. / Křehké. Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Εύθραστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Őm, kásitsege ettevaatlíult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сығыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkites atsargiai. / Trausls; rikoties uzmanīgi / Breekbaar, voorzichtig behandelen. / Ømtålí, håndter forsigtig. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşıyin. / Тендітна, звертатися з обережністю / 易碎, 小心轻放

R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."



bd.com/e-labeling

KEY-CODE: 8008200

Europe, CH, GB, NO:		+800 135 79 135
International:		+31 20 794 7071
AR	+800 135 79 135	LT 8800 30728
AU	+800 135 79 135	MT +31 20 796 5693
BR	0800 591 1055	NZ +800 135 79 135
CA	+1 855 805 8539	RO 0800 895 084
CO	+800 135 79 135	RU +800 135 79 135
EE	0800 0100567	SG 800 101 3366
GR	00800 161 22015 7799	SK 0800 606 287
HR	0800 804 804	TR 00800 142 064 866
IL	+800 135 79 135	US +1 855 236 0910
IS	800 8996	UY +800 135 79 135
LI	+31 20 796 5692	VN 122 80297



Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA



Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC® is a trademark of the American Type Culture Collection.

BD, the BD Logo, BACTEC, BBL, MGIT and Trypticase are trademarks of Becton, Dickinson and Company or its affiliates.

© 2019 BD. All rights reserved.