

# BD BACTEC™ MGIT™ 960 SIRE Kit

For the Antimycobacterial Susceptibility Testing of *Mycobacterium tuberculosis*



88-2041-1JAA(04)  
2019-09

## INTENDED USE

The BD BACTEC™ MGIT™ 960 SIRE Kit is used as a rapid qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*, from culture, to streptomycin, isoniazid, rifampin and ethambutol using the BD BACTEC MGIT 960 and BD BACTEC MGIT 320 Systems.

## SUMMARY AND EXPLANATION

Antimycobacterial susceptibility testing is necessary for the proper treatment of patients with tuberculosis. The treatment of tuberculosis is commonly through a multiple drug regimen which includes the antimycobacterial drugs streptomycin, isoniazid, rifampin and/or ethambutol. It is important that the antimycobacterial drugs prescribed show appropriate activity against *Mycobacterium tuberculosis*, i.e., susceptibility of the isolate to the drug.

Multidrug resistant *Mycobacterium tuberculosis* (MDR-TB) has recently become a serious public health problem.<sup>1</sup> Resistance to any of the four primary drugs, streptomycin (STR), isoniazid (INH), rifampin (RIF), and ethambutol (EMB), makes the disease more difficult and expensive to treat. The rapid detection of these strains is critical to the effective treatment of the patient.

Two methods have been widely used for antimycobacterial susceptibility testing. The first method, known as the Method of Proportion,<sup>2</sup> uses Middlebrook and Cohn 7H10 Agar. It compares colony counts on drug-containing and drug-free media. Resistance to a drug is detected when 1% or more of the bacterial population is resistant to the drug concentration under test. Results are generally available after 21 days of incubation. The second method, known as the BD BACTEC 460TB radiometric susceptibility method,<sup>3</sup> generally takes from 4 to 12 days. It is based on the production of radioactive <sup>14</sup>C-labeled carbon dioxide by the growing mycobacteria, manifested by a growth index increase in the system.

Historically, the Method of Proportion (MOP) procedure has included susceptibility testing of *Mycobacterium tuberculosis* using two concentrations of antimicrobials. The Clinical and Laboratory Standards Institute (CLSI) continues to recommend that the MOP test procedure include two concentrations of the primary drugs for testing except rifampin. The recommended low concentrations for the MOP procedure are generally considered to be the critical concentrations for these drugs. The critical concentration is defined as the drug concentration which inhibits the wild type sub-population while allowing sufficient growth of the mutant resistant sub-population in order to determine resistance at the critical proportion of 1%. The high drug concentration is not considered to be the critical concentration. However, resistance at the high concentration is evidence that resistance is generally distributed among the population of the test *M. tuberculosis* strain. Some clinicians use the susceptibility test results at the high concentrations to profile the degree of resistance of the test strain.

The BD BACTEC MGIT 960 SIRE test provides the susceptibility result in approximately the same time frame as the BD BACTEC 460TB system. In addition, this method is non-radiometric and allows appropriate antibiotic susceptibility results to be reported earlier, in most cases, than with the MOP procedure.

The BD BACTEC MGIT system has been developed to allow susceptibility testing at the critical concentration for streptomycin, isoniazid, rifampin and ethambutol, and at a higher concentration for streptomycin, isoniazid, and ethambutol. These concentrations correlate with the two concentrations used in the MOP procedure. A susceptible result at the critical concentration can be reported and no other tests are necessary. However, any strain determined to be resistant with the BD BACTEC MGIT 960 SIRE kit to streptomycin, isoniazid, and ethambutol at the critical concentration may be, at a minimum, tested at the high concentration. In this case, a final result of resistant at the critical concentration may be reported, with notification that an additional test at the high concentration is being performed.

## PRINCIPLES OF THE PROCEDURE

The BD BBL™ MGIT 7 mL Mycobacteria Growth Indicator Tube is a tube containing a modified Middlebrook 7H9 Broth which supports the growth and detection of mycobacteria (see BD BBL MGIT 7 mL package insert). The BD MGIT tube contains a fluorescent compound embedded in silicone on the bottom of a 16 x 100 mm round-bottom tube. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. The initial concentration of dissolved oxygen quenches the emission from the compound, and little fluorescence can be detected. Later, actively respiring microorganisms consume the oxygen which allows the compound to fluoresce.

The BD BACTEC MGIT 960 SIRE Kit is a 4–13 day qualitative test. The test is based on growth of the *Mycobacterium tuberculosis* isolate in a drug-containing tube compared to a drug-free tube (Growth Control). The BD BACTEC MGIT instrument monitors tubes for increased fluorescence. Analysis of fluorescence in the drug-containing tube compared to the fluorescence of the Growth Control tube is used by the instrument to determine susceptibility results.

The BD BACTEC MGIT instrument automatically interprets these results using predefined algorithms (which compare growth in the drug containing tube to that in the Growth Control tube) and the susceptible or resistant result is reported by the instrument.

## REAGENTS

BD BACTEC MGIT 960 SIRE Kit contains one each lyophilized vials of streptomycin, isoniazid, rifampin and ethambutol, and eight vials of SIRE Supplement.

Approximate Formula\* Per Vial Lyophilized drug: Streptomycin (STR)..... 332 µg

Approximate Formula\* Per Vial Lyophilized drug: Isoniazid (INH) ..... 33.2 µg

Approximate Formula\* Per Vial Lyophilized drug: Rifampin (RIF)..... 332 µg

Approximate Formula\* Per Vial Lyophilized drug: Ethambutol (EMB)..... 1,660 µg

BD BACTEC MGIT STR 4.0 Kit contains one vial lyophilized streptomycin and two vials of SIRE Supplement.

Approximate Formula\* Per Vial Lyophilized drug: Streptomycin ..... 664 µg

BD BACTEC MGIT INH 0.4 Kit contains one vial lyophilized isoniazid and two vials of SIRE Supplement.

Approximate Formula\* Per Vial Lyophilized drug: Isoniazid ..... 66.4 µg

BD BACTEC MGIT EMB 7.5 Kit contains one vial lyophilized ethambutol and two vials of SIRE Supplement.

Approximate Formula\* Per Vial Lyophilized drug: Ethambutol ..... 1,245 µg

BD BACTEC MGIT 960 SIRE Supplement contains 20 mL Middlebrook OADC enrichment

Approximate Formula\* Per L Purified Water

Bovine albumin ..... 50.0 g      Catalase ..... 0.03 g

Dextrose ..... 20.0 g      Oleic Acid ..... 0.6 g

*\*Adjusted and/or supplemented as required to meet performance criteria.*

**Storage and reconstitution of reagents:** BD BACTEC MGIT 960 SIRE Drug vials – on receipt, store the lyophilized drug vials at 2–8 °C. Once reconstituted, the antibiotic solutions may be frozen and stored at -20 °C or colder up to six months, not to exceed the original expiration date. Once thawed, use immediately. Discard unused portions.

BD BACTEC MGIT SIRE Supplement – On receipt, store in dark at 2–8 °C. Avoid freezing or overheating. Open and use prior to the expiration date. Minimize exposure to light.

### Directions For Use:

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Streptomycin lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Isoniazid lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 8.3 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Rifampin lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each BD BACTEC MGIT 960 SIRE Kit Ethambutol lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 415 µg/mL.

**NOTE: The following are reconstituted with a different volume. Failure to use the appropriate volume of sterile distilled water for reconstitution of the higher drug concentrations will invalidate those test results.**

Reconstitute each BD BACTEC MGIT 960 STR 4.0 Kit streptomycin lyophilized drug vial with **2 mL** of sterile distilled/deionized water to make a stock solution of 332 µg/mL.

Reconstitute each BD BACTEC MGIT 960 INH 0.4 Kit isoniazid lyophilized drug vial with **2 mL** of sterile distilled/deionized water to make a stock solution of 33.2 µg/mL.

Reconstitute each BD BACTEC MGIT 960 EMB 7.5 Kit ethambutol lyophilized drug vial with **2 mL** of sterile distilled/deionized water to make a stock solution of 622.5 µg/mL.

### WARNINGS:

For *in vitro* Diagnostic Use.

**POTENTIALLY INFECTIOUS TEST SPECIMEN: Observe “Universal Precautions”<sup>4</sup> and institutional guidelines when handling and disposing of infectious materials.**

BD BACTEC MGIT 960 – Catalog number 245127

BD BACTEC MGIT 960 Ethambutol

### Danger



**H360** May damage fertility or the unborn child.

**P201** Obtain special instructions before use. **P202** Do not handle until all safety precautions have been read and understood.

**P280** Wear protective gloves/protective clothing/eye protection/face protection. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Working with *M. tuberculosis* growth in culture requires Biosafety Level (BSL) 3 practices, containment equipment and facilities. Read and follow directions contained in all appropriate package inserts including the BD BACTEC MGIT 7 mL Mycobacteria Growth Indicator Tube.

Prior to use, the user should examine the tubes and vials for evidence of contamination or damage. Discard any tubes or vials if they appear unsuitable or if BD MGIT tubes exhibit fluorescence prior to use.

In the event of tube breakage: 1) Close the instrument drawers; 2) Turn off the instrument; 3) Vacate the area immediately; 4) Consult your facility/CDC guidelines. An inoculated leaking or broken tube may produce an aerosol of mycobacteria; appropriate handling should be observed.

Autoclave all inoculated BD MGIT tubes prior to disposal.

### SPECIMEN PREPARATION

All preparations detailed below are from cultures of *Mycobacterium tuberculosis*. The laboratory should confirm, by appropriate identification techniques, that the isolate to be tested is a pure culture.

#### Preparation of the Isolate from Solid Media:

1. Add 4 mL of BD BBL Middlebrook 7H9 Broth (or BD BACTEC MGIT broth) to a 16.5 x 128 mm sterile tube with cap containing 8–10 glass beads.
2. Scrape with a sterile loop as many colonies as possible from growth no more than 14 days old, trying not to remove any solid medium. Suspend the colonies in the Middlebrook 7H9 Broth. The suspension should exceed in turbidity a 1.0 McFarland standard.
3. Vortex the suspension for 2–3 min to break up the larger clumps.
4. Let the suspension sit for 20 min without disturbing.
5. Transfer the supernatant fluid to another 16.5 x 128 mm sterile tube with cap (avoid transferring any of the sediment) and let sit for another 15 min.
6. Transfer the supernatant fluid (it should be smooth, free of any clumps) to a third 16.5 x 128 mm sterile tube.
7. Adjust the suspension to a 0.5 McFarland standard by visual comparison to a 0.5 McFarland turbidity standard.
8. Dilute 1 mL of the adjusted suspension in 4 mL of sterile saline (1:5 dilution).

#### Preparation from a Positive BD BACTEC MGIT Tube:

1. For the preparation of the test inoculum, a positive 7 mL BD MGIT tube should be used the day **after** it first becomes positive on the BD BACTEC MGIT instrument (Day 1), up to and including the fifth day (Day 5) after instrument positivity. A tube which has been positive longer than five days should be subcultured to a fresh 7 mL BD MGIT tube containing BD BACTEC MGIT Growth Supplement and tested on the BD BACTEC MGIT instrument until positive, and used from one to five days following positivity.
2. If the tube is a Day 1 or Day 2 positive, proceed to "Inoculation Procedure for Susceptibility Test."
3. If the tube is a Day 3, Day 4, or Day 5 positive, then dilute 1 mL of the positive broth in 4 mL of sterile saline (1:5 dilution). Use the diluted suspension for the inoculation procedures. Proceed to "Inoculation Procedure for Susceptibility Test."

### PROCEDURE

**Materials Provided:** BD BACTEC MGIT 960 SIRE Kit containing one vial each lyophilized drug and eight vials of SIRE Supplement (approximately 40 tests per drug per kit). BD BACTEC MGIT 960 STR 4.0 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit), BD BACTEC MGIT 960 INH 0.4 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit), and BD BACTEC MGIT 960 EMB 7.5 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit).

**Materials Required But Not Provided:** BD BACTEC MGIT 7 mL Mycobacteria Growth Indicator Tubes, ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

#### IMPORTANT PROCEDURAL NOTE:

Be sure to tightly recap the BD MGIT tubes after all additions are made. Mix tubes thoroughly by gentle inversion three to four times. Thorough mixing of inoculated tubes is important. Failure to mix the tubes adequately can lead to false resistant results.

#### Inoculation Procedure for BD BACTEC MGIT 960 SIRE Kit Susceptibility Test:

1. Label five 7 mL BD MGIT tubes for each test isolate. Label one as GC (Growth Control), one as STR, one as INH, one as RIF, and the last as EMB. Place the tubes in the correct sequence in the appropriate size AST set carrier (see BD BACTEC MGIT Instrument User's Manual).
2. Aseptically add 0.8 mL of BD BACTEC MGIT SIRE Supplement to each tube. NOTE: It is important to use the correct supplement.
3. Aseptically pipet, using a micropipet, 100 µL of 83 µg/mL BD BACTEC MGIT STR solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 8.3 µg/mL BD BACTEC MGIT INH solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 83 µg/mL BD BACTEC MGIT RIF solution to the appropriately labeled BD MGIT tube. Aseptically pipet 100 µL of 415 µg/mL BD BACTEC MGIT EMB solution to the appropriately labeled BD MGIT tube. No drug solutions should be added to the BD MGIT GC tube.

| Drug               | Concentration of Drug after Reconstitution** | Volume Added to BD MGIT Tubes for Test | Final Concentration in BD MGIT Tubes |
|--------------------|----------------------------------------------|----------------------------------------|--------------------------------------|
| BD BACTEC MGIT STR | 83 µg/mL                                     | 100 µL                                 | 1.0* µg/mL                           |
| BD BACTEC MGIT INH | 8.3 µg/mL                                    | 100 µL                                 | 0.1* µg/mL                           |
| BD BACTEC MGIT RIF | 83 µg/mL                                     | 100 µL                                 | 1.0* µg/mL                           |
| BD BACTEC MGIT EMB | 415 µg/mL                                    | 100 µL                                 | 5.0* µg/mL                           |

\*Equivalent to CDC<sup>4</sup> recommended critical drug concentrations.

\*\* These drugs must be reconstituted using 4 mL sterile/deionized water to achieve concentrations indicated.

- Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see "Specimen Preparation") into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate **0.5** mL of the 1:100 Growth Control suspension into the BD MGIT tube labeled "GC."
- Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see "Specimen Preparation") into each of the FOUR remaining drug tubes (STR, INH, RIF, EMB).
- Tightly recap the tubes and mix well.
- Enter the AST set into the BD BACTEC MGIT 960 SIRE Kit using the AST set entry feature (refer to the BD BACTEC MGIT Instrument User's Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the Set Carrier definitions selected when performing the AST set entry feature.
- Streak 0.1 mL of the organism suspension to a BD Trypticase™ Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35–37 °C.
- Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the BD BACTEC MGIT Instrument User's Manual) and repeat testing with pure culture.

#### Inoculation Procedure for BD BACTEC MGIT 960 STR 4.0, INH 0.4 and EMB 7.5 Kits Susceptibility Test:

If resistance occurs at the critical concentration, an optional profile test may be performed which, at minimum, tests the high concentration of the drug to which it was originally resistant:

**Isolate Source** – The isolate used for this testing must have been prepared as described in "Specimen Preparation." A seed tube may be prepared from the drug free Growth Control tube from the previously tested AST set of the isolate, from which 0.5 mL may be inoculated to a fresh BD MGIT 7 mL tube containing BD BACTEC MGIT Growth Supplement. Once the seed tube is instrument positive, proceed as described in "Specimen Preparation: Preparation from a Positive BD MGIT tube."

- Label enough BD MGIT 7 mL tubes for the test isolate to have a BD MGIT GC (Growth Control) and a BD MGIT drug tube for each antimicrobial tested.
- Aseptically add 0.8 mL of BD BACTEC MGIT SIRE Supplement to each tube. NOTE: it is important to use the correct supplement.
- Aseptically pipet, using a micropipet, 100 µL of the drug solution to the appropriately labeled BD MGIT tube. No antibiotics should be added to the BD MGIT GC tube.

| Drug                   | Concentration of Drug after Reconstitution** | Volume Added to BD MGIT Tubes for Test | Final Concentration in BD MGIT Tubes |
|------------------------|----------------------------------------------|----------------------------------------|--------------------------------------|
| BD BACTEC MGIT STR 4.0 | 332 µg/mL                                    | 100 µL                                 | 4.0* µg/mL                           |
| BD BACTEC MGIT INH 0.4 | 33.2 µg/mL                                   | 100 µL                                 | 0.4* µg/mL                           |
| BD BACTEC MGIT EMB 7.5 | 622.5 µg/mL                                  | 100 µL                                 | 7.5* µg/mL                           |

\* Equivalent to MOP high drug concentrations.<sup>4</sup>

\*\* These drugs must be reconstituted using 2 mL sterile/deionized water to achieve concentrations indicated.

- Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see "Specimen Preparation") into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate **0.5** mL of the 1:100 Growth Control suspension into the BD MGIT tube labeled "GC."
- Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see "Specimen Preparation") into each of the drug tubes.
- Tightly recap the tubes and mix well.
- Enter the AST set into the BD BACTEC MGIT instrument using the AST set entry feature (refer to the BD BACTEC MGIT Instrument User's Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the Set Carrier definitions selected when performing the AST set entry feature.
- Streak 0.1 mL of the organism suspension to a BD Trypticase Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35–37 °C.
- Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the BD BACTEC MGIT Instrument User's Manual) and repeat testing with pure culture.

**NOTE** – The susceptibility test may be configured in a variety of formats. For example, a five tube carrier set containing only the critical concentrations can be configured into the system. A variety of other tube carrier sets can be configured depending on the optional profile tests being run (refer to BD BACTEC MGIT Instrument User's manual).

**User Quality Control** – Upon receipt of a new shipment or lot number of BD BACTEC MGIT 960 SIRE Kit vials, it is suggested that the control organism shown below be tested (see “Inoculation Procedure for Susceptibility Test”). Observation of the proper results, as shown below, within 4–13 days indicates that the BD BACTEC MGIT 960 SIRE Kit is ready for use in testing patient isolates. If the proper results are not observed, repeat the test. If, after repeating the test, the proper results are still not observed, do not use the product until you have contacted Technical Services at 1.800.638.8663 (United States only).

The same control organism should be run as batch QC once each week when susceptibility testing is performed. If the batch QC fails, do not report patient results for the drug (s) that failed for that testing period. Repeat the QC for the drug(s) and patient isolates affected by the initial QC failure. If the repeat QC does not perform as expected, do not report patient results. Do not use the product until you have contacted Technical Services at 1.800.638.8663 (United States only).

| Strains                            | GC       | BD BACTEC<br>MGIT STR | BD BACTEC<br>MGIT INH | BD BACTEC<br>MGIT RIF | BD BACTEC<br>MGIT EMB |
|------------------------------------|----------|-----------------------|-----------------------|-----------------------|-----------------------|
| <i>M. tuberculosis</i> ATCC® 27294 | Positive | Susceptible           | Susceptible           | Susceptible           | Susceptible           |

| Strains                           | GC       | BD BACTEC<br>MGIT STR 4.0 | BD BACTEC<br>MGIT INH 0.4 | BD BACTEC<br>MGIT EMB 7.5 |
|-----------------------------------|----------|---------------------------|---------------------------|---------------------------|
| <i>M. tuberculosis</i> ATCC 27294 | Positive | Susceptible               | Susceptible               | Susceptible               |

## REPORTING OF RESULTS

The BD BACTEC MGIT instrument will monitor susceptibility test sets until a susceptible or resistant determination is made. Once the set testing is completed, the results are reported by the BD BACTEC MGIT instrument (refer to the BD BACTEC MGIT Instrument User's Manual).

## LIMITATIONS OF THE PROCEDURE

Suspensions made from solid media must be allowed to settle for the prescribed times prior to standardization. Inoculum preparations made from solid media without the use of a 0.5 McFarland turbidity standard and the appropriate dilutions may give inaccurate results.

Use only pure cultures of *M. tuberculosis*. Cultures which are contaminated or which contain multiple strains of mycobacteria may give erroneous results.

Failure to reconstitute the drugs with the appropriate volume of sterile deionized water may give inaccurate results.

Failure to use a 1:100 dilution of the isolate for the inoculation of the Growth Control tube may give inaccurate results.

Failure to load the tubes of the AST set into the AST Set Carrier in the proper sequence may give inaccurate results.

Failure to use the SIRE Supplement in the AST set may give inaccurate results. Do not add **BD BACTEC MGIT** Growth Supplement to the AST set.

## EXPECTED RESULTS

The average time-to-result for susceptibility testing, using the critical concentrations, is 8.0 days (range 4 to 13).

## PERFORMANCE CHARACTERISTICS

The performance of the BACTEC MGIT 960 SIRE test was internally evaluated using a panel of challenge strains obtained from the Centers for Disease Control and Prevention (CDC), GA, USA. The panel consisted of thirty strains of *M. tuberculosis* with known susceptibility patterns (using agar proportion method). Table 1 shows the BACTEC MGIT 960 SIRE test results compared to the expected results.

Evaluation of the BACTEC MGIT 960 SIRE test is being performed at five geographically diverse clinical sites, to include regional reference centers, university hospital-based laboratories, and one internal site. The BACTEC MGIT 960 SIRE test is being compared to the modified Method of Proportion (MOP).<sup>6</sup> Table 2 shows the initial BACTEC MGIT 960 SIRE test results compared to the equivalent drug concentration in MOP. Table 3 shows the second tier BACTEC MGIT 960 SIRE test results compared to the equivalent drug concentration in MOP. This data includes testing of both fresh and/or reference clinical isolates.

Reproducibility of the BACTEC MGIT 960 SIRE test was evaluated at the clinical sites using a panel of ten qualified strains, including several strains resistant to each of the drugs, and comparing the BACTEC MGIT 960 SIRE test results to the expected result. The reproducibility results are: 99.6% for STR, 97.7% for INH, 99.4% for RIF and 94.0% for EMB. Individual site reproducibility results ranged from 94.5% to 99.5% for combined drug results.

Table 1

| CDC Challenge Panel |     | No. Tested | Category Agreement |       |
|---------------------|-----|------------|--------------------|-------|
| Drug                |     |            | No.                | %     |
| STR                 | All | 29         | 26                 | 90.0  |
|                     | S   | 24         | 21                 | 88.0  |
|                     | R   | 5          | 5                  | 100.0 |
| INH                 | All | 29         | 29                 | 100.0 |
|                     | S   | 8          | 8                  | 100.0 |
|                     | R   | 21         | 21                 | 100.0 |
| RIF                 | All | 29         | 29                 | 100.0 |
|                     | S   | 19         | 19                 | 100.0 |
|                     | R   | 10         | 10                 | 100.0 |
| EMB                 | All | 29         | 29                 | 100.0 |
|                     | S   | 25         | 25                 | 100.0 |
|                     | R   | 4          | 4                  | 100.0 |
| STR 4.0             | All | 8          | 8                  | 100.0 |
|                     | S   | 3          | 3                  | 100.0 |
|                     | R   | 5          | 5                  | 100.0 |
| INH 0.4             | All | 21         | 21                 | 100.0 |
|                     | S   | 8          | 8                  | 100.0 |
|                     | R   | 13         | 13                 | 100.0 |

Note: Results of MOP test for EMB 7.5 not available for comparison.

Table 2

| Clinical Strains |     | No. Tested | Category Agreement |       |
|------------------|-----|------------|--------------------|-------|
| Drug             |     |            | No.                | %     |
| STR              | All | 113        | 104                | 92.0  |
|                  | S   | 80         | 72                 | 90.0  |
|                  | R   | 33         | 32                 | 97.0  |
| INH              | All | 117        | 114                | 97.4  |
|                  | S   | 70         | 67                 | 95.7  |
|                  | R   | 47         | 47                 | 100.0 |
| RIF              | All | 116        | 115                | 99.1  |
|                  | S   | 82         | 82                 | 100.0 |
|                  | R   | 34         | 33                 | 97.1  |
| EMB              | All | 75         | 73                 | 97.0  |
|                  | S   | 67         | 65                 | 97.0  |
|                  | R   | 8          | 8                  | 100.0 |

Table 3

| Clinical Strains |     | No. Tested | Category Agreement |       |
|------------------|-----|------------|--------------------|-------|
| Drug             |     |            | No.                | %     |
| STR 4.0          | All | 31         | 22                 | 71.0  |
|                  | S   | 16         | 10                 | 62.5  |
|                  | R   | 15         | 12                 | 80.0  |
| INH 0.4          | All | 39         | 36                 | 92.3  |
|                  | S   | 6          | 5                  | 83.3  |
|                  | R   | 33         | 31                 | 93.9  |
| EMB 7.5          | All | 74         | 71                 | 96.0  |
|                  | S   | 66         | 66                 | 100.0 |
|                  | R   | 8          | 5                  | 63.0  |

## AVAILABILITY

### Cat. No. Description

245123 BD BACTEC™ MGIT™ 960 SIRE Kit, carton of 4 lyophilized drug vials, and 8 SIRE Supplements.  
 245125 BD BACTEC™ MGIT™ 960 STR 4.0 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.  
 245126 BD BACTEC™ MGIT™ 960 INH 0.4 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.  
 245127 BD BACTEC™ MGIT™ 960 EMB 7.5 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.

## REFERENCES

1. Barenfanger, J. 1993. Making your lab safe against multi-drug resistant *Mycobacterium tuberculosis*. Clin. Microbiol. News. 15: 76–80.
2. Clinical and Laboratory Standards Institute. 2003. Approved Standard: M24-A. Susceptibility testing of mycobacteria, nocardiae, and other aerobic actinomycetes. CLSI, Wayne, Pa.
3. BD Diagnostic Systems. BD BACTEC 460TB System Product and Procedure Manual.
4. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
5. Data on file at BD Diagnostic Systems.
6. Clinical and Laboratory Standards Institute. 1994. Tentative Standard: M24-T. Antimycobacterial Susceptibility Testing for *Mycobacterium tuberculosis*. CLSI, Wayne, PA.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or bd.com.

## Change History

| Revision | Date    | Change Summary                                                                                                                                                                                                                                                                                                                                    |
|----------|---------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (04)     | 2019-09 | <p>Converted printed instructions for use to electronic format and added access information to obtain the document from BD.com/e-labeling.</p> <p>Per Safety Data Sheet for catalog number 245127 added health hazard pictogram, signal word "Danger", all hazard &amp; precautionary codes and statements for BD BACTEC MGIT 960 Ethambutol.</p> |

US Customers only: For symbol glossary, refer to [bd.com/symbols-glossary](http://bd.com/symbols-glossary)



Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производитель / Výrobca / Proizvodač / Tillverkare / Üretici / Виробник / 生产厂商



Use by / Исползвайте до / Spotřebujte do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použít do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)  
 ГТГГ-ММ-ДД / ГТГГ-ММ (ММ = края на месеца)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)  
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)  
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)  
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)  
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)  
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)  
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)  
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)  
 AAAA-MM-GG / AAAA-MM (MM = fine mese)  
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы)  
 YYYY-MM-DD/YYYY-MM (MM = 월말)  
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)  
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)  
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)  
 AAAA-MM-DD / AAAA-MM (MM = slutten av måneden)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)  
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)  
 ГТГГ-ММ-ДД / ГТГГ-ММ (ММ = конец месяца)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)  
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)  
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)  
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)  
 YYYY-MM-DD / YYYY-MM (MM = 月末)



Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номер / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号



Authorized Representative in the European Community / Оторизан представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizuirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유럽 공동체의 위임 대표 / Įgaliotasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Reprezentante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Europskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Автура Топилулуҗу Yetkilil Temsilcisi / Уповноважений представник у країнах ЄС / 欧洲共同体授权代表



In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiniaparatuur / Dispositif médical de diagnostic in vitro / Medicinska romagala za In Vitro Dijagnostiku / In vitro diagnosztikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жүргізілетін медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietaisai / Medicinas ierices, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики ин витро / Medicinska romôcka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Dijagnostik Tibbi Cihaz / Медицинский прибор для диагностики ин витро / 体外诊断医疗设备



Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrensning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūros ierėbožojumi / Temperaturulimiet / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatura / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制



Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(로트) / Partijas numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)  
Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenuto sufficiente per <n> test / <n> тесттегі үшін жеткілікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder tilstrækkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достаточно для <n> тестов(а) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli malmeme içerir / Вистачить для аналізу: <n> / 足够进行 <n> 次检测



Consult Instructions for Use / Направете справка в инструктите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používání / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明



Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolosiți / Не использовать повторно / Nepoužívejte opakovane / Ne upotrebľajte ponovo / Får ej återanvändas / Tekrar kullannayin / Не використовувати повторно / 请勿重复使用



Serial number / Серийн номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / N° de serie / Seerianumber / Numéro de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Sērijas numurs / Serie nummer / Numer serijny / Numéro de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号  
For IVD Performance evaluation only / Само за оценка качеството на работа на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Μόνο για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réservé à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnozu / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасаңды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietaisų veikimo charakteristikoms tikrinti / Vienīgi IVD darbības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluering av IVD-ytelse / Tylko do oceny wydajności IVD / Uso exclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики in vitro / Ürçenē iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans deęerlendirilmesi için / Tільки для оцінювання якості діагностики in vitro / 仅限 IVD 性能评估



For US: "For Investigational Use Only"  
Lower limit of temperature / Долен лимит на температурата / Dolni hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температурный нижний порог / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurilimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限



Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / Контроль / 对照



Positive control / Положительный контроль / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitivá kontrola / Positiieve controle / Kontrola dodatna / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂



Negative control / Отрицательный контроль / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontrol / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativá kontrola / Negatieve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂



Method of sterilization: ethylene oxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: ethylenoxid / Steriliseringmethode: ethylenoxid / Sterilisationsmethode: ethylene oxide / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Sterilisierungsmethode: ethylenoxid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Sterilizacija azici – этилен тотыры / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksis / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringmethode: ethylenoxid / Metodă sterilizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metodă sterilizácie: etylénoxid / Metoda sterilizacije: etilen oksit / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: етиленоксид / 灭菌方法: 环氧乙烷  
Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmethode: bestraling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisierungsmethode: kiurguss / Méthode de stérilisation: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Sterilizacija azici – сәуле түсірі / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmethode: bestraling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metodă sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringsmethod: stråling / Sterilizasyon yöntemi: irradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射



Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Bioloogilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологичный тауекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологічна небезпека / 生物学风险



Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatus! Lugeda kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Uprozorenje, koristi prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, žiūrėkite priedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vystřa, pozri spřievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увага: див. супутню документацію / 小心, 请参阅附带文档。



Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температурный верхний порог / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurilimiet / Øvre temperaturgrænse / Gornja granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sicaklık üst sınırı / Максимальна температура / 温度上限



Keep dry / Пазете сухо / Skladujte v suchém prostředí / Orbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Százaz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausu / Uzglabāt sausu / Droog houden / Hoides tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávajte v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Берегти від вологи / 请保持干燥



Collection time / Време на събиране / Čas odběru / Opsamlingsstidspunkt / Entnahmehszeit / Ώρα συλλογής / Hora de recogida / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жынау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colecta / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingsstid / Toplama zamanı / Час забору / 采集时间  
Peel / Obenere / Otevřete zde / Åbn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстиңи қабатын алып таста / 벗기다 / Plești čia / Atīmēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odrhňte / Oljuštiti / Dra isär / Ayırma / Відклеїти / 撕下



Perforation / Перфорация / Perforace / Perforering / Διότρηση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Текік тесу / 절취선 / Perforacija / Perforăcia / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforăcia / Perforasyon / Перфорация / 穿孔



Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιό. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Егер пакет бұзылған болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuočių pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не използovať при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用



Keep away from heat / Пазете от топлина / Nevystavujte přilišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Óvja a melegtől / Tenere lontano dal calore / Салкын жерде сакта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не награвать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тепла / 请远离热源



Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Rezi / Vágja ki / Tagliare / Keciңiz / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciać / Cortar / Decupați / Отрезать / Odstrihnite / Iseći / Klipp / Kesme / Розрізати / 剪下



Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Pæmimo data / Savākšanas datums / Verzameldatum / Date prøvetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Datum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期 / μL/test / μL/тест / μL/Test / μL/εξέταση / μL/prueba / μL/teszt / μL/테스트 / мкл/тест / μL/tyrimas / μL/pārbaude / μL/teste / мкл/анализ / μL/검測



Keep away from light / Пазете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қазақыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Feriti de lumină / Хранить в темноте / Uchovávať mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Берегти від дії світла / 请远离光线 / Hydrogen gas generated / Образуван е водород газ / Možnost úniku plynného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadrží hydrogen vodík / Hidrogén gázt fejeleszt / Produzione di gas idrogeno / Газтөкрес сүтєгі пайда болды / 수소 가스 생성됨 / Išskiria vandenilio dujas / Rodas udegradis / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Oslobada se vodonik / Genererad vätgas / Açığa çıkan hidrojen gazı / Реакция з виділенням водню / 会产生氢气



Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттін идентификациялық нөмірі / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificatienummer van de patiënt / Pasientens ID-nummer / Numer ID pacienta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентификатор пациента / 患者标识号 / Fragile, Handle with Care / Чупливо, Работете с необходимото внимание. / Křehké. Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Εύθραστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Öm, käsitlege ettevaatlikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сығыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkités atsargiai. / Trausis; rikoties uzmanīgi / Breekbaar, voorzichtig behandelen. / Ømtålig, håndter forsigtig. / Vncha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşıyın. / Тендітна, звертатися з обережністю. / 易碎, 小心轻放



bd.com/e-labeling

KEY-CODE: 88-2041-1JAA

| Europe, CH, GB, NO: |                      | +800 135 79 135      |
|---------------------|----------------------|----------------------|
| International:      |                      | +31 20 794 7071      |
| AR                  | +800 135 79 135      | LT 8800 30728        |
| AU                  | +800 135 79 135      | MT +31 20 796 5693   |
| BR                  | 0800 591 1055        | NZ +800 135 79 135   |
| CA                  | +1 855 805 8539      | RO 0800 895 084      |
| CO                  | +800 135 79 135      | RU +800 135 79 135   |
| EE                  | 0800 0100567         | SG 800 101 3366      |
| GR                  | 00800 161 22015 7799 | SK 0800 606 287      |
| HR                  | 0800 804 804         | TR 00800 142 064 866 |
| IL                  | +800 135 79 135      | US +1 855 236 0910   |
| IS                  | 800 8996             | UY +800 135 79 135   |
| LI                  | +31 20 796 5692      | VN 122 80297         |



Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA



Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

#### Australian Sponsor:

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC® is a trademark of the American Type Culture Collection.

BD, the BD Logo, BACTEC, BBL, MGIT, and Trypticase are trademarks of Becton, Dickinson and Company or its affiliates.

All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved.