Mycobacteria Growth Indicator Tube 7 mL
(Epruveta od 7 mL sa indikatorom rasta mikobakterija)
Sa BD BACTEC MGIT 960 Supplement Kit (dodatni pribor)

L000180JAA(05)
2019-09
Srpski

NAMENA

BD BBL MGIT Mycobacteria Growth Indicator Tube (epruveta sa indikatorom rasta mikobakterija) pojačana supstancama BD BACTEC MGIT Growth Supplement (dodatak rastu) i BD BBL MGIT PANTA antibiotic mixture (antibiotska mešavina) namenjena je za otkrivanje i izolovanje mikobakterija pomoću sistema BD BACTEC MGIT 960 i BD BACTEC MGIT 320. Prihvatljivi tipovi uzorka su rastvoreni i dekontaminirani klinički uzorci (osim urina) i sterilne telesne tečnosti (osim krvi).

REZIME I OBJAŠNJENJE

Od 1985. do 1992. broj prijavljenih slučajeva infekcije vrstom *Mycobacterium tuberculosis* (MTB) porastao je za 18%. Tuberkuloza i dalje odnosi oko 3 miliona života godišnje u celom svetu, što je čini najsmrtonosnijom zaraznom bolešću.¹ Između 1981. i 1987. praćenje obolelih od SIDE ukazivalo je na to da je 5,5% pacijenata koji boluju od SIDE prethodno raširilo netuberkulozne mikobakterijske infekcije; npr. u granicama maksimalne dozvoljene koncentracije. Do 1990. porast raširenih netuberkuloznih mikobakterijskih infekcija doveo je do kumulativne učestalosti od 7,6%.² Pored novih žarišta infekcije vrstom MTB i MTB otporna na više lekova (multidrug resistant MTB, MDR-TB) postala je rastući problem. Kašnjenje laboratorija pri gajenju, identifikaciji i prijavljivanju ovih slučajeva infekcije vrstom MDR-TB donekle je doprinela širenju ove bolesti.³

Američki centri za kontrolu i prevenciju bolesti (U.S. Centers for Disease Control and Prevention, CDC), preporučili su da se učini sve kako bi laboratorije koristile najbrže metode na raspolaganju za dijagnostičko testiranje mikobakterija. Ove preporuke obuhvataju upotrebu kako tečnih, tako i čvrstih podloga za mikobakterijsku kulturu.^{3,4}

MGIT Mycobacteria Growth Indicator Tube (epruveta sa indikatorom rasta mikobakterija) sadrži 7 mL modifikovane baze Middlebrook 7H9 bujona.^{5,6} Kompletni medijum, zajedno sa OADC obogaćenjem i PANTA antibiotskom mešavinom, jedan je od najčešće korišćenih tečnih podloga za kultivaciju mikobakterija.

Svi tipovi kliničkih uzoraka, plućni kao i vanplućni (osim krvi i urina), mogu se obrađivati u cilju primarne izolacije u MGIT epruveti pomoću konvencionalnih metoda.⁴ Obradeni uzorak se inokuliše u MGIT epruveti, koja se stavlja u aparat BD BACTEC MGIT u cilju stalnog praćenja dok ne bude pozitivna ili dok se ne završi postupak testiranja.

PRINCIPI PROCEDURE

Fluorescentno jedinjenje je umetnuto u silikon na dnu epruveta sa okruglim dnom dimenzija 16 x 100 mm. Ovo fluorescentno jedinjenje je osetljivo na prisustvo kiseonika rastvorenog u bujonu. U početku, velika količina rastvorenog kiseonika gasi emisiju jedinjenja, pa se može otkriti malo fluorescencije. Kasnije, mikroorganizmi sa aktivnom respiracijom troše kiseonik i omogućuju da se otkrije fluorescencija.

Epruvete koje su unete u aparat BD BACTEC MGIT stalno se inkubiraju na 37 °C i svakih 60 minuta se pregleda povećanje fluorescencije. Analiza fluorescencije koristi se da bi se odredilo da li je epruveta pozitivna u instrumentu; tj. test uzorak sadrži odgovarajuće organizme. Epruveta koja je pozitivna u instrumentu sadrži oko 10⁵ do 10⁶ jedinica za formiranje kolonije po mililitru (Colony Forming Units, CFU/mL). Bočice sa kulturom koje ostaju negativne najmanje 42 dana (do 56 dana) i koje ne pokazuju vidljive znake pozitivnosti vade se iz instrumenta kao negativne i sterilišu pre bacanja.

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
REFERENCE

1. Bloom, B.R., and C.J.L. Murray. 1992. Tuberculosis: commentary on a reemergent killer. *Science* 257:1055–1064.
2. Horsburg, C.R., Jr., 1991. *Mycobacterium avium* complex infection in the acquired immunodeficiency syndrome. *N. Engl. J. Med.* 324:1332–1338.
3. Tenover, F.C., et al., 1993. The resurgence of tuberculosis: is your laboratory ready? *J. Clin. Microbiol.* 31:767–770.
4. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS, Centers for Disease Control, Atlanta.
5. Cohn, M.L., R.F. Waggoner and J.K. McClatchy. 1968. The 7H11 medium for the cultivation of mycobacteria. *Am. Rev. Respir. Dis.* 98:295–296.
6. Youmans, G.P. 1979. Cultivation of mycobacteria, the morphology and metabolism of mycobacteria, p. 25–35. *Tuberculosis*. W.B. Saunders Co., Philadelphia.
7. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
9. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.
11. Isenberg, Henry D. (ed.) 1992. *Clinical microbiology procedures handbook*. vol. 1. American Society for Microbiology, Washington, D.C.
12. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control for commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
13. Lindeboom, J. A., et al. 2011. Clinical Manifestations, Diagnosis, and Treatment of *Mycobacterium haemophilum* Infections. *Clinical Microbiology Reviews* 24, 701–717.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC® is a trademark of the American Type Culture Collection.

BD, the BD Logo, BACTEC, BBL, MGIT, MycoPrep, and PANTA are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved.