Part Number:	8361716	BALTSO0191 Version 12.0 Template 4 Inserts	
Category and Description:	Package Insert, Sputum Collection System	Rev from: 02	Rev to: 03
		Job Number:	311-17

Catalog Number: 290020

Blank (Sheet) Size: Length: 11" Width: 8.5"


Number of Pages: 8 Number of Sheets: N/A

Page Size: Length: 11" Width: 8.5" Final Folded Size: N/A

Ink Colors: Number of Colors: 1 PMS #: Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 1

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

OEM Specification Controlled by Vendor.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Sonia Thompson at 10:41 am, Feb 05, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 2:26 pm, Feb 05, 2018

Third Eye By: **THIRD EYE BY**
By Nichole Graham at 12:15 pm, Feb 15, 2018

BD Sputum Collection System

English: page 1 Italiano: page 4
Français : page 2 Español: page 5
Deutsch: page 3

R_x Only

8361716(03)
2018-02

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жерінікіті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The plastic BD Sputum Collection System is for the collection, transport and processing of sputum specimens.

SUMMARY AND EXPLANATION

The BD Sputum Collection System consists of an individually wrapped funnel with a hinged lid, protective amber body, a detachable graduated 50 mL conical tube and removable base. The system as a unit is packaged in a plastic bag and is processed to a sterility assurance level (SAL) of 10^{-3} . The screw cap on the 50 mL conical tube is specially designed to accomplish the sealing of the tube by an interference fit between the plug inside the cap and the inside of the tube. To screw the cap into the sealed position, rotate it clockwise until it is securely seated. The wall of the centrifuge tube is etched to provide an area for marking the identity of the patient from whom the sample was obtained.

PRINCIPLES OF THE PROCEDURE

The Sputum Collection System is designed to allow sputum to be collected and transported to the laboratory without personnel contact with contaminated parts. The patient raises the lid, delivers the specimen into the funnel and closes the lid. The sputum goes directly into a detachable conical centrifuge tube. The collection system should be covered with the unit bag after the patient has deposited the specimen. **Gloves should be worn when handling the system for protection from any contamination on the outside of the system occurring at the time the specimen is collected.** In the laboratory, the base is detached from the body of the collection system for removal of the cap and the tube containing the specimen.

Warnings and Precautions

For *in vitro* Diagnostic Use.

Warning: The system will not withstand autoclaving prior to use.

Use only if package is intact.

SPECIMEN COLLECTION AND TRANSPORT

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving. Whenever possible, the laboratory should provide suitable collection containers and instructions for their proper use.

Specimens should be obtained before antibiotics or other antimicrobial agents have been administered. If therapy was initiated prior to collection of the specimen, this should be noted on the forms accompanying the specimen. Material should be collected from the site where the suspected organism is most likely to be found with as little external contamination as possible, and in a clinical situation at the stage of disease most likely to yield a positive culture. Specimens should be of sufficient quantity to permit complete examination. To prevent overgrowth by contaminants and to ensure that pathogens are viable upon receipt, provision must be made for prompt delivery to the laboratory. Specimen containers should be properly labeled with all pertinent information such as name of patient, time specimen was taken, etc. The laboratory must be furnished sufficient clinical information to enable the microbiologist to select the most suitable media and appropriate techniques.

If at all possible, specimens should be transported direct to the laboratory.

Use a Double Mailing Container for shipping these tubes. Current postal and other transportation regulatory agencies' requirements for labeling and shipping of hazardous or etiologic specimens must be observed.

For more detailed information, appropriate references should be consulted.^{5,6}

PROCEDURE

Materials Provided: BD Sputum Collection System.

Materials Required But Not Provided: Sealable transport bag, solutions, culture media, reagents, quality control organisms and laboratory equipment as required for the procedures and tests to be performed on the specimens.

Instructions: Use protective barriers to prevent exposure to blood, body fluids containing visible blood, and other fluids to which universal precautions apply. The type of protective barrier(s) should be appropriate for the procedure being performed and the type of exposure anticipated.

Immediately and thoroughly wash hands and other skin surfaces that are contaminated with blood, body fluids containing visible blood, or other body fluids to which these universal precautions apply.

Handle all biologic samples and materials contaminated with blood in accordance with the policies and procedures of your facility.

Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since the samples may transmit HBV (hepatitis), HIV (AIDS), or other infectious diseases.⁷

Sample Collection

Instructions are printed on the unit bag.

Nursing personnel should open the unit bag by cutting along dotted line with scissors. The patient should be given the opened bag containing the system. It should be held upright.

The patient should remove the bag, open the lid, and expectorate into the container. The lid should then be closed and the bag replaced over the top of the system. This process should be repeated until approximately 20 mL of sputum is in the container. **NOTE: Collect no more than 20 mL of sputum into the container.**

Sample Processing

Instructions are printed on the unit bag.

1. Wearing gloves, remove the system from the patient's bedside by grasping the unit bag and system. Ensure the lid has been firmly closed. Label the system with the appropriate patient information and place the system and the unit bag into a sealable bag for transportation to the laboratory. The system should be kept in an upright position until the tube containing the specimen is removed.
2. A laboratory staff member will receive the specimen. Wearing gloves, hold the system in an upright position and remove the base for access to the tube and cap.
3. Set the base aside for later disposal.
4. Holding the unit with transport bag in place over the system, pull the tube down to separate it from the funnel. Turn the funnel assembly on its side; use the U-shaped access area to put the cap on the tube.
5. Tighten the cap slightly by turning the tube clockwise, then remove the tube and tighten the cap completely.
6. Label the tube with laboratory identification or with specific patient information writing on the frosted area.
7. Replace the base on the bottom of the funnel assembly and place in the transport bag.
8. Place the transport bag with contents into a biohazard container for disposal.

Directions for Processing the Tubes: In the laboratory, all procedures involving hazardous material, such as sputum, should be performed in a suitable biological containment cabinet while utilizing proper personal protective equipment.

1. Prior to handling the centrifuge tube, check for proper positioning and sealing of the cap.
2. When removing the cap, always turn slowly to avoid a sudden release of pressure. Before removing the cap, tap it gently several times to cause excess beads of moisture to fall inside the tube rather than on the threads. Allow two sec drain time.
3. Before repositioning the cap, remove any liquid droplets inside the cap or on the top and outside thread area of the tube by using methods and procedures of universal precautions.
4. Before performing any procedures, retighten the cap on the tube, particularly if the tube is subject to variations in temperature.
5. Use the proper centrifuge equipment and recommended RCF. Make sure there is clearance between the cap of the tube and the shield of the centrifuge, so there is no contact during centrifugation.
6. Extreme CAUTION should be used in disposing of contaminated material. Dispose of the sputum system by autoclaving or incinerating. **NOTE:** Protect autoclave from molten material.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of Sputum Collection Systems are tested to verify specific product characteristics. Certificate of Processing and Dosimetry Records for the sterilization of Sputum Collection Systems are reviewed for acceptable dose (kGy).

An extensive series of toxicity tests as specified in *The United States Pharmacopoeia (USP)* were conducted to assure that the 50 mL conical tube does not contain concentrations of bioactive substances which could potentially interfere with biological experiments.⁸ All tests verified that there was no biological interference from the material used in the manufacturing of this device.

AVAILABILITY

Cat. No.	Description
290020	BD™ Sputum Collection System, Case of 72 units.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.

3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
5. Isenberg, H.D. (ed.) 1992. *Clinical Microbiology Procedures Handbook*. American Society for Microbiology, Washington, D.C.
6. Federal Register, Vol 45, No. 141, Monday, July 21, 1980, 42 CFR part 72 - Interstate Shipment of Etiologic Agents.
7. OSHA Final Standard for Occupational Exposure to Bloodborne Pathogens. 56 Fed. Reg. 64175, Dec. 6, 1991; 29 CFR Part 1910.1030.
8. United States Pharmacopoeial Convention, Inc. 1994. *The United States pharmacopoeia 23/The national formulary 18-1995*. United States Pharmacopoeial Convention, Inc., Rockville, Md.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Made in Mexico.

© 2018 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.