POSTUPCI KONTROLE KVALITETE

I UVOD

Medij se koristi za uzgoj dermatofita i ostalih patogenih i nepatogenih gljivica iz kliničkih i nekliničkih uzoraka.

II POSTUPAK ISPITIVANJA UČINKOVITOSTI

- Inokulirajte reprezentativne uzorke s dolje navedenim kulturama.
 - Za izravni inokulat sojevima *Trichophyton mentagrophytes* ATCC 9533 i *Aspergillus brasiliensis* ATCC 16404 pomoću ušice od 0,01 mL kulture razvijene u glukoznom agaru **BBL Sabouraud Dextrose Agar**.
 - Za inokulat sojevima *Candida albicans* ATCC 10231 i *Escherichia coli* ATCC 25922 pomoću 0,01 mL razrijeđene fiziološke otopine za 10^3 do 10^4 CFU-a.
- Inkubirajte epruvete s otpuštenim poklopcima pri temperaturi od 25 do 30 °C do 7 dana u aerobnoj atmosferi.
- Očekivani rezultati

Organizmi	ATCC	Izolacija
* <i>Candida albicans</i>	10231	Prosječan do velik rast
* <i>Trichophyton mentagrophytes</i>	9533	Prosječan do velik rast
* <i>Aspergillus brasiliensis</i>	16404	Inhibicija (djelomična do potpuna)
* <i>Escherichia coli</i>	25922	Inhibicija (djelomična do potpuna)

*Preporučeni soj organizama za korisničku kontrolu kvalitete.

III DODATNA KONTROLA KVALITETE

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XIV REFERENCE

1. Sabouraud, R. 1910. Les teignes, p. 553. Masson et Cie, Paris.
2. Ajello, L., L.K. Georg, W. Kaplan, and L. Kaufman. 1963. CDC Laboratory manual for medical mycology. PHS Publication No. 994, U.S. Government Printing Office, Washington, D.C.
3. Weitzman, I., J. Kane, and R.C. Summerbell. 1995. *Trichophyton*, *Microsporum*, *Epidermophyton*, and agents of superficial mycoses, p. 791-808. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
4. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol 17:53-80.
6. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from the risks related to exposure to the biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
8. Haley, L.D., H. Trandel, and M.B. Coyle. 1980. Cumitech 11, Practical methods for culture and identification of fungi in the clinical mycology laboratory. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.
9. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
10. Merz, W.G., and G.D. Roberts. 1995. Detection and recovery of fungi from clinical specimens, p. 709-722. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
11. Kwon-Chung, K.J., and J.E. Bennett. 1992. Medical mycology. Lea & Febiger, Philadelphia.
12. Larone, D.H. 1995. Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
13. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1992. Mycology, p. 791-877. Color atlas and textbook of diagnostic microbiology, 4th ed., J.B. Lippincott Co., Philadelphia.

Tehnički servis i podrška na BD Diagnostics: obratite se lokalnom predstavniku tvrtke BD ili posjetite www.bd.com/ds.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2015 BD.