






## Revisions

BALTSO0191 Version 8.0 Template 4  
Inserts

| Rev from | Rev to | JOB #   |
|----------|--------|---------|
| 01       | 02     | 7043-14 |

## NOTES:

1. BD Catalog Number: 261201
2. Blank (Sheet) Size: Length: 11" Width: 22.5"
3. Number of Pages: 10 Number of Sheets: 1
4. Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
5. Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
6. Printed two sides: Yes ☒ No ☐
7. Style (see illustrations below): # 4

8. See specification control no. N/A for material information.
9. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

|              |                                                         |                                                                                                                                                            |  |                                                                                                                                                                          |
|--------------|---------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Label Design | REVISED BY<br>By Nancy Carlsen at 3:15 pm, Apr 30, 2015 | COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION. |  |  <b>BD</b> Becton, Dickinson and Company<br>7 Loveton Circle<br>Sparks, MD 21152 USA |
| Proofer      | <i>[Signature]</i> 05-04-15                             |                                                                                                                                                            |  |                                                                                                                                                                          |
| Checked By   | <i>[Signature]</i> 05-04-15                             | Category and Description                                                                                                                                   |  | Sheet: 1 of 11                                                                                                                                                           |
| Part Number: | L001184                                                 | Package Insert,<br>Ninhydrin Reagent Droppers                                                                                                              |  | Scale: N/A                                                                                                                                                               |
|              |                                                         |                                                                                                                                                            |  | <b>A</b>                                                                                                                                                                 |

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## INTENDED USE

**BD BBL™** Ninhydrin Reagent Droppers are used in the determination of the hippurate reaction to aid in the identification of certain bacteria.

## SUMMARY AND EXPLANATION

The hydrolysis of sodium hippurate is a qualitative test useful in the identification of *Streptococcus agalactiae* (group B  $\beta$ -hemolytic streptococci), *Campylobacter jejuni*, *Gardnerella vaginalis*, *Listeria* sp. and other aerobic bacteria.<sup>1-4</sup> The presence of the enzyme hippurate hydrolase is visualized by the purple color formed when ninhydrin reagent oxidizes the amino acid produced during hippurate hydrolysis.

## PRINCIPLES OF THE PROCEDURE

The end products of hydrolysis of sodium hippurate include glycine and benzoic acid. Glycine is deaminated by the ninhydrin reagent, which becomes reduced during the process. The end products of the ninhydrin oxidation react to form a purple colored dye.<sup>5</sup>

## REAGENTS

Ninhydrin Reagent Droppers contain 0.5 mL of 3.5 % ninhydrin in 50 % aqueous ethanol.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

### Danger

**H225** Highly flammable liquid and vapour. **H302** Harmful if swallowed. **H350** May cause cancer. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P233** Keep container tightly closed. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P370+P378** In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray. **P403+P235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Storage Instructions:** Store at room temperature 15 – 30 °C (59 – 86 °F). Each dropper is good for one day's use after breaking ampule.

**Product Deterioration:** This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. Do not use if reagent is other than yellow.

## PROCEDURE

**Material Provided:** BD BBL Ninhydrin Reagent Droppers.

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

### Test Procedure

To use Ninhydrin Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom of the dropper on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent.

**Sodium Hippurate Broth Preparation:** Make 1% sodium hippurate solution in distilled water. Dispense into 12 x 75 mm plastic tubes with caps in 0.4 mL aliquots. Store frozen at -20 °C for up to 6 months.<sup>5</sup>

1. Thaw the needed number of tubes of 1% sodium hippurate solution at room temperature. Sodium hippurate may be made fresh and used immediately.
2. Inoculate solution with suspect colonies. Make a heavy suspension.
3. Cover tubes and incubate tubes in a 35 ± 2 °C waterbath for 2 h.
4. Add 5 drops of ninhydrin reagent to each tube.
5. Cover tubes and return to the waterbath for 10 min. Do not exceed 30 min as false reactions may occur.
6. Read tubes immediately for purple coloration.

### User Quality Control

Examine the reagent for signs of deterioration (see "Product Deterioration").

A reagent control should be performed daily or with each batch of tests and when new 1% sodium hippurate solution is made. The reagent control is tested using the procedure above except that no organism is added. There should be no purple coloration after the 10 min incubation period.

When tested according to the recommended procedure, *Streptococcus agalactiae* ATCC™ 27956 should give a positive reaction and *Enterococcus faecalis* ATCC 19433 should give a negative reaction.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control

procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

In positive reactions a purple color will form within the 10 min incubation period. Negative reactions are colorless or pale yellow.

Consult appropriate texts for expected results.<sup>1,3-5</sup>

LIMITATIONS OF THE PROCEDURE

Use only pure cultures. Mixed cultures may produce misleading results.

Do not read reactions after 30 min.<sup>2</sup>

Hippurate hydrolysis should be used in conjunction with other biochemical tests to ensure accurate identification.

PERFORMANCE CHARACTERISTICS

An internal study compared **BD BBL™** Ninhydrin Reagent and another commercially available Ninhydrin Reagent using ATCC and stock organisms.<sup>6</sup> As noted in the table, all reactions were identical for 25 strains tested using both reagents.

| Organism                                   | Expected Results | BD BBL Reagent | Commercial Reagent |
|--------------------------------------------|------------------|----------------|--------------------|
| <i>Streptococcus agalactiae</i> ATCC 12386 | +                | +              | +                  |
| <i>Streptococcus pyogenes</i> ATCC 19615   | -                | -              | -                  |
| <i>Campylobacter coli</i> 802              | -                | -              | -                  |
| <i>Campylobacter coli</i> 834              | -                | -              | -                  |
| <i>Campylobacter jejuni</i> 1118           | +                | +              | +                  |
| <i>Campylobacter jejuni</i> ATCC 29428     | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> ATCC 14028    | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> JH3           | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD12          | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> BD16          | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD11          | +                | W+             | W+                 |
| <i>Gardnerella vaginalis</i> BD10          | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> BD8           | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD7           | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD6           | +                | W+             | W+                 |
| <i>Gardnerella vaginalis</i> BD3           | +                | W+             | W+                 |
| <i>Gardnerella vaginalis</i> BD2           | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> BD1           | +                | W+             | +                  |
| <i>Gardnerella vaginalis</i> BD5           | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> JH4           | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> JH7           | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> BD9           | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD14          | -                | -              | -                  |
| <i>Gardnerella vaginalis</i> BD15          | +                | +              | +                  |
| <i>Gardnerella vaginalis</i> AMMS17        | +                | +              | +                  |
| Control ( no organism)                     | -                | -              | -                  |

| Key |               |
|-----|---------------|
| -   | Negative      |
| +   | Positive      |
| W+  | Weak Positive |

AVAILABILITY

**Cat. No.**      **Description**  
261201      **BD BBL™** Ninhydrin Reagent Droppers, 50.

REFERENCES

1. Balows, A., W.J. Hausler, Jr., K.L. Hermann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
2. MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
3. Baron, S. 1986. Medical microbiology, 2nd ed. Addison-Wesley, Menlo Park, Calif.
4. Washington, J.A. (ed.). 1985. Laboratory procedures in clinical microbiology, 2nd ed. Springer-Verlag, New York.
5. Baron, E.J., and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Co., St. Louis.
6. Data on file, BD Diagnostic Systems.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

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Becton Dickinson Pty Ltd.
Benex Limited