
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	8548-16

NOTES:

- BD Catalog Number: 261197, 261198, 261207
- Blank (Sheet) Size: Length: 11" Width: 31.5"
- Number of Pages: 14 Number of Sheets: 1
- Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue, 032 Red, Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. D11835RM for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Mary Schmidt at 2:59 pm, Jun 01, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 3:44 pm, Jun 01, 2016		
Checked By	THIRD EYE BY By Kim Stone at 8:12 am, Jun 21, 2016		
Part Number:	L001190	Category and Description Package Insert, Nitrate A, B & C Reagent Droppers	Sheet: 1 of 15 Scale: N/A
		A	

BD BBL™ Nitrate A, Nitrate B and Nitrate C Reagent Droppers

English: pages 1 – 3 Italiano: pagine 8 – 10
Français: pages 3 – 5 Español: páginas 10 – 12
Deutsch: Seiten 5 – 8

 L001190(03)
2016-06

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyks lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жерінікі BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL™ Nitrate Reagent A, BD BBL Nitrate Reagent B and BD BBL Nitrate Reagent C Droppers are for use in determining the ability of bacteria to reduce nitrate to nitrite or free nitrogen gas.

SUMMARY AND EXPLANATION

The *Enterobacteriaceae*, many other gram-negative bacteria, mycobacteria and fungi reduce nitrate to nitrite. The nitrate-reducing characteristic of a species is constant for certain genera and species. The Nitrate A and B Reagents, when added in equal parts, indicate the presence of a catabolic end product or the absence of nitrate in the medium.^{1,2} In the event of a negative reaction, the addition of Nitrate C Reagent shows the presence of nitrate in the medium and confirms the negative reaction.

PRINCIPLES OF THE PROCEDURE

The reaction indicates the ability of the organism to reduce nitrate (NO₃) to nitrite (NO₂), nitrous oxide (N₂O), nitric oxide (NO) or ammonia (NH₃). In the presence of any of these end products, the Nitrate A and B Reagents, added in equal proportions, will turn the medium pink or red if the organism is nitrate positive. The formation of other end products that do not turn the medium pink or red is also possible; therefore, the addition of Nitrate C Reagent (zinc dust) is used to detect unreduced nitrate. In the presence of unreduced nitrate, zinc dust will turn the medium pink or red. The absence of color after the addition of zinc dust indicates that nitrate has been reduced to nitrogen gas (N₂) and the organism is nitrate positive.

REAGENTS

Nitrate A Reagent Droppers contain 0.5 mL of 0.8% Sulfanilic Acid in 5N Acetic Acid. Nitrate B Reagent Droppers contain 0.5 mL of 0.6% N,N-Dimethyl- α -naphthylamine in 5N Acetic Acid. Nitrate C Reagent Droppers contain 1 g of zinc dust.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Nitrate B Reagent, although not classified as carcinogenic by the U.S. Occupational Safety and Health Administration (OSHA), is closely related to α -naphthylamine, which is considered a carcinogen.³ Nitrate C Reagent in contact with water liberates extremely flammable gases. Never use water.

Follow proper laboratory procedures in handling and disposing of infectious materials.

Cat. No.	Description
261197	BD BBL™ Nitrate A Reagent Droppers, 50.

Danger

H314 Causes severe skin burns and eye damage. **H317** May cause an allergic skin reaction. **P101** If medical advice is needed, have product container or label at hand. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Cat. No.	Description
261198	BD BBL™ Nitrate B Reagent Droppers, 50.

Danger

H314 Causes severe skin burns and eye damage. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P301+P330+P331** IF SWALLOWED: rinse mouth. Do NOT induce vomiting. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Cat. No.	Description
261207	BD BBL™ Nitrate C Reagent Droppers, 50.

Danger

H250 Catches fire spontaneously if exposed to air. **H260** In contact with water releases flammable gases which may ignite spontaneously. **H400** Very toxic to aquatic life. **H410** Very toxic to aquatic life with long lasting effects. **P222** Do not allow contact with air. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P370+P378** In case of fire: Use for extinction: CO₂, sand, extinguishing powder. **P402** Store contents under inert gas. **P402+P404** Store in a dry place. Store in a closed container. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store Reagent Droppers at room temperature 15 – 30 °C (59 – 86 °F). Protect from light. Each dropper is good for one day's use after breaking ampule.

Product Deterioration: These reagents are hermetically sealed in ampules, which afford protection from chemical instability until the expiration date. A change in color of reagents A and B from clear or light yellow to brown indicates improper storage, which may cause weaker reactions.

PROCEDURE

Materials Provided: Nitrate A Reagent Droppers, Nitrate B Reagent Droppers or Nitrate C Reagent Droppers, as ordered.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedures:

To use Nitrate A or B Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom of dropper on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent. To use Nitrate C Reagent Droppers, remove the cap and invert to dispense the zinc dust.

Liquid Medium Method: Method for *Enterobacteriaceae* and nonfermenting gram-negative rods.

1. Inoculate Nitrate Broth with Durham tube with the organism to be tested using an 18 – 24 h culture.
2. Incubate inoculated tubes at 35 – 37 °C for 18 – 24 h. Occasionally, prolonged incubation is needed to achieve a positive reaction for some organisms.¹
3. Remove from incubator and examine Durham tube for presence of gas. If gas is present, the organism is nitrate positive. If no gas is present, continue with step 4.
4. Add 10 drops of Nitrate A and 10 drops of Nitrate B to tube. (Alternative procedure: Remove approximately 2 mL of broth and place in clean 12 x 75 mm test tube. Add 5 drops of Nitrate A and 5 drops of Nitrate B to tube. Shake gently to mix reagents. If negative, remaining broth may be reincubated and tested the next day.)
5. Observe broth for pink/red color development within 2 min.
6. If pink/red color develops, the organism is nitrate positive.
7. If no color develops within 2 min, add a small amount of Nitrate C (one inversion) to the tube.
8. If the medium turns pink/red within 2 min, nitrate is present (i.e., not reduced by the organism) and the organism is nitrate negative.
9. If there is no color change noted within 2 min, nitrate is not present in the broth (i.e., the organism reduced the nitrate to nitrogen gas) and the organism is nitrate positive.

Semisolid Medium Method: Method for fastidious gram-negative rods, anaerobes and gram-positive organisms.

1. Inoculate Indole Nitrite Medium with the organism to be tested using an 18 – 24 h culture.
2. Incubate inoculated tubes at 35 – 37 °C for 18 – 24 h. Occasionally, prolonged incubation is needed to achieve a positive reaction for some organisms.¹
3. Add 10 drops of Nitrate A and 10 drops of Nitrate B to tube. (Alternative procedure: Remove approximately 2 – 3 mL of medium and place in clean 12 x 75 mm test tube. Add 5 drops of Nitrate A and 5 drops of Nitrate B to tube. Shake gently to mix reagents. If negative, remaining broth may be reincubated and tested the next day.)
4. Observe broth for pink/red color development within 2 min.
5. If pink/red color develops, the organism is nitrate positive.
6. If no color develops within 2 min, add a small amount of Nitrate C (one inversion) to the tube.
7. If the medium turns pink/red within 2 min, nitrate is present (i.e., not reduced by the organism) and the organism is nitrate negative.
8. If there is no color change within 2 min, nitrate is not present in the medium (i.e., the organism reduced the nitrate to nitrogen gas) and the organism is nitrate positive.

Other Methods: These reagents can be used in any variation of the nitrate reduction test that requires Nitrate A, Nitrate B and Nitrate C.

User Quality Control: Examine the reagents for signs of deterioration (see "Product Deterioration"). *Escherichia coli* ATCC™ 25922 will yield a positive result, while *Acinetobacter baumannii* ATCC 17961 will yield a negative result.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Nitrate Positive:

- a. Gas in Durham tube *or*,
- b. After addition of reagents, pink/red color develops *or*,
- c. After addition of reagents no color develops and after addition of zinc dust no color develops.

Nitrate Negative: After addition of reagents no color develops *and* after addition of zinc dust pink/red color develops.

LIMITATIONS OF THE PROCEDURE

Allow at least 2 min for the color to develop before considering the nitrate test negative.

Nitrate C (zinc dust) must be added to any negative reaction to eliminate the possibility of a false negative reaction due to the organism's ability to reduce the nitrate to nitrogen gas.

Use only pure cultures. Mixed cultures may produce misleading results.
Use only the recommended amount of zinc dust. Excess zinc dust may cause false reactions.⁴

PERFORMANCE CHARACTERISTICS

An internal study compared **BD BBL™** Nitrate A, B and C Reagents to another set of commercially available Nitrate Reagents.⁵ As noted in the table, all reactions were identified for the 7 ATCC strains tested using both sets of Nitrate Reagents.

Organism	Expected Results	BD BBL™ Reagents	Commercial Reagents
<i>Acinetobacter anitratus</i> ATCC™ 7961	–	–	–
<i>Escherichia coli</i> ATCC 25922	+	+	+
<i>Pasteurella multocida</i> ATCC 11096	+	W+	W+
<i>Branhamella cattharalis</i> ATCC 38098	+	W+	W+
<i>Eikenella corrodens</i> ATCC 7622	+	+	+
<i>Branhamella cattharalis</i> ATCC 25238	+	W+	W+
<i>Enterobacter cloacae</i> ATCC 23355	+	+	+

Key	
–	Negative
+	Positive
W+	Weak Positive

AVAILABILITY

Cat. No.	Description
261197	BD BBL™ Nitrate A Reagent Droppers, 50.
261198	BD BBL™ Nitrate B Reagent Droppers, 50.
261207	BD BBL™ Nitrate C Reagent Droppers, 50.

REFERENCES

1. MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
2. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey and Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St Louis.
3. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.) 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
4. Porres, J.M., and V. Porter. 1974. Rapid nitrate reduction test. Am. J. Med. Technol. 40: 257-259.
5. Data on file, BD Diagnostic Systems.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

 BD BBL Nitrate A, Nitrate B and Nitrate C Reagent Droppers

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Benex Limited
Becton Dickinson Pty Ltd.