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REFERENSER

1. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
2. Baron, E.J., and S.M. Finegold. 1990. Bailey and Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.
3. Engelkirk, P.G., J. Duben-Engelkirk, and V.R. Dowell, Jr. (ed.). 1992. Principles and practice of clinical anaerobic bacteriology. Star Publishing Company, Belmont, Calif.
4. Holdeman, L.V., E.P. Cato and W.E.C. Moore. 1977. Anaerobe laboratory manual, 4th edition. Virginia Polytechnic Institute and State University, Blacksburg.
5. Holdeman, L.V., E.P. Cato and W.E.C. Moore. 1987. Anaerobe laboratory manual update. Supplement to the 4th edition. Virginia Polytechnic Institute and State University, Blacksburg.
6. Holdeman, L.V., E.P. Cato and W.E.C. Moore. 1993. Anaerobe laboratory manual update. Supplement to the 4th edition. Virginia Polytechnic Institute and State University, Blacksburg.
7. Mandell, G.L., R.G. Douglas, Jr. and J.E. Bennett. 1990. Principles and practice of infectious diseases, 3rd ed. Churchill Livingstone Inc., New York.
8. Rodloff, A.C., P.C. Appelbaum, and R.J. Zabransky. 1991. Cumitech 5A, Practical anaerobic bacteriology, Coordinating ed., A.C. Rodloff. American Society for Microbiology, Washington, D.C.
9. Summanen, P., E.J. Barron, D.M. Citron, C.A. Strong; H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Company, Belmont, Calif.
10. Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. Am. J. Public Health. 8:922-923.
11. Cowan, S.T., and K.J. Steel. 1974. Manual for the identification of medical bacteria. 2nd ed. Cambridge University Press, Cambridge.
12. Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
13. Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
14. Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Public Health. Trop. Med. :96-100.
15. Edberg, S.C., and C.M. Konnick. 1986. Comparison of α -glucuronidase-based substrate systems for identification of *Escherichia coli*. J. Clin. Microbiol. 24:368-371.
16. Kämpfer, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877-2879.
17. Kilian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae* 1: detection of bacterial glycosidases. Acta Pathol. Microbiol. Scand. Sect. B. 84:245-251.
18. Maddocks, J.L., and M. Greenan. 1975. Rapid method for identifying bacterial enzymes. J. Clin. Pathol. 28:686-687.
19. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
20. Mangels, J., I. Edvalson, and M. Cox. 1993. Rapid Identification of *Bacteroides fragilis* group organisms with the use of 4-methylumbelliferone derivative substrates. Clin. Infect. Dis. 16(54):5319-5321.
21. Moncla, B.J., P. Braham, L.K. Rabe, and S. L. Hiller. 1991. Rapid presumptive identification of black-pigmented gram-negative anaerobic bacteria by using 4-methylumbelliferone derivatives. J. Clin. Microbiol. 29:1955-1958.
22. Qadri, S.M., and S. Johnson. 1981. Rapid test for esculin hydrolysis by anaerobic bacteria. Antonie van Leeuwenhoek 47:371-379.
23. Sneath, P.H.A. 1957. The application of computers to taxonomy. J. Gen. Microbiol. 17:201-221.
24. Hansen, S.L., and B.J. Stewart. 1978. Slide catalase. A reliable test for differentiation and presumptive identification of certain clinically significant anaerobes. Am. J. Clin. Microbiol. 13:444-448.
25. Data on file at BD Diagnostics.

BD Diagnostics teknisk service: Kontakta närmaste BD-representant.

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Atқарушы / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvođač / Tilverkare / Üretici / Виробник

Use by / Исполняйте до / Spotføjbt do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / Upotrebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použít do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до ліне
YYYY-MM-DD / YYYY-MM (MM = end of month)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
AAAA-MM-DD / AAAA-MM (MM = fin del mes)
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)
AAAA-MM-GG / AAAA-MM (MM = fine mese)
ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)
MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
GGGG-MM-DD / GGGG-MM (MM = mēneša beigas)
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
AAAA-MM-DD / AAAA-MM (MM = sluten av måneden)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = fim do mês)
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
PPPP-MM-DD / PPPP-MM (MM = кінець місяця)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номери / Katalogo numeris / Kataloga numurs / Catalogus number / Numer katalogowy / Număr de catalog / Номер по каталог / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталог

Authorized Representative in the European Community / Оторизиран представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatus esindaja Euroopa Nõukogu / Représentant autorisé pour la Communauté européenne / Autorizirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közöségben / Rappresentante autorizzata nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / Įgaliotasis atstovas Europos Bendrijos / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupce v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Аврпа Topluluğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro diagnostickeho prietlaku / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiniaparatu / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnostikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жүргізілетін медициналық диагностика аспабы / In vitro diagnostikos prietaisais / Medicināas ierīces, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro / Medicinska pomočka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Dijagnostik Tibbi Cihaz / Медицинский прибор для диагностики ин витро

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dovoljena temperatura / Hömersékleti határ / Limiti di temperatura / Temperatuuri shtreke / Laikymo temperatūra / Temperatūras ierobežojumi / Temperaturlimit / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури

Batch Code (Lot) / Код на партида / Kód (číslo) šarže / Batch-code (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kod / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serie / Partnummer (Lot) / Parti Kodu (Lot) / Код партії

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koistuti upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / Skatytite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanni Talimatları'na başvurun / Див. інструкції за використання

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road Laoighaire
Co. Dublin, Ireland

Tween is a trademark of ICI Americas, Inc.

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, BBL and BBL Crystal are trademarks of Becton, Dickinson and Company. ©2014 BD