
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	8368-16

NOTES:

- BD Catalog Number: 231401
- Blank (Sheet) Size: Length: 8.5" Width: 25.5"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 4.25" Final Folded Size: 2 1/8" x 4 1/4"
- Ink Colors: No. of Colors: 1 PMS#: Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. BALTL0001610 for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 5:06 pm, May 25, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Terrence Means at 5:34 pm, May 25, 2016		
Checked By	THIRD EYE BY By Helaine Jones at 1:54 pm, May 26, 2016		
Part Number:	L0001610JAA	Category and Description Package Insert, BBL Gram Slide	Sheet: 1 of 13 Scale: N/A A

For evaluating Gram stain reagents and staining techniques

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INTENDED USE

The **BD BBL™** Gram Slide is used for evaluating and controlling the quality of Gram stain reagents and techniques.

SUMMARY AND EXPLANATION

Quality control procedures are performed to help ensure that the information reported by laboratories is accurate, reliable and reproducible. Reagents and personnel are monitored at established intervals to document the validity of the test method.¹⁻⁸ The **BD BBL** Gram Slide offers a standardized, pretested stable control for use when testing Gram stain reagents. The prepared slide eliminates the need for maintaining stock cultures to prepare slides.

PRINCIPLES OF THE PROCEDURE

Prepared slides containing known quality control cultures (*Staphylococcus aureus* ATCC® 25923 and *Escherichia coli* ATCC 25922) are used to test Gram stain reagents and staining technique. As a result, variables that could cause incorrect test results may be detected and corrected.


A crystal violet-iodine complex forms in the protoplast (not the cell wall) of all organisms stained by this procedure. Organisms able to retain this dye complex after decolorization are classified as gram-positive while those that can be decolorized and counterstained are classified as gram-negative.

Upon disruption or removal of the cell wall, the protoplast of gram-positive (as well as gram-negative) cells can be decolorized and the gram-positive attribute lost. Thus, the mechanism of the Gram stain appears to be related to the presence of an intact cell wall able to act as a barrier to decolorization of the primary stain.

Generally, the cell wall is nonselectively permeable. It is theorized that during the Gram stain procedure, the cell wall of gram-positive cells is dehydrated by the alcohol in the decolorizer and loses permeability, thus retaining the primary stain. However, the cell wall of the gram-negative cells has a higher lipid content and becomes more permeable when treated with alcohol, resulting in loss of the primary stain.

REAGENTS

The **BD BBL** Gram Slide is a conventional 1" x 3" microscope slide imprinted with 10 squares. One square contains unstained control organisms. The remaining nine squares are available for staining test isolates. The control square, labeled (C±), consists of a mixture of gram-positive cocci and gram-negative bacilli, *Staphylococcus aureus* ATCC 25923 and *Escherichia coli* ATCC 25922, respectively. The **BD BBL** Gram Slide must be heat-fixed prior to staining.

WARNINGS AND PRECAUTIONS:

For *in vitro* Diagnostic Use.

The control organisms have been chemically inactivated and air-dried on the slide. However, all slides should be handled as if they have infectious material. Follow proper established laboratory procedure in handling and disposing of infectious materials.

Due to possible carryover of organisms during the staining procedure, this product is not recommended for use with clinical specimens.

Storage: Store the **BD BBL Gram Slide** below 30 °C. Do not expose the slides to extreme temperatures.

The expiration date applies to the product in its intact container when stored as directed.

Stained slides are stable indefinitely and may be kept as a permanent record.

Product Deterioration: Do not use a product if it fails to meet performance specifications for identity and Gram reaction.

PROCEDURES

Material Provided: **BD BBL Gram Slide**.

Materials Required But Not Provided: Gram stain reagents, slide warmer or Bunsen burner, staining rack, forceps and conventional microscope with oil-immersion lens.

TEST PROCEDURE

Heat-fix the **BD BBL Gram Slide** by passing the slide two to three times through the flame of a Bunsen burner. Alternatively, hold the slide in front of a micro-incinerator for 5–10 sec. Do not overheat.

1. Stain the **BD BBL Gram Slide** along with test slides, using Gram stain reagents and following your laboratory's recommended procedures.
2. Keep slides well separated during the staining procedure to avoid cross-contaminating carry-over of stain reagents from slide to slide.
3. Read the stained slide microscopically under an oil-immersion lens and record results.

User Quality Control: Stain slides by the Gram stain technique and examine microscopically. The "C±" square should contain a mixture of gram-positive cocci and gram-negative rods.

INTERPRETATION OF RESULTS

Read Gram stained slides microscopically under an oil-immersion lens. Record appearance of organisms observed (i.e., morphology and color).

LIMITATIONS OF THE PROCEDURE

Abnormal staining or partial loss of the smear on the control section of the slide may result from overheating during fixing of the test smear, improper decolorization, an overly forceful washing procedure or deterioration of the staining reagents used.¹

Antimicrobial agents may cause increased susceptibility of the test isolate to decolorization in the Gram stain procedure.¹

As with any procedure that involves the addition of multiple isolates to a single slide, the possibility exists that some organisms may loosen and float on the slide during the staining process. Isolates exhibiting questionable or unexpected patterns of staining should be reexamined using a single smear per slide.¹

The Gram stain reaction is altered by physical disruption of the bacterial cell wall or protoplast. The cell walls of gram-positive bacteria interpose a barrier which prevents leaching of the dye complex from the cytoplasm. Cell walls of gram-negative bacteria contain lipids soluble in organic solvents, which are then free to decolorize the cytoplasm. Therefore, a microorganism that is physically disrupted by excess heating will not react to Gram staining as expected.

"Careful adherence to procedure and interpretive criteria is required for accurate results. Accuracy is highly dependent on the training and skill of the microbiologist."⁹

Gram stain results, including organism morphology, can be affected by the age of the isolate, bacteria containing autolytic enzyme systems, cultures transferred from antibiotic-containing media, as well as specimens collected from patients on antibiotics.¹⁰ "Background material and artifacts can also interfere with interpretation. Precipitated gram-positive stain generally appears as irregular coccoid shapes or asters resembling fungal hyphae."¹⁰

PERFORMANCE CHARACTERISTICS

Gram-positive organisms will appear blue to purple. Gram-negative organisms will appear pink to red.

AVAILABILITY

Cat. No.	Description
231401	BD BBL™ Gram Slide , 50.

REFERENCES

1. Sewell, D.L. 1994. Laboratory records, p. 13.2.1-13.2.35. *In* H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 2. American Society for Microbiology, Washington, D.C.
2. August, M.J., J.A. Hindler, T.W. Huber, and D.L. Sewell. 1990. Cumitech 3A, Quality control and quality assurance practices in clinical microbiology. Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
3. Commission on Laboratory Accreditation. 1990. Inspection checklist. Diagnostic immunology and syphilis serology. College of American Pathologists, Northfield, Ill.
4. Miller, J.M. 1987. Quality control in microbiology. Centers for Disease Control, Atlanta.
5. Miller, J.M., and B.B. Wentworth (ed.). 1985. Methods for quality control in diagnostic microbiology. American Public Health Association, Washington, D.C.

- National Committee for Clinical Laboratory Standards. 1996. Approved guideline: GP2-A3. Clinical laboratory technical procedure manuals, 3rd ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
- Weissfeld, A.S., and R.C. Bartlett. 1987. Quality control, p. 35-65. In B.J. Howard, J. Klass II, S.J. Rubin, A.S. Weissfeld, and R.C. Tilton (ed.), Clinical and pathogenic microbiology. The C.V. Mosby Co., St. Louis.
- Health Care Financing Administration. 1988. Medicare, Medicaid, and CLIA programs; revision of the clinical laboratory regulations for Medicare, Medicaid, and Clinical Laboratories Improvement Act of 1967 programs. Fed. Regist. 53: 29590-29632.
- Kruczak-Filipov, P., and R.G. Shively. 1992. Gram stain procedure, p. 1.5.1-1.5.18. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol.1. American Society for Microbiology, Washington, D.C.
- Chapin, K. 1995. Clinical microscopy, p. 33-51, In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed), Manual of clinical microbiology, 6th ed., American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

© 2016 BD. BD, the BD Logo and BD BBL are trademarks of Becton, Dickinson and Company.