
Revisions

BALTSO0191 Version 11.0 Template 4
Inserts

Rev from	Rev to	JOB #
03	04	9060-17

NOTES:

- BD Catalog Number: 231391
- Blank (Sheet) Size: Length: 17.00" Width: 12.75"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 4.25" Final Folded Size: 2.125" L x 4.25" W
- Ink Colors: No. of Colors: 1 PMS#: Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 7

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. BALTL0001620 for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 3:24 pm, Apr 14, 2017	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Natalie Morio at 7:42 am, Apr 17, 2017		
Checked By	THIRD EYE BY By Terrence Means at 4:38 pm, Apr 17, 2017		
Part Number:	L0001620JAA	Category and Description Package Insert, BD BBL™ AFB Slide	Sheet: 1 of 13 Scale: N/A A

BD BBL™ AFB Slide

For evaluating the quality of acid-fast and fluorescent stain reagents and staining techniques

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INTENDED USE

The **BD BBL™ AFB Slide** is used for evaluating and controlling the quality of acid-fast and fluorescent stain reagents and staining techniques for the detection of acid-fast bacilli (AFB).

SUMMARY AND EXPLANATION

Quality control procedures are performed to help ensure that the information reported by laboratories is accurate, reliable and reproducible. Reagents and methods are required to be monitored at established intervals to document the validity of the test.¹⁻⁸ The **BD BBL AFB Slide** offers a standardized, stable preparation for use as control when testing acid-fast and fluorescent stain reagents. The prepared slide eliminates the need for maintaining stock cultures to prepare slides.

PRINCIPLES OF THE PROCEDURE

Prepared slides contain an area negative for acid-fast bacilli with a mixture of inactivated non-acid-fast cocci (*Staphylococcus aureus* ATCC® 25923) and bacilli (*Klebsiella pneumoniae* ATCC 13883). A positive area contains these organisms plus inactivated *Mycobacterium tuberculosis* H37Ra. The slide is used to check the quality of staining reagents and techniques.

Staining applications used for evaluating AFB include both acid-fast and fluorescent stains. Examples of these stains include Ziehl-Neelsen, Kinyoun and Auramine-Rhodamine. In the Ziehl-Neelsen stain the primary stain is basic carbol fuchsin and methylene blue is the counterstain.⁹ The Kinyoun stain differs from the Ziehl-Neelsen stain in that phenol (instead of steam) is used as the mordant to allow penetration of the primary stain into the organism cell walls.⁹ Auramine and rhodamine are nonspecific fluorochromes that bind to mycolic acids in the mycobacterial cell wall.⁹ With the Auramine-Rhodamine stain, potassium permanganate serves as a counterstain which suppresses nonspecific background fluorescence.

REAGENTS

The **BD BBL AFB Slide** is a conventional 1" x 3" microscope slide imprinted with two circles. The positive circle contains a mixture of unstained *Mycobacterium tuberculosis*, H37Ra, *S. aureus* ATCC 25923 and *K. pneumoniae* ATCC 13883. This circle serves as the positive control. The negative circle contains an unstained mixture of *S. aureus* ATCC 25923 and *K. pneumoniae* ATCC 13883. This circle serves as the negative control. The **BD BBL AFB Slide** must be heat-fixed prior to staining.

Warnings and Precautions:

Contains mercury  

For *in vitro* Diagnostic Use.

The control organisms have been chemically inactivated and air-dried on the slide. However, all slides should be handled as if they have infectious material. Follow proper established laboratory procedure in handling and disposing of infectious materials.

Storage: Store the **BD BBL AFB Slide** below 30 °C. Do not expose the slides to extreme temperatures. The expiration date applies to the product in its intact container when stored as directed.

Stained slides are stable indefinitely and may be kept as a permanent record.

Product Deterioration: Do not use a product if it fails to meet performance specifications for identity and stain reaction.

PROCEDURES

Material Provided: BD BBL AFB Slide.

Materials Required But Not Provided: Acid-fast stain reagents and/or fluorescent stain reagents, slide warmer or Bunsen burner, staining rack, forceps, conventional microscope with oil-immersion lens and/or fluorescence microscope assembly.

Test Procedure:

Reagent Preparation

Heat-fix the **BD BBL AFB Slide** by passing slide two to three times through the flame of a Bunsen burner. Alternatively, hold the slide in front of a micro-incinerator for 5–10 sec. Do not overheat.

Test Procedure

1. Stain the **BD BBL AFB Slide** along with test slides, using a fluorescent or acid-fast stain kit according to your laboratory's recommended procedures.
2. Keep slides well separated during the staining procedure to avoid cross-contaminating carry-over of stain reagents from slide to slide.
3. Read the stained slide under an oil-immersion lens with a microscope appropriate for the type of stain used and record results.

User Quality Control: After staining with an appropriate AFB stain, the positive circle should contain clearly visible acid-fast rods mixed with non-acid-fast cells and debris. The negative circle should contain non-acid-fast cells and debris with not more than an occasional acid-fast rod.

RESULTS

Read acid-fast-stained slides microscopically under an oil-immersion lens. Observe slides stained with one of the fluorescent stains using a fluorescence microscope. Record organism morphology and color observed for the appropriate staining procedure used.⁴

LIMITATIONS OF THE PROCEDURE

Abnormal staining or partial loss of the smear on the slide may result from improper decolorization, an overly forceful washing procedure or deterioration of the staining reagents used.

Improper heat-fixing of the smear on the slide may result in removal of cells from the slide.

PERFORMANCE CHARACTERISTICS

Kinyoun (cold acid-fast)^{10,11}: Acid-fast bacilli: dark red to pink. Non-acid-fast bacilli: green. **Ziehl-Neelsen (hot acid-fast)¹²:** Acid-fast bacilli: dark red to pink. Non-acid-fast bacilli: blue. **Auramine-rhodamine⁹:** Acid-fast bacilli: fluorescent orange-yellow with UV light.

AVAILABILITY

Cat. No.	Description
231391	BD BBL™ AFB Slide, 50.

REFERENCES

1. Sewell, D.L. 1994. Laboratory records, p. 13.2.1-13.2.35. *In* H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
2. August, M.J., J.A. Hindler, T.W. Huber, and D.L. Sewell. 1990. Cumitech 3A, Quality control and quality assurance practices in clinical microbiology. Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
3. Commission on Laboratory Accreditation. 1990. Inspection checklist. Diagnostic immunology and syphilis serology. College of American Pathologists, Northfield, Ill.
4. Miller, J.M. 1987. Quality control in microbiology. Centers for Disease Control, Atlanta.
5. Miller, J.M., and B.B. Wentworth (ed.). 1985. Methods for quality control in diagnostic microbiology. American Public Health Association, Washington, D.C.
6. National Committee for Clinical Laboratory Standards. 1996. Approved guideline: GP2-A3. Clinical laboratory technical procedure manuals, 3rd ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
7. Weissfeld, A.S., and R.C. Bartlett. 1987. Quality control, p. 35-65. *In* B.J. Howard, J. Klass II, S.J. Rubin, A.S. Weissfeld, and R.C. Tilton (ed.), Clinical and pathogenic microbiology. The C.V. Mosby Co., St. Louis.
8. Health Care Financing Administration. 1988. Medicare, Medicaid, and CLIA programs; revision of the clinical laboratory regulations for Medicare, Medicaid, and Clinical Laboratories Improvement Act of 1967 programs. Fed. Regist. 53: 29590-29632.
9. Chapin, K. 1995 Clinical microscopy, pp. 40-46. *In* Murray, P. R., Baron, E. J., Pfaller, M.A., Tenover, F.C. and R. H. Tenover (ed). Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
10. Vet. Adm. Arm. Forces Coop. Study on the Chemo. of Tuber. Nov. 1962. VA Dept. of Med, & Sur., Wash., D.C.
11. Strong, B.E., and G.P. Kubica. Isolation and identification of *Mycobacterium tuberculosis*, a guide for the Level II Laboratory, U.S. Public Health Service Publication. 81-8390, Center for Disease Control, Atlanta.
12. Lab. Methods Clin, and Publ. Hlth. Mycobact., USPHS. 1968

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuriipiir / Limite inférieure de température / Najnižza dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite mínimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контрольно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrole / Controle / Controllo / Контроль / kontroll / Контроль / 对照

Positive control / Положительный контроль / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozityvā kontrole / Positive controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂

Negative control / Отрицательный контроль / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negatiivne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontrol / Controllo negativo / 음성 컨트롤 / Neigiama kontrolė / Negativā kontrole / Negative controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

© 2017 BD. BD, the BD Logo and BD BBL are trademarks of Becton, Dickinson and Company.