
Revisions

SO 0191-5

Rev from	Rev to	ECO #
0604	2010/07	5391-10

Notes:

- BD Cat. Number 231109, 231108, 231107, 231106, 231105, 231104
- Blank (Sheet) Size : Length: 11" Width: 25.5"
 Number of Pages: 12 Number of Sheets: 1
 Page Size: Length 11" Width 4.25" Final Folded Size: 4.25" x 2.25"
- Style (see illustrations below): # 7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8810091JAA		Category and Description Package Insert, BBL Strips for Differentiation of Haemophilus Species	Sheet: 1 of 13	A
			Scale: 1:1	

BD BBL™ Strips for Differentiation of *Haemophilus* Species

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INTENDED USE – Taxo™ Strips impregnated with X, V and X and V factors are used in a qualitative procedure for the isolation and differentiation of *Haemophilus* species. The growth of these organisms is partially dependent upon the presence of the X or V factor or both.¹ These strips provide a simple method for determining the growth requirements of organisms in this group.

Taxo X and V Factor Strips may be substituted for the *Staphylococcus* streak used in some laboratories for isolation of *Haemophilus* when the base medium is **Trypticase™** Soy Agar with or without sheep blood.

SUMMARY AND EXPLANATION – The growth requirements (“factors”) of *Haemophilus* are used to speciate the genus.²

Haemophilus influenzae, which may be a serious primary or secondary pathogen when isolated from patients with respiratory infections, does not grow well, if at all, on blood agar containing sheep blood. This organism requires the X-factor or heat-stable substance of hemoglobin (hemin) and the V-factor or heat labile substance (nicotinamide adenine dinucleotide, NAD).³ Both of these components are contained within intact red blood cells and are not readily available to *Haemophilus* on unheated sheep blood agar medium, although some X factor diffuses into the medium.⁴ Colonies may become discernible if they happen to grow near colonies of “feeder” organisms that produce large amounts of V factor during growth; e.g., staphylococci, *Neisseria* and *Pseudomonas*.⁴

In addition, sheep blood is inhibitory to *Haemophilus*, due partly to an anti-V factor enzyme (NADase) present in the blood.⁵

Rabbit and horse blood agar plates support growth of *Haemophilus* species. However, the “nonpathogenic” *H. haemolyticus* also grows and may be confused with β-hemolytic streptococci, although a Gram stain readily differentiates them.⁶

For isolation from throat cultures, a **Sensi-Disc™** bacitracin 10 unit disc may be placed on the area of heaviest inoculation, usually where the inoculation is initiated. After incubation, colonies of *Haemophilus*, e.g., *H. influenzae*, may be present in the zone around the disc.

PRINCIPLES OF THE PROCEDURE – Taxo X and V Factor Strips provide X and V factors, permitting growth of *Haemophilus* species from specimens and differentiation of isolates in pure culture based on their X and V factor requirements.

In differentiating species, growth only in the area of the strip indicates dependence on the factor(s) provided.

REAGENTS

Taxo™ X Factor Strips are impregnated with hemin.

Taxo™ V Factor Strips are impregnated with nicotinamide adenine dinucleotide (NAD).

Taxo™ X and V Factor Strips are impregnated with both hemin and NAD.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. “Standard Precautions”⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

Directions should be read and followed carefully.

Storage Instructions: On receipt, store at -20 to 8°C. After opening, store at 2 to 8°C to protect product integrity.

Product Deterioration: Do not use strips if they show evidence of discoloration or other signs of deterioration.

SPECIMENS – These strips may be used with either clinical specimens or pure cultures of organisms suspected to be *Haemophilus* species. Consult appropriate texts for more information.^{4,11}

PROCEDURE

Material Provided: Taxo V, X and / or VX Factor Strips.

Materials Required But Not Provided: Prepared plated **Trypticase** Soy Agar; **Trypticase** Soy Broth; inoculating loop; swabs; sterile forceps; incubator, 35°C, with source of supplemental CO₂.


Test Procedure:

A. Isolation from Primary Culture


1. Streak a culture plate with the specimen, e.g., throat culture, in a conventional manner.
2. With sterile forceps, aseptically place a **Taxo** XV Factor Strip on a less-heavily inoculated area.
3. Incubate for 24 h at 35°C in a 5 to 10% carbon dioxide atmosphere.
4. After incubation, “satellite” growth of *Haemophilus* organisms will be observed around the XV strip.
5. Pick well-isolated colonies and transfer to a second plate for differentiation (see procedure B). A pure culture must be used for determination of growth factor requirements. Contaminated or mixed cultures may yield erroneous results. If the inoculum is obtained from a medium containing any of the essential growth factors, dilute the inoculum in 5 mL sterile **Trypticase** Soy Broth prior to inoculating the plate. This procedure will help avoid carry-over of the X factor from the primary to the differentiation plate.⁴

- B. Differentiation from Pure Cultures
- Using a loop or swab, inoculate a **Trypticase** Soy Agar plate over its entire surface with the organism.
 - With sterile forceps, aseptically place **Taxo** X, V and XV strips on the inoculated agar so that there is about a 20 mm distance between strips, as shown.

Methods for Applying Strips

ISOLATION

DIFFERENTIATION

- Incubate at 35°C for 18 to 24 h in a 5 to 10% carbon dioxide atmosphere.
- Observe plates for growth around the strips.

User Quality Control: Examine strips as described under “Product Deterioration.” Check performance by testing a representative sample of strips with pure cultures of stable control organisms that give known, desired reactions. The following test strains are recommended:

Test Strain	Growth Along Strips		
	XV	V	X
<i>H. influenzae</i> , ATCC™ 10211	+	-*	-
<i>H. parainfluenzae</i> , ATCC 9796	+	+	-

+ indicates growth - indicates no growth

* Very light, diffuse growth (possibly covering the entire plate) may be observed. This phenomenon may be attributed to trace amounts of hemin in the medium (**Trypticase** Soy Agar) or to carry over of hemin from the primary (chocolate) isolation plate. Growth between the X and V strips should be much greater than any background growth.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory’s standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS¹¹ – On primary cultures, observe plates for the presence of “satellite” colonies around the XV strip.

In differentiating species, growth only in the area of the strip indicates dependence on the factor(s) provided. Some typical reactions are tabulated below:

Species	Growth Along Strips		
	XV	V	X
<i>H. influenzae</i>	+	-	-
<i>H. aegyptius</i>	+	-	-
<i>H. parainfluenzae</i>	+	+	-
<i>H. parahaemolyticus</i>	+	+	-
<i>H. haemolyticus</i>	+	-	-
<i>H. aphrophilus</i> *	+	-/+	+
<i>H. paraphrophilus</i>	+	+	-

+ indicates growth - indicates no growth

* X factor requirement may be observed on primary isolation but is often lost on subculture.¹¹ If X factor is required, growth will not occur around the V factor strip.

LIMITATIONS OF THE PROCEDURE – Tests using these strips are presumptive and should be followed with serological or biochemical identification tests, as appropriate.^{4,11,12}

PERFORMANCE CHARACTERISTICS – Taubitz et al. studied methods for identification and serotyping of encapsulated strains of *Haemophilus influenzae*. Two hundred fifty-two (252) strains of *H. influenzae* and 200 strains of *H. parainfluenzae* were tested. Growth factor requirements were determined using V, X and XV factor strips. Opalescent growth, characteristic of encapsulated strains of *H. influenzae*, could be observed at the same time that growth factor requirements were determined.¹³ In another study by Parker and Hoeprich, 255 nasopharyngeal swabs and 1,317 throat swabs were examined using X and V factor disks. Of the nasopharyngeal swabs 62 were positive for *H. influenzae* and 12 were positive for *H. parainfluenzae*. Of the 1,317 throat swabs, 129 were positive for *H. influenzae* and 302 were positive for *H. parainfluenzae*.¹

AVAILABILITY

Cat. No.	Description
231109	Taxo™ V Factor Strips, Pkg. of 6 Vials.
231108	Taxo™ V Factor Strips, 1 vial of 50 strips.
231107	Taxo™ X Factor Strips, Pkg. of 6 Vials.
231106	Taxo™ X Factor Strips, 1 vial of 50 strips.
231105	Taxo™ XV Factor Strips, Pkg. of 6 Vials.
231104	Taxo™ XV Factor Strips, 1 vial of 50 strips.

REFERENCES

- Parker, R.H., and P.D. Hoeprich. 1962. Disk method for rapid identification of *Haemophilus* species. Am. J. Clin. Pathol. 37:319-327.
- Kilian, M., and E.L. Biberstein. 1984. Genus II. *Haemophilus*, p. 558-569. In N.R. Krieg and J.G. Holt (ed.), Bergey’s manual of systematic bacteriology, vol. 1. Williams & Wilkins, Baltimore.
- Holt, L.B. 1962. The growth-factor requirements of *Haemophilus influenzae*. J. Gen. Microbiol. 27:317-322.
- Howard, B.J. 1987. *Haemophilus*, p. 279-288. In B.J. Howard (ed.), Clinical and pathogenic microbiology. The C.V. Mosby Company, St. Louis.
- Waterworth, P.M. 1955. The stimulation and inhibition of the growth of *Haemophilus influenzae* in media containing blood. Br. J. Exp. Pathol. 36:186-194.
- Krumwiede, E., and A.G. Kuttner. 1938. A growth inhibitory substance for the influenzae group of organisms in the blood of various animal species. J. Exp. Med. 67:429-441.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

11. Campos, J.M. 1995. *Haemophilus*, p. 556-565. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
12. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1997. *Haemophilus*, p. 363-394. Color atlas and textbook of diagnostic microbiology, 5th ed. Lippincott-Raven Publishers, Philadelphia.
13. Taubitz, I.S., and H. Brandis, 1988. A comparison between methods of identification and serotyping of encapsulated strains of *Haemophilus influenzae*. Zbl. Bakt. Hyg. A 270, 83-97.

BD BBL Strips for Differentiation of *Haemophilus* Species

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© 2010 BD.
ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL, Sensi-Disc, Taxo and Trypticase are trademarks of Becton, Dickinson and Company.
Benex Limited