Part Number:	0216446JAA	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, BD Directigen™ Neisseria Meningitidis Test	Rev from: 03	Rev to: 04
		Job Number:	730-18

Catalog Number: 250160

Blank (Sheet) Size: Length: 8.5" Width: 11"


Number of Pages: 24 Number of Sheets: 6

Page Size: Length: 8.5" Width: 5.5" Final Folded Size: N/A

Ink Colors: Number of Colors: 3 PMS #: 2755 Blue; 032 Red; Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 5

Vendor Printed: ☒

Online / In House Printed: ☐

Web Printed: ☐

See Specification control no. BALT0216446 for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Tori Pagani at 10:17 am, Dec 06, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Daria Kuznetsova at 7:26 am, Dec 07, 2018

Third Eye By: **THIRD EYE BY**
By Terrence Means at 4:58 pm, Dec 10, 2018

BD Directigen™ *Neisseria meningitidis* Test

For the detection of *Neisseria meningitidis* groups A, C, Y and W135

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INTENDED USE

The BD Directigen™ *N. meningitidis* Groups A, C, Y and W135 Test is a presumptive latex slide agglutination test for the qualitative detection of antigens to *Neisseria meningitidis* Groups A, C, Y and W135 directly in cerebrospinal fluid (CSF), serum or urine. In addition, the test kit provides confirmation and serogrouping capabilities from suspected colonies of *N. meningitidis* Groups A/Y or C/W135. Visible agglutination occurs when a sample containing any of these bacterial antigens is reacted with its respective antibody-coated latex beads.

SUMMARY AND EXPLANATION

The diagnosis of bacteremia and meningitis, especially in young children, can be difficult. As many as 55% of children are seen by a physician and started on antibiotics before meningitis is detected.¹ Detection of microbial antigens in the CSF is a rapid and helpful method for diagnostic microbiology. It may be the single most important test in cases of partially treated meningitis since Gram stain and culture may be negative. Detection of a specific antigen is a clinically significant finding and a valuable aid when choosing antimicrobial therapy.²

Haemophilus influenzae type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.³

Immunological methods for detecting characteristic exoantigens of pathogenic microorganisms in patient fluids (CSF, serum, urine) are typically faster than traditional methods such as culture. These techniques include counterimmunoelectrophoresis (CIE) and latex agglutination.⁴⁻⁷ The latex agglutination procedure has been found to be more rapid and sensitive than CIE in the detection of purified antigen.^{8,9}

PRINCIPLES OF THE PROCEDURE

Specific antibodies are bound to the surface of latex beads. Latex particle aggregation becomes large enough to allow rapid visualization of positive agglutination in the presence of specific antigens. These specific soluble polysaccharide antigens accumulate in CSF, serum or urine as a result of infection by *N. meningitidis* Groups A, B, C, Y and W135, *H. influenzae* type b, *S. pneumoniae*, and Group B *Streptococcus*. These other antigens can be detected with the BD Directigen Meningitis Combo Test Kit, or other BD Directigen Test Kits (see "Availability").

REAGENTS

BD Directigen *N. meningitidis* Groups A, C, Y and W135:

Reagent 1	(1.0 mL),	Anti- <i>N. meningitidis</i> Groups C and W135, Rabbit Antibody-Coated Latex Suspension, with 0.2% sodium azide (preservative),
Reagent 2	(1.0 mL),	Anti- <i>N. meningitidis</i> Groups A and Y, Rabbit Polyclonal Antibody-Coated Latex Suspension, with 0.2% sodium azide (preservative),
Reagent A	(0.5 mL),	Control Latex, Rabbit Immunoglobulin-Coated Latex Suspension, with 0.2% sodium azide (preservative),
Control +	(3.0 mL),	Polyvalent Positive Antigen Control, <i>N. meningitidis</i> Groups A/Y and C/W135, <i>H. influenzae</i> type b, <i>S. pneumoniae</i> , and Group B <i>Streptococcus</i> Antigens, with 0.2% sodium azide (preservative),
Control –	(3.0 mL),	Negative Antigen Control, glycine buffered saline, with 0.2% sodium azide (preservative).

WARNING

H302 Harmful if swallowed.

P264 Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Precautions: For *in vitro* Diagnostic Use.

Warning: Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁰⁻¹³ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

Reagents: Do not use beyond the expiration date. Upon removal from the refrigerator, allow the reagents to warm to room temperature (15–30 °C) before use.

To assure proper drop delivery, reagent dispensing bottles must be held vertically, dispensing one free-falling drop at a time.

Warning: Reagents contain sodium azide, which is very toxic by inhalation, in contact with skin and if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Controls: Do not use the kit if **Control +** and **Control –** do not yield appropriate results.

Test Cards: Cards must be flat for proper reactions. If necessary, flatten cards by bowing back in a direction opposite to that of the curl. Care should be taken not to finger-mark the test areas, since this may result in an oily deposit and improper test results. Use each card once and discard. Store cards in the original package in a dry area at room temperature. When spreading within confines of test areas, avoid scratching the card surface with the plastic stirrers. If the specimen does not spread to the outer perimeter of test area, use another test area of the card.

Test Slide (Glass): If the glass slide is used, disinfect the slide after each use, and wash before reuse. A phenolic disinfectant is suggested. Be sure all detergent and/or disinfectants are thoroughly rinsed from the slide before reuse.

Rotation: The recommended speed for mechanical rotation is 100 ± 2 rpm, but rotation between 95 and 110 rpm does not significantly affect the results obtained. The rotator should circumscribe a circle approximately 2 cm in diameter in the horizontal plane. A moistened humidifying cover should be used to prevent drying of test specimens during rotation.

Storage of Reagents: Upon receipt, remove the carton containing reagents and refrigerate at 2–8 °C. DO NOT FREEZE. Reagents should be recapped and returned to refrigeration when not in use, taking care not to mix color-coded caps.

SPECIMEN COLLECTION AND HANDLING

Refer to appropriate texts for details of specimen collection and handling procedures. Specimens should be tested as soon as possible; however, if the sample cannot be tested immediately, it should be stored at 2–8 °C (for up to 48 h), or at -20 °C.

Serum must be separated from whole blood prior to testing and storage.

SPECIMEN PRETREATMENT

The use of covered glass test tubes (borosilicate) in a 100 °C (boiling) water bath is the most effective means of specimen treatment. Although glass test tubes in a heat block may be used, greater variation in the thorough heating of a sample may be noted as a result of variations in tube size and sample volume.

PROCEDURES

The testing area, all reagents, test specimens and test components should be at room temperature (15–30 °C) when used.

Materials Provided: All materials as listed under "Reagents," work station and accessories.

Materials Required But Not Provided: BD Directigen *Meningitis* Combo Test Cards, BD Directigen Specimen Buffer, Rotator, 100 ± 2 rpm, Humidifying cover, Micropipettor, 50 µL delivery, Pipette tips. } (see "Availability")

Also required are the necessary equipment and labware used for preparation, storage and handling of CSF, serum and urine specimens.

Specimen Preparation (CSF):

1. Heat specimens for 3 min at 100 °C (e.g., water bath or heat block) and allow to cool to room temperature before use.
2. For CSF specimens showing turbidity, centrifuge after heating for 10 min at 1400 x g prior to testing. The supernatant fluid is to be used as the test specimen.
3. Test specimens as described in "Test Procedure."

Specimen Preparation (Serum):

1. Dilute serum specimens of at least 0.6 mL 1:1 with BD Directigen Specimen Buffer and mix.
2. Heat specimens for 5 min at 100 °C (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Using a wooden applicator stick, break up the protein "clot" formed, and vortex vigorously (approximately 5 s).
4. Centrifuge at a minimum of 1400 x g for 15 min.
5. Test supernatant fluid as described in "Test Procedure."

Specimen Preparation (Unconcentrated Urine):

1. Dilute urine specimens of at least 0.4 mL 1:1 with BD Directigen Specimen Buffer and mix.
2. Heat specimens for 5 min at 100 °C (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Centrifuge at a minimum of 1400 x g for 10 min.
4. Test supernatant fluid as described in "Test Procedure."

Specimen Preparation (Concentrated Urine):

1. Urine samples that are turbid or have particulate material should be centrifuged at 1400 x g for 10 min before concentrating.
2. Urine samples may be concentrated 25-fold with a Minicon B-15 concentrator (Amicon Corporation, Danvers, Massachusetts).
3. Dilute at least 200 µL of urine concentrate 1:1 with BD Directigen Specimen Buffer and mix.
4. Heat specimens for 5 min at 100 °C (e.g., water bath or heat block) and allow to cool to room temperature before use.
5. Centrifuge at a minimum of 1400 x g for 10 min.
6. Test supernatant fluid as described in "Test Procedure."

Preparation for Confirmation of Colonies from Culture:

1. Locate suspected colonies on the agar surface from 18–24 h cultures that meet morphological and Gram stain characteristics of organisms that are appropriate for testing with BD Directigen meningitis latex reagents. **Note:** A Gram stain should be performed prior to testing to ensure that organisms are appropriate for testing with the BD Directigen meningitis latex reagents.
2. Pipette 0.5 mL (approximately 10 drops) of Control – reagent into a small glass test tube (10 x 75 mm or equivalent).
3. Select several (2–3) isolated colonies of similar morphology from the original or subculture plate using a sterile loop and suspend into the above tube to achieve a suspension equal to a McFarland #1 turbidity standard. **Note:** Over inoculating will yield an excessively heavy suspension which may cause erroneous results.
4. Heat the suspension for 3 min at 100 °C (e.g., boiling water bath or heat block) and allow to cool to room temperature before testing.
5. Centrifuge at a minimum of 1400 x g for 10 min.
6. Test supernatant as described under “Test Procedure.” **Note:** If atypical agglutination patterns are observed, refer to “Limitations of the Procedure.”

User Quality Control: Include **Control +** and **Control –** testing with each batch of specimens tested as described in step 1, “Test Procedure.”

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Test Procedure (CSF, serum, urine, concentrated urine, and colony confirmation):

Remove the reagents from refrigerated storage and slide the reagent tray into the slot in the kit Work Station.

For testing a specimen with Reagents 1–2, use the following procedure:

Refer to Procedural Chart Illustrations, page 22.

The Meningitis Combo test card or glass slide may be used to test samples. If using the Meningitis Combo card for more than one sample, label circles appropriately. Before use of the glass slide, thoroughly clean with a lint-free tissue.

1. Dispense one drop of **Control +** onto circles “+” in the NmC/W and the NmA/Y columns. Place one drop of **Control –** onto circles “–” in the NmC/W and the NmA/Y columns.
2. Micropipette 50 µL of test sample onto circle(s) of sample row S in the NmC/W and NmA/Y columns and circle A.
3. Holding the dispensing bottle by the cap, vigorously swing (without inverting) to thoroughly mix each reagent. Before uncapping the bottle, gently tap the base on the counter top to assure that no latex remains in the tip.
4. Dispense one drop of **Reagent A** onto circle A.
5. Dispense one drop of **Reagent 1** onto the circles in S, “+” and “–” in the NmC/W column. Dispense one drop of **Reagent 2** onto circles in S, “+” and “–” in the NmA/Y column.
6. Mix the samples and the latex reagents in each circle with a plastic stirrer, alternately using first one end of the stirrer and then the opposite end for the next circle. Discard the stirrer.
7. Place the test card or glass slide on a mechanical rotator and rotate at a speed of 100 ± 2 rpm for 10 min. Use a moistened humidifying cover to prevent evaporation.
8. Immediately at the end of 10 min, read the test results macroscopically under a high intensity incandescent light.

INTERPRETATION OF TEST RESULTS

Record **Control +** and **Control –** test results first.

The **Control +** should yield strong agglutination within 10 min. The **Control –** should show no agglutination.

Agglutination in any of the circles containing **Reagent A** or **Control –** renders the reaction uninterpretable.

Record patient test results:

Positive Test – Should show agglutination. Any degree of agglutination present in one of the latex reagents indicates the presence of the corresponding antigen. *Agglutination in two or more latex reagents or the corresponding BD Directigen™ Reagent-A renders the reaction uninterpretable.*

Negative Test – Should show no agglutination.

LIMITATIONS OF THE PROCEDURE

These latex agglutination tests are not intended as a substitute for bacterial culture. Confirmatory diagnosis of bacterial meningitis infection is only possible with appropriate culture procedures and should be routinely performed.

Samples with extremely low levels of antigen, for instance early in the course of infection, may yield negative results. Furthermore, samples with exceedingly high antigen concentrations may exhibit prozone effects producing inappropriately negative results. Although not extensively studied, prozone phenomena have only been observed in specimens seeded with extremely high antigen levels, and not in clinical specimens.

This assay is for the qualitative detection of antigens to *N. meningitidis* groups A, C, Y or W135 in CSF, serum or urine. In addition, the test kit can also be used with suspected colonies of *N. meningitidis* Groups A/Y or C/W135. Performance with other specimen types has not been assessed. The BD Directigen *Neisseria meningitidis* test cannot be used with blood culture media.

N. meningitidis Groups C and W135 reagent is known to cross react in the presence of *Escherichia coli* K92. *Other cross reactions, uninterpretable and false positives may occur which are reduced or eliminated by proper specimen preparation (see “Specimen Pretreatment” and “Procedures-Specimen Preparation”).*

Untreated serum and urines may yield uninterpretable results with the latex reagents and should be treated as in “Procedures-Specimen Preparation.” CSF specimens or treated urine specimens showing turbidity should be centrifuged for 10 min at 1400 x g after heating, and prior to testing. The supernatant fluid is to be used as the test specimen.

Urine specimens (Processed according to "Procedures-Specimen Preparation") yielding uninterpretable or suspected false positive results may be filtered using a MILLIPORE™ Millex-HA 0.45 µm filter (#SLHA-0250S) and retested with the reactive test latex(es) and the corresponding control latex(es).

Colony confirmation samples that yield uninterpretable or atypical agglutination reactions may require a 1:10 dilution of the supernatant fluid with Control – reagent. Care should be taken to avoid bacterial suspensions greater than a McFarland #1 turbidity standard.

EXPECTED VALUES

Haemophilus influenzae type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.³

PERFORMANCE CHARACTERISTICS

Sensitivity: BD Directigen *N. meningitidis* Groups A and Y, C and W135 latex reagents were compared with counter-immunoelectrophoresis for the ability to detect a partially purified *N. meningitidis* Group polysaccharide antigen in CSF, serum and urine. Commercial antisera specific for each *N. meningitidis* group were used in the CIE test. The BD Directigen *N. meningitidis* latex reagents were equal to or more sensitive than CIE in detecting *N. meningitidis* Group antigens.

A series of retrospective clinical evaluations of the *N. meningitidis* A and Y, C and W135 latex reagents were performed with specimens originally derived from patients with culture, CIE, or commercial latex reagent positive confirmed cases of *N. meningitidis* Groups A, C, and W135 infection. The data presented in **Table 1 (see page 20)** are an accumulation of these evaluations. In some cases CIE testing was not included in some of the evaluation.

Specificity: The specificity of the BD Directigen *N. meningitidis* Groups A and Y, C and W135 latex reagents was determined by testing retrospective and prospective CSF, serum and urine during a series of multiple site clinical evaluations (**Table 2, see page 21**). The latex reagents were tested against culture negative and nonindexed culture positive specimens. The data indicates specificity ranges of 97–100% depending on the latex reagent and specimen tested.

Colony Confirmation: Latex reagents were tested using suspensions of isolated colonies that met morphological characteristics of the organisms. Performance characteristics are listed in **Table 3 (see page 21)**.

AVAILABILITY

Cat. No.	Description
250160	BD Directigen™ <i>Neisseria meningitidis</i> Test, 30 tests.
255460	BD Directigen™ Group B Strep Test, 30 tests.
252260	BD Directigen™ <i>H. influenzae</i> type b Test, 30 tests.
255560	BD Directigen™ <i>Neisseria meningitidis</i> Group B and <i>E. coli</i> K1 Test, 30 tests.
251960	BD Directigen™ <i>S. pneumoniae</i> Test, 30 tests.
252360	BD Directigen™ Meningitis Combo Test, 30 tests.
252480	BD Directigen™ Meningitis Combo Test Card (single use), Box 30.
256391	BD Directigen™ Specimen Buffer Solution, 8 mL.
278051	BD Macro-Vue™ Card Test Rotator, with humidifying cover.

REFERENCES

1. Rothrock, S.G., S.M. Green, J. Wren, D. Letai, L. Daniel-Underwood, E. Pillar: Pediatric bacterial meningitis: is prior antibiotic therapy associated with unaltered clinical presentation?, *Annals Emerg. Med.* 21:146–152, 1992.
2. Fasola, E. and P. Ferrieri: Laboratory diagnostic methods for central nervous system infections, *Neuro. Clinics N. Amer.* 3:279–290, 1992.
3. Morbidity and Mortality Weekly Report (CDC), Surveillance summary: bacterial meningitis and meningococcemia, United States, 28:277–279, 1978.
4. Edwards, E.A.: Immunologic investigations of meningococcal disease, group-specific *Neisseria meningitidis* antigens present in the serum of patients with fulminant meningococcemia., *J. Immunol.* 106:314–317, 1971.
5. Whittle, H.C., Egler, J.J., Tugwell, P. and Greenwood, B.M.: Rapid bacteriological diagnosis of pyogenic meningitis by latex agglutination., *Lancet* 2:619–621, 1974.
6. Ingram, D.L., Anderson, P. and Smith, D.H.: Counterimmunoelectrophoresis in the diagnosis of disease caused by *Haemophilus influenzae* type b., *J. Ped.*, 81:1156–1159, 1972.
7. Newman, R.B., Stevens, R. W. and Gaafar, H.H.: Latex agglutination tests for the diagnosis of *Haemophilus influenzae* Meningitis., *J. Lab. Clin. Med.*, 76:107–113, 1970.
8. Data on file at BD Diagnostic Systems.
9. Tilton, R. C., Dias, F. and Ryan, R.W.: Comparative evaluation of three commercial products and counterimmunoelectrophoresis for the detection of antigens in cerebrospinal fluid, *J. Clin. Microbiol.*, 20:231–234, 1984.
10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
11. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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6

7

Directigen™ MENINGITIS TEST			
	Anti-Hib Latex	Anti-Sp Latex	Anti-GBS Latex
+			
-			
S			

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant /
 Proizvođač / Gyártó / Fabricante / Атқарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirkir / Producent / Producător /
 Производител / Výrobca / Proizvođač / Tilverkare / Üretici / Виробник / 生产厂商

Use by / Исползвайте до / Spøtfebuje do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date
 de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Деяін пайдалануға / Naudokite iki /
 Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použite
 do / Upotrebiti do / Använd före / Son kullanna tarilni / Використати до/line / 使用截止日期
 YYYY-MM-DD / YYYY-MM (MM = end of month)
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)
 AAAA-MM-GG / AAAA-MM (MM = fine mese)
 ЖОЖОЖ-АА-КК / ЖОЖОЖ-АА / (АА = айдың соңы)
 YYYY-MM-DD/YYYY-MM (MM = 월말)
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
 AAAA-MM-DD / AAAA-MM (MM = sluttan av måneden)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
 YYYY-MM-DD / YYYY-MM (MM = 月末)

REF Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo /
 Katalognummer / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Καταλογ αριθμός / Каталогно
 번호 / Katalogo / numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по
 каталогу / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号

EC REP Authorized Representative in the European Community / Оторизирани представител в Европейската общност /
 Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber /
 Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogu / Représentant autorisé
 pour la Communauté européenne / Autorizuirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai
 Közöségben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл /
 유럽 공동체의 위임 대표 / Įgaliotasis atstovas Europos Bendrijose / Pilnvarotais Europos Kopienā / Bevoegde
 vertegenwoordiger in de Europese Gemeenschap / Autorisert representant i EU / Autoryzowane przedstawicielstwo
 we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru
 Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca
 v Európskom spoločenstve / Autorizovano predstavništvo v Evropskoj uniji / Auktoriserad representant i Europeiska
 gemenskapen / Автура Топлупулу Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 歐洲共同體授權代表

IVD In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro
 diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική
 συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiniparatuur / Dispositif médical
 de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnosztikai orvosi eszköz / Dispositivo
 medicale per diagnostica in vitro / Жасанды жағдайда жүзгізетін медициналық диагностика аспабы / In Vitro Diagnostik
 의료 기기 / In vitro diagnostikos prietaisai / Medicinas ierices, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-
 vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico
 para diagnóstico in vitro / Dispozitiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro /
 Medicinska pomôcka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-
 diagnostik / In Vitro Diagnostik Tibbi Cihaz / Медицинский прибор для диагностики in vitro / 体外诊断医疗设备

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung /
 Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena
 temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra
 / Temperatūras ierobežojumi / Temperaturilimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de
 temperatura / Limite de temperatură / Ограничение температуры / Ograničenje teploty / Ograničenje temperature /
 Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

LOT Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας
 (partida) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтам
 коды / 배치 코드(로트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-kode (parti) / Kod partii
 (serie) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serie / Partinummer (Lot) / Partii
 Kodu (Lot) / Код партії / 批号 (亚批)

Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů /
 Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido
 suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova /
 <n> teszthez elegendő / Contenuto sufficiente per <n> test / <n> тесттеpi үшiн жеткінікті / <n> 테스트가 충분히 포함됨
 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder
 tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut
 suficient pentru <n> teste / Достаточно для <n> тестов(а) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova
 / Innehåller tillräckligt för <n> analyser / <n> test için yeterli malzeme içerir / Вистачить для аналізів: <n> / 足够进行 <n>
 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se
 brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones
 de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati
 utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite
 naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz
 instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство
 по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanım
 Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

CONTROL+

Positive control / Положителен контрол / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας /
 Control positivo / Positiivne kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontroll / Controllo positivo / Оң бақылай
 / 양성 컨트롤 / Teigiama kontrolė / Pozitívá kontrolé / Positieve controle / Kontrola dodatnia / Controllo positivo / Control
 pozitiv / Положительный контроль / Pozitif kontrol / Позитивний контроль / 阳性对照试剂

CONTROL-

Negative control / Отрицателен контрол / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας /
 Control negativo / Negatiivne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo /
 Негативтік бақылай / 음성 컨트롤 / Neigiama kontrolė / Negatívá kontrolé / Negatieve controle / Kontrola ujemna / Controllo
 negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативний контроль / 阴性对照试剂

Rx Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." /
 S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed
 practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a
 licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a
 licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of
 a licensed practitioner."

Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA

Benex Limited
 Pottery Road, Dun Laoghaire
 Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
 4 Research Park Drive
 Macquarie University Research Park
 North Ryde, NSW 2113
 Australia

MILLIPORE is a trademark of Millipore Corporation.

© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.