|                                  |                                               |                                                   |                   |
|----------------------------------|-----------------------------------------------|---------------------------------------------------|-------------------|
| <b>Part Number:</b>              | 8809701JAA                                    | <b>BALTSO0191 Version 13.0 Template 4 Inserts</b> |                   |
| <b>Category and Description:</b> | Package Insert,<br>BD BBL™ Crystal™ GP ID Kit | <b>Rev from:</b> 03                               | <b>Rev to:</b> 04 |
|                                  |                                               | <b>Job Number:</b>                                | 558-18            |

Catalog Number: 245140

Blank (Sheet) Size: Length: 11" Width: 8.5"


Number of Pages: 40 Number of Sheets: 10

Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 8.5" X 5.5"

Ink Colors: Number of Colors: 1 PMS #: 513 Purple

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 5

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. BALT8809701 for material information.

**BD** Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**  
By Nancy Carlsen at 2:09 pm, Oct 10, 2018

Proofing Approved By: **PROOFING APPROVED BY**  
By Daria Kuznetsova at 12:07 pm, Oct 11, 2018

Third Eye By: **THIRD EYE BY**  
By Terrence Means at 12:22 pm, Oct 19, 2018

## Gram-Positive ID Kit

English: pages 1 – 6 Italiano: pagine 15 – 19  
 Français: pages 6 – 10 Português: páginas 19 – 23  
 Deutsch: Seiten 10 – 14 Español: páginas 24 – 28

8809701JAA(04)  
 2018-10

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

### INTENDED USE

The BD BBL™ Crystal™ Gram-Positive (GP) Identification (ID) system is a miniaturized identification method employing modified conventional, fluorogenic and chromogenic substrates. It is intended for the identification of frequently isolated aerobic gram-positive bacteria.<sup>1,2,13,16</sup>

### SUMMARY AND EXPLANATION

Micromethods for the biochemical identification of microorganisms were reported as early as 1918.<sup>3</sup> Several publications reported on the use of the reagent-impregnated paper discs and micro-tube methods for differentiating enteric bacteria.<sup>3,4,7,17,19</sup> The interest in miniaturized identification systems led to the introduction of several commercial systems in the late 1960s, and they provided advantages in requiring little storage space, extended shelf life, standardized quality control and ease of use.

In general, many of the tests used in the BD BBL Crystal ID Systems are modifications of classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition, there are chromogen and fluorogen linked substrates, as in the BD BBL Crystal GP ID panel, to detect enzymes that microbes use to metabolize various substrates.<sup>5,7,8,9,11,12,14,15</sup>

The BD BBL Crystal GP ID kit is comprised of (i) BD BBL Crystal GP ID panel lids, (ii) BD BBL Crystal bases and (iii) BD BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid (IF) tubes. The lid contains 29 dehydrated substrates and a fluorescence control on tips of plastic prongs. The base has 30 reaction wells. Test inoculum is prepared with the inoculum fluid and is used to fill all 30 wells in the base. When the lid is aligned with the base and snapped in place, the test inoculum rehydrates the dried substrates and initiates test reactions.

Following an incubation period, the wells are examined for color changes or presence of fluorescence that result from metabolic activities of the microorganisms. The resulting pattern of the 29 reactions is converted into a ten-digit profile number that is used as the basis for identification.<sup>18</sup> Biochemical and enzymatic reaction patterns for the 29 BD BBL Crystal GP ID substrates for a wide variety of microorganisms are stored in the BD BBL Crystal GP ID data base. Identification is derived from a comparative analysis of the reaction pattern of the test isolate to those held in the database. A complete list of taxa that comprises the current database is provided in Table 1.

### PRINCIPLES OF THE PROCEDURE

The BD BBL Crystal GP ID panels contain 29 dried biochemical and enzymatic substrates. A bacterial suspension in the inoculum fluid is used for rehydration of the substrates. The tests used in the system are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Enzymatic hydrolysis of fluorogenic substrates containing coumarin derivatives of 4-methylumbelliferone (4MU) or 7-amino-4-methylcoumarin (7-AMC), results in increased fluorescence that is easily detected visually with a UV light source.<sup>11,12,14,15</sup> Chromogenic substrates upon hydrolysis produce color changes that can be detected visually. In addition, there are tests that detect the ability of an organism to hydrolyze, degrade, reduce or otherwise utilize a substrate in the BD BBL Crystal ID Systems.

Reactions employed by various substrates and a brief explanation of the principles employed in the system are described in Table 2. Panel location in referred tables indicates the row and column where the well is located (example: 1J refers to Row 1 in column J).

## REAGENTS

The **BD BBL™ Crystal™** GP ID panel contains 29 enzymatic and biochemical substrates. Refer to Table 3 for a list of active ingredients.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

After use, all infectious materials including plates, cotton swabs, inoculum fluid tubes, and panels must be autoclaved prior to disposal or incineration.

## STORAGE AND HANDLING/SHELF LIFE

**Lids:** Lids are individually packaged and must be stored unopened in a refrigerator at 2–8 °C. **DO NOT FREEZE.** Visually inspect the package for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Lids in the original packaging, if stored as recommended, will retain expected reactivity until the date of expiration.

**Bases:** Bases are packaged in two sets of ten, in BD BBL Crystal incubation trays. The bases are stacked facing down to minimize air contamination. Store in a dust-free environment at 2–30 °C, until ready to use. Store unused bases in the tray, in plastic bag. Empty trays should be used to incubate inoculated panels.

**Inoculum Fluid:** BD BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid (IF) is packaged in two sets of ten tubes. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store tubes at 2–25 °C. Expiration dating is shown on the tube label. Only BBL Crystal ANR, GP, RGP, N/H Inoculum Fluid should be used with BD BBL Crystal GP ID panels.

On receipt, store the BD BBL Crystal GP ID kit at 2–8 °C. Once opened, only the lids need to be stored at 2–8 °C. The remaining components of the kit may be stored at 2–25 °C. If the kit or any of the components are stored refrigerated, each should be brought to room temperature prior to use.

## SPECIMEN COLLECTION AND PROCESSING

BD BBL Crystal ID Systems are **not** for use directly with clinical specimens. Use isolates from media such as Trypticase™ Soy Agar with 5% Sheep Blood (TSA II) or Columbia Agar with 5% Sheep Blood (Columbia). Use of selective media such as Phenylethyl Alcohol Agar with 5% Sheep Blood (PEA) or Columbia CNA Agar with 5% Sheep Blood (CNA) is also acceptable. Media containing esculin should not be used. The test isolate must be a pure culture, no more than 18–24 h old for most genera; for some slow growing organisms up to 48 h may be acceptable. When swabs are utilized, only cotton-tipped applicators should be used to prepare the inoculum suspensions. Some polyester swabs may cause problems with inoculation of the panels. (See “Limitations of the Procedure.”) Once lids are removed from the sealed pouches, they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator used should be humidified to prevent evaporation of fluid from the wells during incubation. The recommended humidity level is 40–60%. The usefulness of BD BBL Crystal ID Systems or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the *Manual of Clinical Microbiology* for specimen collection, transport and inoculation onto primary isolation media.<sup>1,16</sup>

## TEST PROCEDURE

**Materials Provided:** BD BBL Crystal GP ID Kit –


- 20 BD BBL Crystal GP ID Panel Lids,
- 20 BD BBL Crystal Bases,
- 20 BD BBL Crystal ANR, GP, RGP, N/H ID IF Tubes. Each tube has approximately 2.3 ± 0.15 mL of Inoculum Fluid containing: KCl 7.5 g, CaCl<sub>2</sub> 0.5 g, Tricine N-[2-Hydroxy-1, 1-bis (hydroxymethyl)methyl] glycine 0.895 g, purified water to 1000 mL.
- 2 incubation trays,
- 1 BD BBL Crystal GP ID Report Pad.

**Materials Required But Not Provided:** Sterile cotton swabs (do not use polyester swabs), incubator (35–37 °C) non-CO<sub>2</sub> (40–60% humidity), McFarland No. 0.5 standard, BD BBL Crystal Panel Viewer, BD BBL Crystal ID System Electronic Codebook or BD BBL Crystal GP Manual Codebook, and appropriate culture media.


Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.

**Test Procedure:** BD BBL Crystal GP ID System requires a Gram stain.


1. Remove lids from pouch. Discard desiccant. Once removed from the pouch, covered lids should be used within 1 h. Do not use the panel if there is no desiccant in the pouch.
2. Take an inoculum fluid tube and label with patient's specimen number. Using aseptic technique, with the tip of a sterile cotton swab (do not use a polyester swab) or a wooden applicator stick or disposable plastic loop, pick colonies of the same morphology from one of the recommended media (see section “Specimen Collection and Processing”).
3. Suspend colonies in a tube of BD BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid.


4. Recap tube and vortex for approximately 10–15 s. The turbidity should be equivalent to a McFarland No. 0.5 standard. If the inoculum suspension concentration is in excess of the recommended McFarland standard, one of the following steps is recommended:
  - a. Use a fresh tube of inoculum fluid to prepare a new inoculum suspension equivalent to a McFarland No. 0.5 standard.
  - b. If additional colonies are unavailable for preparation of a new inoculum suspension, using aseptic techniques, dilute the inoculum by adding the minimum required volume (not to exceed 1.0 mL) of 0.85% sterile saline or inoculum fluid to bring down the turbidity equivalent to a McFarland No. 0.5 standard. Remove the excess amount added to the tube with a sterile pipette so that the final volume of inoculum fluid is approximately equivalent to that of the original volume in the tube ( $2.3 \pm 0.15$  mL). Failure to make this adjustment in volume will result in spilling of the inoculum suspension over the black portion of the base rendering the panel unusable.
5. Take a base, and mark the patient's specimen number on the side wall.
6. Pour entire contents of the inoculum fluid tube into target area of the base.


7. Hold base in both hands and roll inoculum gently along the tracks until all of the wells are filled. Roll back any excess fluid to the target area and place the base on a bench top.

8. Align the lid so that the labeled end of the lid is on top of the target area of the base.


9. Push down until a slight resistance is felt. Place thumb on edge of lid towards middle of panel on each side and push downwards simultaneously until the lid snaps into place (listen for two "clicks").


**Purity Plate:** Using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the base and inoculate an agar slant or plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the slant or plate for 24–48 h at 35–37 °C under appropriate conditions. The purity plate or slant may also be used for any supplementary tests or serology, if required.

**Incubation:** Place inoculated panels in incubation trays. Ten panels can fit in one tray (5 rows of 2 panels). All panels should be incubated **face down** (larger windows facing up; label facing down) in a non-CO<sub>2</sub> incubator with 40–60% **humidity**. Trays should not be stacked more than two high during incubation. The incubation time for panels is 18–24 h at 35–37 °C. If panels are incubated for 24 h, they should be read within 30 min after removing from incubator.

**Reading:** After the recommended period of incubation, remove the panels from the incubator. All panels should be read **face down** (larger windows up; label facing down) using the BD BBL™ Crystal™ Panel Viewer. Refer to the color reaction chart and/or Table 3 for an interpretation of the reactions. Use the BD BBL Crystal GP Report Pad to record reactions. Alternatively, the BD BBL Crystal AutoReader may be used to read the panels.

- Read columns E thru J first, using the regular (white) light source.
- Read columns A thru D (fluorescent substrates) using the UV light source in the panel viewer. A fluorescent substrate well is considered positive *only if* the intensity of the fluorescence observed in the well is *greater* than the Negative Control well (4A).

**Calculation of BD BBL Crystal Profile Number:** Each test result (except 4A, which is used as a fluorescence negative control) scored positive is given a value of 4, 2, or 1, corresponding to the row where the test is located. A value of 0 (zero) is given to any negative result. The values resulting from each positive reaction in each column are then added together. A 10-digit number is generated; this is the profile number.

| Example:       | A        | B        | C        | D        | E        | F        | G        | H        | I        | J        |
|----------------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|
| 4              | +        | +        | —        | —        | +        | +        | +        | —        | +        | —        |
| 2              | —        | +        | +        | +        | —        | +        | —        | +        | +        | —        |
| 1              | +        | —        | +        | —        | +        | —        | —        | +        | +        | —        |
| <b>Profile</b> | <b>1</b> | <b>6</b> | <b>3</b> | <b>2</b> | <b>5</b> | <b>6</b> | <b>4</b> | <b>3</b> | <b>7</b> | <b>0</b> |

\* (4A) = fluorescent negative control

The resulting profile number and cell morphology, if known, should be entered on a PC in which the BD BBL Crystal mind software has been installed to obtain the identification. A manual codebook is also available. If a PC is not available contact BD Technical Services for assistance with the identification. If using the BD BBL Crystal AutoReader, organisms are automatically identified by the PC.

**User Quality Control:** Quality control testing is recommended for each lot of panels as follows —

- Inoculate a panel with *Streptococcus pyogenes* ATCC® 19615 per recommended procedure (refer to "Test Procedure").
- Incubate panel for 18–20 h at 35–37 °C.
- Read panel with the panel viewer and color reaction chart; record reactions using the report pad. Alternatively, read the panel on the BD BBL Crystal AutoReader.
- Compare recorded reactions with those listed in Table 4. If discrepant results are obtained, confirm purity of quality control strain before contacting BD Technical Services.

Expected test results for additional quality control test strains are listed in Table 5.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

## LIMITATIONS OF THE PROCEDURE

The BD BBL Crystal GP ID System is designed for the taxa provided. Taxa other than those listed in Table 1 are not intended for use in this system.

The BD BBL Crystal GP ID database was developed with BBL™ brand media. Reactivity of some substrates in miniaturized identification systems may be dependent upon the source media used in inoculum preparations. We recommend the use of the following media for use with the BD BBL Crystal GP ID System: TSA II and Columbia Blood Agar. Use of selective media, such as PEA or CNA is also acceptable. Media containing esculin should not be used.

BD BBL Crystal Identification Systems use a modified microenvironment; therefore, expected values for its individual tests may differ from information previously established with conventional test reactions. The accuracy of the BD BBL Crystal GP ID System is based on statistical use of specially designed tests and an exclusive database.

While BD BBL Crystal GP ID System aids in microbial differentiation, it should be recognized that minor variations may exist in strains within species. Use of panels and interpretation of results require a competent microbiologist. The final identification of the isolate should take into consideration the source of the specimen, aerotolerance, cell morphology, colonial characteristics on various media as well as metabolic end products as determined by gas-liquid chromatography, when warranted.

While the majority of *Enterococcus faecium* isolates identify correctly in the BD BBL Crystal GP System, some strains of Vancomycin Resistant *Enterococcus faecium* produce atypical substrate reactions that can lead to an identification of *Enterococcus durans* or, less frequently, *Helcococcus kunzii*. Therefore, confirmatory testing is recommended when either *Enterococcus durans* or *Helcococcus kunzii* is reported as the identification.

Only cotton-tipped applicator swabs should be used to prepare the inoculum suspension as some polyester swabs may cause the inoculum fluid to become viscous. This may result in insufficient inoculum fluid to fill the wells. Once lids are removed from the sealed pouches they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator where panels are placed should be humidified to prevent evaporation of inoculum fluid from the wells during incubation. The recommended humidity level is 40–60%.

The panels, after inoculation, should only be incubated **face down** (larger windows facing up; label facing down) to maximize the effectiveness of substrates.

If the BD BBL Crystal test profile yields a “No identification” result and culture purity has been confirmed, then it is likely that (i) the test isolate is producing *atypical* BD BBL Crystal *reactions* (which may also be caused by procedural errors), (ii) the test species is not part of the intended taxa or (iii) the system is unable to identify the test isolate with the required level of confidence. Conventional test methods are recommended when user error has been ruled out.

## PERFORMANCE CHARACTERISTICS

**Reproducibility:** In an external study involving four clinical laboratories, (total of four evaluations), the reproducibility of BD BBL Crystal GP ID substrate (29) reactions was studied by replicate testing. The reproducibility of the individual substrate reactions ranged from 79.2% to 100%. The overall reproducibility of BD BBL Crystal GP ID panel was determined to be 96.7%.<sup>20</sup>

**Accuracy of Identification:** The performance of BD BBL Crystal GP ID System was compared to currently available commercial systems using **clinical isolates and stock cultures**. A total of four studies were conducted in four independent laboratories. Fresh, routine isolates arriving in the clinical laboratory, as well as previously identified isolates of the clinical trial sites’ choice were utilized to establish performance characteristics.

Out of 735 total isolates tested from the studies, 668 (90.9%) were correctly identified (including isolates that required supplemental testing) by the BD BBL Crystal GP Identification System. A total of 56 (7.6%) isolates were incorrectly identified, and a message of “No Identification” was obtained for 11 (1.5%) isolates.<sup>20</sup>

## AVAILABILITY

| Cat. No. | Description                                                            | Cat. No. | Description                                                             |
|----------|------------------------------------------------------------------------|----------|-------------------------------------------------------------------------|
| 245140   | BD BBL™ Crystal™ GP Gram-Positive ID System, 1 Kit.                    | 221165   | BD BBL™ Columbia Agar with 5% Sheep Blood, pkg. of 20.                  |
| 245038   | BD BBL™ Crystal™ ANR, GP, RGP, N/H ID Inoculum Fluid.                  | 221263   | BD BBL™ Columbia Agar with 5% Sheep Blood, ctn. of 100.                 |
| 245031   | BD BBL™ Crystal™ Panel Viewer, Domestic model, 110 V, 60 Hz.           | 221352   | BD BBL™ Columbia CNA Agar with 5% Sheep Blood, pkg. of 20.              |
| 245032   | BD BBL™ Crystal™ Panel Viewer, European model, 220 V, 50 Hz.           | 221353   | BD BBL™ Columbia CNA Agar with 5% Sheep Blood, ctn. of 100.             |
| 245033   | BD BBL™ Crystal™ Panel Viewer, Japanese model, 100 V, 50/60 Hz.        | 221179   | BD BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, pkg. of 20.       |
| 245034   | BD BBL™ Crystal™ Panel Viewer, Longwave UV Tube.                       | 221277   | BD BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, ctn. of 100.      |
| 245036   | BD BBL™ Crystal™ Panel Viewer, White Light Tube.                       | 221239   | BD BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), pkg. of 20.  |
| 245037   | BD BBL™ Crystal™ Identification Systems Gram-Positive Manual Codebook. | 221261   | BD BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), ctn. of 100. |
| 441010   | BD BBL™ Crystal™ Mind Software                                         | 212539   | BD BBL™ Gram Stain Kit, pkg. of 4 x 250 mL bottles.                     |

## REFERENCES

- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott’s diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
- Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. Am. J. Public Health. 8:922–923.
- Cowan, S.T., and K.J. Steel. 1974. Manual for the identification of medical bacteria. 2nd ed. Cambridge University Press, Cambridge.
- Edberg, S.C., and C.M. Kontrick. 1986. Comparison of b-glucuronidase-based substrate systems for identification of *Escherichia coli*. J. Clin. Microbiol. 24:368–371.
- Ferguson, W.W., and A.E. Hook. 1943. Urease activity of *Proteus* and *Salmonella* organisms. J. Lab. Clin. Med. 28:1715–1720.
- Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
- Kampfer, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877–2879.
- Killian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae* 1: detection of bacterial glycosidases. Acta Pathol. Microbiol. Scand. Sect. B. 84:245–251.
- MacFaddin, J.F. 1980. Biochemical tests for identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
- Maddocks, J.L., and M. Greenan. 1975. Rapid method for identifying bacterial enzymes. J. Clin. Pathol. 28:686–687.

12. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. *Microbiol. Rev.* 55:335–348.
13. Mandell, G.L., R.G. Douglas, Jr. and J.E. Bennett. 1990. Principles and practice of infectious diseases, 3rd ed. Churchill Livingstone Inc., New York.
14. Mangels, J., I. Edvalson, and M. Cox. 1993. Rapid identification of *Bacteroides fragilis* group organisms with the use of 4-methylumbelliferone derivative substrates. *Clin. Infect. Dis.* 16(54):5319–5321.
15. Moncla, B.J., P. Braham, L.K. Rabe, and S.L. Hiller. 1991. Rapid presumptive identification of black-pigmented gram-negative anaerobic bacteria by using 4-methylumbelliferone derivatives. *J. Clin. Microbiol.* 29:1955–1958.
16. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
17. Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. *Appl. Microbiol.* 5:36–40.
18. Sneath, P.H.A. 1957. The application of computers to taxonomy. *J. Gen. Microbiol.* 17:201–221.
19. Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. *Puerto Rican J. Public Health. Trop. Med.* 25:96–100.
20. Data on file at BD Diagnostics.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or [www.bd.com](http://www.bd.com).

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(...continued...)

Table 3 / Tableau 3 / Tabelle 3 / Tabella 3 / Quadro 3 / Tabla 3

| Panel Location /<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate<br>Substrat<br>Substrato<br>Localização no Painel /<br>Tampão                                                                                                                                 | Code<br>Código<br>Pos.<br>Código | Neg. / Nég.                                                                                                                                                                                           | Active Ingredients<br>Réactifs<br>Reaktive Bestandteile<br>Ingredienti attivi<br>Principios activos<br>Ingredientes activos                                                                             | Approx Amt. (g/L)<br>Quantité approx. (g/L)<br>Unggef. Menge (g/L)<br>Qda approx. (g/L)<br>Qtd. Aprox. (g/L)<br>Cantidad aprox. (g/L) |
|-----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| 4C                                                                                                                                      | L-tryptophan-AMC<br>L-tryptophane-AMC<br>L-Tryptophan-AMC<br>L-triptofano-AMC<br>L-triptófano-AMC                                                                                                       | FTF                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | L-tryptophan-AMC<br>L-tryptophane-AMC<br>L-Tryptophan-AMC<br>L-triptofano-AMC<br>L-triptófano-AMC                                                                                                       | ≤1                                                                                                                                    |
| 2C                                                                                                                                      | L-arginine / Arginin / arginina-AMC                                                                                                                                                                     | FAR                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | L-arginine / Arginin / arginina-AMC                                                                                                                                                                     | ≤1                                                                                                                                    |
| 1C                                                                                                                                      | 4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-Azetyl-β-D-Glukosaminid<br>4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-acetyl-β-D-glucosaminida<br>4MU-N-acetyl-β-D-glucosaminida | FGA                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | 4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-Azetyl-β-D-Glukosaminid<br>4MU-N-acetyl-β-D-glucosaminide<br>4MU-N-acetyl-β-D-glucosaminida<br>4MU-N-acetyl-β-D-glucosaminida | ≤1                                                                                                                                    |
| 4D                                                                                                                                      | 4MU-phosphate / Phosphat / fosfato                                                                                                                                                                      | FHO                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | 4MU-phosphate / Phosphat / fosfato                                                                                                                                                                      | ≤1                                                                                                                                    |
| 2D                                                                                                                                      | 4MU-β-D-glucuronide<br>4MU-β-D-Glukuronid<br>4MU-β-D-glucuronide<br>4MU-β-D-glucurónido                                                                                                                 | FGN                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | 4MU-β-D-glucuronide<br>4MU-β-D-Glukuronid<br>4MU-β-D-glucuronide<br>4MU-β-D-glucurónido                                                                                                                 | ≤1                                                                                                                                    |
| 1D                                                                                                                                      | L-isoleucine / Isoleucin / isoleucina-AMC                                                                                                                                                               | FIS                              | blue fluorescence >FCT well<br>fluorescence blue >puits FCT<br>blaue Fluoreszenz >FCT Vertiefung<br>fluorescência blu >pozzetto FCT<br>fluorescência azul >pogo FCT<br>Fluorescência azul >poçilo FCT | L-isoleucine / Isoleucin / isoleucina-AMC                                                                                                                                                               | ≤1                                                                                                                                    |

(...continued...)

(...continued...)

Table 3 / Tableau 3 / Tabelle 3 / Tabella 3 / Quadro 3 / Tabla 3

| Panel Location /<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate<br>Substrat<br>Substrato<br>Tampão                                                                                       | Code<br>Codice<br>Código | Pos.                                                                                         | Neg. / Nég.                                                                                    | Active Ingredients<br>Réactifs<br>Reaktive Bestandteile<br>Ingredienti attivi<br>Principios activos<br>Ingredientes activos        | Approx Amt. (g/L)<br>Quantité approx. (g/L)<br>Ungef. Menge (g/L)<br>Qta approx. (g/L)<br>Qtd. Approx. (g/L)<br>Cantidad aprox. (g/L) |
|-----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|--------------------------|----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| 4E                                                                                                                                      | Trehalose / Trehalose /<br>Trealosio / Trehalosa                                                                                   | TRE                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo /<br>Dorado/amarillo | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Trehalose / Trehalose /<br>Trealosio / Trehalosa                                                                                   | ≤300                                                                                                                                  |
| 2E                                                                                                                                      | Lactose<br>Lattosio<br>Lactosa                                                                                                     | LAC                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo / Dorado/amarillo    | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Lactose<br>Lattosio<br>Lactosa                                                                                                     | ≤300                                                                                                                                  |
| 1E                                                                                                                                      | Methyl-α & β-glucoside / Méthyl-α<br>& β-glucoside / Methyl-α & β-<br>Gukosid / Methyl-α & β-glucoside /<br>Methyl-α y β-glucósido | MAB                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo /<br>Dorado/amarillo | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Methyl-α & β-glucoside / Méthyl-α<br>& β-glucoside / Methyl-α & β-<br>Gukosid / Methyl-α & β-glucoside /<br>Methyl-α y β-glucósido | ≤300                                                                                                                                  |
| 4F                                                                                                                                      | Sucrose / Saccharose /<br>Saccarosio / Sacarose /<br>Sacarosa                                                                      | SUC                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo / Dorado/amarillo    | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Sucrose / Saccharose /<br>Saccarosio / Sacarose /<br>Sacarosa                                                                      | ≤300                                                                                                                                  |
| 2F                                                                                                                                      | Mannitol<br>Mannitolo<br>Manitol                                                                                                   | MNT                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo / Dorado/amarillo    | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Mannitol<br>Mannitolo<br>Manitol                                                                                                   | ≤300                                                                                                                                  |
| 1F                                                                                                                                      | Maltotriose<br>Maltotrioso<br>Maltotriose<br>Maltotriosa                                                                           | MTT                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo /<br>Dorado/amarillo | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Maltotriose<br>Maltotrioso<br>Maltotriose<br>Maltotriosa                                                                           | ≤300                                                                                                                                  |
| 4G                                                                                                                                      | Arabinose<br>Arabinosio<br>Arabinosa                                                                                               | ARA                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo / Dorado/amarillo    | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Arabinose<br>Arabinosio<br>Arabinosa                                                                                               | ≤300                                                                                                                                  |
| 2G                                                                                                                                      | Glycerol / Glycérol /<br>Glyzerin / Glicerolo / Glicerol                                                                           | GLR                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo /<br>Dorado/amarillo | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Glycerol / Glycérol /<br>Glyzerin / Glicerolo / Glicerol                                                                           | ≤300                                                                                                                                  |
| 1G                                                                                                                                      | Fructose / Fruktose /<br>Fruttosio / Fucose /<br>Fructosa                                                                          | FRU                      | Gold/Yellow / Or/Jaune /<br>Gold/gelb / oro/giallo /<br>Dourado/Amarelo /<br>Dorado/amarillo | Orange/Red / Orange/Rouge /<br>orange/rot / Arancio/rosso /<br>Laranja/Vermelho / Naranja/rojo | Fructose / Fruktose /<br>Fruttosio / Fucose /<br>Fructosa                                                                          | ≤300                                                                                                                                  |

(...continued...)

(...continued...)

Table 3 / Tableau 3 / Tabelle 3 / Tabella 3 / Quadro 3 / Tabla 3

| Panel Location /<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate<br>Substrat<br>Substrato                                                                                                                                                                               | Code<br>Codice<br>Código | Pos.                                                                                                                          | Neg. / Nég.                                                                                                   | Active Ingredients<br>Réactifs<br>Reaktive Bestandteile<br>Ingredienti attivi<br>Principios activos<br>Ingredientes activos                                                                                       | Approx Amt. (g/L)<br>Quantité approx. (g/L)<br>Ungef. Menge (g/L)<br>Qtà approx. (g/L)<br>Qtd. Approx. (g/L)<br>Cantidad aprox. (g/L) |
|-----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| 4H                                                                                                                                      | p-n-p-β-D-glucoside /<br>Glukosid /<br>glucosido                                                                                                                                                                 | BGL                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | p-n-p-β-D-glucoside /<br>Glukosid /<br>glucosido                                                                                                                                                                  | ≤10                                                                                                                                   |
| 2H                                                                                                                                      | p-n-p-β-D-cellobioside /<br>Cellobiosid /<br>cellobiosido                                                                                                                                                        | PCE                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | p-n-p-β-D-cellobioside /<br>Cellobiosid /<br>cellobiosido                                                                                                                                                         | ≤10                                                                                                                                   |
| 1H                                                                                                                                      | Prolin & Leucine-p-nitroanilide<br>Nitroanilide p-proline & leucine<br>Prolin & Leucin-p-Nitroanilid<br>Prolina & Leucina-p-nitroanilido<br>Prolina & Leucina-p-nitroanilido<br>Prolina y Leucina-p-nitroanilida | PLN                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | Proline & Leucine-p-nitroanilide<br>Nitroanilide p-proline & leucine<br>Prolin & Leucin-p-Nitroanilid<br>Prolina & Leucina-p-nitroanilido<br>Prolina & Leucina-p-nitroanilido<br>Prolina y Leucina-p-nitroanilida | ≤10                                                                                                                                   |
| 4I                                                                                                                                      | p-n-p-phosphate /<br>Phosphat /<br>fosfato                                                                                                                                                                       | PHO                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | p-n-p-phosphate /<br>Phosphat /<br>fosfato                                                                                                                                                                        | ≤10                                                                                                                                   |
| 2I                                                                                                                                      | p-n-p-α-D maltoside /<br>maltosid /<br>maltósido                                                                                                                                                                 | PAM                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | p-n-p-α-D maltoside /<br>maltosid /<br>maltósido                                                                                                                                                                  | ≤10                                                                                                                                   |
| 1I                                                                                                                                      | ONPG & p-n-p-α-D galactoside /<br>Galaktosid /<br>galactósido                                                                                                                                                    | PGO                      | Yellow / Jaune / gelb /<br>Giallo / Amarelo / Amarillo                                                                        | Colorless / Incolore / farblos /<br>Incolore / Incolor                                                        | ONPG & p-n-p-α-D galactoside /<br>Galaktosid /<br>galactósido                                                                                                                                                     | ≤10                                                                                                                                   |
| 4J                                                                                                                                      | Urea<br>Urée<br>Harnstoff<br>Ureia                                                                                                                                                                               | URE                      | Aqua/Blue / Bleu/Bleu-vert /<br>aquamarin/bleu /<br>Verde acqua/blu /<br>Verde acqua/azul /<br>Azul/Azul / Verde azulado/azul | Yellow/Green / Jaune/Vert /<br>gelb/grün / Giallo/verde /<br>Amarelo/Verde /<br>Amarillo/verde                | Urea<br>Urée<br>Harnstoff<br>Ureia                                                                                                                                                                                | ≤50                                                                                                                                   |
| 2J                                                                                                                                      | Esculin /<br>Esculina                                                                                                                                                                                            | ESC                      | Brown/Maroon / Brun/Bordeaux /<br>kastanien/braun /<br>Bruno/Marone /<br>Castanho/Bordeaux<br>Pardor/rojo parduzco            | Clear/Tan / Clair/Ocre /<br>klar/gelbbraun / Chiaro/tanino /<br>Transparente/Bronzeado<br>Incoloro/café claro | Esculin /<br>Esculina                                                                                                                                                                                             | ≤25                                                                                                                                   |
| 1J                                                                                                                                      | Arginine /<br>Arginina                                                                                                                                                                                           | ARG                      | Purple / Violet / Illa / Porpora /<br>Roxo / Púrpura                                                                          | Yellow/Gray / Jaune/Gris /<br>gelb/grau / Giallo/grigio /<br>Amarelo/Cincento / Amarillo/gris                 | Arginine /<br>Arginina                                                                                                                                                                                            | ≤200                                                                                                                                  |

**Table 4 / Tableau 4 / Tabelle 4 / Tabella 4 / Quadro 4 / Tabla 4**

Quality Control Chart for BD BBL Crystal GP ID System After 18–20 Hours Incubation from TSA II or Columbia Blood Agar / Tableau de contrôle de qualité du système BD BBL Crystal GP ID Après 18–20 heures d'incubation à partir des géloses TSA II ou des géloses de sang Columbia / Qualitätskontroll-Tabelle für das BD BBL Crystal GP-ID System Nach 18- bis 20-stündiger Inkubation auf TSA II oder Columbia Blutagar / Prospetto del controllo qualità per il sistema BD BBL Crystal GP ID dopo 18–20 ore d'incubazione su agar TSA II o agar sangue Columbia / Gráfico de Control de Qualidade para o Sistema BD BBL Crystal GP ID Após 18 a 20 Horas de Incubação de um TSA II ou Agar de Sangue Columbia / Plantilla del control de calidad para el sistema BD BBL Crystal GP ID después de 18–20 horas de incubación en TSA II o agar sangre Columbia

| Panel Location/<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate / Substrat / Substrato / Tampão                                                                                                                                                                          | Code<br>Codice<br>Código | <i>Streptococcus<br/>pyogenes</i><br>ATCC 19615 |
|----------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|-------------------------------------------------|
| 4A                                                                                                                                     | Fluorescent negative control / Contrôle négatif fluorescent /<br>Fluoreszierende negative Kontrolle / Controllo negativo fluorescente /<br>Controlo negativo da fluorescência / Control fluorescente negativo      | FCT                      | –                                               |
| 2A                                                                                                                                     | 4MU-β-D-glucoside / Glukosid / glucósido                                                                                                                                                                           | FGC                      | –                                               |
| 1A                                                                                                                                     | L-valine / Valin / valina-AMC                                                                                                                                                                                      | FVA                      | +                                               |
| 4B                                                                                                                                     | L-phenylalanine / phénylalanine / Phenylalanin / fenilalanina-AMC                                                                                                                                                  | FPH                      | +                                               |
| 2B                                                                                                                                     | 4MU-α-D-glucoside / Glukosid / glucósido                                                                                                                                                                           | FGS                      | +                                               |
| 1B                                                                                                                                     | L-pyroglutamic acid-AMC / Acide L-pyroglutamique-AMC /<br>L-Pyroglutaminsäure-AMC / L-acido piroglutamico-AMC /<br>L-ácido piroglutâmico-AMC / L-ácido piroglutámico-AMC                                           | FPY                      | +                                               |
| 4C                                                                                                                                     | L-tryptophan / tryptophane / Tryptophan / triptofano / triptófano-AMC                                                                                                                                              | FTR                      | +                                               |
| 2C                                                                                                                                     | L-arginine / Arginin / arginina-AMC                                                                                                                                                                                | FAR                      | +                                               |
| 1C                                                                                                                                     | 4MU-N-acetyl-β-D-glucosaminide / 4MU-N-acétyl-β-D-glucosaminide /<br>4MU-N-Azetyl-β-D-Glucosaminid / 4MU-N-acetil-β-D-glucosaminide /<br>4MU-N-acetyl-β-D-glicosaminida / 4MU-N-acetil-β-D-glucosaminida           | FGA                      | –                                               |
| 4D                                                                                                                                     | 4MU-phosphate / Phosphat / fosfato                                                                                                                                                                                 | FHO                      | +                                               |
| 2D                                                                                                                                     | 4MU-β-D-glucuronide / Glukuronid / glucurónido                                                                                                                                                                     | FGN                      | –                                               |
| 1D                                                                                                                                     | L-isoleucine / Isoleucin / isoleucina-AMC                                                                                                                                                                          | FIS                      | +                                               |
| 4E                                                                                                                                     | Trehalose / Tréhalose / Trealosio / Trehalosa                                                                                                                                                                      | TRE                      | +                                               |
| 2E                                                                                                                                     | Lactose / Lattosio / Lactosa                                                                                                                                                                                       | LAC                      | +                                               |
| 1E                                                                                                                                     | Methyl-α & β-glucoside / Méthyl-α & β-glucoside / Methyl-α & β-Glukosid /<br>Metil-α & β-glucoside / Metil-α y β-glucósido                                                                                         | MAB                      | +                                               |
| 4F                                                                                                                                     | Sucrose / Saccharose / Saccarosio / Sacarose / Sacarosa                                                                                                                                                            | SUC                      | +                                               |
| 2F                                                                                                                                     | Mannitol / Mannitolo / Manitol                                                                                                                                                                                     | MNT                      | –                                               |
| 1F                                                                                                                                     | Maltotriose / Maltotriosio / Maltotriosia                                                                                                                                                                          | MTT                      | +                                               |
| 4G                                                                                                                                     | Arabinose / Arabinosio / Arabinosa                                                                                                                                                                                 | ARA                      | –                                               |
| 2G                                                                                                                                     | Glycerol / Glycérol / Glyzerin / Glicerolo / Glicerol                                                                                                                                                              | GLR                      | +                                               |
| 1G                                                                                                                                     | Fructose / Fruktose / Fruttosio / Frutose / Fructosa                                                                                                                                                               | FRU                      | +                                               |
| 4H                                                                                                                                     | p-n-p-β-D-glucoside / Glukosid / glucósido                                                                                                                                                                         | BGL                      | V                                               |
| 2H                                                                                                                                     | p-n-p-β-D-cellobioside / Cellobiosid / celobiósido                                                                                                                                                                 | PCE                      | –                                               |
| 1H                                                                                                                                     | Proline & Leucine-p-nitroanilide / Nitroanilide p-proline & leucine /<br>Prolin & Leucin-p-Nitroanilid / Prolina & Leucina-p-nitroanilide /<br>Prolina & Leucina-p-nitroanilido / Prolina y Leucina-p-nitroanilida | PLN                      | +                                               |
| 4I                                                                                                                                     | p-n-p-phosphate / Phosphat / fosfato                                                                                                                                                                               | PHO                      | V                                               |
| 2I                                                                                                                                     | p-n-p-α-D-maltoside / Maltosid / maltósido                                                                                                                                                                         | PAM                      | –*                                              |
| 1I                                                                                                                                     | ONPG & p-n-p-α D-galactoside / Galaktosid / galattoside / galactósido                                                                                                                                              | PGO                      | –                                               |
| 4J                                                                                                                                     | Urea / Urée / Harnstoff / Ureia                                                                                                                                                                                    | URE                      | –                                               |
| 2J                                                                                                                                     | Esculin / Esculine / Eskulin / Esculina                                                                                                                                                                            | ESC                      | –                                               |
| 1J                                                                                                                                     | Arginine / Arginin / Arginina                                                                                                                                                                                      | ARG                      | V                                               |

\* = variable when tested from Columbia Blood Agar / variable à partir de la gélose de sang Columbia / variabel auf Columbia-Blutagar /  
variable quando testato su agar sangue Columbia / variável quando testado de Agar de Sangue Columbia / variable cuando se prueba en  
agar sangre Columbia.

Table 5 / Tableau 5 / Tabelle 5 / Tabella 5 / Quadro 5 / Tabla 5

Additional Quality Control Strains for BD BBL Crystal GP ID System After 18–20 Hours Incubation from TSA II or Columbia Blood Agar / Tableau supplémentaire de contrôle de qualité du système BD BBL Crystal GP ID après 18–20 heures d'incubation à partir des géloses TSA II ou des géloses de sang Columbia / Zusätzliche Qualitätskontroll-Stämme für das BD BBL Crystal GP-ID System nach 18– bis 20-stündiger Inkubation auf TSA II oder Columbia Blutagar / Ceppi supplementari di controllo qualità per il sistema BD BBL Crystal GP ID dopo 18–20 ore d'incubazione su agar TSA II o agar sangue Columbia / Estripes de Control de Qualidade Adicionais para o Sistema BD BBL Crystal GP ID Após 18 a 20 Horas de Incubação de um TSA II ou Agar de Sangue Columbia / Cepas adicionales del control de calidad para el sistema BD BBL Crystal GP ID después de 18–20 horas de incubación en TSA II o el agar sangre Columbia

| Panel Location /<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate<br>Substrat<br>Substrato<br>Tampão                                                                                                                                                            | Code<br>Codice<br>Código | Staphylococcus<br>epidermidis<br>ATCC 12228 | Bacillus<br>brevis<br>ATCC 8246 | Enterococcus<br>faecalis<br>ATCC 19433 | Staphylococcus<br>xylosum<br>ATCC 35033 |
|-----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|---------------------------------------------|---------------------------------|----------------------------------------|-----------------------------------------|
| 4A                                                                                                                                      | Fluorescent negative control / Contrôle négatif fluorescent / Fluoreszierende negative Kontrolle / Controllo negativo fluorescente / Controlo negativo da fluorescência / Control fluorescente negativo | FCT                      | –                                           | –                               | –                                      | –                                       |
| 2A                                                                                                                                      | 4MU-β-D-glucoside / Glucosid / glucósido                                                                                                                                                                | FGC                      | –                                           | +                               | +                                      | –                                       |
| 1A                                                                                                                                      | L-valine / Valin / valina-AMC                                                                                                                                                                           | FVA                      | –                                           | V                               | –                                      | –                                       |
| 4B                                                                                                                                      | L-phenylalanine / phénylalanine / Fenilalanina-AMC                                                                                                                                                      | FPH                      | –                                           | +                               | +                                      | –                                       |
| 2B                                                                                                                                      | 4MU-α-D-glucoside / Glucosid / glucósido                                                                                                                                                                | FGS                      | –*                                          | +                               | +                                      | –                                       |
| 1B                                                                                                                                      | L-pyroglutamic acid-AMC / Acide L-pyroglutamique-AMC / L-Pyroglutaminsäure-AMC / L-pyroglutámico-AMC / L-ácido piroglutámico-AMC / L-ácido piroglutámico-AMC                                            | FPY                      | –                                           | +                               | +                                      | V                                       |
| 4C                                                                                                                                      | L-tryptophan-AMC / L-tryptophane-AMC / L-Tryptophan-AMC / L-triptófano-AMC                                                                                                                              | FTR                      | –                                           | +                               | +                                      | V                                       |
| 2C                                                                                                                                      | L-arginine / Arginin / arginina-AMC                                                                                                                                                                     | FAR                      | V                                           | +                               | –                                      | –                                       |
| 1C                                                                                                                                      | 4MU-N-acetyl-β-D-glucosamine / 4MU-N-acetyl-β-D-glucosaminide / 4MU-N-Azetyl-β-D-Glukosaminid / 4MU-N-acetil-β-D-glucosaminide / 4MU-N-acetyl-β-D-glucosaminida / 4MU-N-acetil-β-D-glucosaminida        | FGA                      | –                                           | +                               | +                                      | –                                       |
| 4D                                                                                                                                      | 4MU-phosphate / Phosphat / fosfato                                                                                                                                                                      | FHO                      | +                                           | V                               | V                                      | +                                       |
| 2D                                                                                                                                      | 4MU-β-D-glucuronide / Glukuronid / glucurónido                                                                                                                                                          | FGN                      | –                                           | –                               | –                                      | +                                       |
| 1D                                                                                                                                      | L-isoleucine / Isoleucin / isoleucina-AMC                                                                                                                                                               | FIS                      | –                                           | V                               | –                                      | –                                       |
| 4E                                                                                                                                      | Trehalose / Trehalose / Trealosio / Trehalosa                                                                                                                                                           | TRE                      | –                                           | –                               | +                                      | +                                       |
| 2E                                                                                                                                      | Lactose / Lattosio / Lactosa                                                                                                                                                                            | LAC                      | +                                           | –                               | +                                      | +                                       |
| 1E                                                                                                                                      | Methyl-α & β-glucoside / Méthyl-α & β-glucoside / Methyl-α & β-glucoside / Metil-α & β-glucosido                                                                                                        | MAB                      | –                                           | –                               | +                                      | +                                       |
| 4F                                                                                                                                      | Sucrose / Saccharose / Saccarosio / Sacarose / Sacarosa                                                                                                                                                 | SUC                      | +                                           | –                               | +                                      | +                                       |

\* = variable when tested from Columbia Blood Agar / variable à partir de la gélose de sang Columbia / variabile quando testato su agar sangue Columbia / variavel quando testado de Agar de Sangue Columbia / variable cuando se prueba en agar sangre Columbia.

(...continued...)

(...continued...)

Table 5 / Tableau 5 / Tabelle 5 / Tabella 5 / Quadro 5 / Tabla 5

| Panel Location /<br>Pos. sur panel /<br>Panel-Position /<br>Posizione nel pannello /<br>Localização no Painel /<br>Posición en el panel | Substrate<br>Substrat<br>Substrato<br>Tampão                                                                                                                                                                    | Code<br>Codice<br>Código | Staphylococcus<br>epidermidis<br>ATCC 12228 | Bacillus<br>brevis<br>ATCC 8246 | Enterococcus<br>faecalis<br>ATCC 19433 | Staphylococcus<br>xylosum<br>ATCC 35033 |
|-----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|---------------------------------------------|---------------------------------|----------------------------------------|-----------------------------------------|
| 2F                                                                                                                                      | Mannitol / Mannitolo / Manitol                                                                                                                                                                                  | MNT                      | –                                           | –                               | +                                      | +                                       |
| 1F                                                                                                                                      | Maltotriose / Maltotrioso / Maltotriosa                                                                                                                                                                         | MTT                      | +                                           | –                               | +                                      | V                                       |
| 4G                                                                                                                                      | Arabinose / Arabinosio / Arabinosa                                                                                                                                                                              | ARA                      | –                                           | –                               | –                                      | V                                       |
| 2G                                                                                                                                      | Glycerol / Glicérol / Glicerolo / Glicerol                                                                                                                                                                      | GLR                      | +                                           | –                               | +                                      | +                                       |
| 1G                                                                                                                                      | Fructose / Fructose / Fruttosio / Fructose / Fructosa                                                                                                                                                           | FRU                      | +                                           | –                               | +                                      | +                                       |
| 4H                                                                                                                                      | p-n-p-β-D-glucoside / Glucosid / glucósido                                                                                                                                                                      | BGL                      | –                                           | V                               | +                                      | +                                       |
| 2H                                                                                                                                      | p-n-p-β-D-cellobioside / Cellobiosid / cellobiósid                                                                                                                                                              | PCE                      | –                                           | –                               | +                                      | –                                       |
| 1H                                                                                                                                      | Proline & Leucine-p-nitroanilide / Nitroanilide p-proline & leucine / Prolina y Leucin-p-Nitroanilid /<br>Prolin & Leucine-p-nitroanilide / Prolina & Leucina-p-nitroanilido / Prolina y Leucina-p-nitroanilide | PLN                      | V                                           | V                               | –                                      | –                                       |
| 4I                                                                                                                                      | p-n-p-phosphate / Phosphat / fosfato                                                                                                                                                                            | PHO                      | V                                           | V                               | V                                      | +                                       |
| 2I                                                                                                                                      | p-n-p-α-D-maltoside / Maltosid / maltósido / maltoside                                                                                                                                                          | PAM                      | –*                                          | V                               | +                                      | –*                                      |
| 1I                                                                                                                                      | ONPG & p-n-p-α-D-galactoside / Galaktosid / galactósido / galactósido                                                                                                                                           | PGO                      | V                                           | –                               | –                                      | V                                       |
| 4J                                                                                                                                      | Urea / Urée / Harnstoff / Ureia                                                                                                                                                                                 | URE                      | +                                           | V                               | V                                      | +                                       |
| 2J                                                                                                                                      | Esculin / Esculine / Eskulin / Esculina                                                                                                                                                                         | ESC                      | –                                           | V                               | +                                      | –                                       |
| 1J                                                                                                                                      | Arginine / Arginin / Arginina                                                                                                                                                                                   | ARG                      | V                                           | +                               | +                                      | V                                       |

\* = variable when tested from Columbia Blood Agar / variable à partir de la gélose de sang Columbia / variabile quando testato su agar sangue Columbia / variável quando testado de Agar de Sangue Columbia / variable cuando se prueba en agar sangre Columbia.

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvodač / Gyártó / Fabricante / Атқарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvodač / Tilverkare / Üretici / Виробник / 生产厂商

Use by / Използвайте до / Spottebujtje do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebi do / Felhasználhatósági dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Использовать до / Pouzítie do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до/line / 使用截止日期  
 YYYY-MM-DD / YYYY-MM (MM = end of month)  
 ГТТТ-ММ-ДД / ГТТТ-ММ (ММ = края на месеца)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)  
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)  
 EEEE-MM-NH / EEEE-MM (MM = τέλος του μήνα)  
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)  
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)  
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)  
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)  
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)  
 AAAA-MM-GG / AAAA-MM (MM = fine mese)  
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)  
 YYYY-MM-DD/YYYY-MM (MM = 월말)  
 MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga)  
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)  
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)  
 AAAA-MM-DD / AAAA-MM (MM = sluten av måneden)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)  
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)  
 ГТТТ-ММ-ДД / ГТТТ-ММ (ММ = конец месяца)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)  
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)  
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)  
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)  
 YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номери / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus number / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalógové číslo / Kataloški broj / Katalog numaras / Номер за каталогом / 目录号

Authorized Representative in the European Community / Авторизирани представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret representant i De Europæiske Fællesskaber / Autoriserter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizad en la Comunidad Europea / Volititud esindaja Euroopa Liidukogus / Représentant autorisé pour la Communauté européenne / Autorizirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Rappresentante autorizzato nella Comunità Europea / Европа крауымдастындагы уәкәлеттi өкiл / 유럽 공동체의 위임 대표 / Įgaliojasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autorisert representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Représentant autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupce v Evropském společenství / Autorizovanio predstavništvo u Evropskoj uniji / Autoriserad representant i Europeiska gemenskapen / Автура Топлулугу Yetkilii Temsilcisi / Упономочений представник у крaїнах ЄС / 歐洲共同體授權代表

In Vitro Diagnostic Medical Device / Медицински уред за диагностика in vitro / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiiniaparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Diagnostiku / In vitro diagnostikaitali orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жагдайда жүргізілетін медициналық диагностика аспабы / In Vitro Diagnostico 의료 기기 / In vitro diagnostikos prietaisais / Medicinas ierices, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro / Medicinska pomočka na diagnostiku in vitro / Medicinski uređaj za in vitro diagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Diyagnostik Tibbi Cihaz / Медицинский прибор для диагностики in vitro / 体外诊断医疗设备

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hörséskleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperatuurlimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kod / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топлама коды / 배치 코드(로트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-code (parti) / Kod partii (serie) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kod série (šarža) / Kod serie / Partinummer (Lot) / Parti Kodu (Lot) / Код партии / 批号 / 亚批

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitlykite naujojoimo instrukcijas / Skatiti lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i brugsanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultati instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se brugsanvisningen / Kullanni Talimatları'na başvurun / Див. інструкції за використання / 请参阅使用说明

Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

## Australian Sponsor:

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.