8808191JAA
2010/07
Italiano

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REFERENCES

1. Trepeta, R.W., and S.C. Edberg. 1984. Methylumbelliferyl- β -D-glucuronide- based medium for rapid isolation and identification of *Escherichia coli*. J. Clin. Microbiol. 19:172-174.
2. Killian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae*. I. Detection of bacterial glycosidases. Acta Pathol. Microbiol. Scand. Sec.B 84:245-251.
3. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
4. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol 17:53-80.
5. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
6. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
7. Isenberg, H.D., F.D. Schoenknecht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
8. Miller, J.M., and H.T. Holmes. 1999. Specimen collection, transport, and storage, p. 33-63. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
9. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
10. Feng, P.C.S., and P.A. Hartman. 1982. Fluorogenic assays for immediate confirmation of *Escherichia coli*. Appl. Environ. Microbiol. 43:1320-1329.
11. Robison, B.J. 1984. Evaluation of a fluorogenic assay for detection of *Escherichia coli* in foods. Appl. Environ. Microbiol. 48:285-288.
12. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual™ of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.
13. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1999. Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.

-  Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Uretici / Proizvođač / Производител / Аткарушы
-  Use by / Spotřebuje do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeikäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Pouzítě do / Usar antes de / Använd före / Исполняйте до / A se utiliza până la / Son kullanna tarihi / Upotrebiti do / Исползовать до / дейи пайдаланура / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = konec měsíce) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / ÁÁÁÁ-KK-PP / ÁÁÁÁ-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / ÁÁÁÁ-MM-JJ / ÁÁÁÁ-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / ÁÁÁÁ-MM-GG / ÁÁÁÁ-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mêsio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = fim do mês) / RRRR-MM-DD / RRRR-MM (MM = koniec mesiac) / áááá-mm-dd / áááá-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) / ГТТТ-ММ-ДД / ГТТТ-ММ (ММ = края на месеца) / ÁÁÁÁ-LL-ZZ / ÁÁÁÁ-LL (LL = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayın sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГТТТ-ММ-ДД / ГТТТ-ММ (ММ = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
- REF** Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalogszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numerasi / Kataloški broj / Номер по каталогу / Каталог номери
- EC REP** Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatus esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирани представител в EU / Reprezentant autorizat în Uniunea Europeană / Автура Төрлүлүгү Yetkilii Төмсүлсү / Овлашчени представник у Европској заједници / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizuirani predstavnik u EU
- IVD** In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro / In vitro diagnostikos prietaisas / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska romčicka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики ин витро / Жасанды жагдайда жүргізітін медициналық диагностика аспабы / Medicinska romagala za In Vitro Dijagnostiku
-  Temperature limitation / Teplotni omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturbereich / Οριο θερμοκρασίας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohrančenie teploty / Limitación de temperatura / Temperaturbegrænsning / Температурни ограничения / Limitare de temperatură / Sicaklik sinirlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dozvoljena temperatura
- LOT** Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (lotul) / Parti Kodu (Lot) / Kod serije / Код партии (лот) / Топтама коды / Lot (kod)
-  Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направете справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na bakurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайданалуу нускаулығымән танысын алыңыз / Koristi upute za upotrebu
-  Keep away from light / Nevystavujte světlu / Må ikke udsættes for lys / Weghouden van licht / Hoida eemal valgusest / Suojattava valolta / Conserver à l'abri de la lumière / Vor Licht schützen / Φυλάξτε το μακριά από φως / Fény nem érheti / Tenere al riparo dalla luce / Laikyti atokiau nuo šilumos šaltinių / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Uchovávať mimo dosahu svetla / Mantener alejado de la luz / Får ej utsättas för ljus / Παετέ от светлина / A se feri de lumină / Işıktan uzak tutun / Držite dalje od svetlosti / Хранить в темноте / Қараңғыланған жерде ұста / Držati dalje od svetla

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
800-638-8663
www.bd.com/ds

EC REP Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare, Ireland