


## Revisions

SO 0191-5

Rev from	Rev to	ECO #
0703	2010/07	5393-10

### Notes:


- BD Cat. Number 297941, 297841
- Blank (Sheet) Size: Length: 11"      Width: 8.5"  
 Number of Pages: 2      Number of Sheets: 1  
 Page Size: Length 11"      Width 8.5"      Final Folded Size: 4.5" x 1.5"
- Style (see illustrations below): # 1




#1




#2




#3




#4



#5




#6



#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No  
 No. of Colors: 1      PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 <b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8808511JAA		Category and Description  Package Insert, TSA II Deeps	Sheet: 1 of 3	A
			Scale: N/A	

## BD BBL™ Prepared Tubed General Purpose Medium

## Trypticase™ Soy Agar, Modified (TSA II) Deeps



8808511JAA  
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD-representant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítás kérje a BD helyi képviselőtől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakt lokal Becton Dickinson-representant for anvisningar. / **Свяжитесь с местным представителем на BD за инструкциями.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / **Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз.** / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

**Trypticase™ Soy Agar, Modified (TSA II)** supplemented with blood is used for cultivating fastidious microorganisms and for the visualization of hemolytic reactions produced by many bacterial species.

## SUMMARY AND EXPLANATION

The nutritional composition of **Trypticase** Soy Agar has made it a popular medium, both unsupplemented and as a base for media containing blood. TSA II is an improved version of the original **Trypticase** Soy Agar formulation for use with animal blood supplements. With 5 or 10% sheep blood, it is extensively used for the recovery and cultivation of fastidious microbial species and for the determination of hemolytic reactions that are important differentiating characteristics for bacteria, especially *Streptococcus* species.

## PRINCIPLES OF THE PROCEDURE

The combination of casein and soy peptones renders the medium highly nutritious by supplying organic nitrogen, particularly amino acids and longer-chained peptides. The sodium chloride maintains osmotic equilibrium.

Defibrinated sheep blood is the most widely used blood for enriching agar base media.<sup>1</sup> Hemolytic reactions of streptococci are proper and growth of *Haemophilus hemolyticus*, a nonpathogen whose hemolytic colonies are indistinguishable from those of beta-hemolytic streptococci, is inhibited.

## REAGENTS

**Trypticase™ Soy Agar, Modified (TSA II)**

Approximate Formula\* Per Liter Purified Water

Pancreatic Digest of Casein	14.5 g
Papaic Digest of Soybean Meal	5.0 g
Sodium Chloride	5.0 g
Agar	14.0 g
Growth Factors	1.5 g

\* Adjusted and/or supplemented as required to meet performance criteria.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Tubes with tight caps should be opened carefully to avoid injury due to breakage of glass.

**Storage Instructions:** On receipt, store tubes in the dark at 2 to 25°C. Avoid freezing and overheating. Do not open until ready to use. Tubed media stored as labeled until just prior to use may be inoculated up to the expiration date and incubated for the recommended incubation times. Minimize exposure to light.

**Product Deterioration:** Do not use tubes if they show evidence of microbial contamination, discoloration, drying or other signs of deterioration.

## SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be obtained using various techniques. For detailed information, consult appropriate texts.<sup>2,3</sup> Specimens should be obtained before antimicrobial agents have been administered. Provision must be made for prompt delivery to the laboratory.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>4-7</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

## PROCEDURE

**Material Provided: Trypticase Soy Agar, Modified (TSA II)**

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

**Test Procedure:** Observe aseptic techniques. To prepare plated medium, place agar deeps with loosened caps in a boiling water bath until the medium becomes liquid (clear). Cool to 45 to 50°C, add blood, if desired, and pour into sterile Petri dishes. Allow the medium to solidify and dry before use. The agar surface should be smooth and moist, but without excessive moisture.

Inoculate the medium as soon as possible after the specimen arrives at the laboratory. The streak plate is used primarily to isolate pure cultures from specimens containing mixed flora. Alternatively, if material is being cultured directly from a swab, roll the swab over a small area of the surface at the edge; then streak from this inoculated area. Since many pathogens require carbon dioxide on primary isolation, plates may be incubated in an atmosphere containing approximately 3 to 10% CO<sub>2</sub>. Incubate plates at 35 ± 2°C for 18 to 24 h.

**User Quality Control:**

1. Examine the tubes for signs of deterioration as described under "Product Deterioration".
2. Check performance by inoculating a representative sample of plates with pure cultures of stable control organisms that give known, desired reactions. The following test strains are recommended:

Test Strain	Expected Results
Medium without the addition of blood.	
<i>Shigella flexneri</i> ATCC™ 12022	Growth
<i>Staphylococcus aureus</i> ATCC 25923	Growth
Medium with the addition of 5% defibrinated sheep blood.	
<i>Streptococcus pneumoniae</i> ATCC 6305	Growth. Colonies surrounded by zones of alpha hemolysis (green).
<i>Streptococcus pyogenes</i> ATCC 19615	Growth. Colonies surrounded by zones of beta hemolysis.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

After sufficient incubation, the plates should show isolated colonies in streaked areas and confluent growth in areas of heavy inoculation.

Hemolytic reactions should be noted for organisms inoculated on the medium containing blood.

### LIMITATIONS OF THE PROCEDURE

For identification, organisms must be in pure culture. Morphological, biochemical, and/or serological tests should be performed for final identification. Consult appropriate texts for detailed information and recommended procedures.<sup>2,8</sup>

## PERFORMANCE CHARACTERISTICS

**Trypticase Soy Agar (ISA)** with 5% Sheep Blood was used as a control in a study using broth enhanced culture (Todd Hewitt) and Optical Immunoassay method for the diagnosis of  $\beta$ -hemolytic streptococcal infection. Five hundred two (502) specimens were tested. TSA with 5% Sheep Blood had a sensitivity and specificity of 92.5% and 99.4%, respectively.<sup>9</sup> Nguyen et al. used **Trypticase Soy Agar with 5% Sheep Blood** as the 'gold standard' for the detection of group B *Streptococcus* from the lower genital tract of pregnant women.<sup>10</sup> In another study, Rossmann et al. successfully reisolated *Lautropia mirabilis* on **Trypticase Soy Agar** with 5% Sheep Blood from the oral cavities of human immunodeficiency virus infected children.<sup>11</sup> Of the 85 children evaluated in this study, 35 (41.4%) were positive for *L. mirabilis*. Isenberg et al. used **Trypticase Soy Agar with 5% Sheep Blood** as a control to evaluate the recovery of *Enterococcus* from a selective medium under study.<sup>12</sup> Two hundred fifty (250) group D streptococcal strains isolated from clinical material and 8 strains obtained from the National Communicable Disease Center (Atlanta, Ga.) were used. Kantor et al. maintained stock cultures at room temperature using **Trypticase Soy Agar** slants covered with sterile mineral oil for a study on the identification of nonfermentative gram-negative bacteria in the clinical laboratory.<sup>13</sup>

## AVAILABILITY

Cat. No.	Description
297941	<b>BBL™ Trypticase™</b> Soy Agar, Modified (TSA II) Deepes, 20 mL, Ctn. of 100 size A tubes
297841	<b>BBL™ Trypticase™</b> Soy Agar, Modified (TSA II) Deepes, 9 mL, Ctn. of 100 size D tubes

## REFERENCES

1. Vera, H.D. and D.A. Power. 1980. Culture media, p. 969. In E.H. Lennette, A. Balows, W.J. Hausler, Jr., and J.P. Truant (ed.), Manual of clinical microbiology, 3rd ed. American Society for Microbiology, Washington, D.C.
2. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1999. Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
3. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
4. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
6. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
8. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual™ of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.
9. Fries, S.M. 1995. Diagnosis of group A streptococcal pharyngitis in a private clinic: comparative evaluation of an optical immunoassay method and culture. J. of Ped. vol. 126, number 6.
10. Nguyen, T.M. et al. 1998. Detection of group B streptococcus: comparison of an optical immunoassay with direct plating and broth-enhanced culture methods. J. Matern. Fetal. Med. Jul-Aug; 7 (4): 172-176.
11. Rossmann, S.N. et al. 1998. Isolation of *Lautropia mirabilis* from oral cavities of human immunodeficiency virus infected children. J. Clin. Microbiol. 36: 1756-1760.
12. Isenberg, H.D., Goldberg, D. and J. Sampson. 1970. Laboratory studies with a selective medium. Appl. Microbiol. Sept. 1970, p. 443-436.
13. Kantor, L.T., Spyros, D.K. and R.B. Yee. 1975. Identification of nonfermentative gram-negative bacteria in the clinical laboratory. Amer. J. of Med. Tech. vol. 41, number 1.



Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA  
800-638-8663  
www.bd.com/ds



Benex Limited  
Rineanna House  
Shannon Free Zone  
Shannon, County Clare, Ireland

ATCC is a trademark of the American Type Culture Collection.  
BD, BD Logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company.  
© 2010 BD.



Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производител / Аткарушы



Use by / Spotřebujte do / Användes för / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Usar antes de / Använd före / Используйте до / А se utiliza până la / Son kullanna tarihi / Upotrebiti do / Исползовать до / дейін пайдаланыра / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = konec měsíce) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / AAAA-MM-JJ / AAAA-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / AAAA-MM-GG / AAAA-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mêsio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / AAAA-MM-DD / AAAA-MM (MM = fim do mês) / RRRR-MM-DD / RRRR-MM (MM = koniec mesiacu) / aaaa-mm-dd / aaaa-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = sluet på månaden) / ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца) / AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayın sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)



Catalog number / Katalógové číslo / Katalognummer / Catalognummer / Kataloogi number / Tuotenumber / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalógové číslo / Número de catálogo / Kataložen nomer / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу / Каталог нөмірі



Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatus esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgalotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирован представитель в ЕС / Représentant autorizat în Uniunea Europeană / Автура Топлугу Yetkili Temsilcisi / Ovlašteni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizirani predstavnik u EU



In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostikkalalaite / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnostikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska pomůcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uredaj za in vitro dijagnostiku / Медицинский прибор для диагностики in vitro / Жасанды жағдайда жүргізетін медициналық диагностика аспабы / Medicinska pomagala za In Vitro Dijagnostiku



Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturenbereich / Οριο θερμοκρασίας / Hömörsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenie teploty / Limitación de temperatura / Temperaturbegrænsning / Температури ограничения / Limitare de temperatură / Sıcaklık sınırlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dovoljena temperatura



Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod serije / Код партии (лот) / Топтама коды / Lot (kod)



Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направете справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na başvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайдалану нұсқаулығымен танысып алыңыз / Koristi upute za upotrebu