
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	8507-16

NOTES:

- BD Catalog Number: 231040, 231041, 231042, 231552
- Blank (Sheet) Size: Length: 22 1/2" Width: 10"
- Number of Pages: 10 Number of Sheets: 1
- Page Size: Length: 4 1/2" Width: 10" Final Folded Size: 4 1/2"L x 1 1/4"W
- Ink Colors: No. of Colors: 1 PMS#: Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. VS-88-0067-1PR for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VENDOR SPEC HELD BY BD CARIBE, LTD.

Label Design	REVISED BY By Sonia Thompson at 1:37 pm, Jun 22, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Daria Kuznetsova at 3:13 pm, Jun 22, 2016		
Checked By	THIRD EYE BY By Nichole Graham at 12:53 pm, Sep 22, 2016		
Part Number:	8800671JAA	Category and Description Package Insert, BD BBL Taxo A Discs for Goup A Streptococci	Sheet: 1 of 11 Scale: N/A A

BD BBL™ Taxo™ A Discs for Differentiation of Group A Streptococci

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INTENDED USE

BD BBL Taxo™ A discs are for the presumptive identification of group A beta-hemolytic streptococci based on susceptibility to a low level of bacitracin. Discs are intended for use with pure cultures, with the exception noted under "Specimens."

SUMMARY AND EXPLANATION

BD BBL Taxo A discs are impregnated with a low level of bacitracin. According to the work of Maxted, of Levinson and Frank, and others, the group A streptococci may be differentiated from the other Lancefield groups of hemolytic streptococci by the formation of a zone of inhibition around the disc.¹⁻⁶

PRINCIPLES OF PROCEDURE

Group A beta-hemolytic streptococci are sensitive to small amounts of bacitracin, while beta-hemolytic streptococci of other serologic groups are more resistant. The **BD BBL Taxo** A disc on a blood agar plate can be used for presumptive identification of group A beta-hemolytic streptococci after overnight incubation.

REAGENTS

BD BBL Taxo A discs are impregnated with approximately 0.04 unit of bacitracin per disc.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

BD BBL Taxo A discs *are not for susceptibility testing.*

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, test plates and other contaminated materials must be sterilized by autoclaving before discarding. Directions for use should be read and followed carefully.

Storage Instructions: On receipt, store at -20 to +8 °C. After use, store vial or cartridge to protect product integrity at 2 to 8 °C.

Use oldest discs first and discard expired discs. Allow containers to come to room temperature before opening. Return unused discs to the refrigerator when application of discs has been completed. Vials and cartridges from which discs have been frequently removed during one week and discs left out overnight in the laboratory should be discarded, or the discs should be tested for performance with control organisms prior to continued use.

SPECIMENS

BD BBL Taxo A discs are not for use directly with clinical specimens or other sources containing mixed flora. The organism to be presumptively identified must first be isolated as separate colonies by streaking the specimen onto appropriate culture media; e.g., **Trypticase™** Soy Agar with 5% Sheep Blood (TSA II). However, the **BBL™** Group A Selective Strep Agar with 5% Sheep Blood (**ssA™**) plate, which may be used for the primary isolation of group A streptococci from throat specimens, was designed specifically for use with the **BD BBL Taxo** A disc.⁷

PROCEDURE

Material Provided: **BD BBL Taxo** A Discs.

Materials Required But Not Provided: Ancillary culture media, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure

1. Inoculate a **Trypticase** Soy Agar with 5% Sheep Blood plate with the test organism exhibiting beta-hemolysis on the primary isolation plate. If the plate is inoculated with a suspension, it should be adjusted to provide just confluent growth over the surface of the plate. With sterile forceps or single disc dispenser place the **BD BBL Taxo** A disc in the center of the inoculated area. If the organism is streak-inoculated, the disc should be placed in the primary streak area or at the junction of the primary and secondary streak area. For further information regarding the use of the **BD BBL Taxo** A disc with **ssA**, consult **ssA** product literature.
2. Incubate plate(s) in ambient air (or in an atmosphere enriched with 5 to 10% CO₂) at 35 to 37 °C for 18 to 24 hours.
3. Observe plate for presence of a zone of growth inhibition around the bacitracin disc.

User Quality Control: At the time of use, check performance with pure cultures of stable control organisms producing known, desired reactions. The use of *Streptococcus pyogenes* ATCC® 12384 is recommended to demonstrate zone formation. One or more beta-hemolytic streptococcal species belonging to groups B, C, D and/or G may be employed to demonstrate lack of zone formation. Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A zone of inhibition is formed around the **BD BBL Taxo A** disc if the organism is a group A streptococcus. It is recommended that any zone of inhibition, regardless of diameter, be reported as "beta-hemolytic *Streptococcus*, presumptively group A by bacitracin."⁸
No zone of inhibition (growth up to the edge of the disc) is reported as "beta-hemolytic *Streptococcus*, presumptively not group A by bacitracin."⁸

LIMITATIONS OF THE PROCEDURE

The **BD BBL Taxo A** disc test is presumptive, and a positive result should be followed with more specific physiological and/or serological tests.⁹
It should be noted that an excessive inoculum concentration may result in the absence of a zone of inhibition (false-negative) with some group A streptococci. It has been reported that 6% of group B and 7.5% of groups C and G streptococci may produce zones of inhibition (false-positive result).¹⁰ Also, some streptococci may fail to grow in the absence of CO₂.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BD BBL Taxo A** Discs are tested for performance characteristics. Representative samples of the lot are assayed for bacitracin content using a microbiological assay procedure. The acceptable minimal concentration of bacitracin per disc is determined.
Additionally, plates of **Trypticase Soy Agar** with 5% Sheep Blood are swab-inoculated with cultures of *Streptococcus pyogenes* group A (A451), *S. pyogenes* group A (ATCC 12384), *S. pyogenes* group A (ATCC 10389), *S. pyogenes* group A (ATCC 19615), *Streptococcus* group B (ATCC 12400), *Streptococcus* group B (ATCC 12386), *Streptococcus* group C (152), *Streptococcus* group D (1305) and *Streptococcus* group G (R58/2336) diluted to 10⁻¹. Representative samples of the lot are placed on the inoculated plates and incubated at 35 ± 2 °C for 18 to 24 h. The plates are then read for zone size around the **BD BBL Taxo A** disc. All strains of *S. pyogenes* group A show zones of inhibition around the disc. No zones are observed around the **BD BBL Taxo A** disc with *Streptococcus* groups B, C, D and G.

AVAILABILITY

Cat No.	Description
231040	BD BBL Taxo™ A , 50
231041	BD BBL Taxo™ A , 6 x 50
231042	BD BBL Taxo™ A , 50
231552	BD BBL Taxo™ A , 10 x 50

REFERENCES

1. Maxted, W.R. 1953. The use of bacitracin for identifying group A haemolytic streptococci. J. Clin. Pathol. 6:224-226.
2. Levinson, M.L., and P.F. Frank. 1955. Differentiation of group A from other beta-hemolytic streptococci with bacitracin. J. Bacteriol. 69:284-287.
3. Streamer, C.W., P.M. Williams, W.L. Wang, R.S. Johnson, C.D. McGuire, I.J. Abelow, and R.J. Glaser. 1962. Identification of group A streptococci. Bacitracin disc and fluorescent antibody techniques compared with the Lancefield precipitin method. Am. J. Dis. Child. 104:157-160.
4. Petran, E.I. 1964. Comparison of the fluorescent antibody and the bacitracin disk methods for the identification of group A streptococci. Am. J. Clin. Pathol. 41:224-226.
5. Moffet, H.L., H.G. Cramblett, and J.P. Black. 1964. Group A streptococcal infections in a children's home. 1. Evaluation of practical bacteriologic methods. Pediatrics 32:5-10.
6. Chitwood, L.A., M.B. Jennings, and H.D. Riley, Jr. 1969. Time, cost, and efficacy study of identifying group A streptococci with commercially available reagents. Appl. Microbiol. 18:193-197.
7. Evans, G.L., and T.E. O'Neill. 1984. Development of an improved selective medium for the isolation of group A streptococci from throat cultures, abstr. C-136, p. 259. Abstr. 84th Annu. Meet. Am. Soc. Microbiol. 1984.
8. Facklam, R.R. and R.B. Cary. 1985. Streptococci and aerococci, p. 154-175. In E.H. Lennette, A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. Am. Soc. Microbiol., Washington, D.C.
9. Baron, E.J., and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.
10. Facklam, R.R., J.F. Padula, L.G. Thacker, E.C. Wartham, and B.J. Sconyers. 1974. Presumptive identification of group A, B and D streptococci. Appl. Microbiol. 27:107-113.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
© 2016 BD. BD, the BD Logo and all other trademarks are property of
Becton, Dickinson and Company.