PROCEDURES DE CONTROLE DE QUALITE

I INTRODUCTION

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7. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
9. Isenberg, H.D., F.D. Schoenknecht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
10. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby Company, St. Louis.
11. Miller, J.M., and H.T. Holmes. 1999. Specimen collection, transport, and storage, p. 33-63. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
12. Farmer, J.J., III. 1999. *Enterobacteriaceae*: introduction and identification p. 442-458. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
13. Downes and Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods. 4th ed. American Public Health Association, Washington, D.C.
14. U.S. Pharmacopeial Convention, Inc. 2001. The U.S. pharmacopeia 25/The national formulary 20-2002. U.S. Pharmacopeial Convention, Inc. Rockville, Md.
15. Cunniff, P. (ed.). 1995. Official methods of analysis of AOAC International, 16th ed. AOAC International, Arlington, Va.
16. U.S. Food and Drug Administration. 1995. Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, Md.
17. Bopp, C.A., F.W. Brenner, J.G. Wells, and N.A. Strockbine. 1999. *Escherichia*, *Shigella*, and *Salmonella*, p. 459-474. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.

 Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA
 800-638-8663
www.bd.com/ds

 Benex Limited
 Pottery Road, Dun Laoghaire
 Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
 BD, BD Logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company. © 2013 BD