Part Number:	8820281JAA	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, BD BBL™ Taxo™ Taxo Nitrate	Rev from: 02	Rev to: 03
		Job Number:	745-18

Catalog Number: 231734

Blank (Sheet) Size: Length: 13.75" Width: 16.125"


Number of Pages: 6 Number of Sheets: 1

Page Size: Length: 13.75" Width: 5.375" Final Folded Size: 1.25" X 5.375"

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. VS-88-2028-1PR for material information.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Shanthi Sagayamary Ramu at 4:12 pm, Oct 10, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Sonia F. Thompson at 3:13 pm, Oct 10, 2018

Third Eye By: **THIRD EYE BY**
By Natalie Morio at 11:00 am, Oct 16, 2018

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INTENDED USE

BD BBL™ Taxo™ Differentiation Discs Nitrate are used for detecting nitrate reductase activity by anaerobes in the nitrate reduction test.

SUMMARY AND EXPLANATION

The nitrate disc test is a rapid qualitative procedure first developed by Wideman et al.¹ The test determines the ability of anaerobic bacteria to reduce nitrate to nitrite. The nitrate disc test is recommended for use in the differentiation and identification of anaerobic bacteria.^{1–4}

PRINCIPLES OF THE PROCEDURE

The nitrate reduction test detects the ability of bacteria to use nitrate as a final electron acceptor in their respiratory chain.⁴ In this reaction, nitrate is reduced to nitrite, which may be further reduced to nitrogen gas or ammonia.^{1,4} The reduction of nitrate to nitrite is detected by the addition of sulfanilic acid and N,N-dimethyl- α -naphthylamine that combine with nitrite to form a red compound, *p*-sulfobenzene-azo- α -naphthylamine (indicating a positive test). If no color change takes place after the addition of the reagents (negative reaction), zinc dust is added to determine if unreduced nitrate or products other than nitrite (i.e., ammonia, nitrogen, nitric oxide, hydroxylamine) are present. If nitrate has been reduced beyond nitrite, zinc ions will produce no color change (positive result). A positive zinc and a negative nitrite reaction presumptively indicates that nitrate was reduced to end products beyond the nitrite stage.^{1,4}

REAGENTS

BD BBL Taxo Differentiation Discs Nitrate are paper discs imprinted with the letters NIT and impregnated with a solution containing 40% potassium nitrate and 0.1% sodium molybdate.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BD BBL Taxo Differentiation Discs Nitrate at 2–8 °C.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use the product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens in sterile containers or with sterile swabs and transport immediately to the laboratory in accordance with recommended guidelines.^{2,3,5}

Process each specimen as appropriate for that specimen.^{2,3,5}

PROCEDURE

Material Provided: BD BBL Taxo Differentiation Discs Nitrate.

Materials Required But Not Provided: Nitrate Reagent A (Sulfanilic Acid), Nitrate Reagent B (N,N-dimethyl- α -naphthylamine), Nitrate Reagent C (zinc dust), forceps, inoculating loop or needle, anaerobic jar or chamber, Bunsen burner or incinerator, reduced blood agar plates, incubator (35 °C) and quality control organisms.

Test Procedure:

Nitrate Reduction Test

1. Streak a reduced blood agar plate with the test organism to obtain a confluent lawn of growth.
2. Using aseptic technique, place one BD BBL Taxo Nitrate Disc in the center of the inoculated area. Using flamed and cooled forceps, gently press down on the disc to ensure a good contact with the medium.
3. Invert and incubate plates anaerobically at 35 $\pm$ 2 °C for 40–48 h.
4. Remove plates from incubation and add one drop each of Nitrate Reagent A (sulfanilic acid) and Nitrate Reagent B (N,N-dimethyl- α -naphthylamine) to the disc.
5. Examine the disc for development of a red color within 3–5 min.
6. If results are negative (no color development), confirm the negative finding by sprinkling a small amount of Nitrate Reagent C (zinc dust) on the disc.
7. Examine the disc for development of a pink to red color in 5–10 min.

User Quality Control:

Identity Specifications – Round, off-white (ivory) 1/4" paper discs with "NIT" printed on both sides.

Cultural Response – Place BD BBL Taxo Differentiation Discs Nitrate on blood plates heavily inoculated with the test organisms. Incubate the plates under anaerobic conditions at 35 $\pm$ 2 °C for 40–48 h or longer if necessary for good growth. Use Nitrate Reagent A and Nitrate Reagent B to indicate reaction.

Organism	ATCC®	Nitrate Reduction Reaction
<i>Propionibacterium acnes</i>	11827	+
<i>Clostridium difficile</i>	9689	–

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Procedure	Observation	Interpretation of Reaction
Addition of Nitrate Reagents A and B to disc	Development of red color	Positive – Nitrate reduced
	No color change	Presumptive negative – Nitrate not reduced or reduced to products other than nitrite
Addition of Nitrate Reagent C (zinc dust) to negative test disc.	Development of pink to red color	Negative – Nitrate not reduced
	No color change	Positive – Nitrate reduced

LIMITATIONS OF THE PROCEDURE

The nitrate reduction test may be used as a diagnostic aid in the identification of bacteria. Perform additional biochemical and serological testing using pure cultures for complete identification. For further information, consult appropriate references.^{2,3,5}

The nitrate discs may turn a tan color from hemolysis and/or metabolism when testing rapidly growing organisms. Addition of test reagents may produce only a very subtle color change or no color change. If this reaction occurs, use another means of nitrate reduction testing (tube test).¹

Organisms that yield light or non-confluent growth may fail to produce a sufficient quantity of nitrate reductase and may give negative test results.

Use only fresh cultures. Reincubate cultures that produce little or no growth before adding reagents.¹

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BD BBL Taxo Differentiation Discs Nitrate are tested to verify specific product characteristics. Plates of BD Trypticase™ Soy Agar with 5% Sheep Blood are swab-inoculated with cultures of *Clostridium difficile* (ATCC 9689), *Propionibacterium acnes* (ATCC 11827) and *Veillonella parvula* (ATCC 10790). Samples are placed on the inoculated plates and incubated at 35 $\pm$ 2 °C under anaerobic conditions. After 48 h, the plates are removed from incubation and one drop each of Nitrate Reagent A (sulfanilic acid) and Nitrate Reagent B (N,N-dimethyl- α -naphthylamine) are added to the disc. The discs are examined for the development of red color within 3 to 5 min. If the initial results are negative (no color change), a small amount of Nitrate Reagent C (zinc dust) is added directly to the disc. After 10 min the discs are examined again. A positive reaction is indicated by no color change on the disc. *P. acnes* and *V. parvula* reduce nitrate and give a positive reaction. *C. difficile* does not reduce nitrate and gives a negative reaction (pink to red color) after the addition of the zinc dust.

AVAILABILITY

Cat. No.	Description
231734	BD BBL™ Taxo™ Differentiation Discs Nitrate
261197	BD BBL™ Nitrate Reagent A Droppers
261198	BD BBL™ Nitrate Reagent B Droppers
261207	BD BBL™ Nitrate Reagent C

REFERENCES

- Wideman, P.A., D.M. Citronbaum, and V.L. Sutter. 1977. Simple disc technique for detection of nitrate reduction by anaerobic bacteria. J. Clin. Microbiol. 5: 315-319.
- Isenberg, H.D. (ed.) 1994. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
- MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore.
- Baron, E.J., L.R. Peterson and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Benex Limited