PROCEDURES DE CONTROLE DE QUALITE

I INTRODUCTION

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XIV REFERENCES

1. Warren, J.R., and B. Marshall. 1983. Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* *i*:1273-1275.
2. Rauws, E.A.J., W. Langenberg, H.J. Houthoff, H.C. Zanen, and G.N.J. Tytgat. 1988. *Campylobacter pyloridis*-associated chronic active antral gastritis. *Gastroenterology* *94*:33-40.
3. Marshall, B.J., J.A. Armstrong, D.B. McGeachie, and R.J. Glancy. 1985. Attempt to fulfill Koch's postulates for pyloric campylobacter. *Med. J. Aust.* *142*:436-439.
4. Buck, G.E., W.K. Gourley, W.K. Lee, K. Subramanyam, J.M. Latimer, and A.R. DiNuzzo. 1986. Relation of *Campylobacter pyloridis* to gastritis and peptic ulcer. *J. Infect. Dis.* *153*:664-669.
5. Czinn, S.J., and H. Carr. 1987. Rapid diagnosis of *Campylobacter pyloridis*-associated gastritis. *J. Pediatr.* *110*:569-570.
6. McNulty, C.A.M., and R. Wise. 1985. Rapid diagnosis of *Campylobacter*-associated gastritis. *Lancet* *i*:1443-1444.
7. Bates, H.M. 1987. *C. pylori*-induced gastritis and peptic ulcers. *Lab. Management* *25*(11):27-31.
8. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
9. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* *17*:53-80.
10. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
11. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
12. Coudron, P.E., and D.F. Kirby. 1989. Comparison of rapid urease tests, staining techniques, and growth on different solid media for detection of *Campylobacter pylori*. *J. Clin. Microbiol.* *27*:1527-1530.
13. Nachamkin, I. 1995. *Campylobacter* and *Arcobacter*, p. 483-491. *In* P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
14. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
15. Westblom, T.U., E. Madan, J. Kemp, and M. Subik. 1988. Evaluation of a rapid urease test to detect *Campylobacter pylori* infection. *J. Clin. Microbiol.* *26*:1393-1394.

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 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

 Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD