
## Revisions

SO 0191-5


| Rev from | Rev to  | ECO #   |
|----------|---------|---------|
| J        | 2010/07 | 5394-10 |

### Notes:


- BD Cat. Number 444764
- Blank (Sheet) Size : Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): N/A


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number VS-0250 for Material Information
- Ink Colors: Printed two sides ☒ Yes      ☐ No  
 No. of Colors: 1      PMS# 199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

|                              |      |                                                                                                                                                           |                                                                                                                                                                       |   |
|------------------------------|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|
| Label Design                 | Date | COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION |  <b>Becton, Dickinson and Company</b><br>7 Loveton Circle<br>Sparks, MD 21152 USA |   |
| Proofer                      | Date |                                                                                                                                                           |                                                                                                                                                                       |   |
| Checked By                   | Date |                                                                                                                                                           |                                                                                                                                                                       |   |
| Part Number:<br><br>PP115JAA |      | Category and Description<br><br>Package Insert,<br>BACTEC PANTA Plus Kit                                                                                  | Sheet: 1 of 21                                                                                                                                                        | A |
|                              |      |                                                                                                                                                           | Scale: 1:1                                                                                                                                                            |   |

# **BD BACTEC™ PANTA™ PLUS Kit**

|            |        |       |            |         |         |
|------------|--------|-------|------------|---------|---------|
| English:   | pages  | 1 – 3 | Italiano:  | pagine  | 9 – 11  |
| Français : | pages  | 4 – 6 | Português: | páginas | 11 – 13 |
| Deutsch:   | Seiten | 6 – 8 | Español:   | páginas | 14 – 16 |

PP115JAA  
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati útastást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

**BACTEC™ PANTA™ PLUS** is an antimicrobial supplement composed of **BACTEC PANTA Supplement** and **BACTEC Reconstituting Fluid**. Principal use is with the **BACTEC 460TB System** instrument.

**BACTEC PANTA Supplement** is recommended as an addition to **BACTEC 12B Medium** during the primary isolation of mycobacteria and direct drug susceptibility testing from clinical specimens.<sup>1,2</sup> **PANTA Supplement** suppresses the normal flora associated with the digested and decontaminated clinical specimens.

## SUMMARY AND EXPLANATION

Since the **BACTEC** system utilizes a liquid medium, it is important to add an antimicrobial supplement to suppress growth of contaminating microorganisms which may survive the decontamination process of specimens. Contamination may be reduced by supplementing the medium with a mixture of antimicrobials prior to inoculation. **PANTA Supplement** which contains polymyxin B, amphotericin B, nalidixic acid, trimethoprim and azlocillin is available in a lyophilized form. The concentrations of antimicrobials have been selected after a careful evaluation of their effect on the growth of mycobacteria and their ability to reduce non-mycobacterial contamination.<sup>1,3</sup> Reconstituted **PANTA** is added into 12B vials prior to inoculation of the specimen (inoculation of the specimen must be performed within 2 h of the addition of **PANTA** to the 12B vial). The leftover **PANTA** should be frozen at -20°C or colder for up to 6 months.

The addition of **PANTA** should enhance primary isolation of mycobacteria from clinical specimens. Reconstituting Fluid (RF) contains a growth promoting substance POES (polyoxyethylene stearate).<sup>4</sup> Thus, it is recommended that when **PANTA** is not used, in the case of an aseptically collected specimen, 0.1 mL of RF alone be added to 12B Medium prior to inoculation of the specimen.

## PRINCIPLES OF THE PROCEDURE

Most of the non-sterile clinical specimens processed for mycobacterial inoculation contain bacterial flora other than mycobacteria. These bacteria sometimes survive the decontamination procedure which is generally used in a clinical lab. Since mycobacteria grow slower than the contaminating bacteria, these bacteria may overgrow the culture, especially in liquid medium. Carefully selected antimicrobials, which inhibit gram-positive bacteria and gram-negative bacteria and fungi but do not significantly inhibit mycobacteria, can be used to selectively inhibit contamination. At the same time a growth promoting substance, POES (polyoxyethylene stearate), which is present in Reconstituting Fluid, is added to the medium which minimizes the inhibitory effect of the antimicrobials on mycobacteria as well as promotes growth of mycobacteria, especially of slowly growing *Mycobacterium tuberculosis*.

## REAGENTS

The **BACTEC PANTA PLUS** Kit vials contain the following active ingredients prior to processing:

### **BACTEC PANTA Supplement**

Lyophilized antimicrobial supplement: each vial is formulated with:

|                      |              |
|----------------------|--------------|
| Polymyxin B .....    | 10,000 units |
| Amphotericin B ..... | 1,000 µg     |
| Nalidixic acid ..... | 4,000 µg     |
| Trimethoprim .....   | 1,000 µg     |
| Azlocillin .....     | 1,000 µg     |

### **BACTEC Reconstituting Fluid, 10 mL each**

Per mL of Processed Water:

|                                       |      |
|---------------------------------------|------|
| Polyoxyethylene Stearate (POES) ..... | 4 mg |
|---------------------------------------|------|

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>5-8</sup> and institutional guidelines should be followed in handling all items contaminated with blood or other body fluids.

**DO NOT USE** if lyophilized product shows evidence of moisture.

Prior to use, each vial should be examined for evidence of damage, deterioration or contamination. **DO NOT USE** any vial showing evidence of contamination. Vials that exhibit evidence of damage should be discarded. On rare occasions, the glass bottle may be cracked and the neck may break during removal of the flip-off cap or in handling.

Positive 12B Mycobacteria Medium culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which may build up due to microbial metabolism. Aseptic techniques and established precautions against microbial hazards should be used at all times. Procedures should be performed in a suitable biological safety cabinet in a room with an appropriate ventilation system as recommended by the CDC.<sup>9</sup> Use a proper protective gown, mask and gloves while handling specimens and cultures of potential pathogens. Follow the CDC and OSHA recommendations.<sup>9</sup>

To minimize the potential of leakage during inoculation, use syringes with permanently attached needles or securely fastened **Luer-Lok™** brand tips.

## Storage Instructions

Prior to reconstitution, **BACTEC PANTA** Supplement vials should be stored in a cool, 2 – 8°C, dry place out of direct sunlight. Once reconstituted, they may be stored at 2 – 8°C for a maximum of one week. For low volume users, small aliquots may be frozen at -20°C or colder so that a tube may be thawed for each use (up to 6 months). Do not refreeze after thawing. Protect from light.

The expiration date applies to the unopened lyophilized bottle when stored as directed.

## SPECIMEN COLLECTION

Follow the procedures described in the **BACTEC 12B Mycobacteria Medium** package insert (PP-116JAA) and in the **BACTEC 460TB System Product and Procedure Manual** (MA-0029).

## PROCEDURE

### Material Provided:

- BACTEC PANTA** Supplement
- BACTEC** Reconstituting Fluid

### Materials Required But Not Provided:

- BACTEC 12B Medium**
- BACTEC 460TB Instrument**
- 5 mL syringe
- Tuberculin syringe with permanently attached needle
- Disinfectant solution
- 70% isopropyl alcohol swab

Add 5 mL of the **BACTEC** Reconstituting Fluid to the lyophilized **PANTA** Supplement. Mix until completely dissolved. Add 0.1 mL of the reconstituted **PANTA** to each 12B Mycobacteria Medium vial before inoculation of the specimen. Inoculation of the specimen must be performed within 2 h of the addition of **PANTA** to the vial. Aseptically collected specimens may be inoculated without the addition of **PANTA** into 12B Mycobacteria Medium. Some laboratories prefer to inoculate such specimens into two 12B Medium vials, one with **PANTA**, and one without **PANTA**.

Prior to inoculation add 0.1 mL of **BACTEC** Reconstituting Fluid only if a specimen is inoculated without **PANTA**, or if a subculture is made into a 12B Mycobacteria Medium vial.

### Quality Control

**DO NOT USE** vials past their expiration date.

**DO NOT USE** vials that exhibit any cracks or defects; discard the vial in an appropriate manner.

Quality Control Certificates are provided with each **PANTA PLUS** Kit. Quality Control Certificates show test organisms utilized for this type of product.

**BACTEC PANTA** Supplement vials displaying deterioration or moisture and **BACTEC** Reconstituting Fluid vials displaying turbidity, contamination, or discoloration (darkening) should not be used.

Follow the procedures described in the **BACTEC 12B Mycobacteria Medium** package insert (PP-116JAA) and in the **BACTEC 460TB System Product and Procedure Manual** (MA-0029).

As part of the daily maintenance schedule, Performance Test Vials (**BACTEC** Performance Test Kit) should be tested on the **BACTEC 460TB** System instrument. If low readings are obtained, yet Performance Test Vials were tested properly, possible instrument problems may be indicated (see **BACTEC** PTK package insert PP-046JAA).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulation or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

After the addition of reconstituted **PANTA** Supplement into 12B Medium, the contamination is suppressed while viable mycobacteria, if present in the specimen, grow and can be detected by the **BACTEC** 460TB System instrument.

## LIMITATIONS OF THE PROCEDURE

**BACTEC PANTA** Supplement has been designed to minimize overgrowth of contaminating flora; however, when sputum specimens have been insufficiently or improperly decontaminated, overgrowth can occur. In addition, extreme transport times or inadequate specimen transport conditions may significantly impact specimen contamination rates. Contaminants resistant to the antimicrobials contained in **PANTA** Supplement may be encountered. Occasionally, some mycobacterial strains (*M. kansasii*) may show inhibition by **PANTA**.

Due to the nature of biological materials in media products and inherent organism variability, the user shall be cognizant of potential variable results in the recovery of certain organisms.

## EXPECTED VALUES

The contamination rate in the **BACTEC** 12B Medium for the isolation of mycobacteria from clinical specimens may be equal to that experienced in conventional solid media. If specimen processing is performed according to CDC recommendations,<sup>9</sup> the contamination rate is usually within 5% of the total vials inoculated as reported in several publications.<sup>10-12</sup>

If a higher contamination rate is observed, review the decontamination procedure of the specimen. Increasing the sodium hydroxide concentration (from 4% to 5% with NALC) may decrease the contamination rate. Increasing the concentration of **PANTA** by rehydrating it with 4 mL (instead of 5 mL) Reconstituting Fluid may also decrease the high contamination rate. Consult CDC recommendations<sup>9</sup> for changes to the process.

## PERFORMANCE CHARACTERISTICS

See the **BACTEC** 12B Mycobacteria Medium package insert (PP-166JAA) for details.

## AVAILABILITY

| Cat No. | Description                                                                                 |
|---------|---------------------------------------------------------------------------------------------|
| 444764  | <b>BATEC PANTA PLUS</b> Kit, lyophilized antimicrobial supplement ( <b>PANTA</b> ), 5 vials |
| 442004  | <b>BACTEC</b> 12B Medium, 4 mL, 100 per carton                                              |

## REFERENCES

1. Libonati, J.P. et al. Effect of increased inoculum size on recovery rate of mycobacteria in **BACTEC** 7H12 Middlebrook Medium. Abstract C223, A.S.M. Annual Meeting, St. Louis, MO, 1984.
2. Libonati, J.P. et al. Identification and direct drug susceptibility testing of mycobacteria by the radiometric method. Abstract C300, A.S.M. Meeting, St. Louis, MO 1984.
3. Siddiqi, S.H. et al. A new antimicrobial mixture (polymyxin B, amphotericin B, nalidixic acid, trimethoprim and azlocillin) for effective suppression of non-mycobacterial contamination during primary isolation of mycobacteria. Abstract U-35. A.S.M. Meeting, Washington, D.C., 1986.
4. Siddiqi, S.H. Studies to further enhance mycobacterial growth in radiometric 7H12 Medium. Abstract U-46. A.S.M. Annual Meeting, Las Vegas, NV, 1985.
5. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.
6. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
7. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4<sup>th</sup> ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
9. Kent, P.T. et al. Public health mycobacteriology, a guide for the level III laboratory. Centers for Disease Control, Division of Laboratory Training and Consultation. Atlanta, GA, 1985.
10. Libonati, J.P. et al. Evaluation of the recovery and contamination rates of mycobacteria in the **BACTEC** system using increased media volume, increased inoculum size and new antimicrobial supplement (**PANTA**). Abstract U37, A.S.M. Meeting, Washington, D.C. 1986.
11. Roberts, G.D. et al. Evaluation of **BACTEC** radiometric method for recovery of mycobacteria and drug susceptibility testing of *Mycobacterium tuberculosis* from acid-fast smear-positive specimens. J. Clin. Microbiol. 1983, 18:689-696.
12. Takahashi, H. et al. Detection and recovery of mycobacteria by a radiometric procedure. J. Clin. Microbiol. 1983, 17:380-381.

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Benex Limited