
Revisions

SO 0191-5

Rev from	Rev to	ECO #
0607	2010/07	5394-10

Notes:

- BD Cat. Number 442004
- Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number VS-0250 for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# 199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP116JAA		Category and Description Package Insert, BACTEC 12B Mycobacteria Culture Vials Middlebrook 7H12	Sheet: 1 of 31	A
			Scale: 1:1	

BD BACTEC™ 12B Mycobacteria Culture Vials Middlebrook 7H12

English: pages 1 – 5 Italiano: pagine 14 – 18
Français: pages 5 – 9 Português: páginas 18 – 22
Deutsch: Seiten 10 – 14 Español: páginas 23 – 27

PP116JAA
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratíte se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The qualitative **BACTEC™ 12B Mycobacteria Medium** is recommended for the culture and recovery of mycobacteria from clinical specimens, sputum, gastric, urine, tissue, mucopurulent specimens, other body fluids and other respiratory secretions, differentiation of the *Mycobacterium tuberculosis* complex from other mycobacteria and drug susceptibility testing of *M. tuberculosis*.

SUMMARY AND EXPLANATION

The **BACTEC Radiometric Technique** has been widely used for the rapid recovery of mycobacteria from sputum and other clinical specimens.¹⁻⁴ All types of clinical specimens, pulmonary as well as extra-pulmonary, can be processed for **BACTEC** primary isolation in a manner similar to the conventional procedures.^{2,5-7} Details are provided in the CDC Manual, Procedure for the Isolation and Identification of Mycobacteria.³

The sample to be tested is inoculated into one or more **12B Mycobacteria Medium** vials, with a syringe through the rubber septum, and incubated. The culture vial is periodically placed into the **BACTEC 460TB System** instrument for testing, which consists of aspiration of the head space gas and assay of its radioactive content. A positive reading indicates the presence of viable microorganisms in the vial.

PRINCIPLES OF THE PROCEDURE

BACTEC 12B Mycobacteria Medium is an enriched Middlebrook 7H9 broth base. If microorganisms are present in the test sample inoculated into the **BACTEC** vial, they utilize a ¹⁴C labeled substrate (fatty acid) present in the medium and release ¹⁴CO₂ into the atmosphere above the medium. When the vials are tested on the **BACTEC 460TB System** instrument, the gas is aspirated from the vial and the ¹⁴CO₂ radioactivity is determined quantitatively in terms of numbers on a scale from 0 to 999. These numbers are designated as the Growth Index (GI). The GI numbers are displayed by the **BACTEC 460TB System** instrument and are also printed along with the identifying rack and bottle numbers (100 GI units are approximately equal to 0.025 µCi). The daily increase in the GI is directly proportional to the rate and amount of growth in the medium.

If an inhibitory agent is introduced into the medium, inhibition of metabolism is indicated by reduced production of ¹⁴CO₂ when compared to a control having no inhibitory agent. This basic principle is applied for drug susceptibility testing and in differentiating *M. tuberculosis* complex from other mycobacteria.

The **BACTEC 460TB System** instrument must be used with a special TB hood when employed for Mycobacteriology. The TB hood provides HEPA filtered exhaust air and negative pressure in the test area. In addition, the TB hood is equipped with an ultraviolet light source in the test area. The unit is designed for automatic testing of vials and must not be used for inoculation or subculturing in place of a biological cabinet.

REAGENTS

The **BACTEC 12B Medium** culture vials contain the following reactive ingredients prior to processing:

Amount added per liter of Processed Water	
7H9 Broth	4.7 g
Casein Hydrolysate	1.0 g
Bovine Serum Albumin	5.0 g
Catalase	48,000 units
¹⁴ C-substrate	1,000 µCi

Composition may have been adjusted to meet specific performance requirements.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Pathogenic microorganisms including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁸⁻¹¹ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Perform all procedures, including processing of specimens, smear preparation, inoculum preparation, making dilutions, inoculations of media, subculturing, etc., in a suitable biological safety cabinet in a room with an appropriate ventilation system as recommended by the CDC.³ Use a proper protective gown, mask, and gloves while handling specimens and cultures of potential pathogens. Follow the CDC and OSHA recommendations.

Prior to use, each vial should be examined for evidence of damage, deterioration or contamination such as cloudiness, bulging or depressed septum. **DO NOT USE** any vial showing evidence of contamination. A contaminated vial could contain positive pressure. Vial contamination may not be readily apparent.

Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. Vials that exhibit evidence of damage should be discarded prior to inoculation. On rare occasions, the glass bottle may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases, the contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which may build up due to microbial metabolism. Inoculation and sampling should be performed in a biological cabinet, and appropriate protective clothing, including gloves and masks should be worn. See **Procedure** section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or securely fastened **Luer-Lok™** brand tips.

The rubber septum of a vial should be disinfected each time the vial is entered. Clean the top of each vial with disinfectant, followed by 70% alcohol. Use an appropriate disinfectant for cleaning the work area and swabbing the vial septum and other surfaces. The CDC states: "With so many disinfectants available, it is important to consult the product brochures to make certain that the disinfectant is bactericidal for mycobacteria."¹²

Working with a syringe and needle requires special attention and safety precautions. Use syringes with permanently attached needles or securely fastened **Luer-Lok** brand tips. Handle the syringe with appropriate caution. Do not put the cover of the needle back after use. Discard syringes in sharps containers and handle the disposal containers carefully. Refer to the NCCLS Approved Guideline M29, "Protection of Laboratory Workers from Occupationally Acquired Infections."⁸

Storage Instructions

BACTEC 12B Mycobacteria vials should be stored at 2 – 25°C in dry place.

The medium vials are stable at the recommended storage conditions until the labeled expiration date.

It is good practice to segregate radioactive material by storage in a cabinet accessible only to authorized personnel. The cabinet should be labeled with the words "Caution: Radioactive Materials."

Do not expose the medium to direct bright light. Always initially test each 12B vial on the **BACTEC** brand 460TB System instrument to establish 5 – 10% CO₂ atmosphere in the vial and to screen out vials with high background readings. If a vial reads GI 20 or more on initial screening, discard it.

SPECIMEN COLLECTION

Specimen collection and handling should be carried out according to the standard methods which are adequately reviewed in the CDC manual³ and in the **BACTEC 460TB** System Product and Procedure Manual (MA-0029).

PROCEDURE

Materials Provided:

BACTEC 12B Medium

Materials Required But Not Provided:

BACTEC PANTA™ PLUS Kit

BACTEC NAP Differentiation Kit

BACTEC 460TB Instrument

Tuberculin syringe with permanently attached needle

Disinfectant solution

70% isopropyl alcohol swab

Sterile needle with filter

Processing of the Specimen: Numerous procedures for digestion and decontamination have been recommended and are being followed in various laboratories. Some procedures are known to be compatible only with egg-based media and are not suitable for use with the **BACTEC** method. Sodium Hydroxide (NaOH) and N-acetyl-L-cysteine-NaOH (NALC) methods are recommended for the **BACTEC 460TB** System. Other methods may be compatible with the **BACTEC** procedure, but have insufficient data available to demonstrate this conclusively.

For specimen collected aseptically, no decontamination procedure is required. For details see the **BACTEC 460TB System Product and Procedure Manual**.

Inoculation of 12B Medium: Prior to inoculation, remove the flip-off cap from the **BACTEC** vial top and inspect for cracks, contamination, excessive cloudiness, and bulging or indented septum. **DO NOT USE** if any defect is noted.

Test all 12B Medium vials on the **BACTEC 460TB System** instrument to eliminate vials with high background (GI 20) readings and to establish a CO₂ enriched atmosphere in the vial.

12B Medium may be used as a stand-alone medium for the recovery of mycobacteria. However, for maximum recovery of mycobacteria, along with the 12B vial one tube of Lowenstein-Jensen (LJ), other modifications of LJ medium, or 7H10/7H11 (regular plate or a bi-plate) may also be inoculated.¹³ The decision to use 12B medium alone or with additional conventional media should be made after reviewing each institution's own experience and requirements, as well as after review of published data.^{4-7,14,15}

Contamination may be reduced by supplementing the medium with a mixture of antimicrobials prior to inoculation.¹⁶ **PANTA** Supplement which contains polymyxin B, amphotericin B, nalidixic acid, trimethoprim and azlocillin is available in a lyophilized form (see the **BACTEC 460TB System Product and Procedure Manual** for details). Inoculation of the specimen must be performed within 2 h of the addition of **PANTA** to the 12B Medium vial. The Reconstituting Fluid (RF) contains a growth promoting substance POES (polyoxyethylene stearate).¹⁷ Thus, it is recommended that when **PANTA** is not used, 0.1 mL of RF alone be added to 12B Medium prior to inoculation of the specimen.

Using a new sterile syringe with a permanently attached needle or **Luer-Lok** brand tip for each specimen, inoculate 0.5 mL of the well-mixed specimen concentrate into each 12B Medium vial. A higher volume will disturb the pH of the medium due to the high pH of the inoculum. However, specimens collected aseptically and not decontaminated can be inoculated with as much as 1.0 mL. Volumes greater than 1.0 mL will dilute the medium and may affect growth.

Disinfect the rubber septum of each vial with a proper disinfectant after inoculation or after sampling aliquots from an inoculated vial followed by cleaning with a 70% alcohol swab.

Inoculated vials should be incubated at 37° ± 1°C without shaking and should be tested on a **BACTEC 460TB System** instrument using 5 – 10% CO₂ in air. Test the vials every two to three days for the first two weeks and weekly thereafter, for a total of six weeks. The testing schedule may vary depending upon the work load of a laboratory.

If a specimen is suspected of containing mycobacteria which require an optimum temperature other than 37°C, (i.e., *M. marinum* and *M. chelonae* require 30°C), two sets of media should be inoculated, one set for incubation at 37°C and the second for a more optimum temperature.

Prepare a smear from a **BACTEC 12B Medium** vial at GI 50 – 100. If it is positive for AFB, report as culture positive pending identification. Subculture to a solid medium. Incubate the subculture at 37° ± 1°C.

TB Differentiation: When a 12B Medium vial GI reaches 50, perform the **BACTEC NAP** test following the recommended procedure. This procedure will differentiate TB Complex from mycobacteria other than tuberculosis (MOTT).

Speciation of mycobacteria should be carried out by the conventional techniques as recommended by the Centers for Disease Control and Prevention.³

Subculturing: Prior to sampling, subculturing, or staining, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

Drug Susceptibility Testing: Any drug susceptibility test for *M. tuberculosis* may be set up at GI 300 + 1 additional day of incubation or GI ≥500 (see the **BACTEC 460TB System Product and Procedure Manual** for details).

QUALITY CONTROL

DO NOT USE vials past their expiration date.

DO NOT USE vials that exhibit any cracks or defects; discard the vial in an appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates show test organisms, including ATCC[™] cultures specified in the NCCLS standard, *Quality Assurance for Commercially Prepared Microbiological Culture Media*.¹⁸

As part of the daily maintenance schedule, Performance Test Vials (**BACTEC Performance Test Kit**) should be tested on the **BACTEC** instrument. If low readings are obtained, yet Performance Test Vials were tested properly, possible instrument problems may be indicated (see **BACTEC PTK** package insert PP-046JAA).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

The instrument reads out directly on a scale of 0 – 999 called a Growth Index. No calculations are required. The Growth Index is a measure of the ¹⁴CO₂ radioactivity aspirated into the ion chamber and is directly proportional to the amount of active growth in the vial.

Current Federal Regulations (CLIA '88) require that instrument printouts of test results must be associated with the appropriate patient sample. These printouts must be retained for a period of at least two years. Local or state regulations may be more restrictive than the Federal requirements: please consult those regulations for additional compliance requirements.

Once the GI reaches 10 or more, the culture should be considered “presumptive H positive” and the vial should be tested daily. The AFB smear is important for the confirmation of presence of mycobacteria. Once AFB are observed on the smear made from a positive vial, the specimen may be reported as culture positive for AFB (identification pending). Negative cultures should be reported after six weeks of incubation.

The detection time for a positive culture varies considerably due to such factors as the type of specimen, number of organisms present in the specimen, the treatment status of the patient, and the incubation temperature. Procedural differences, especially the decontamination method, may play an important role in detection time. The average detection time for *M. tuberculosis* from clinical specimens (smear-negative and smear-positive) is reported to be 7 to 12 days. MOTT bacilli are generally detected two to three days earlier than *M. tuberculosis*.^{4,6,7,14,15}

The incidence of contamination varies from laboratory to laboratory, depending upon the decontamination procedure used. Contamination in the 12B isolation vial is usually lower than routinely observed in conventional media where a contamination rate of up to 5% is acceptable. Maintain a contamination record and calculate the contamination rate periodically. Review the processing protocol and adjust if needed. A contaminated 12B Medium vial may be reprocessed, refer to the **BACTEC 460TB System Product and Procedure Manual** for the procedure.

LIMITATIONS OF THE PROCEDURE

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BACTEC** vial. A contaminated sample will give a positive reading but will not indicate a clinically relevant positive sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, appropriate equipment checks, etc. Contamination may also result from improper maintenance of **BACTEC** needles and tubing, poor septum cleansing or malfunction of the **BACTEC** needle heating unit (see the **BACTEC 460TB Operation and Maintenance Manual**).

The GI level reflects the rate and amount of growth occurring in the medium but cannot be as accurately quantified as by colony counts on solid media. Colonial morphology and pigmentation characteristics can only be determined on solid media.

Due to the nature of biological materials in media products and inherent organism variability, the user shall be cognizant of potential variable results in the recovery of certain organisms.

PERFORMANCE CHARACTERISTICS

In a study by Tortoli et al., recovery of mycobacteria was compared in two automatic systems and one manual system. The **BACTEC 460** system using 12B Mycobacteria Medium yielded 201 isolates, in comparison with 190 isolates with the **BACTEC MGIT 960** system and 167 isolates with Lowenstein-Jensen medium. The average time to detection using the **BACTEC 460** and 12B medium was 14.8 days.¹⁹

TABLE 1. Recovery of mycobacteria by individual systems and combinations of systems^a

Mycobacterium or specimen	Total no. (%) recovered					
	All media	BACTEC MGIT 960 system	BACTEC 460 system	Lowenstein-Jensen medium	BACTEC MGIT 960 system + Lowenstein-Jensen medium	BACTEC 460 system + Lowenstein-Jensen medium
All	236	190 (80)	201 (85)	167 (71)	212 (90)	225 (95)
<i>M. tuberculosis</i> complex	169	149 (88)	153 (92)	124 (74)	158 (94)	160 (95)
<i>M. xenopi</i>	24	13 (54)	11 (46)	17 (71)	21 (87)	22 (92)
<i>M. avium</i> complex	22	21 (95)	22 (100)	16 (73)	21 (95)	22 (100)
Other MOTT ^b	21	7 (33)	15 (71)	10 (48)	12 (57)	21 (100)
Smear positive	228	122 (53)	121 (53)	103 (45)	125 (55)	127 (56)
Smear negative	108	68 (63)	80 (74)	64 (59)	87 (80)	98 (91)

^a More frequently occurring taxa are considered separately, as were smear-positive and smear-negative specimens.

^b MOTT, mycobacteria other than *M. tuberculosis*.

AVAILABILITY

Cat. No.	Description
442004	BACTEC™ 12B Medium, 4 mL, 100 per carton
444764	BACTEC™ PANTA™ Plus Kit, lyophilized antimicrobial supplement (PANTA), 5 vials; Reconstituting Fluid, 10 mL, 5 vials
442103	BACTEC™ NAP Differentiation Kit, 5 µg NAP per vial, 10 vials per carton
445203	BACTEC™ 460TB System Instrument
442102	BACTEC™ SIRE Drug Kit (100 tests)
442104	BACTEC™ Diluting Fluid (10 tests)
442146	Isoniazid, 1 vial
442139	PZA Test Medium (5 Tests)
442143	PZA Drug/Reconstituting Fluid (50 tests)

REFERENCES

1. Middlebrook, G. et al. Automated radiometric detection of growth of *Mycobacterium tuberculosis* in selective media. *Am. Rev. Respir. Dis.* 115:1066-1069, 1977.
2. Bannister, E.R. et al. Comparison of detection, identification and drug susceptibility testing of mycobacteria by using the BACTEC radiometric method and conventional methods. Abstract C260, A.S.M. Annual Meeting, Las Vegas, 1985.
3. Kent, P.T. et al. Public health mycobacteriology, a guide for the level III laboratory. Centers for Disease Control, Division of Laboratory Training and Consultation. Atlanta, GA, 1985.
4. Roberts, G.D. et al. Evaluation of BACTEC radiometric method for recovery of mycobacteria and drug susceptibility testing of *Mycobacterium tuberculosis* from acid-fast smear-positive specimens. *J. Clin. Microbiol.* 18:689-696, 1983.
5. Fadda, G. et al. Recovery and susceptibility testing of *Mycobacterium tuberculosis* from extrapulmonary specimens by the BACTEC radiometric method. *J. Clin. Microbiol.* 19:720-721, 1984.
6. Morgan, M.A. et al. Comparison of a radiometric method (BACTEC) and conventional culture media for recovery of mycobacteria from smear-negative specimens. *J. Clin. Microbiol.* 18:384-388, 1983.
7. Takahashi, H. et al. Detection and recovery of mycobacteria by a radiometric procedure. *J. Clin. Microbiol.* 17:380-381, 1983.
8. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.
9. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
10. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
11. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
12. Kent, P.T. et al. Public health mycobacteriology: Centrifugal efficiency and digestant toxicity. Abstract U-2. A.S.M. Annual Meeting, St. Louis, MO, 1984.
13. Libonati, J.P. et al. The role of solid media in the BACTEC system for mycobacterial isolation and identification. Abstract U-71. A.S.M. Meeting, Atlanta, GA, 1987.
14. Kirihaara, J. M. et al. Improved detection times of *Mycobacterium avium* complex and *M. tuberculosis* with the BACTEC radiometric system. *J. Clin. Microbiol.* 22:841-845, 1985.
15. Park, C.H. et al. Rapid recovery of mycobacteria from clinical specimens using automated radiometric technique. *Am. J. Clin. Path.* 81:341-345, 1984.
16. Siddiqi, S.H. et al. A new antimicrobial mixture (polymyxin B, amphotericin B, nalidixic acid, trimethoprim and azlocillin) for effective suppression of non-mycobacterial contamination during primary isolation of mycobacteria. Abstract U-35. A.S.M. Meeting, Washington, D.C., 1986.
17. Siddiqi, S.H. Studies to further enhance mycobacterial growth in radiometric 7H12 Medium. Abstract U-46. A.S.M. Annual Meeting, Las Vegas, NV, 1985.
18. National Committee for Clinical Laboratory Standards. 1996. Approved Standard M22-A2, Quality assurance for commercially prepared microbiological culture media, 2nd ed., NCCLS, Wayne, PA.
19. Tortoli, E., P. Cinchero, C. Piersimoni, M.T. Simonetti, G. Gesu, and D. Nista. 1999. Use of BACTEC MGIT 960 for recovery of mycobacteria for clinical specimens: multicenter study. *J. Clin Microbiol.* 37:3578-3582.

BD BACTEC 12B Mycobacteria Culture Vials Middlebrook 7H12

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BACTEC, Luer-Lok, MGIT and PANTA are trademarks of Becton, Dickinson and Company. ©2010 BD.
Benex Limited