
Revisions

SO 0191-5


Rev from	Rev to	ECO #
0607	2010/07	5394-10

Notes:


- BD Cat. Number 442103
- Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): N/A


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number VS-0250 for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# 199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP119JAA		Category and Description Package Insert, BACTEC NAP TB Differentiation Test Kit	Sheet: 1 of 21	A
			Scale: 1:1	

BD BACTEC™ NAP TB **Differentiation Test Kit**

English:	pages	1 – 4	Italiano:	pagine	10 – 12
Français :	pages	4 – 7	Português:	páginas	13 – 15
Deutsch:	Seiten	7 – 10	Español:	páginas	16 – 18

PP119JAA
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свяжитесь с местным представителем компании BD. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratíte se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті БД өкіліне жүгініп нұсқа алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BACTEC NAP TB Differentiation Test is a test to differentiate the *Mycobacterium tuberculosis* complex from other mycobacteria. Principal use is with the **BACTEC 460TB System** instrument. **BACTEC NAP** (p-nitro- α -acetyl-amino- β -hydroxy-propionophenone) TB Differentiation Test is recommended for use with **BACTEC 12B Medium** for rapid differentiation of the TB complex (*M. tuberculosis*, *M. bovis*, *M. africanum* and *M. microti*) from mycobacteria other than *M. tuberculosis* (MOTT). This test should be used in conjunction with bacterial morphology on smear and growth characteristics in **BACTEC 12B Medium**.

SUMMARY AND EXPLANATION

BACTEC 12B Medium does not allow observation of colonial morphology, however, the GI output can be characteristic for certain species of mycobacteria. *M. tuberculosis* and *M. bovis* exhibit slow growth with a 2 to 3 fold increase in the daily GI output while other mycobacteria (MOTT bacilli) may show daily increases as high as 10 fold. Visually, *M. tuberculosis* growth in **BACTEC 12B Medium** displays small specks or clumps with a clear medium. Other mycobacteria grow in very small particles and may give slight turbidity to the medium.

Smears made from a positive vial (GI = 100 or more) may give an indication of the type of mycobacteria growing in the medium. Presence of serpentine cords and clumps is a characteristic of *M. tuberculosis* (*M. bovis* and *M. kansasii* may also form cords). Other mycobacteria have some other characteristic morphology, such as *M. avium*/*M. intracellulare* which show branching with long strands of bacilli.

The **BACTEC NAP Test** further helps in a rapid differentiation of *M. tuberculosis*, the most important species of mycobacteria. This test is independent of the Niacin test. A sample of actively growing culture is inoculated into the NAP vial which contains a disc containing 5 μ g of NAP. The vial, along with the control, is incubated and growth in the presence of NAP is monitored on the **BACTEC 460TB System** instrument. Analysis of the daily GI readings within the next 3 to 5 days helps in the differentiation of the TB complex from other mycobacteria.

PRINCIPLES OF THE PROCEDURE

NAP (p-nitro- α -acetyl-amino- β -hydroxy-propionophenone), an intermediate compound in the synthesis of chloramphenicol, inhibits mycobacteria belonging to the TB Complex (*M. tuberculosis*, *M. bovis*, *M. africanum* and *M. microti*) almost completely, while other mycobacteria show either slight or no inhibition. When mycobacterial growth is inhibited in the presence of NAP, the evolution of $^{14}\text{CO}_2$ is also inhibited, as indicated either by no increase or a decrease in the GI. This effect on the GI is used as a tool for identification.

REAGENTS

The **BACTEC NAP Kit** vials contain a paper disk with the following active ingredient:

p-nitro- α -acetyl-amino- β -hydroxy-propionophenone 5 μ g

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood or other body fluids.

DO NOT USE a NAP Differentiation Test vial if it is obviously contaminated, discolored, or does not contain a disk. Moisture, especially at higher temperatures, causes rapid deterioration of NAP.

Aseptic techniques and established procedures against microbial hazards should be used at all times.⁵ Always use within a biological safety cabinet. Avoid any contact with the disc. After use, specimen containers and other contaminated material should be sterilized by autoclaving.

Prior to use, each vial should be examined for evidence of damage.

Vials that exhibit evidence of damage should be discarded prior to inoculation. On rare occasions, the glass bottle may be cracked and the neck may break during removal of the flip-off cap or in handling. Also on rare occasions a vial may not be sealed sufficiently. In both cases, the contents of the vials may leak or spill, especially if the vial

is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining: Before sampling it is necessary to release gas which may build up due to microbial metabolism. Aseptic techniques and established precautions against microbial hazards should be used at all times. Procedures should be performed in a suitable biological safety cabinet in a room with an appropriate ventilation system as recommended by the CDC.⁵ Use a proper protective gown, mask, and gloves while handling specimens and cultures of potential pathogens. Follow the CDC and OSHA recommendations.

To minimize the potential of leakage during inoculation, use syringes with permanently attached needles or securely fastened **Luer-Lok** brand tips.

The toxicity of this compound has not been studied. Its use should be accompanied by appropriate procedures to avoid exposure.

Storage Instructions

BACTEC NAP vials should be kept refrigerated at 2 – 8°C out of direct sunlight. Do not expose to heat or high temperature, as NAP is heat labile.

The expiration date applies to the unopened vial when stored as directed.

SPECIMEN COLLECTION

Actively growing cultures in **BACTEC** 12B Medium may be used to perform the NAP test. The test culture should be pure and free from non-mycobacterial contamination.

PROCEDURES

Materials Provided: Ten vials, each containing one NAP disk.

Materials Required But Not Provided: **BACTEC** 460TB Instrument, **BACTEC** 12B Medium, Tuberculin syringe with permanently attached needle, Disinfectant solution, Alcohol swab, **BACTEC** Diluting Fluid.

For details see the **BACTEC** 460TB System Product and Procedure Manual (MA-0029).

When a **BACTEC** 12B Medium vial shows a GI 10 or more, test the vial daily until the GI reaches 50 – 100. Ensure that the culture is not contaminated by performing a Gram stain or subculture. Homogenize the culture with a tuberculin syringe (permanently attached needle) and, using the same syringe, aseptically transfer 1 mL to the NAP vial. Swab the top of the NAP vial and the control (original culture) vial with an appropriate disinfectant, followed by cleaning with an alcohol swab. Test both vials on the **BACTEC** instrument immediately to purge with 5 – 10% CO₂. Disregard the readings. Incubate at 37°C ± 1°C or at the optimum temperature and test daily for the next 2 to 6 days. Record the daily GI of the NAP vial and the original vial which acts as a control.

If the GI is found to be more than 100, the culture should be diluted as follows prior to inoculating the NAP vial:

GI 50 – 100	no dilution
GI 101 – 200	0.8 mL into a fresh BACTEC 12B Medium vial
GI 201 – 400	0.6 mL into a fresh BACTEC 12B Medium vial
GI 401 – 600	0.4 mL into a fresh BACTEC 12B Medium vial
GI 601 – 800	0.3 mL into a fresh BACTEC 12B Medium vial
GI 801 – 999	0.2 mL into a fresh BACTEC 12B Medium vial
GI 999 more than 1 day	0.1 mL into a fresh BACTEC 12B Medium vial

Do not let the culture grow beyond the peak GI value (999). If a culture is overgrown, subculture in a **BACTEC** 12B Medium vial, incubate and test daily until the desired GI is obtained.

If growth on a solid medium is to be tested, make a suspension by scraping a few representative colonies and emulsifying them in a tube of sterile water or **BACTEC** diluting fluid. Adjust the turbidity of the suspension so that it is comparable to the McFarland #1 standard. Inoculate 0.1 mL of this suspension into a **BACTEC** 12B Medium vial. Incubate and test daily on the **BACTEC** 460TB System instrument. Once the GI reaches 50 – 100, follow the procedure described above.

QUALITY CONTROL

The following positive and negative controls should be included in each batch of testing to ensure the accuracy of the test procedure.

<i>M. tuberculosis</i> , H37Rv	ATCC™ No. 27294
<i>M. kansasii</i>	ATCC No. 35775
<i>M. intracellulare</i>	ATCC No. 13950
<i>M. avium</i>	ATCC No. 35717

*For details, see the **BACTEC** 460TB System Product and Procedure Manual (MA-0029).

DO NOT USE vials past their expiration date.

DO NOT USE vials that exhibit any cracks or defects; discard the vial in an appropriate manner.

Quality Control Certificates are provided with each NAP Kit. Quality Control Certificates show test organisms utilized for this type of product.

As part of the daily maintenance schedule, Performance Test Vials (**BACTEC** Performance Test Kit) should be tested on the **BACTEC** instrument. If low readings are obtained, yet Performance Test Vials were tested properly, possible instrument problems may be indicated (see the **BACTEC** PTK package insert, PP-046JAA).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is

recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

INTERPRETATION OF RESULTS

The daily GI of the control (original culture vial) will continue to increase. In the NAP test vial, the decrease or increase depends upon the species of mycobacteria present. A decrease or unchanging GI indicates TB Complex, whereas an increase in GI indicates MOTT bacilli. The following criteria would apply to almost all mycobacteria in the NAP vial:

TB Complex (*M. tuberculosis*, *M. bovis*, *M. microti*, *M. africanum*)

- two consecutive significant decreases in GI after the inoculation
- slight, but not significant increase, in the first two days and then a decrease or no increase in GI

MOTT bacilli

- daily GI reading increases to over 400 within four days
- a slight decrease or no increase in the first one to three days after inoculation, and then two consecutive daily significant increases following day two

A “significant” increase or decrease means a 20% or more change from one day to the next.

LIMITATIONS OF THE PROCEDURE

Separation of TB Complex from MOTT by the BACTEC NAP TB Differentiation Test is possible with isolated, pure cultures only. A mixed culture of TB Complex and MOTT bacilli will result in an increase of GI, indicating MOTT bacilli only. Contamination with bacteria other than mycobacteria can also result in a false increase in the GI when the contaminant utilizes the substrate. NAP vials should be checked visually for any obvious contamination and/or turbidity. A smear for AFB from the NAP vial in the event of increasing GI is helpful in confirming the presence of MOTT bacilli and TB, or contaminating bacteria.

Occasional strains of *M. kansasii* may be inhibited by NAP to the extent that the results may be incorrectly reported. Further speciation is to be carried out by other methods.

Due to the nature of biological materials in media products and inherent organism variability, the user shall be cognizant of potential variable results in the recovery of certain organisms.

EXPECTED VALUES

The reporting time of results varies from 2 to 6 days depending on the mycobacterial species. Results should be reported only when a clear trend is achieved. If the GI is unchanged, identification as TB Complex should not be interpreted in less than 4 days. Average time for the completion of the NAP test has been reported to be 4 days.^{6,7-11}

NAP test results should be supported by growth characteristics and morphology of AFB on an acid-fast stained smear. The smear could also indicate the presence of mixed cultures or contamination.

PERFORMANCE CHARACTERISTICS

Although most NAP tests are easy to interpret, some tests require additional incubation time in order to demonstrate a clear trend. This is most important with the mycobacteria which are partially inhibited by NAP and show a slow or a delayed increase in GI or a longer lag phase in the presence of NAP (e.g., some strains of *M. kansasii*, *M. gastri*, *M. szulgai*, *M. terrae*, and *M. triviale*). All MOTT bacilli studied during the NAP evaluation demonstrated an increase in GI within 5 to 6 days.^{5,7} However, occasional strains of *M. kansasii* may be inhibited more than usual by NAP and could be misidentified by the test. Incubation at the optimal temperature is necessary to obtain correct results. Mycobacteria such as *M. chelonae* and *M. marinum*, which grow at lower than 37°C, and *M. xenopi* which grows at a higher temperature, may give incorrect results if tested at 37°C. In cases where the optimum incubation temperature is not certain, the test should be performed at more than one temperature. In previous studies, overall sensitivity and specificity of the test were estimated to be 98 – 99%.^{7,10}

AVAILABILITY

Cat. No.	Description
----------	-------------

442103	BACTEC NAP TB Differentiation Kit, 5 µg NAP per vial, 10 vials per carton
--------	---

442004	BACTEC 12B Medium, 4 mL, 100 per carton
--------	---

REFERENCES

1. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Kent, P.T. et al. Public health mycobacteriology, a guide for the level III laboratory. Centers for Disease Control, Division of Laboratory Training and Consultation. Atlanta, GA, 1985.
6. Gross, W.M. et al. Radiometric selective inhibition tests for differentiation of *Mycobacterium tuberculosis*, *Mycobacterium bovis*, and other mycobacteria. J. Clin. Microbiol., 21:565-568, 1985.

- Laszlo, A. et al. Evaluation of a rapid radiometric differentiation test for the *Mycobacterium tuberculosis* complex by selective inhibition with *p*-nitro- α -acetyl amino- β -hydroxy-propionophenone. J. Clin. Microbiol., 19:694-695, 1984.
- Libonati, J.P. et al. Identification and direct drug susceptibility testing of mycobacteria by the radiometric method. Abstract C300, A.S.M. Meeting, St. Louis, MO 1984.
- Morgan, M.A. et al. Evaluation of the *p*-nitro- α -acetyl amino- β -hydroxy-propionophenone differential test for identification of *Mycobacterium tuberculosis* complex. J. Clin. Microbiol., 21:634-635, 1985.
- Siddiqi, S.H. et al. Rapid radiometric methods to detect and differentiate *Mycobacterium tuberculosis* / *M. bovis* from other mycobacterial species. Amer. Rev. Respir. Dis., 130:634-640, 1984.
- Stager, C.E. et al. Identification of *Mycobacterium tuberculosis* from acid-fast smear-positive specimens by direct inoculation of the **BACTEC** NAP Test. Abstract C267, A.S.M. Annual Meeting, Las Vegas, NV, 1985.

BACTEC NAP TB Differentiation Test Kit

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BACTEC and Luer-Lok are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited