## Revisions

SO 0191-5

| Rev from | Rev to  | ECO #   |
|----------|---------|---------|
| 01       | 2011/01 | 5632-10 |

### Notes:


- BD Cat. Number 448057
- Blank (Sheet) Size: Length: 8.5"      Width: 14.0"  
 Number of Pages: 2      Number of Sheets: 1  
 Page Size: Length 8.5"      Width 14.0"      Final Folded Size: No Fold.
- Style (see illustrations below): # 1


#1

#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes      ☐ No  
 No. of Colors: 1      PMS# Standard Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

PC-Generated in-house

|                             |      |                                                                                                                                                           |                                                                                                                                                                       |   |
|-----------------------------|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|
| Label Design                | Date | COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION |  <b>Becton, Dickinson and Company</b><br>7 Loveton Circle<br>Sparks, MD 21152 USA |   |
| Proofer                     | Date |                                                                                                                                                           |                                                                                                                                                                       |   |
| Checked By                  | Date |                                                                                                                                                           |                                                                                                                                                                       |   |
| Part Number:<br><br>L009644 |      | Category and Description<br><br>Package Insert,<br>Phoenix PMIC/ID-60                                                                                     | Sheet: 1 of 3                                                                                                                                                         | A |
|                             |      |                                                                                                                                                           | Scale: N/A                                                                                                                                                            |   |

 **BD Phoenix™**  
PMIC/ID – 60

**English**

**INTENDED USE**

The **Phoenix™** PMIC/ID panel is used for the rapid identification and susceptibility testing of most aerobic and facultative anaerobic gram-positive bacteria of human origin with selected antimicrobial agents (for a complete listing of taxa, refer to the **Phoenix** System User's Manual). This panel is only for use with the **Phoenix** Automated Microbiology System instrument.

**PRINCIPLES OF THE PROCEDURE**

Refer to the **Phoenix** System User's Manual.

**PANELS**

*Identification*: Refer to the **Phoenix** System User's Manual.

*Susceptibility*: The PMIC/ID panel contains the antimicrobial agents and concentrations in doubling dilutions found in Table 1.

**STORAGE**

Store at 15-25°C. Do not use panel if the pouch is punctured or opened. Do not use the panel if desiccant is missing or if the desiccant pouch is torn.

Panels must be used within 2 h of being removed from the pouch.

**Precautions:**

For *in vitro* Diagnostic Use.

**WARNINGS**

Observe established precautions against microbiological hazards throughout all procedures. "Standard Precautions"<sup>1,2</sup> and institutional guidelines should be followed in handling all items contaminated with specimens and microorganisms. Prior to discarding, sterilize specimen containers and other inoculated materials by autoclaving.

**QUALITY CONTROL**

*Identification*: See expected results below.

| Test Organism                               | Expected Result              |
|---------------------------------------------|------------------------------|
| <i>Staphylococcus aureus</i><br>ATCC™ 29213 | <i>Staphylococcus aureus</i> |
| <i>Enterococcus faecalis</i><br>ATCC 29212  | <i>Enterococcus faecalis</i> |

*Susceptibility*: See expected results in Table 1.

NOTE: QC may be reported as less than or equal to the lowest, or greater than the highest concentration of the antimicrobic.

**LIMITATIONS OF THE PROCEDURE**

The clinical relevance of a specific antimicrobial and organism combination is determined by the activity of that antimicrobial against the organism and whether the antimicrobic is indicated for treatment of a disease state associated with that organism. The **Phoenix** System provides results for combinations, whether they are clinically relevant or not relevant. Overall, the results obtained using the **Phoenix** System compare favorably with the CLSI (formerly NCCLS) broth microdilution reference methodology;<sup>3</sup> however, Table 2 contains clinically relevant combinations that did not meet the strict standards of BD during clinical trials and will either not be reported or an alternate method is recommended for confirmation of the result. Similarly, Table 3 contains combinations which are not clinically relevant and are either not reported or an alternate method is recommended for confirmation of the result.

**PERFORMANCE CHARACTERISTICS**

Refer to the **Phoenix** System User's Manual.

**REFERENCES**

- Clinical and Laboratory Standards Institute. 2005. Approved guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, Pa.
- U.S. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. 1991. Occupational exposure to bloodborne pathogens. Federal Register 56:64175-64182.
- Clinical and Laboratory Standards Institute. 2009. Approved standard M7-A8. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 8th ed., CLSI, Wayne, Pa.

ATCC is a trademark of the American Type Culture Collection. BD, BD Logo and Phoenix are trademarks of Becton, Dickinson and Company. © 2011 BD.

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