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BIBLIOGRAFIA

1. Barton, A.L., and G.V. Doern. 1995. Selective media for detecting gastrointestinal carriage of vancomycin-resistant enterococci. *Diagn. Microbiol. Infect. Dis.* 23:119-122.
2. Jett, B.D., M.M. Huycke, and M.S. Gilmore. 1994. Virulence of enterococci. *Clin. Microbiol. Rev.* 7:462-478.
3. Moellering, R.C., Jr. 1992. Emergence of *Enterococcus* as a significant pathogen. *Clin. Infect. Dis.* 14:1173-1178.
4. Emori, T.G., and R.P. Gaynes. 1993. An overview of nosocomial infections, including the role of the microbiology laboratory. *Clin. Microbiol. Rev.* 6:428-442.
5. Landry, S.L., D.L. Kaiser, and R.P. Wenzel. 1989. Hospital stay and mortality attributed to nosocomial enterococcal bacteremia: a controlled study. *Am. J. Infect. Control* 17:323-329.
6. Subcommittee on Prevention and Control of Antimicrobial-Resistant Microorganisms in Hospitals. 1994. Preventing the spread of vancomycin resistance: a report from the hospital infection control practices advisory committee. *Fed. Regist.* 59:25758-25763.
7. MacFaddin, J.F. 1980. *Biochemical tests for identification of medical bacteria*, 2nd ed. Williams & Wilkins, Baltimore.
8. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
9. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
10. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
11. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000*, p. 0021-0045.
12. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
13. Isenberg, H.D., F.D. Schoenknecht, and A. von Graevenitz, 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Ruben. American Society for Microbiology, Washington, D.C.
14. Woodford, N., A.P. Johnson, D. Morrison, and D.C.E. Speller. 1995. Current perspectives on glycopeptide resistance. *Clin. Microbiol. Rev.* 8:585-615.

Manufacturer
Producent
Fabrikant
Valmistaja
Fabricant
Hersteller
Κατασκευαστής
Ditta produttrice
Fabrikant
Fabricante
Fabricante
Tillverkare

Use by
Anvendes før
Houdbaar tot
Viimeinkäyttöpäivä
A utiliser avant
Verwendbar bis
Ημερομηνία λήξης
Usare entro
Brukes for
Utilizar em
Usar antes de
Använd före

YYYY-MM-DD / YYYY-MM (MM = end of month)
AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä)
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
AAAA-MM-GG / AAAA-MM (MM = fine mese)
AAAA-MM-DD / AAAA-MM (MM = slutten av måneden)
AAAA-MM-DD / AAAA-MM (MM = fim do mês)
aaaa-mm-dd / aaaa-mm (mm = fin del mes)
AAAA-MM-DD / AAAA-MM (MM = slutet på månaden)

Catalog number
Katalognummer
Catalogusnummer
Tuotenumero
Numéro catalogue
Bestellnummer
Αριθμός καταλόγου
Numero di catalogo
Katalognummer
Número do catálogo
Número de catálogo
Katalognummer

Authorized Representative in the
European Community
Autoriseret repræsentant i EU
Erkend vertegenwoordiger in de Europese Unie
Valtuutettu edustaja Euroopan yhteisössä
Représentant agréé pour la C.E.E.
Autorisierte EG-Vertretung
Εξουσιοδοτημένος αντιπρόσωπος στην
Ευρωπαϊκή Κοινότητα
Rappresentante autorizzato nella Comunità
europea
Autorisert representant i EU
Representante autorizado na União Europeia
Representante autorizado en la Comunidad
Europea
Auktoriserad representant i EU

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
800-638-8663
www.bd.com/ds

In Vitro Diagnostic Medical Device
In vitro diagnostisk medicinsk anordning
Medisch hulpmiddel voor in vitro
diagnose
Lääkinnällinen in vitro -diagnostiikkalaitte
Dispositif médical de diagnostic in vitro
Medizinisches In-vitro-Diagnostikum
In vitro διαγνωστική ιατρική συσκευή
Dispositivo medico diagnostico in vitro.
In vitro diagnostisk medisinsk utstyr
Dispositivo médico para diagnóstico
in vitro
Dispositivo médico de diagnóstico
in vitro
Medicinsk anordning för in
vitro-diagnostik

Temperature limitation
Temperaturbegrensning
Temperatuurlimiet
Lämpötilarajoitus
Température limite
Zulässiger Temperaturbereich
Όριο θερμοκρασίας
Temperatura limite
Temperaturbegrensning
Limitação da temperatura
Limitación de temperatura
Temperaturbegrensning

Batch Code (Lot)
Batch kode (Lot)
Chargenummer (lot)
Eräkoodi (LOT)
Code de lot (Lot)
Chargencode (Chargenbezeichnung)
Κωδικός παρτίδας (Παρτίδα)
Codice del lotto (partita)
Batch-kode (Serie)
Código do lote (Lote)
Código de lote (Lote)
Satskod (parti)

Consult Instructions for Use
Læs brugsanvisningen
Raadpleeg gebruiksaanwijzing
Tarkista käyttöohjeista
Consulter la notice d'emploi
Gebrauchsanweisung beachten
Συμβουλευτείτε τις οδηγίες χρήσης
Consultare le istruzioni per l'uso
Se i brugsanvisningen
Consulte as instruções de utilização
Consultar las instrucciones de uso
Se brugsanvisningen

Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare, Ireland