
Revisions

SO 0191-5


Rev from	Rev to	ECO #
0107	2011/01	5662-11

Notes:


- BD Cat. Number 448077
- Blank (Sheet) Size: Length: 14" Width: 8.5"
 Number of Pages: 2 Number of Sheets: 1
 Page Size: Length 14" Width 8.5" Final Folded Size: No Fold.
- Style (see illustrations below): # 1


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS Standard Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

PC-Generated in-house

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L009736		Category and Description Package Insert, Phoenix PMIC-64	Sheet: 1 of 3	A
			Scale: N/A	

Phoenix™

PMIC-64

English

INTENDED USE

The **Phoenix™** PMIC panel is used for the susceptibility testing of most aerobic and facultative anaerobic gram-positive bacteria of human origin with selected antimicrobial agents (for a complete listing of taxa, refer to the **Phoenix** System User's Manual). This panel is only for use with the **Phoenix** Automated Microbiology System instrument.

PRINCIPLES OF THE PROCEDURE

Refer to the **Phoenix** System User's Manual.

PANEL

Susceptibility: The PMIC panel contains the antimicrobial agents and concentrations in doubling dilutions found in Table 1.

STORAGE

Store at 15-25°C. Do not use panel if the pouch is punctured or opened. Do not use the panel if desiccant is missing or if the desiccant pouch is torn.

Panels must be used within 2 h of being removed from the pouch.

Precautions:

For *in vitro* Diagnostic Use.

WARNINGS

Observe established precautions against microbiological hazards throughout all procedures. "Standard Precautions"^{1,2} and institutional guidelines should be followed in handling all items contaminated with specimens and microorganisms. Prior to discarding, sterilize specimen containers and other inoculated materials by autoclaving.

QUALITY CONTROL

Susceptibility: See expected results in Table 1.

NOTE: QC may be reported as less than or equal to the lowest, or greater than the highest concentration of the antimicrobic.

LIMITATIONS OF THE PROCEDURE

The clinical relevance of a specific antimicrobic and organism combination is determined by the activity of that antimicrobic against the organism and whether the antimicrobic is indicated for treatment of a disease state associated with that organism. The **Phoenix** System provides results for combinations, whether they are clinically relevant or not relevant. Overall, the results obtained using the **Phoenix** System compare favorably with the CLSI (formerly NCCLS) broth microdilution reference methodology;³ however, Table 2 contains clinically relevant combinations that did not meet the strict standards of BD during clinical trials and will either not be reported or an alternate method is recommended for confirmation of the result. Similarly, Table 3 contains combinations which are not clinically relevant and are either not reported or an alternate method is recommended for confirmation of the result.

PERFORMANCE CHARACTERISTICS

Refer to the **Phoenix** System User's Manual.

REFERENCES

- Clinical and Laboratory Standards Institute. 2005. Approved guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- U.S. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. 1991. Occupational exposure to bloodborne pathogens. Federal Register 56:64175-64182.
- Clinical and Laboratory Standards Institute. 2009. Approved standard M7-A8. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 8th ed. CLSI, Wayne, Pa.

ATCC is a trademark of the American Type Culture Collection. BD, BD Logo and Phoenix are trademarks of Becton, Dickinson and Company. © 2011 BD.

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Benex Limited