

# **BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials** **Soybean-Casein Digest Broth in a Plastic Vial**



R<sub>x</sub> Only



8089974(08)

2019-09

English

## **INTENDED USE**

BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials (prereduced enriched Soybean-Casein Digest broth with CO<sub>2</sub>) are for anaerobic blood cultures. Principal use is with the BD BACTEC fluorescent series instruments for the qualitative culture and recovery of anaerobic microorganisms from blood.

## **SUMMARY AND EXPLANATION**

The sample to be tested is inoculated into one or more vials which are inserted into the BD BACTEC fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO<sub>2</sub> produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO<sub>2</sub> present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

## **PRINCIPLES OF THE PROCEDURE**

If microorganisms are present in the test sample inoculated into the BD BACTEC vial, CO<sub>2</sub> will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO<sub>2</sub> are monitored by the BD BACTEC fluorescent series instrument. Analysis of the rate and amount of CO<sub>2</sub> increase enables the BD BACTEC fluorescent series instrument to determine if the vial is positive; i.e., that the test sample contains viable organisms.

## **REAGENTS**

The BD BACTEC Lytic/10 Anaerobic/F Culture Vials contain the following active ingredients prior to processing:

### **List of Ingredients**

|                                         |              |
|-----------------------------------------|--------------|
| Processed Water .....                   | 40 mL        |
| Soybean-Casein Digest Broth .....       | 2.75% w/v    |
| Yeast Extract .....                     | 0.2% w/v     |
| Animal Tissue Digest .....              | 0.05% w/v    |
| Dextrose .....                          | 0.2% w/v     |
| Hemin .....                             | 0.0005% w/v  |
| Menadione .....                         | 0.00005% w/v |
| Sodium Citrate .....                    | 0.02% w/v    |
| Thiols .....                            | 0.1% w/v     |
| Sodium Pyruvate .....                   | 0.1% w/v     |
| Saponin .....                           | 0.26% w/v    |
| Antifoaming Agent .....                 | 0.01% w/v    |
| Sodium Polyanetholsulfonate (SPS) ..... | 0.035% w/v   |

All BD BACTEC media are dispensed with added CO<sub>2</sub>. Anaerobic media are prereduced and dispensed with added CO<sub>2</sub> and N<sub>2</sub>. Composition may have been adjusted to meet specific performance requirements.

## **Warnings and Precautions**

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.**

Prior to use, each vial should be examined for evidence of damage, contamination or deterioration. Vials displaying evidence of damage or contamination such as leakage, cloudiness, discoloration (darkening), bulging or depressed septum should not be used.

A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. When using direct draw procedures, monitor the process closely to avoid refluxing materials into the patient.

On rare occasions, a vial may not be sealed sufficiently, which may result in the contents of the vial leaking or spilling. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or BD Luer-Lok™ brand tips.

Molecular tests performed on positive blood cultures will detect both viable and non-viable organisms commonly found in culture media. Therefore, Molecular test results should be evaluated in conjunction with Gram Stain results in accordance with standard-of-care practices as well as manufacturer's instructions for use.

### Storage Instructions

The BD BACTEC vials are ready for use as received and require no reconstitution or dilution. Store at 2–25 °C in a dry location, **out of direct light**.

### SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. Published studies have shown that the recommended specimen volume is 8–10 mL.<sup>5,6</sup> It is recommended that the specimen be inoculated into the BD BACTEC vials at bedside. Most commonly, a 10 cc or 20 cc syringe with a BD Luer-Lok brand tip is used to draw the sample. If appropriate, a BD Vacutainer® brand Needle Holder and a BD Vacutainer brand Blood Collection Set, BD Vacutainer Safety-Lok™ Blood Collection Set or other tubing “butterfly” set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 8–10 mL has been drawn, the flow should be stopped by crimping the tubing and removing the tubing set from the BD BACTEC vial. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. **The inoculated BD BACTEC vial should be transported as quickly as possible to the laboratory.**

### PROCEDURE

Remove the flip-off cap from BD BACTEC vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented stoppers. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **not** recommended). Aseptically inject or draw directly 8–10 mL of specimen per vial. If sample volumes of 3–4 mL are used, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated anaerobic vials should be placed in the BD BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the BD BACTEC fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period.

Positive vials will be determined by the BD BACTEC fluorescent series instrument and identified as such (see the appropriate BD BACTEC fluorescent series instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the BD BACTEC fluorescent series instrument can determine a difference in fluorescence.

Blood will lyse immediately upon addition to the BD BACTEC Lytic/10 Anaerobic/F Medium. The blood will appear chocolatized or very dark initially. If at the end of the testing period a BD BACTEC Lytic/10 Anaerobic/F vial is observed to have a bulging septum, it should be subcultured, Gram-stained or treated as presumptive positive.

Positive vials should be subcultured and a Gram-stained slide prepared. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to selective media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the BD BACTEC vials.

**Subculturing:** Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

For maximum yield of isolates, negative cultures may be checked by stain and/or subcultured at some point prior to discarding as negative.

### QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

**DO NOT USE** culture vials past their expiration date.

**DO NOT USE** culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Standard, *Quality Control for Commercially Prepared Microbiological Culture Media*.<sup>7</sup>

The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

*Clostridium perfringens* ATCC 13124

*Bacteroides fragilis*\* ATCC 25285

*Bacteroides vulgatus* ATCC 8482

*Streptococcus pneumoniae* ATCC 6305

*Escherichia coli* ATCC 25922

*Staphylococcus aureus* ATCC 25923

*Clostridium histolyticum* ATCC 19401

\*CLSI strain

For information on Quality Control for the BD BACTEC fluorescent series instrument, refer to the appropriate BD BACTEC fluorescent series instrument User's Manual.

## LIMITATIONS OF THE PROCEDURE

### Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the BD BACTEC vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical result. Such a determination must be made by the user based on such factors as type of organisms recovered, occurrence of the same organism in multiple cultures, patient history, etc.

### Recovery of SPS Sensitive Organisms from Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS (such as *P. anaerobius*), the presence of maximum volumes of blood (i.e., up to 10 mL) can help to optimize recovery of these organisms. To enhance the growth of SPS sensitive organisms when less than 8 mL of blood is inoculated, additional whole human blood may be added.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is 3.0 mL or less, an appropriate supplement may be required for recovery of these organisms. BD BACTEC FOS™ Fastidious Organism Supplement may be used as a nutritional supplement.

### Nonviable Organisms

A Gram-stained smear from culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.<sup>8</sup>

### General Considerations

Optimum recovery of isolates will be achieved by adding 8–10 mL of blood.<sup>5,6</sup> Use of lower or higher volumes may adversely affect recovery and/or detection times. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO<sub>2</sub> to be detected by the system or significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high. The default 5 day protocol was utilized for all analytical testing with this device and protocols longer than 5 days have not been evaluated.

## EXPECTED RESULTS

Performance of the BD BACTEC Lytic/10 Anaerobic/F medium in glass vials has been established by a number of published external clinical studies.<sup>9,10</sup> Seeded laboratory studies performed by BD have shown equivalent performance of the BD BACTEC Lytic/10 Anaerobic/F medium in plastic vials compared to BD BACTEC Lytic/10 Anaerobic/F medium in glass vials.<sup>11</sup>

A total of 342 paired sets at 10 to 100 CFU per vial were evaluated for recovery, with 100% recovering in both the BD BACTEC Lytic/10 Anaerobic/F medium contained in a plastic vial and the BD BACTEC Lytic/10 Anaerobic/F medium contained in a glass vial. This study included a diverse set of anaerobic and aerobic microorganisms frequently isolated in blood. The median time to detection (TTD) difference between the paired sets was 10 minutes, in favor of the BD BACTEC Lytic/10 Anaerobic/F medium contained in a plastic vial. Ninety-five percent of the TTD differences between the paired sets were between -1.68 hours faster for the glass vial and 3 hours faster for the plastic vial.

The following anaerobes were evaluated in the analytical studies: *Bacteroides fragilis*, *B. ovatus*, *B. thetaiotaomicron*, *B. vulgatus*, *Clostridium histolyticum*, *C. novyi*, *C. perfringens*, *Fusobacterium nucleatum*, *Porphyromonas asaccharolytica* (formerly *Bacteroides melaninogenicus* subsp. *asaccharolyticus*) and *Veillonella parvula*. The facultative anaerobe *S. pneumoniae* was also tested.

A subset of organisms, including *Fingoldia magna* (formerly *Peptostreptococcus magnus*) and *Peptoniphilus asaccharolyticus* (formerly *Peptostreptococcus asaccharolyticus*) were evaluated on the BD BACTEC FX instrument at 10 to 100 CFU per vial and demonstrated 100% recovery in both the BD BACTEC Lytic/10 Anaerobic/F medium contained in a plastic vial and the BD BACTEC Lytic/10 Anaerobic/F medium contained in a glass vial.

In microbial detection limit testing, a total of 312 paired sets at inoculum levels of 0 to 1 and 1 to 10 CFU per vial were evaluated. This study was designed to assess the capability of the BD BACTEC blood culture media tested to detect one CFU, when present. Of the 312 paired sets tested, 191 grew and detected in both devices and 44 did not detect in either. Twenty-nine (29) cultures grew and detected only in the BD BACTEC Lytic/10 Anaerobic/F medium contained in a glass vial. Forty-eight (48) cultures grew and detected only in BD BACTEC Lytic/10 Anaerobic/F medium contained in a plastic vial. One of 12 replicates of *Porphyromonas asaccharolytica* (ATCC 25260, 4 CFU per bottle) failed to detect in the BD BACTEC Lytic/10 Anaerobic/F medium contained in a plastic vial. Signal analysis demonstrated no evidence of growth in the replicate and a terminal subculture yielded no growth; indicating that there were likely no viable organisms inoculated into the vial.

## AVAILABILITY

### Cat. No. Description

442021 BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials

## REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.
5. Lin, H-H, *et al.* 2012. Evaluation of the blood volume effect on the diagnosis of bacteremia in automated blood culture systems. *Doi:10.1016/j.jmii. 2012. 03.2012*.
6. Reimer, L.G. *et al.* 1997. Update on detection of bacteremia and fungemia. *Clin. Micro. Rev.* 10:444–465.
7. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control of commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
8. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.). 2007. *Manual of clinical microbiology*, 9th ed. American Society for Microbiology, Washington, D.C.
9. Hollick, G.E., *et al.* 1996. *Diagnostic Microbiology and Infectious Disease Journal.* 24:191–196.
10. Rohner, P., *et al.* 1997. Advantage of combining resin with Lytic BACTEC blood culture media. *J. Clin. Micro.* 35:2634–2638.
11. Data available from BD Life Sciences.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or [bd.com](http://bd.com).

## Change History

| Revision | Date    | Change Summary                                                                                                                                                                                                                                                                                                                                                                             |
|----------|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (08)     | 2019-09 | Converted printed instructions for use to electronic format and added access information to obtain the document from <a href="http://BD.com/e-labeling">BD.com/e-labeling</a> .<br><br>In Warnings and Precautions section, added recommendation to perform molecular testings on positive blood cultures according to standard-of-care practices and manufacturer's instructions for use. |

US Customers only: For symbol glossary, refer to [bd.com/symbols-glossary](http://bd.com/symbols-glossary)



Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbicante / Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvođač / Tillverkare / Üretici / Виробник / 生产厂商



Use by / Използвайте до / Spotføjbejite do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейин пайдалануға / Naudokite iki / Izljetot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Използovať до / Použite do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)  
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)  
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce)  
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)  
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)  
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)  
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)  
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)  
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)  
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)  
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)  
 AAAA-MM-GG / AAAA-MM (MM = fine mese)  
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)  
 YYYY-MM-DD/YYYY-MM (MM = 월말)  
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)  
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)  
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)  
 AAAA-MM-DD / AAAA-MM (MM = slutten av månaden)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 AAAA-MM-DD / AAAA-MM (MM = fim do mês)  
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)  
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)  
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)  
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)  
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu)  
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)  
 YYYY-MM-DD / YYYY-MM (MM = 月末)



Catalog number / Каталоген номер / Katalógové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номер / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalógové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号



Authorized Representative in the European Community / Оторизирани представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret representant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizuirani predstavnik u Evropskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유럽 공동체의 위임 대표 / Įgaliotasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Avrupa Topluluğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 欧洲共同体授权代表



In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro diagnostisk iatrický suokceř / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiiniparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Diagnostiku / In vitro diagnosztikai orvosi eszköz / Dispositivo medical per diagnostica in vitro / Жасанды жағдайда жүргізетін медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietais / Medicinas ierīces, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispozitiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики in vitro / Medicinska pomůcka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Diagnostik Tibbi Cihaz / Медицинский прибор для диагностики in vitro / 体外诊断医疗设备



Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperaturlimit / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制



Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)



Contains sufficient for <n> tests / Съдържанието е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenu suffisant pour <n> test / <n> тесттері үшін жеткілікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Continut suficient pentru <n> teste / Достаточное для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli malmaze içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测



Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaityti naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i brugsanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se brugsanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明



Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Nepoužívejte opakovaně / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kulanmayın / Не використовувати повторно / 请勿重复使用



Serial number / Серийен номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / N° de serie / Seerianumber / Numéro de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Sėrijas numurs / Serie nummer / Numer seryjny / Número de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号



|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                   | For IVD Performance evaluation only / Само за оценка качество на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Μόνο για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réserve à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietaisų veikimo charakteristikoms tikrinti / Vienīgi IVD darbības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluering av IVD-ytelse / Tylko do oceny wydajności IVD / Uso exclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики in vitro / Určenie iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirmesi için / Тільки для оцінювання якості діагностики in vitro / 仅限 IVD 性能评估 |
|                   | For US: "For Investigational Use Only"                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|                   | Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūros žemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrense / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Мінімальна температура / 温度下限                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>CONTROL</b>    | Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / 对照                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>CONTROL +</b>  | Positive control / Положителен контрол / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitivna kontrola / Pozitív kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívá kontrolé / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивний контроль / 阳性对照试剂                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>CONTROL -</b>  | Negative control / Отрицателен контрол / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negatiivne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo / Негативті бақылау / 음성 컨트롤 / Neigiama kontrolė / Negatívá kontrolé / Negatieve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативний контроль / 阴性对照试剂                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>STERILE EO</b> | Method of sterilization: ethylene oxide / Метод на стерилизация: етиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringensmetode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseringensmetode: etüleenoksiid / Méthode de stérilisation : oxyde d'éthylène / Metoda sterilizácie: etilén-oxid / Sterilizálás módszere: etilén-oxid / Metoda di sterilizzazione: ossido di etilene / Стерилизация әдісі – этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksisds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringensmetode: etylenoksid / Metoda sterylizacji: tienek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metodá sterilizácie: etylenoxid / Metoda sterilizacije: etilen oksid / Steriliseringensmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷                                                                                                                                                         |
| <b>STERILE R</b>  | Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringensmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseringensmetode: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizácie: zračenie / Sterilizálás módszere: besugárzás / Metoda di sterilizzazione: irradiazione / Стерилизация әдісі – сәуле тисы / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringensmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metodá sterilizácie: ožiarenie / Metoda sterilizacije: ozračevanje / Steriliseringensmetod: strålning / Sterilizasyon yöntemi: ırradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射                                                                                                                                                                                                                                      |
|                   | Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologilised riskid / Risques biologiques / Biološki rizik / Biológiai veszélyes / Rischio biologico / Биологический тауекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisches Risiko / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Biologische опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Біологічна небезпека / 生物学风险                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                   | Caution, consult accompanying documents / Внимание, направете справка в придружаващите документи / Pozor! Pristudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatust! Lugeka kaasnevat dokumentatsiooni! / Attention, consulter les documents joints / Urozorenje, koristi prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, žiūrėkite priedamus dokumentus / Piesardziba, skatit pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vystraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увага: див. супутню документацію / 小心, 请参阅附带文档。                     |
|                   | Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Ausgējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Górná granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sıcaklık üst sınırı / Мінімальна температура / 温度上限                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|                   | Keep dry / Пазете сухо / Skladujte v suchém prostředí / Opbevarer tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausiai / Uzglabāt sausu / Droog houden / Holdes tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Бергerti від вологості / 请保持干燥                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|                   | Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid föravtagning / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamanı / Час забору / 采集时间                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|                   | Peel / Обелете / Otevfete zde / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skin / Húzza le / Staccare / Ўстиғи қабатын алып таста / 벗기 / Plęsti čia / Atfimt / Schillen / Trek av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odrhните / Oljuštiti / Dra isär / Ayırma / Відкрити / 撕下                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|                   | Perforation / Перфорация / Perforace / Perforering / Διάρθρωση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 찢취선 / Perforacija / Perforácia / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|                   | Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Inne kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használni, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Ереп пакет бұзылған болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用                                                                     |
|                   | Keep away from heat / Пазете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Сапқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не награвать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Fär ej utsättas för värme / Isidan uzak tutun / Бергerti від дії тепла / 请远离热源                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|                   | Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Reži / Vágja ki / Tagliare / Keciңiz / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciać / Cortar / Decupați / Отрезать / Odstřihnite / Iseći / Klipp / Kesme / Pozpizati / 剪下                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |



Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paemimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Dátum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期



µL/test / µL/тест / µL/Test / µL/εξέταση / µL/prueba / µL/teszt / µL/테스트 / мкл/тест / µL/tyrimas / µL/pårbaude / µL/teste / мкл/анализ / µL/检测



Keep away from light / Παзете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svjetla / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródła światła / Manter ao abrigo da luz / Feriți de lumină / Хранить в темноте / Uchovávajte mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Беретти від дії світла / 请远离光线



Hydrogen gas generated / Образуван е водород газ / Možnost úniku plynného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadrží hydrogen vodik / Hidrogén gázt fejleszt / Produzione di gas idrogeno / Газтөктес сутері пайда болды / 수소 가스 생성됨 / Išskiria vandenilio dujas / Rodas ūdenradis / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Osloбаda se vodonik / Genererad vätgas / Açığa çıkan hidrojen gazı / Реакция з виділенням водню / 会产生氢气



Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттің идентификациялық нөмірі / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificationnummer van de patiënt / Pasientens ID-nummer / Numer ID pacjenta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентифікатор пацієнта / 患者标识号



Fragile, Handle with Care / Чупливо, Работете с необходимото внимание. / Křehké, Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Εύθραυστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Őm, kásitsege ettevaatlukult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сыныш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkites atsargiai. / Trausls; rikoties uzmanīgi / Breekbaar, voorzichtig behandelen. / Ømtål, håndter forsigtig. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşıyın. / Тендітна, звертатися з обережністю / 易碎, 小心轻放

R<sub>x</sub> Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."



bd.com/e-labeling

KEY-CODE: 8089974

|                            |                      |                        |
|----------------------------|----------------------|------------------------|
| <b>Europe, CH, GB, NO:</b> |                      | <b>+800 135 79 135</b> |
| <b>International:</b>      |                      | <b>+31 20 794 7071</b> |
| AR                         | +800 135 79 135      | LT 8800 30728          |
| AU                         | +800 135 79 135      | MT +31 20 796 5693     |
| BR                         | 0800 591 1055        | NZ +800 135 79 135     |
| CA                         | +1 855 805 8539      | RO 0800 895 084        |
| CO                         | +800 135 79 135      | RU +800 135 79 135     |
| EE                         | 0800 0100567         | SG 800 101 3366        |
| GR                         | 00800 161 22015 7799 | SK 0800 606 287        |
| HR                         | 0800 804 804         | TR 00800 142 064 866   |
| IL                         | +800 135 79 135      | US +1 855 236 0910     |
| IS                         | 800 8996             | UY +800 135 79 135     |
| LI                         | +31 20 796 5692      | VN 122 80297           |



Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA



Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

#### Australian Sponsor:

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC® is a trademark of American Type Culture Collection.

BD, the BD Logo, BACTEC, FOS, Luer-Lok, Safety-Lok, and Vacutainer® are trademarks of Becton, Dickinson and Company or its affiliates. © 2019 BD. All rights reserved.