BD CLED Agar/MacConkey II Agar (Biplate)

NAMJENA

BD CLED Agar/MacConkey II Agar (Biplate) upotrebljava se za mikrobiološku analizu urina. CLED Agar diferencijalna je hranjiva podloga koja se koristi za izolaciju i brojenje bakterija u urinu. MacConkey II Agar selektivna je i diferencijalna podloga za izolaciju i diferencijaciju *enterobakterija* i mnoštvo drugih gram-negativnih sojeva iz kliničkih i nekliničkih uzoraka.

NAČELA I OBJAŠNJENJE POSTUPKA

Mikrobiološka metoda.

CLED Agar: 1960. godine Sandys je izvijestio o razvoju nove metode sprječavanja množenja bakterije *Proteus* na čvrstoj podlozi ograničavanjem elektrolita na hranjivoj podlozi koja je naknadno nekoliko puta izmijenjena za upotrebu u kulturi urina.¹⁻³ Označena je kao podloga s nedostatkom cistina, laktoze i elektrolita (CLED) te se preporučuje kao idealna za upotrebu u urinarnoj bakteriologiji.

Peptoni i goveđi ekstrakt snabdijevaju podlogu CLED Agar hranjivim tvarima. Laktoza je izvor energije za organizme koji je mogu iskoristiti mehanizmom fermentacije. Bromtimol plavo indikator je pH vrijednosti za diferencijaciju bakterija koje fermentiraju laktozu od onih koje je ne fermentiraju. Organizmi koji fermentiraju laktozu imat će snižene pH vrijednosti i promijenit će boju podloge iz zelene u žutu. Cistin omogućuje rast koliformnih "kolonija nanobakterija".³ Izvori elektrolita smanjeni su da bi se minimaliziralo množenje bakterije vrste *Proteus*.

MacConkey II Agar: jednu od najranijih formula podloge za diferencijaciju *enterobakterija* razvio je MacConkey te je objavio 1900. i 1905. godine.^{4,5} Ta formula temelji se na poznavanju činjenice da se žučne soli talože uslijed djelovanja kiselina, a da određeni enterični mikroorganizmi fermentiraju laktozu, dok drugi nemaju tu sposobnost. Ta je podloga naknadno izmijenjena nekoliko puta.^{6,7}

Podloga MacConkey Agar samo je blago selektivna budući da je koncentracija žučnih soli koje inhibiraju gram-pozitivne mikroorganizme niska u usporedbi s drugim enteričnim podlogama u petrijevim zdjelicama. Ova podloga preporučuje se za upotrebu s kliničkim uzorcima za koje se pretpostavlja da sadrže miješanu mikrobnu floru, kao što su urin i mnogi drugi uzorci, zbog toga što omogućuje preliminarno grupiranje *enterobakterija* i drugih gram-negativnih sojeva u one koji fermentiraju laktozu i u one koji je ne fermentiraju.⁷⁻¹⁰

Formula podloge MacConkey II Agar osmišljena je da poboljša inhibiciju množenja vrsta *Proteus* radi dobivanja bolje diferencijacije na one bakterije koje fermentiraju laktozu i one koje je ne fermentiraju i za izrazit rast.

Peptoni snabdijevaju podlogu MacConkey II Agar hranjivim tvarima. Kristalno ljubičasta boja inhibira gram-pozitivne bakterije, osobito enterokoke i stafilokoke. Diferencijacija enteričnih mikroorganizama postiže se kombinacijom laktoze i neutralnog crvenog pH indikatora. Stvaraju se bezbojne ili ružičaste do crvene kolonije, ovisno o sposobnosti izolata da fermentira ugljikohidrate.

Prisutnost podloga CLED Agar i MacConkey II Agar u toj petrijevoj zdjelici s dva dijela omogućuje određivanje ukupnog broja i izolaciju gram-pozitivnih i gram-negativnih bakterija iz uzoraka urina.

REAGENSI

CLED Agar/MacConkey II Agar (Biplate)

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REFERENCE

1. Sandys, G.H. 1960. A new method of preventing swarming of *Proteus* sp. with a description of a new medium suitable for use in routine laboratory practice. J. Med. Lab. Technol. 17:224-233.
2. Mackey, J.P., and G.H. Sandys. 1965. Laboratory diagnosis of infection of the urinary tract in general practice by means of a dip-inoculum transport medium. Br. Med. J. 2:1286-1288.
3. Mackey, J.P., and G.H. Sandys. 1966. Diagnosis of urinary infections. Br. Med. J. 1:1173.
4. MacConkey, A.T. 1900. Note on a new medium for the growth and differentiation of the *Bacillus coli communis* and the *Bacillus typhi abdominalis*. The Lancet, Part II:20.
5. MacConkey, A. 1905. Lactose-fermenting bacteria in faeces. J. Hyg. 5:333-379.
6. Levine, M., and H.W. Schoenlein. 1930. A compilation of culture media for the cultivation of microorganisms. The Williams & Wilkins Company, Baltimore.
7. MacFaddin, J.F. 1985. Media for isolation-cultivation- identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
8. Thomson, R.B., and J.M. Miller. 2003. Specimen collection, transport, and processing: bacteriology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
9. Gattermann, S., et al. 1997. Harnwegsinfektionen. In: Mauch, H, et al. (ed.). MiQ - Qualitätsstandards in der mikrobiologischen Diagnostik, vol. 2. Fischer Verlag, Stuttgart, Jena.
10. Gallien, R. 1988. Mikrobiologische Diagnostik in der ärztlichen Praxis – ein Leitfad. Fischer Verlag, Stuttgart.
11. Clarridge, J.E., M.T. Pezzlo, and K.L. Vosti. 1987. Cumitech 2A, Laboratory diagnosis of urinary tract infections. Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
12. Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
13. Mazura-Reetz, G., T.R. Neblett, and J.M. Galperin. 1979. MacConkey agar: CO₂ vs. ambient incubation, abstr. C 179, p. 339. Abstr. 79th Annu. Meet. Am. Soc. Microbiol. 1979.

AMBALAŽA/DOSTUPNOST

CLED Agar/MacConkey II Agar (Biplate)

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BD, logotip BD i svi drugi zaštitni znakovi u vlasništvu su tvrtke Becton, Dickinson and Company. © 2016 BD
Becton Dickinson GmbH