Part Number:	L011104	BALTSO0191 Version 12.0 Template 4 Inserts	
Category and Description:	Package Insert, BD Phoenix™ NMIC/ID-503 Panel	Rev from: 01	Rev to: 02
		Job Number:	362-18

Catalog Number: 449026

Blank (Sheet) Size: Length: 8.5" Width: 14.0"


Number of Pages: 2 Number of Sheets: 1

Page Size: Length: 8.5" Width: 14.0" Final Folded Size: No Fold

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 1

Vendor Printed: ☐ Online / In House Printed: ☒ Web Printed: ☐

See Specification control no. N/A for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Sonia Thompson at 1:47 pm, Feb 02, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 2:38 pm, Feb 02, 2018

Third Eye By: **THIRD EYE BY**
By Terrence Means at 2:37 pm, Feb 05, 2018

INTENDED USE

The **BD Phoenix™** NMIC/ID panel is used for the rapid identification and susceptibility testing of most aerobic and facultative anaerobic gram-negative bacteria of human origin with selected antimicrobial agents (for a complete listing of taxa, refer to the **BD Phoenix** System User's Manual). This panel is only for use with the **BD Phoenix** Automated Microbiology System instrument.

PRINCIPLES OF THE PROCEDURE

Refer to the **BD Phoenix** System User's Manual.

PANEL CONTENTS

Identification: Refer to the **BD Phoenix** System User's Manual.

Susceptibility: The NMIC/ID panel contains the antimicrobial agents and concentrations in doubling dilutions found in Table 1.

Confirmatory CPO Detect Test: Utilizes meropenem, doripenem, and cloxacillin in amounts required for the detection of carbapenemase-producing organisms.

STORAGE

Store at 15–25 °C. Do not use panel if the pouch is punctured or opened. Do not use the panel if desiccant is missing or if the desiccant pouch is torn.

Panels must be used within 2 h of being removed from the pouch.

Precautions:

For *in vitro* Diagnostic Use.

WARNINGS

Observe established precautions against microbiological hazards throughout all procedures. "Standard Precautions"^{1,2} and institutional guidelines should be followed in handling all items contaminated with specimens and microorganisms. Prior to discarding, sterilize specimen containers and other inoculated materials by autoclaving.

QUALITY CONTROL

Identification: See expected results below.

Test Organism	Expected Result
<i>Escherichia coli</i> ATCC® 25922	<i>Escherichia coli</i>
<i>Pseudomonas aeruginosa</i> ATCC 27853	<i>Pseudomonas aeruginosa</i>

Susceptibility: See expected results in Table 1.

NOTE: QC may be reported as less than or equal to the lowest, or greater than the highest concentration of the antimicrobic.

LIMITATIONS OF THE PROCEDURE

The clinical relevance of a specific antimicrobial and organism combination is determined by the activity of that antimicrobic against the organism and whether the antimicrobic is indicated for treatment of a disease state associated with that organism. The **BD Phoenix** System provides results for combinations, whether they are clinically relevant or not relevant. Overall, the results obtained using the **BD Phoenix** System compare favorably with the CLSI broth microdilution reference methodology;³ however, Table 2 contains clinically relevant combinations that did not meet the strict standards of BD during clinical trials and will either not be reported or an alternate method is recommended for confirmation of the result. Similarly, Table 3 contains combinations which are not clinically relevant and are either not reported or an alternate method is recommended for confirmation of the result.

PERFORMANCE CHARACTERISTICS

Refer to the **BD Phoenix** System User's Manual.

REFERENCES

- Clinical and Laboratory Standards Institute. 2014. Approved guideline M29-A4. Protection of laboratory workers from occupationally acquired infections, 4th ed., CLSI, Wayne, Pa.
- U.S. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. 1991. Occupational exposure to bloodborne pathogens. Federal Register 56:64175-64182.
- Clinical and Laboratory Standards Institute. 2015. Approved standard M07-A10. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 10th ed., CLSI, Wayne, Pa.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

ATCC is a trademark of the American Type Culture Collection.

© 2018 BD, BD, the BD Logo and BD Phoenix are trademarks of Becton, Dickinson and Company.

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