Part Number:	L000021	BALTSO0191 Version 13.0 Template 4 Inserts
Category and Description: Insert, BD BBL™ CultureSwab™	Rev from: 04	Rev to: 05
	Job Number:	246-17

Catalog Number: 220093, 220097, 220099, 220105, 220109, 220129, 220130, 220131, 220132, 220133, 220134

Blank (Sheet) Size: Length: 297 mm Width: 210 mm


Number of Pages: 52 Number of Sheets: 13

Page Size: Length: 210 mm Width: 148 mm Final Folded Size: 210 X 148 mm

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 5

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. COPAN for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By:

REVISED BY
By Tori Pagani at 11:32 am, Jan 16, 2019

Proofing
Approved By:

PROOFING APPROVED BY
By Nichole Graham at 1:33 pm, Jan 16, 2019

Third Eye By:

THIRD EYE BY
By Sonia F. Thompson at 2:08 pm, Jan 16, 2019

Liquid Stuart Medium, Liquid Amies Medium, and Cary-Blair Agar Gel Medium

For In Vitro
Diagnostic Use

R_x Only

STERILE R

0123

L000021(05) / HPC131
2019-01

English:	pages	1–3	Български:	16–18	Polski:	strony	31–33	Türkçe:	sayfalar	46–48
Français :	pages	4–6	Dansk:	side	Português:	páginas	34–36	Українська:	сторінки	49–51
Deutsch:	Seiten	7–9	Ελληνικά:	σελίδες	Русский:	страницы	37–39			
Italiano:	pagina	10–12	Magyar:	oldal	Slovenčina:	stránky	40–42			
Español:	páginas	13–15	Қазақша:	беттер	Svenska:	sidorna	43–45			

BD Swab Applicators – Technical Notes

Technical Note: BD swab applicators are manufactured using natural fibers that have not be treated with chemical additives, whitening agents or bleaches as these substances can compromise the viability of microorganisms and performance of the product. Because BD use natural fibers the appearance of the swab tip can be slightly yellow, this is perfectly normal and has no affect, whatsoever, on product performance or patient safety.

Sodium Thioglycollate – Technical Note

Gel medium formula contains Sodium Thioglycollate, an important component for the performance of the product and the maintenance of organism viability. Sodium Thioglycollate has a natural sulfur-like odor.

It may be possible to detect this sulfur odor momentarily when first opening the swab peel pouch. This odor is a perfectly normal and completely harmless characteristic.

Catalog Number	Transport Medium	Applicator Swab Type	Intended Use/Sampling Site*
220093	Amies Liquid	Regular Single Plastic Applicator	Mouth, Throat, Vagina, Wounds
220097	Cary Blair Agar Gel	Regular Single Plastic Applicator	Mouth, Throat, Vagina, Wounds
220099	Stuart Liquid	Regular Single Plastic Applicator	Mouth, Throat, Vagina, Wounds
220105	Amies Liquid	Two Regular Plastic Applicators	Mouth, Throat, Vagina, Wounds
220109	Stuart Liquid	Two Regular Plastic Applicators	Mouth, Throat, Vagina, Wounds
220129	Amies Liquid	Minitip Aluminum Wire	Eye, ENT, Urogenital, Pediatric
220130	Amies Liquid	Minitip Soft Aluminum Wire	Eye, ENT, Urogenital, Pediatric
220131	Amies Liquid	Minitip Flexible Wire	Eye, ENT, Urogenital, Pediatric
220132	Stuart Liquid	Minitip Aluminum Wire	Eye, ENT, Urogenital, Pediatric
220133	Stuart Liquid	Minitip Soft Aluminum Wire	Eye, ENT, Urogenital, Pediatric
220134	Stuart Liquid	Minitip Flexible Wire	Eye, ENT, Urogenital, Pediatric

ENT = Ear, Nose, Throat

*These are suggested sampling sites. Please refer to your GLP procedures to choose the most appropriate device for the specific sampling site.

INTENDED USE

BD BBL™ CultureSwab™ are sterile ready-to-use systems intended for the collection, transport and preservation of clinical specimens for bacteriological examinations.

Summary and Principles

One of the routine procedures for the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This can be accomplished using the BD BBL CultureSwab. Each culture swab unit comprises of a sterile peel pouch containing a swab applicator used to collect the sample, and a tube containing transport medium into which the swab applicator is placed after sampling. BD BBL CultureSwab is available with a range of different transport media: Liquid Stuart, Liquid Amies, and Cary-Blair Agar Gel. These transport medium are non-nutritious, buffered with phosphate, and provide a reduced environment, due to their formulation with sodium thioglycollate or thioglycollic acid (mercaptoacetic acid).^{1,2,3,4,5} Cary-Blair Transport Medium is specifically recommended for the collection and transport of fecal and rectal swab samples for the investigation of enteric pathogenic bacteria.^{6,7,8} Organisms in the sample material are protected from drying by moisture in the transport medium. The medium is designed to maintain the viability of organisms during transit to the laboratory. Once a swab sample is collected it should be placed in the tube of medium and transported to the laboratory as soon as possible and cultured onto an appropriate primary isolation medium (Blood Agar, Laked Blood Agar, MacConkey etc.).

The CultureSwab is available with different applicator shafts which facilitate the collection of specimens from various sites of the patient as described in the table above. For specific recommendations about collection of specimens for microbiology analysis and primary isolation techniques, consult the following references; Cumitech 9th, Manual of Clinical Microbiology¹⁰ and the Clinical Microbiology Procedures Handbook.¹¹

The transport tube has an hour glass shaped construction which serves two functions; in the case of agar gel medium it helps to keep the 6 cm deep agar gel medium column intact, and in the case of liquid medium products it helps to contain the sponge reservoir in position at the bottom of the tube. BD BBL CultureSwabs are individually packed in plastic film pouches in an outer metallic foil Vi-Pak. The plastic film pouch and metallic foil pack help to minimize the oxidation of medium and evaporation of water from the product during its shelf life.

Reagents - The nominal formula for each medium is as follows:

Liquid Stuart Transport Medium

Sodium Glycero-phosphate 10.0 g
Calcium Chloride 0.1 g
Mercaptoacetic Acid 1.0 ml
Distilled Water 1 liter

Liquid Amies Transport Medium

Sodium Chloride 3.0 g
Potassium Chloride 0.2 g
Calcium Chloride 0.1 g
Magnesium Chloride 0.1 g
Monopotassium Phosphate 0.2 g
Disodium Phosphate 1.15 g
Sodium Thioglycollate 1.0 g
Distilled Water 1 liter

Cary-Blair Agar Gel Transport Medium

Disodium Hydrogen Phosphate 1.1 g
Sodium Thioglycollate 1.5 g
Sodium Chloride 5.0 g
Calcium Chloride 0.09 g
Bacteriological Agar 5.6 g
Distilled Water 1 liter

Technical Notes

BD swab media contains Sodium Thioglycollate, an important component for the performance of the product and the maintenance of organism viability. Sodium thioglycollate has a natural sulfur-like odor. It may be possible to detect this sulfur odor momentarily when first opening the swab peel pouch. This odor is perfectly normal and is completely harmless. From time to time the medium containing tube may demonstrate some yellow coloration to varying degrees. This coloration is natural and a well-known phenomenon which is associated with the medical grade polypropylene used and the process of ionizing irradiation. The coloration has no adverse effect on the quality or performance of the product.

BD swab applicators are manufactured using natural fibers that have not been treated with chemical additives, whitening agents or bleaches as these substances can compromise the viability of microorganisms and performance of the product. Because BD uses natural fibers, the appearance of the swab tip can be slightly yellow. This is perfectly normal and has no effect, whatsoever, on product performance or patient safety.

Precautions

1. ⚠ This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
2. For *In vitro* Diagnostic Use.
3. BD BBL CultureSwab is certified as a Class IIa device under the classification terms of the European Medical Device Directive EC 93/42. In particular, the swab applicator is qualified for short transient contact with the patient in order to collect a specimen. This short contact is made with the patient's external or internal surfaces via body orifices such as the nose, throat, vagina or surgical wounds.
4. When collecting swab samples from patients, care should be taken not to use excessive force or pressure which may result in breakage of the swab shaft.
5. The fiber attachment to the applicator stick is qualified to withstand short transient contact with the patient in order to collect the sample. Prolonged contact must be avoided as this may result in the detachment of fiber.
6. Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
7. When the swab sample is cultured in the laboratory, and the procedure necessitates that the applicator(s) be placed in a tube of culture broth, great care should be taken in detaching the applicator stick from the cap to eliminate any risk of splashes or aerosols. If it is necessary to cut the applicator stick, sterile scissors should be used to facilitate a safe and clean break.
8. Observe aseptic techniques when using the product.
9. It must be assumed that all specimens contain infectious micro-organisms; therefore, all specimens should be handled with the appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.
10. Swab sample processing should be performed inside a protective safety cabinet or under a protective hood. Protective laboratory clothing and eyeglasses should be worn at all times when processing culture swab samples.
11. The product must be used as directed. It must not be subjected to any additional chemical or physical sterilization or micro-cidal or micro-static processes prior to use as this will compromise the performance and function of the product.
12. Certain fiber swabs and transport medium are known to interfere or be incompatible with certain diagnostic test kits and assays. If the intention is to use any part of a BD BBL CultureSwab product with a third party test kit or assay, then the user or manufacturer of such third party test kits or assays must verify acceptability of the BD product or independently validate and qualify the use of the BD BBL CultureSwab with said test kit or assay.

Storage and Stability

Store CultureSwab at 5–25 °C. Do not freeze or overheat. Do not use after the expiration date which is clearly printed on the outer box, each pack of swabs, each individual sterile swab pouch and the specimen transport tube label. If product is stored incorrectly, it can compromise the performance and invalidate the product specifications and performance claims.

Product Deterioration

The contents of unopened or undamaged units are guaranteed sterile. Do not use if they show evidence of damage, dehydration or contamination. Do not use if the expiration date has passed.

Materials Supplied

50 units of sterile CultureSwab are contained in a metallic foil park. Each individual swab pouch contains an applicator and a plastic tube containing transport medium.

Materials Required But Not Supplied

Appropriate materials for isolating, differentiating and culturing bacteria. These materials include culture media plates or tubes, and incubation systems, gas jars or workstations.

Directions for Use

The directions for use are printed on each BD BBL CultureSwab unit along with descriptive diagrams. Directions for use are summarized as follows:

- a. Peel open the BD BBL CultureSwab sterile pouch at the point marked "Peel Here".
- b. Remove cap from transport tube.
- c. Remove applicator swab and collect specimen. During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.
- d. Place applicator swab in transport tube and replace cap firmly to completely seal.
- e. Record patient's name and information on tube label.
- f. Send specimen to the laboratory for immediate analysis.

Precaution - When collecting swab samples from patients, care should be taken not to use excessive force or pressure which may result in breakage of the swab shaft.

Quality Assurance

All raw materials, swab components and batches of finished product are subjected to rigorous quality control. As part of these test procedures, a panel of control organisms is used to test the performance of BD BBL CultureSwab. Certificates of sterility and quality assurance, which describe some of the QC procedures, are available on request. For those laboratories wishing to test the performance of transport swabs, a simple test protocol is described in the Quality Control section in Clinical Microbiology Procedures Handbook.¹¹

Results

The survival of bacteria in a transport medium depends on many factors. These include the type of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium. Specimens should be transported directly to the laboratory and cultured within 24 hours. Published studies have demonstrated that BD BBL CultureSwab with Liquid Stuart and Liquid Amies Transport Medium will maintain the viability of a range of aerobic bacteria for 24 hours.¹³⁻²⁰

Cary and Blair reported in field studies with clinical specimens that *Salmonella* and *Shigella* could be recovered after storage at room temperature for as long as 49 days. In a similar field trial, nonagglutinable *Vibrio* comma strains, Heiberg group III and IV, were isolated after 22 days.⁶ In survey of 162 routine fecal specimens collected in Cary-Blair Medium, strains of *Shigella* were recoverable for as long as 49 days at room temperature.⁷ *Salmonella* continued to be isolated in spite of the presences of *Proteus* and *Pseudomonas aeruginosa* for at least 45 days. In other specimens, not containing these contaminants, *Salmonella* were isolated for as long as 62 days. *Shigella sonnei* I was recovered for up to 34 days from a portion of the bowel placed in the transport medium. Using Cary-Blair Medium, Neuman et. al., reported that *Vibrio parahaemolyticus* survives for 35 days.⁸ The medium of Cary-Blair is also recommended for transporting specimens suspected of containing *Campylobacter jejuni*.¹² Optimal product performances are achieved by BD BBL CultureSwab packed in a barrier plastic peel pouch plus aluminium foil envelope.

Limitations

BD BBL CultureSwab Liquid Stuart Medium, Liquid Amies Medium, and Cary-Blair Agar Gel Medium are intended for the collection and transport of aerobic bacteriological samples only. Samples containing viruses, chlamydia or anaerobic bacteria should be collected and transported using alternative specific transport systems.


Transport media, staining reagents, immersion oil, glass slides and specimens themselves sometimes contain dead organisms visible upon Gram staining.

BD BBL CultureSwab is not validated for environmental sampling and sterility test.

References

1. Stuart R.D. The diagnosis and control of gonorrhea by bacteriological cultures. Glasgow M. J. 27:131-142, May 1946.
2. Moffett M. Young J.L. and Stuart R.D. Centralized gonococcus culture for dispersed clinics. Brit. M. J.: 2: 421-424. Aug. 28, 1948.
3. Stuart R.D., Toshach S.R. and Patsula T.M. The problem of transport of specimens for culture of gonococci. Cand. J. Health, 45: 73-83, February, 1954.
4. Stuart R. D. Transport Medium for specimens in Public Health Bacteriology. Public Health Reports 74: No. 5, 431-438, May, 1959.
5. Amies C.R. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296-300.
6. Cary S. G. and Blair E.B.: New transport medium for shipment of clinical specimens. J. Bact. 88: No. 1,96-98, July 1964.
7. Cary S. G., Matthew M.S.: Fussillo M.S. and Hawkins C. Survival of *Shigella* and *Salmonella* in a New Transport Medium. Am J. Clin. Path. 43: No.3,294-296. March 1965.
8. Neuman D. A., Benenson M. W., Hubster E and Thi Nhu Tuan Am. J. Clin Path. 57: 33-34, Feb.1971.
9. Isenberg H. D., Schoenkench F.D. and Von Graeventz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
10. Balows A., Hausler, Jr. W. J., Herrmann K.L., Isenberg H. D., Shadomy H.J. Manual of Clinical Microbiology. Fifth Edition. American Society for Microbiology, Washington DC, 1991.
11. Isenberg H. D. (Editor in Chief). Clinical Microbiology Procedures Handbook. American Society for Microbiology, Washington DC, 1992.
12. Kaplan R. L. *Campylobacter* p. 236. In Lennette E. H., Balows A., Hausler W. J. Jr. and Truant J. P. (ed) .Manual of Clinical Microbiology, Third Edition. American Society for Microbiology, Washington, DC, 1980.
13. Van Horn K., Toth C. and Wegienek J. Viability of aerobic microorganisms in four swab systems. Poster Session 249/C Abstract C-436.98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
14. Warshauer D. M. and J.T. Paloucek. Comparison of Bacterial Recovery Rates of Four Commercial Swab Transport Systems. Poster Session 249/C. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
15. Perry J. L. Inhibitory properties of a commercially available swab transport device. 97th General Meeting of American Society for Microbiology, Miami Beach, Florida, May 1997.
16. Perry J. L. and Ballou D. R. Inhibitory properties of a swab transport device. J. Clin. Micro, Dec 1997, 3367-3368.
17. Krepl J. The effect of different sponge (pledget) material on the survival of Group A Streptococcus (GAS) during swab transport. Abstracts from Conjoint Meeting on Infectious Diseases, Canadian Association for Clinical Microbiology and Infectious Diseases, October 1997.
18. Perry J. L. Evaluation of fastidious organism survival in swab transport systems. Abstracts of 39th Interscience Conference on Antimicrobial Agents and Chemotherapy (ICAAC), American Society for Microbiology, San Francisco, CA. September 1999.
19. Campos J. M., Ruthman L. and Tshimanga M. Survival of fastidious bacteria on specimen collection swabs stored at room temperature. Abstract C-150.100th General Meeting of the American Society for Microbiology, Los Angeles, CA. May 2000.
20. Toth C. and Van Horn K. Comparison of four swab transport systems for the recovery of aerobic microorganisms. Abstract C-153.100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.

CULTURE SWAB TRANSPORT SYSTEM HOW TO USE SWAB GUIDE

BD BBL CultureSwab

Liquid Stuart Medium, Liquid Amies Medium, and Cary-Blair Agar Gel Medium

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Becton Dickinson Pty Ltd.