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LITERATURA

1. Amies CR. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296–300.
2. Miller JM. A Guide to Specimen Management in Clinical Microbiology. Second Edition. American Society for Microbiology. Washington, DC. 1999.
3. Miller JM, Holmes HT. Specimen collection, transport, and storage. In: Manual of Clinical Microbiology. 6th ed. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. Washington, DC: ASM; 1995:19–20.
4. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2003. Quality Control of Microbiological Transport Systems; Approved Standard. M40-A Vol. 23 No. 34.
5. Sng E-H, Rajan VS, Teo K-L, Goh A-J. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. Sex Trans Dis. 1982; 9:74–78.
6. Sun Y, Taylor T, Williams L, Sautter RL. Comparison of bacterial viability using both the EZ brand collection and transport system with the Difco swab transport pack. Presented at: 96th ASM General Meeting. 1996; Washington DC. Abstract C35.
7. Arbiq JC, Forward KR, LeBlanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. Diagnostic Microbiology and Infectious Disease. 2000; 36:163–168.
8. Perry JL. Effects of temperature on fastidious organism viability during swab transport. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-55.
9. Wilson DA, Tuohy MS, Procop GW, Hall GS. Effects of storage on the recovery of bacteria from three swab transport systems: BD CulturESwab, BD Culturette and Starplex StarSwab II. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-61.
10. Arbiq J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies described in NCCLS document M40-P Quality Control of Microbiology Transport Devices. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-40.
11. Mitchell E, Berman M, Ginocchio CC. Evaluation of two new Liquid Stuart transport systems: Platinum StarSwab II (Starplex Scientific) and BBL CulturESwab (Becton Dickinson). 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-74.
12. Perry JL, Matthews JS. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-42.

13. Robinson A, Gruver ML. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-69.
14. Human RP, Jones GA. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology. 2004; New Orleans, LA. Abstract C-161.
15. Human RP, Jones GA. Evaluation of swab transport systems against a published standard. J Clin Pathol 2004; 57:762–763.
16. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P, Quality Control of Microbiology Transport Devices. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID). 2003; Glasgow, UK. Abstract P-652.
17. Isenberg HD, Schoenkencht FD, Von Graeventiz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, SJ. Rubin. American Society for Microbiology, Washington, DC, 1979.
18. Koneman EW, Allen SD, Janda WM, Schreckenberger PC and Winn, Jr. WC. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
19. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
20. Forbes BA, Sahm DF, Weissfeld AS. 1998. Bailey and Scott's Diagnostic Microbiology. 10th ed. Mosby, St. Louis, MO.
21. Isenberg HD. 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
22. Isenberg HD. 1998. Essential Procedures for Clinical Microbiology. Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances. ASM, Washington, DC.
23. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents; Approved Standard. H5-A3.
24. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH, eds. Manual of Clinical Microbiology. 7th edition. Washington, DC: ASM; 1999.
25. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. (1993). Wadsworth Anaerobic Bacteriology Manual, 5th ed. Star Publishing Company, Belmont, CA.
26. Marler LM, Siders JA, Allen SD. Direct Smear Atlas, A Monograph of Gram-Stained Preparations of Clinical Specimens. Lippincott Williams and Wilkins, 2001.
27. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. Journal of Medical Microbiology, 1991 Vol 35, Issue 2 103–106.
28. Spiegel CA, Amsel R, Holmes KK. Diagnosis of bacterial vaginosis by direct gram stain of vaginal fluid J. Clin Microbiol. 1983 Jul;18 (1):170–177.
29. Benavides MI, Moncada X, Rodriguez B, Castillo C. Gonococcal urethritis in men: clinical experience in 1978–1988. Rev Med Chil. 1992 Oct;120(10):1140–3.
30. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyian G, Grosskurth H, Mabey D. Rapid assessment in Rwandan refugee camps in Tanzania. Genitourin Med, 1997 Feb;73 (1):33–8.
31. Deceuninck G, Asamoah-Adu C, Khonde N, Pepin J, Frost EH, Deslandes S, Asamoah-Abu A, Bekoe V, Alary M. Improvement of clinical algorithms for the diagnosis of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* by the use of Gram-stained smears among female sex workers in Accra, Ghana. Sex Transm Dis. 2000 Aug;27 (7):401–10.
32. Isenberg HD. 1998. Essential Procedures for Clinical Microbiology. Chapter 2.1, Page 41. Gram Stain. ASM, Washington, DC.
33. Isenberg HD. 1998. Essential Procedures for Clinical Microbiology. Chapter 1.1, Page 27. Collection, Transport and Manipulation of Clinical Specimens. Procedure for streaking plates for primary isolation. ASM, Washington, DC.
34. Fleming D. Biological Safety: Principles and Practices. January 2000. ASM, Washington DC.
35. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jyrtext.htm>.
36. Richardson JH. Biosafety in Microbiological and Biomedical Laboratories. December 1994. Diane Publishing Company.
37. Hansen DJ. Healthcare, Laboratories and Biosafety. Vol 2., 1992. CRC Press.
38. Greenberg AE, Clesceri LS, and Eaton AD. 9215 heterotrophic plate count. In: Standard Methods for the Examination of Water and Waste Water. 18th ed. Washington, DC APHA; 1992: 9-33-9-34.
39. Washington JA. 1986. Rapid diagnosis by microscopy. Clin. Microbiol. Newsl. 8:135–137.
40. Van Horn KG, Rankin I. Evaluation and comparison of two Stuart's Liquid Swab transport systems tested by the NCCLS M40 method. 105th General Meeting of the American Society for Microbiology. 2005; Atlanta, Georgia. Abstract C-292.
41. Bourbeau PP, Heiter BJ. Validation of QC standard for bacteriological transport devices as specified in the NCCLS Proposed Standard M40: Quality Control of Microbiological Transport Systems. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-46.

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabbricante / Atқарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производитель / Výrobca / Proizvodač / Tillverkare / Üretici / Виробник / 生产厂商

Use by / Използвайте до / Spotřebujte do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použite do / Upotrebiti do / Använd före / Son kulanma tarihi / Використати до / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
AAAA-MM-DD / AAAA-MM (MM = fin del mes)
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)
AAAA-MM-GG / AAAA-MM (MM = fine mese)
ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА / (АА = айдың соңы)
YYYY-MM-DD/YYYY-MM (MM = 월말)
MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
AAAA-MM-DD / AAAA-MM (MM = sluten av måneden)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = fim do mês)
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталог номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumbr / Numéro catalogue / Kataloški broj / Katalogusszám / Numero di catalogo / Каталог номер / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus number / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperatuurlimiet / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топтама коды / 배치 코드(코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-kode (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (лот) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партии / 批号 (亚批)

Contains sufficient for <n> tests / Съдържа достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkelig til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> teszthez elegendő / Contenido suficiente per <n> test / <n> тесттері үшін жеткілікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder tilstrekkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conținut suficient pentru <n> teste / Достаточо для <n> тестов(а) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli miktarda içerir / Вистачить для аналізу: <n> / 足够进行 <n> 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanım Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Nepoužívejte opakované / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kulanmayın / Не використовувати повторно / 请勿重复使用

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisieringsmethode: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metoda di sterilizzazione: irradiazione / Стерилизация әдісі – сәулелі түрсі / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизації: облучение / Metoda sterilizácie: ožiarenie / Metoda sterilizacije: ozračevanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: irradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射

Identification number of notified body / Identification de l'organisme notifié / Identifizierung der benannten Stelle / Identificación del organismo notificado / Identificazione dell'organismo notificato / Identificação do organismo notificado / Identifikační číslo notifikovaného subjektu / Identifikasjonsnummer ved påvist organisme / Identifieringsnummer av certifierad myndighet

Peel / Обелете / Otefete zde / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстиғи қабатын алып таста / 벗기 / Plăști ăla / Atimēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odtrhnite / Oljušiti / Dra isår / Ayırma / Відклеїти / 撕下

Manufactured by:
Copan Italia SpA
Via F. Perotti, 10
25125 Brescia Italy

Distributed by
Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ESwab® is a trademark of Copan Group.

ATCC is a trademark of the American Type Culture Collection.

© 2019 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.