tekuća podloga Stuart, tekuća podloga Amies i podloga s agarom u gelu Cary-Blair

Za *in vitro*
dijagnostiku

STERILE R

0123

L000021(05) / HPC131

2019-01

Hrvatski

Aplikatori za bris BD – tehničke napomene

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
Reference

1. Stuart R.D. The diagnosis and control of gonorrhea by bacteriological cultures. Glasgow M. J. 27:131–142, May 1946.
2. Moffett M. Young J.L. and Stuart R.D. Centralized gonococcus culture for dispersed clinics. Brit. M. J.. 2: 421–424. Aug. 28, 1948.
3. Stuart R.D., Toshach S.R. and Patsula T.M. The problem of transport of specimens for culture of gonococci. Cand. J. Health, 45: 73–83, February, 1954.
4. Stuart R. D. Transport Medium for specimens in Public Health Bacteriology. Public Health Reports 74: No. 5, 431–438, May, 1959.
5. Amies C.R. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296–300.
6. Cary S. G. and Blair E.B.: New transport medium for shipment of clinical specimens. J. Bact. 88: No. 1,96–98, July 1964.
7. Cary S. G., Matthew M.S., Fussillo M.S. and Hawkins C. Survival of Shigella and Salmonella in a New Transport Medium. Am J. Clin. Path. 43: No.3,294–296. March 1965.
8. Neuman D. A., Benenson M. W., Hubster E and Thi Nhu Tuan Am. J. Clin Path. 57: 33–34, Feb.1971.
9. Isenberg H. D., Schoenkencht F.D. and Von Graeventiz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
10. Balows A., Hausler, Jr. W. J., Herrmann K.L, Isenberg H. D., Shadomy H.J. Manual of Clinical Microbiology. Fifth Edition. American Society for Microbiology, Washington DC, 1991.
11. Isenberg H. D. (Editor in Chief). Clinical Microbiology Procedures Handbook. American Society for Microbiology, Washington DC, 1992.
12. Kaplan R. L. Campylobacter p. 236. In Lennette E. H., Balows A., Hausler W. J. Jr. and Truant J. P. (ed) .Manual of Clinical Microbiology, Third Edition. American Society for Microbiology, Washington, DC, 1980.
13. Van Horn K., Toth C. and Wegienek J. Viability of aerobic microorganisms in four swab systems. Poster Session 249/C Abstract C-436.98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
14. Warshauer D. M. and J.T. Paloucek. Comparison of Bacterial Recovery Rates of Four Commercial Swab Transport Systems. Poster Session 249/C. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
15. Perry J. L. Inhibitory properties of a commercially available swab transport device. 97th General Meeting of American Society for Microbiology, Miami Beach, Florida, May 1997.
16. Perry J. L. and Ballou D. R. Inhibitory properties of a swab transport device. J. Clin. Micro, Dec 1997, 3367–3368.
17. Krepel J. The effect of different sponge (pledget) material on the survival of Group A Streptococcus (GAS) during swab transport. Abstracts from Conjoint Meeting on Infectious Diseases, Canadian Association for Clinical Microbiology and Infectious Diseases, October 1997.
18. Perry J. L. Evaluation of fastidious organism survival in swab transport systems. Abstracts of 39th Interscience Conference on Antimicrobial Agents and Chemotherapy (ICAAC), American Society for Microbiology, San Francisco, CA. September 1999.
19. Campos J. M., Ruthman L, and Tshimanga M. Survival of fastidious bacteria on specimen collection swabs stored at room temperature. Abstract C-150.100th General Meeting of the American Society for Microbiology, Los Angeles, CA. May 2000.
20. Toth C, and Van Horn K. Comparison of four swab transport systems for the recovery of aerobic microorganisms. Abstract C-153.100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.

SUSTAV ZA TRANSPORT BRISA KULTURE VODIČ ZA UPOTREBU ŠTAPIĆA ZA BRIS

Identification number of notified body / Identification de l'organisme notifié / Identifizierung der benannten Stelle / Identificación del organismo notificado / Identificazione dell'organismo notificato / Identificação do organismo notificado / Identifikační číslo notifikovaného subjektu / Identifikasjonsnummer ved påvist organisme / Identifieringsnummer av certifierad myndighet / Identifikacijski broj prijavljenog tijela

Method of sterilization: irradiation / Метод на стерилизација: ирадиација / Způsob sterilizace: záření / Steriliseringmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismeetod: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metoda di sterilizzazione: irradiazione / Стерилизација әдісі – сәуле түсіру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestråling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metóda sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: iradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射

Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Nepoužívajte opakovane / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullanmayın / Не використовувати повторно / 请勿重复使用

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperatuurlimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ochraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Peel / Обелете / Otevřete zde / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Үстіңгі қабатын алып таста / 벗기기 / Pléști čia / Atlîmēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odtrhnite / Oljuštiti / Dra isär / Аyrıрма / Відклеїти / 撕下

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvodač / Gyártó / Fabbicante / Атқарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производитель / Výrobca / Proizvodač / Tillverkare / Üretici / Виробник / 生产厂商

Manufactured by:
Copan Italia SpA
Via F. Perotti, 10
25125 Brescia Italy

Distributed by:
Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia