
Revisions

BALTSO0191 Version 11.0 Template 4
Inserts

Rev from	Rev to	JOB #
03	04	9457-17

Notes:

- BD Catalog Number: 220220, 220221, 220222, 220239, 220240, 220244, 220526, 220527, 220528, 220529, 220531
- Blank (Sheet) Size: Length: 8.5" Width: 11"
- Number of Pages: 32 Number of Sheets: 8
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 8.5" x 5.5"
- Ink Colors: No. of Colors: 1 PMS#: 2755 Blue
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 5

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. OEM for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 9:29 am, Jul 26, 2017	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Natalie Morio at 9:47 am, Jul 26, 2017		
Checked By	THIRD EYE BY By Nichole Graham at 1:57 pm, Jul 31, 2017		
Part Number:	45	Category and Description Package Insert, BD Universal Viral Transport	Sheet: 1 of 33 Scale: N/A A

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жепрінікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinize temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD™ Universal Viral Transport System is intended for the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasmas or ureaplasmas from the collection site to the testing laboratory. This system can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasmal and ureaplasma culture.

SUMMARY AND EXPLANATION

One of the routine procedures in the diagnosis of infections caused by viruses, chlamydiae, mycoplasmas or ureaplasmas involves the collection and safe transportation of biological samples. This can be accomplished using the **BD** Universal Viral Transport System. This system includes a universal transporting medium that is room temperature stable, which can sustain viability (and infectivity) of a plurality of organisms that include clinically important viruses, chlamydiae, mycoplasmas and ureaplasmas during transit to the testing laboratory. The formulation of **BD** Universal Viral Transport medium includes protein for stabilization, antibiotics to minimize bacterial and fungal contamination, and a buffer to maintain a neutral pH.

BD Universal Viral Transport System is provided with labeled capture-cap vials designed for transport of the clinical sample. This system is also supplied as a sample collection kit that comprises a package containing one capture-cap vial of medium and a peel pouch incorporating sterile polyester or nylon flocked specimen collection swabs, with scored shafts for easy breakage. The capture-cap is designed to secure the shaft of the swab sample to the cap, eliminating the use of forceps to remove the swab in the laboratory.

PRINCIPLES OF THE PROCEDURE

BD Universal Viral Transport medium consists of modified Hanks' balanced salt solution supplemented with bovine serum albumin, cysteine, gelatin, sucrose and glutamic acid. The pH is buffered with HEPES buffer. Phenol red is used to indicate pH. Vancomycin, amphotericin B and colistin are incorporated in the medium to inhibit growth of competing bacteria and yeast. The medium is isotonic and non-toxic to mammalian host cells. The presence of sucrose acts as a cryoprotectant which aids in the preservation of viruses and chlamydiae if specimens are frozen (-70 °C) for prolonged storage.

REAGENTS

Universal Viral Transport Medium Components

Hanks' Balanced Salts
Bovine Serum Albumin
L-Cysteine
Gelatin
Sucrose
L-Glutamic Acid
HEPES Buffer
Vancomycin
Amphotericin B
Colistin
Phenol Red
pH 7.3 ± 0.2 @ 25 °C

Warnings and Precautions

For *in vitro* Diagnostic Use.

- Inserting more than one swab in the capture-cap vial may interfere with proper cap closure.
- Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
- Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
- Sterilize all biohazard waste including specimens, containers and media after their use.
- Directions should be read and followed carefully.
- Do not re-sterilize unused swabs.

- Do not re-pack.
- Not suitable to collect and transport microorganisms other than viruses, chlamydiae, mycoplasmas and ureaplasmas.
- Not suitable for any other application than intended use.
- The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
- Do not use if the swab is visibly damaged (i.e., if the swab tip is broken).
- Do not bend flocked swabs prior to specimen collection.
- Do not ingest the medium.
- Do not use the Universal Viral Transport medium for premoistening or prewetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
- Do not use for more than one patient.
- **Ⓢ BD Universal Viral Transport System** is for single use only; reuse may cause a risk of infection and/or inaccurate results.
- Avoid skin contact with medium.

Storage: This product is ready for use and no further preparation is necessary. The product should be transported and stored in its original container at 2–25 °C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual sterile pouch unit and the specimen transport vial label.

Product Deterioration: **BD Universal Viral Transport** should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the color of the medium has changed from light orange-red, (4) the expiration date has passed, (5) the swab pouch is open, or (6) there are other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Specimens for viral, chlamydial, mycoplasmal or ureaplasma investigation should be collected and handled following published manuals and guidelines.⁵⁻¹¹ Once a swab specimen is collected it should be placed immediately into the transport vial where it comes into contact with transport media. To maintain optimum viability, transport the specimen to the laboratory as soon as possible. Best recovery is obtained when specimens are refrigerated at 2–8 °C or kept on wet ice following collection and while in transit. If there will be a long delay before processing, specimens should be frozen at -70 °C or colder and transported on dry ice. Storage at -20 °C is less satisfactory than storage at 4 °C or -70 °C and can result in the loss of infectivity.^{12,13} Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations.^{1,11,14} Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

PROCEDURES

Materials Provided: **BD Universal Viral Transport System** includes a capture-cap vial containing 1 mL or 3 mL of transport medium plus three glass beads. Universal Viral Transport System vials of 1 mL or 3 mL transport medium are supplied with one of the following specimen collection swab options:

Two regular size plastic scored shaft swabs with polyester fiber tips.

One regular size plastic scored shaft swab and one stainless steel wire-plastic scored shaft minitip swab, both with polyester fiber tips.

One regular size plastic scored shaft swab with nylon flocked fiber tip.

One minitip size plastic scored shaft swab with nylon flocked fiber tip.

One flexible minitip size plastic scored shaft swab with nylon flocked fiber tip.

One regular and flexible minitip size plastic scored shaft swab with nylon flocked fiber tip.

These different swab applicator shafts facilitate the collection of specimens from various sites on a patient. Refer to the individual product descriptions for specific information about materials supplied.

Materials Required But Not Provided: Appropriate materials for isolating, differentiating and culturing viruses, chlamydiae, mycoplasmas and ureaplasmas. These materials include tissue culture cell lines, tissue culture medium, incubation systems and reading equipment. Refer to appropriate references for recommended protocols for isolation and identification of viral, chlamydial, mycoplasmal and ureaplasma agents.^{5-8,10}

Test Procedure

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals.⁵⁻¹¹ Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness.

For Universal Viral Transport Medium Vials (1 mL or 3 mL)

1. Aseptically remove cap from vial.
2. Aseptically place vesicle aspirates, corneal or conjunctival scrapings, small pieces of tissue or stool samples into the vial with medium.
3. Replace cap on vial and close tightly.
4. Label with appropriate patient information.
5. Send to the laboratory for immediate analysis.

For Universal Viral Transport Collection Kits

NOTE: Nylon flocked swab should not be bent prior to specimen collection.

1. Collect specimen with one swab.
2. Aseptically remove cap from vial.
3. Insert swab into the vial with medium.
4. Break swab shaft by bending it against the vial wall evenly at the pre-scored line.

5. Replace cap on vial and close tightly.
6. Label with appropriate patient information.
7. Send to the laboratory for immediate analysis.

NOTE: **BD Universal Viral Transport** contains antimicrobial substances intended to inhibit commensal bacteria and fungi. When collecting specimens from sites of the body known to contain high levels of commensal organisms it is good practice to refrigerate specimens and process as soon as possible in order to minimize breakthrough growth of bacteria or fungi. It is also common practice to add an antibiotic mixture to cell culture re-feed medium when the specimen is inoculated. This procedure helps avoid bacterial and fungal contamination of the cell culture. For specific information about process and cultivation techniques for specimens, consult laboratory reference manuals and standards.¹¹

Quality Control

All lots of the Universal Viral Transport medium are tested for microbial contamination, toxicity to host cells and the ability to maintain viability of desired agents. Procedures for quality control of Universal Viral Transport medium and viral culture media are described in a number of publications by the American Society for Microbiology^{6,8,10} and by CLSI.^{15,16} If aberrant quality control results are noted, patient results should not be reported.

RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

NOTE: Lower volumes of **BD Universal Viral Transport** medium will reduce the dilution effect of the specimen introduced into the vial; therefore, the organism(s) under investigation as well as any commensals or normal flora will be more concentrated.

LIMITATIONS OF THE PROCEDURE

1. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.⁵⁻¹¹
2. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
3. Universal Viral Transport Systems are intended for use as a collection and transport medium for viral, chlamydial, mycoplasmal and ureaplasma agents only. The medium can serve as a cryoprotectant for clinical viruses, including cytomegalovirus and varicella-zoster virus.
4. Calcium alginate swabs are toxic for many enveloped viruses and may interfere with fluorescent antibody tests, so they should not be used for specimen collection. Wooden shaft swabs may contain toxins and formaldehydes and should not be used. Polyester-tipped or nylon flocked swabs are suitable when specimen collection by a swab is appropriate.
5. Performance Characteristics of **BD Universal Viral Transport** Systems are validated with **BD Universal Viral Transport Swabs** and nylon flocked swabs. The use of tubes of medium or swabs from any other source has not been validated and could affect the performance of the product.

PERFORMANCE CHARACTERISTICS

Viability studies were performed using **BD Universal Viral Transport** with a variety of viruses, chlamydiae, mycoplasmas and ureaplasmas. Swabs accompanying each transport system were directly inoculated in triplicate with 100 µL of organism suspension. Swabs were then placed in their respective transport medium vials and were held for 0, 24 and 48 h at both 4 °C and room temperature (20–25 °C). At the appropriate time interval, each swab was vortexed, removed from its transport medium vial and then an aliquot of this suspension was inoculated into shell vials or into appropriate culture media. All cultures were processed by standard laboratory culture technique and examined after a specified incubation time. Organism viability was determined by fluorescing foci counts for viral and chlamydial strains and by CFU counts for mycoplasmal and ureaplasma strains. Organisms evaluated were: adenovirus, cytomegalovirus, echovirus type 30, herpes simplex virus type 1, herpes simplex virus type 2, influenza A, parainfluenza 3, respiratory syncytial virus, varicella-zoster virus, *Chlamydomphila pneumoniae*, *Chlamydia trachomatis*, *Mycoplasma hominis*, *Mycoplasma pneumoniae* and *Ureaplasma urealyticum*.

The results for the strains tested using **BD Universal Viral Transport** System are shown in the tables below.

BD Universal Viral Transport System was able to maintain the viability of the following organisms for at least 48 h at both room temperature (20–25 °C) and in the refrigerator (2–8 °C) under the test conditions described above: adenovirus, cytomegalovirus, echovirus type 30, herpes simplex virus type 1, herpes simplex virus type 2, influenza A, parainfluenza 3, respiratory syncytial virus, varicella-zoster virus, *Chlamydomphila pneumoniae*, *Chlamydia trachomatis*, *Mycoplasma hominis*, *Mycoplasma pneumoniae* and *Ureaplasma urealyticum*.

Table 1

Organism	Organism Concentration	Holding Time (hours)	Incubation Time Before Reading (hours)	Viability Challenge at 4 °C Foci of infected cells/200 μL^2	Viability Challenge at RT Foci of infected cells/200 μL^2
Adenovirus	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 70% of cells)	0	24	123	119
		24	24	62	47
		48	24	68	63
	10 ⁻² Neat Virus Stock Suspension* (dilution produces infectivity of 42% of cells)	0	24	17	14
		24	24	5	3
		48	24	5	7
Cytomegalovirus	Neat Virus Stock Suspension* (neat produces infectivity of 3% of cells)	0	24	337	444
		24	24	582	1012
		48	24	394	506
	1:2 Neat Virus Stock Suspension* (dilution produces infectivity of 2% of cells)	0	24	49	195
		24	24	63	80
		48	24	72	228
Echovirus Type 30	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 64% of cells)	0	24	76	79
		24	24	59	75
		48	24	66	60
	10 ⁻² Neat Virus Stock Suspension* (dilution produces infectivity of 35% of cells)	0	24	34	48
		24	24	18	26
		48	24	25	20
Herpes Simplex Virus Type 1	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 100% of cells)	0	24	491	412
		24	24	387	301
		48	24	282	164
	10 ⁻² Neat Virus Stock Suspension* (dilution produces infectivity of 25% of cells)	0	24	98	100
		24	24	68	10
		48	24	21	1
Herpes Simplex Virus Type 2	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 90% of cells)	0	24	TNTC ¹	TNTC ¹
		24	24	615	437
		48	24	525	58
	10 ⁻² Neat Virus Stock Suspension* (dilution produces infectivity of 40% of cells)	0	24	228	315
		24	24	170	73
		48	24	75	7
Influenza A	Neat Virus Stock Suspension* (neat produces infectivity of 59% of cells)	0	16	129	134
		24	16	172	166
		48	16	166	169
	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 47% of cells)	0	16	123	115
		24	16	71	72
		48	16	67	65
Parainfluenza 3	Neat Virus Stock Suspension* (neat produces infectivity of 57% of cells)	0	24	24	32
		24	24	26	28
		48	24	26	19
	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 51% of cells)	0	24	2	8
		24	24	12	10
		48	24	8	4
Respiratory Syncytial Virus	Neat Virus Stock Suspension* (neat produces infectivity of 47% of cells)	0	24	178	248
		24	24	251	208
		48	24	183	232
	10 ⁻¹ Neat Virus Stock Suspension* (dilution produces infectivity of 8% of cells)	0	24	17	13
		24	24	28	21
		48	24	14	16
Varicella-Zoster Virus	Neat Virus Stock Suspension* (neat produces infectivity of 8% of cells)	0	72	TNTC ¹	TNTC ¹
		24	72	TNTC ¹	TNTC ¹
		48	72	283	424
	1:2 Neat Virus Stock Suspension* (dilution produces infectivity of 2% of cells)	0	72	TNTC ¹	TNTC ¹
		24	72	TNTC ¹	TNTC ¹
		48	72	132	159

* 100 μL of suspension dosed onto the swab tip then swab placed in Universal Viral Transport vial containing 3 mL of transport medium

¹ TNTC= Too numerous to count

² Average of triplicate tests performed on 200 μL aliquots of Universal Viral Transport medium at each time point

Table 2

Organism	Organism Concentration	Holding Time (hours)	Incubation Time Before Reading (days)	Viability Challenge at 4 °C Fluorescing cytoplasmic inclusions/200 µL ²	Viability Challenge at RT Fluorescing cytoplasmic inclusions/200 µL ²
<i>Chlamydomophila pneumoniae</i>	Neat <i>Chlamydomophila</i> Stock Suspension* (neat produces TNTC ¹ cytoplasmic inclusions over entire HeLa DHI shell vials coverslip)	0 24 48	3 3 3	TNTC ¹ TNTC ¹ 201	TNTC ¹ TNTC ¹ 136
	10 ⁻¹ Neat <i>Chlamydomophila</i> Stock Suspension* (dilution produces TNTC ¹ cytoplasmic inclusions over entire HeLa DHI shell vials coverslip)	0 24 48	3 3 3	256 175 39	257 276 17
<i>Chlamydia trachomatis</i>	Neat <i>Chlamydia</i> Stock Suspension* (neat produces TNTC ¹ cytoplasmic inclusions over entire BGGMK DHI shell vials coverslip)	0 24 48	3 3 3	TNTC ¹ TNTC ¹ 317	TNTC ¹ TNTC ¹ 50
	10 ⁻¹ Neat <i>Chlamydia</i> Stock Suspension* (dilution produces TNTC ¹ cytoplasmic inclusions over entire BGGMK DHI shell vials coverslip)	0 24 48	3 3 3	216 164 67	171 48 6

* 100 µL of suspension dosed onto the swab tip then swab placed in Universal Viral Transport vial containing 3 mL of transport medium

¹ TNTC= Too numerous to count

² Average of triplicate tests performed on 200 µL aliquots of Universal Viral Transport medium at each time point

Table 3

Organism	Organism Concentration	Holding Time (hours)	Incubation Time Before Reading (days)	Viability Challenge at 4 °C CFU/200 µL ²	Viability Challenge at RT CFU/200 µL ²
<i>Mycoplasma hominis</i>	Neat <i>Mycoplasma</i> Stock Suspension*: Four <i>Mycoplasma hominis</i> Bacti Disks™ reconstituted into 20 mL of PPLO broth and incubated in 5 – 10% CO ₂ at 35–37 °C for 48 h (reference Remel <i>Mycoplasma</i> Bacti Disks™ Pack Insert TI No. 19314)	0 24 48	7 7 7	~ 1000, TNTC ¹ ~ 1000, TNTC ¹ ~ 1000, TNTC ¹	~ 1000, TNTC ¹ ~ 1000, TNTC ¹ ~ 1000, TNTC ¹
	10 ⁻² Neat <i>Mycoplasma</i> Stock Suspension*	0 24 48	7 7 7	17 17 11	16 10 12
<i>Mycoplasma pneumoniae</i>	Neat <i>Mycoplasma</i> Stock Suspension*: Four <i>Mycoplasma pneumoniae</i> Bacti Disks™ reconstituted into 20 mL of SP4 broth with glucose and incubated in ambient air at 35–37 °C for 7 – 14 days until broth becomes yellow (reference Remel <i>Mycoplasma</i> Bacti Disks™ Pack Insert TI No. 19314)	0 24 48	7 7 7	171 219 183	169 238 184
	10 ⁻¹ Neat <i>Mycoplasma</i> Stock Suspension*	0 24 48	7 7 7	17 22 17	18 26 19
<i>Ureaplasma urealyticum</i>	Neat <i>Ureaplasma</i> Stock Suspension*: Ten <i>Ureaplasma urealyticum</i> Bacti Disks™ reconstituted into 18 mL of 10B broth and incubated in ambient air at 35–37 °C for 24 h (reference Remel <i>Ureaplasma</i> Bacti Disks™ Pack Insert TI No. 19315)	0 24 48	3 3 3	1020 1136 1249	1125 1083 1056
	10 ⁻¹ Neat <i>Ureaplasma</i> Stock Suspension*	0 24 48	3 3 3	101 107 116	83 108 103

* 100 µL of suspension dosed onto the swab tip then swab placed in Universal Viral Transport vial containing 3 mL of transport medium

¹ TNTC= Too numerous to count

² Average of triplicate tests performed on 200 µL aliquots of Universal Viral Transport medium at each time point

AVAILABILITY

Cat. No.	Description
220220	BD™ Universal Viral Transport 3 mL Vial, carton of 50.
220221	BD™ Universal Viral Transport Standard Kit; each kit contains 3 mL vial, package of 2 sterile polyester-tip regular swabs with scored plastic shaft, carton of 50.
220222	BD™ Universal Viral Transport Combo Kit; each kit contains 3 mL vial, package of 1 sterile polyester minitip swab with wire/plastic scored shaft and 1 sterile polyester-tip regular swab with scored plastic shaft, carton of 50.
220239	BD™ Universal Viral Transport Regular Swabs, sterile polyester-tip with scored plastic shaft, 2 per pouch, carton of 100.
220240	BD™ Universal Viral Transport Combo Swabs, 1 sterile polyester minitip swab with wire/plastic scored shaft and 1 sterile regular swab with scored plastic shaft, 2 per pouch, carton of 100.
220244	BD™ Universal Viral Transport 1 mL Vial, carton of 50.
220526	BD™ Universal Viral Transport Kit; each kit contains 1 mL vial and 1 sterile nylon flocked flexible minitip swab with scored plastic shaft, carton of 50.
220527	BD™ Universal Viral Transport Kit; each kit contains 3 mL vial and package of 1 sterile nylon flocked regular and 1 sterile nylon flocked flexible minitip swabs with scored plastic shafts, carton of 50.
220528	BD™ Universal Viral Transport Kit; each kit contains 3 mL vial and 1 sterile nylon flocked regular swab with scored plastic shaft, carton of 50.
220529	BD™ Universal Viral Transport Kit; each kit contains 3 mL vial and 1 sterile nylon flocked minitip swab with scored plastic shaft, carton of 50.
220531	BD™ Universal Viral Transport Kit; each kit contains 3 mL vial and 1 sterile nylon flocked flexible minitip swab with scored plastic shaft, carton of 50.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17: 53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller. 2007. Manual of clinical microbiology. 9th ed. American Society for Microbiology, Washington, D.C.
6. Gleaves, C.A., R.L. Hodinka, S.L.G. Johnston, and E.M. Swierkosz. 1994. Cumitech 15A. Laboratory diagnosis of viral infections. ASM, Washington, DC.
7. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's diagnostic microbiology. 11th ed. Mosby, St. Louis, MO.
8. Wardford, A., M. Chernesky, and E. M. Peterson. 1999. Cumitech 19A, Laboratory diagnosis of *Chlamydia trachomatis* infections. ASM, Washington, DC.
9. Miller, J. M. 1999. A guide to specimen management in clinical microbiology, 2nd ed. ASM, Washington, DC.
10. Isenberg, H. D., 2004. Clinical microbiology procedures handbook, 2nd ed. ASM, Washington, DC.
11. Isenberg, H.D., 1998. Essential procedures for clinical microbiology. Chapter 14.12, Page 787. Packaging and shipping infectious substances.
12. Maass, M. and K. Dalhoff. 1995. Transport and storage conditions for cultural recovery of *Chlamydia pneumoniae*. J. Clin. Microbiol. 33:1793-1796.
13. Johnson, F. 2005. Transport of viral specimens. Clin. Microbiol. Rev 3:120-121.
14. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
15. Clinical and Laboratory Standards Institute. 2003. Quality control of microbiological transport systems. Approved Standard M40-A, CLSI, Wayne, PA.
16. Clinical and Laboratory Standards Institute. 2006. Viral culture; Approved Guideline M41-A, CLSI, Wayne, PA.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Bacti Disks is a trademark of Remel, Inc.
© 2017 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.
Benex Limited
Becton Dickinson Pty Ltd.