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1. Parker, R.H., and P.D. Hoerich. 1962. Disk method for rapid identification of *Haemophilus* species. Am. J. Clin. Pathol. 37:319-327.
2. Kilian, M., and E.L. Biberstein. 1984. Genus II. *Haemophilus*, p. 558-569. In N.R. Krieg and J.G. Holt (ed.), *Bergey's manual of systematic bacteriology*, vol. 1. Williams & Wilkins, Baltimore.
3. Holt, L.B. 1962. The growth-factor requirements of *Haemophilus influenzae*. J. Gen. Microbiol. 27:317-322.
4. Howard, B.J. 1987. *Haemophilus*, p. 279-288. In B.J. Howard (ed.), *Clinical and pathogenic microbiology*. The C.V. Mosby Company, St. Louis.
5. Waterworth, P.M. 1955. The stimulation and inhibition of the growth of *Haemophilus influenzae* in media containing blood. Br. J. Exp. Pathol. 36:186-194.
6. Krumwiede, E., and A.G. Kuttner. 1938. A growth inhibitory substance for the influenzae group of organisms in the blood of various animal species. J. Exp. Med. 67:429-441.
7. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
11. Campos, J.M. 1995. *Haemophilus*, p. 556-565. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.

12. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1997. *Haemophilus*, p. 363-394. Color atlas and textbook of diagnostic microbiology, 5th ed. Lippincott-Raven Publishers, Philadelphia.
13. Taubitz, I.S., and H. Brandis, 1988. A comparison between methods of identification and serotyping of encapsulated strains of *Haemophilus influenzae*. Zbl. Bakt. Hyg. A 270, 83-97.

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas /
 Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производител / Аткарушы

Use by / Spottføjte do / Användes för / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης /
 Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Pouzít do / Usar antes de / Använd före / Используйте до /
 A se utiliza până la / Son kulanma taríhi / Uptrebiti do / Исползовать до / дейн пайдаланура /
 YYYY-MM-DD / YYYY-MM (MM = end of month) /
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce) /
 ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) /
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand) /
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) /
 VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) /
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois) /
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) /
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) /
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) /
 AAAA-MM-GG / AAAA-MM (MM = fine mese) /
 MMMM-MM-DD / MMMM-MM (MM = mês do ano) /
 ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) /
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) /
 AAAA-MM-DD / AAAA-MM (MM = fim do mês) /
 RRRR-MM-DD / RRRR-MM (MM = koniec mesiacu) /
 aaaa-mm-dd / aaaa-mm (mm = fin del mes) /
 ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) /
 ТТТТ-ММ-ДД / ТТТТ-ММ (ММ = края на месеца) /
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) /
 YYYY-AA-GG / YYYY-AA (AA = ayın sonu) /
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca) /
 ТТТТ-ММ-ДД / ТТТТ-ММ (ММ = конец месяца) /
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы)

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός
 κατάλογου / Katalogusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo /
 Каталоген номер / Număr de catalog / Katalog numerasi / Kataloški broj / Номер по каталогу / Каталог номері

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend
 vertegenwoordiger in de Europese Unie / Volitatud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé
 pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban /
 Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane
 przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante
 autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизран представител в EU / Reprezentant autorizat în Uniunea Europeană /
 Аутра Төрлөлүгү Yetkilii Төмсилчиси / Ovlašćeni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Еврoпа
 қауымдастығындағы уәкілетті өкіл

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel
 voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro /
 Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnostikai orvosi eszköz / Dispositivo medico diagnostico in vitro /
 In vitro diagnostikos prietais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico
 in vitro / Medicinska pomůcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински
 уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostic Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku /
 Медицинский прибор для диагностики in vitro / Жасанды жагдайда жүргізетін медициналық диагностика аспабы

Temperature limitation / Teplotní omezení / Temperaturbegrensning / Temperaturlimit / Temperaturi piirang / Lämpötilarajoitus / Température
 limite / Zulässiger Temperaturbereich / Όριο θερμοκρασίας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrensning /
 Ograniczenie temperatury / Limitação da temperatura / Ohraničenie teploty / Limitación de temperatura / Temperaturbegrensning / Температури
 ограничения / Limitare de temperatură / Sicaklık sınırlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу

Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode
 (Chargenbezeichnung) / Кодикós партидас (Партида) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii
 (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod serije /
 Код партии (лот) / Топтама коды

Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista
 käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást /
 Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de
 utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направете справка в инструкциите за употреба /
 Consultați instrucțiunile de utilizare / Kullanim Talimatları'na başvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайдалану
 нұсқаулығымен танысыңыз

Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA
 800-638-8663
 www.bd.com/ds

Benex Limited
 Rineanna House
 Shannon Free Zone
 Shannon, County Clare, Ireland

ATCC is a trademark of the American Type Culture Collection.
 BD, BD Logo, BBL, Sensi-Disc, Taxo and Trypticase are trademarks of Becton, Dickinson and Company. © 2010 BD