BD Schaedler Kanamycin-Vancomycin Agar with 5% Sheep Blood (Schaedler-KV Agar)

USO PREVISTO

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BIBLIOGRAFIA

1. Schaedler, R.W., R. Dubos, and R. Costello. 1965. The development of the bacterial flora in the gastrointestinal tract of mice. *J. Exp. Med.* 122:59-66.
2. Mata, L.J., C. Carrillo, and E. Villatoro. 1969. Fecal microflora in healthy persons in a preindustrial region. *Appl. Microbiol.* 17:596-602.
3. Finegold, S.M., A.B. Miller, and D.J. Posnick. 1965. Further studies on selective media for *Bacteroides* and other anaerobes. *Ernährungsforschung* 10:517-528.
4. Gibbons, R.J., and J.B. MacDonald. 1960. Hemin and vitamin K compounds as required factors for the cultivation of certain strains of *Bacteroides melaninogenicus*. *J. Bacteriol.* 80:164-170.
5. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt. 1974. Isolation of anaerobic bacteria. *In*: E.H. Lennette, E.H. Spaulding, and J.P. Truant (ed.), *Manual of clinical microbiology*, 2nd ed. American Society for Microbiology, Washington, D.C.
6. Allen, S.D., J.A. Siders, and L.M. Marler. 1985. Isolation and examination of anaerobic bacteria. *In*: E.H. Lennette, A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.), *Manual of clinical microbiology*, 4th ed. American Society for Microbiology, Washington, D.C.
7. Jousimies-Somer, H.R., et al. 2003. *Bacteroides*, *Porphyromonas*, *Prevotella*, *Fusobacterium*, and other anaerobic gram-negative bacteria. *In*: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). *Manual of clinical microbiology*, 8thed. American Society for Microbiology, Washington, D.C.
8. Miller, J.M., and H.T. Holmes. 1995. Specimen Collection, transport, and storage. *In*: Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
9. Holdeman, L.V., E.P. Cato, and W.E.C. Moore (ed.). 1977. *Anaerobe laboratory manual*, 4th ed. Virginia Polytechnic Institute and State University, Blacksburg.
10. Engelkirk, P.G., J. Duben-Engelkirk, and V.R. Dowell, Jr. 1992. *Principles and practice of clinical anaerobic bacteriology*. Star Publishing Co., Belmont, Calif.

11. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.
12. Murray, P.R., and D.M. Citron. 1991. General processing of specimens for anaerobic bacteria. *In*: A. Balows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
13. van Winkelhoff, A.J., and J. de Graaff. 1983. Vancomycin as a selective agent for isolation of Bacteroides. *J. Clin. Microbiol.* 18:1282-1284.

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection
BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. © 2013 BD

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Becton Dickinson GmbH