
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	8894-16

NOTES:

- BD Catalog Number: 245159
- Blank (Sheet) Size: Length: 11" Width: 8.5"
- Number of Pages: 24 Number of Sheets: 5
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 5.5" x 2.75"
- Ink Colors: No. of Colors: 2 PMS#: # 032 Red, Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 5

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. 8085917 for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 9:23 am, Dec 13, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	PROOFING APPROVED BY By Daria Kuznetsova at 1:59 pm, Dec 13, 2016			
Checked By	THIRD EYE BY By Natalie Morio at 7:35 am, Dec 15, 2016			
Part Number:	8085917	Category and Description	Sheet: 1 of 25	A
			Scale: N/A	
		Package Insert, BD MGIT™ Tbc ID Test		

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INTENDED USE

The **BD MGIT™** Tbc Identification Test (Tbc ID) is a rapid chromatographic immunoassay for the qualitative detection of *Mycobacterium tuberculosis* complex (MTbc) antigen from AFB smear-positive **MGIT** tubes. The device will detect the following species of the MTbc: *M. tuberculosis*, *M. bovis*, *M. africanum*, and *M. microti*.

SUMMARY AND EXPLANATION

Tuberculosis is a chronic disease caused by an infection of *M. tuberculosis*. It is a major global public health problem and one of the leading causes of death from infectious disease in the world.¹ There were an estimated 9.2 million new cases of tuberculosis in 2006 with 1.5 million deaths.² Ninety-five percent of the estimated incidence occurs in developing countries with 80 percent reported from 22 countries.³

Conventional methods for identification of mycobacteria rely on staining specimens for acid-fast bacilli (AFB) followed by culturing on solid media and biochemical testing. It can take as long as two months to speciate an isolate using these standard methods.⁴ The World Health Organization (WHO) has recommended the use of liquid media for the culture of *M. tuberculosis*.⁵ Modern procedures used to identify *M. tuberculosis* from liquid culture include nucleic acid probes and gas-liquid / high performance chromatography.

It is important to distinguish non-tuberculous mycobacteria (NTM) from *Mycobacterium tuberculosis* complex (MTbc).

The **BD MGIT** Tbc ID is a rapid procedure that identifies MTbc from culture grown in liquid media. It is a simple chromatographic immunoassay which requires no sample preparation. The total assay time is 15 minutes with reactivity determined by visual color development.

PRINCIPLES OF THE PROCEDURE

The Tbc ID test is a chromatographic immunoassay for the qualitative detection of MTbc from an AFB smear-positive **MGIT** tube. This product detects MPT64, a mycobacterial protein fraction that is secreted from MTbc cells during culture. When samples are added to the test device, MPT64 antigen binds to anti-MPT64 antibodies conjugated to visualizing particles on the test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by a second specific MPT64 antibody applied to the membrane. If the MPT64 antigen is present in the sample, a color reaction is produced by the labeled colloidal gold particles and is visualized as a pink to red line.

REAGENTS

Materials Provided:

The following components are included:

BD MGIT™ Tbc ID test devices	25 devices	Foil pouched device (containing a single reactive strip) and desiccant. The strip contains a test line of MPT64-specific monoclonal antibody and a control line of anti-species antibody.
-------------------------------------	------------	---

Materials Required But Not Provided: Pipettor (capable of delivering 100 µL), sterile pipette tips and timer.

Warnings and Precautions

For *in vitro* Diagnostic Use

1. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"⁶⁻⁹ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
2. Do not use devices beyond the expiration date.
3. Do not reuse the device.

4. Use a clean, sterile pipette tip for each sample.
5. Do not test clinical specimens directly in this device (e.g., human body fluid, tissue, sputa and bronchial lavage fluid).
6. Use appropriate biosafety precautions for handling mycobacteria. Procedure should be carried out using appropriate biological safety cabinets.^{4,8}
7. Used test devices may contain viable MTbc which could be infectious. Used devices should be discarded according to your institutional guidelines or Standard Precautions requirements.
8. Processed sputum and other non-sterile specimens are typically cultured in the **MGIT** system for the growth and detection of mycobacteria. As a result, primary cultures can contain non-AFB microorganisms. If non-AFB organisms predominate in a positive **MGIT** culture their growth and metabolism may interfere with the **BD MGIT** Tbc ID performance. Care must be taken when using this test with polymicrobial cultures containing both AFB and non-AFB organisms. As a guideline, only use the Tbc ID test on a positive **MGIT** tube if AFB-positive organisms predominate on the smear. Performing the test on a very turbid positive culture (>1.0 McFarland) containing an overgrowth of non-AFB organisms may cause a false positive result.
9. There are reports of rare strains of MTbc that produce MPT64 antigen below detectable limits.¹⁰
10. False positive results may be observed in the presence of protein A-producing strains of bacteria (e.g. *Staphylococcus aureus*).

Danger

H228 Flammable solid.

P280 Wear protective gloves/protective clothing/eye protection/face protection. **P370+P378** In case of fire: Use for extinction: CO₂, powder or water spray.

Storage and Handling: Test devices may be stored at 2–35 °C. DO NOT FREEZE. Devices must be at ambient room temperature at time of testing.

SAMPLE COLLECTION AND PREPARATION

This test is designed to identify MTbc from AFB smear-positive **MGIT** tubes (4 mL and 7 mL). The presence of AFB in a positive **MGIT** tube should be confirmed using an AFB smear prior to conducting the test.

Sample Storage: Positive **MGIT** tubes can be stored at 2–37 °C for up to 10 days after **MGIT** tube positivity and prior to testing with the Tbc ID device. If necessary, positive **MGIT** tubes may be stored and maintained at -20 to 8 °C for up to two months.

PROCEDURE

NOTES:

- AFB smear-positive **MGIT** tubes can be tested in the Tbc ID device within 10 days after **MGIT** tube positivity.
 - If devices are refrigerated, they must be brought to ambient room temperature in the foil pouch prior to testing.
1. Remove the Tbc ID device from its foil pouch immediately before testing. Place the device on a flat surface.
 2. Label one device for each sample to be tested.
 3. Thoroughly mix the sample (AFB smear-positive **MGIT** tube) by inverting or vortexing. Do not centrifuge.
 4. Remove cap from **MGIT** tube and using a sterile pipette tip, pipette 100 µL of the sample into the sample well (as indicated by the teardrop ●) of the appropriately labeled device. Tightly replace cap on **MGIT** tube. Start timer for 15 min.
 5. Read result at 15 min and record test result. Do not interpret test after 60 min.

QUALITY CONTROL

Each device contains both positive and negative internal/procedural controls. The appearance of a control line in the read window at the Control "C" position provides an internal positive control that validates the proper reagent function and assures that the correct test procedure was followed. The membrane area surrounding the test and control lines is the internal negative control for the device. A background area that is white to light pink indicates that the test is performing correctly.

Additional quality control should be performed in accordance with local or country regulations, laboratory accreditation requirements, and your laboratory's standard quality control procedures.

Positive and negative external controls should be tested in the same manner as test samples to provide a means of external quality control. Positive Control: A positive **MGIT** tube prepared by growing a known isolate of MTbc. This should yield a positive result. Negative Control: An uninoculated **MGIT** tube. This should yield a negative result.


External controls should be run, at a minimum, for each new lot or each new shipment received.

If the controls do not perform as expected, do not report sample results. Contact your local BD representative or Technical Services for assistance.


INTERPRETATION OF RESULTS

Positive Test for TBc (MPT64 antigen present) - A pink to red line appears at the Test "T" position and the Control "C" position in the read window. This indicates MPT64 antigen was detected in the sample. The intensity of the C and T lines may vary. The background area should be white to light pink.

Negative Test for TBc (no MPT64 antigen detected) - No pink to red line is visible at the Test "T" position of the read window. This indicates that MPT64 antigen was not detected in the sample. A line at the Control "C" position read window indicates proper performance of the test procedure. The background area should be white to light pink.

Invalid Test - The test is invalid if no pink to red line is visible at the Control "C" position in the read window or if the background area color inhibits test interpretation. If invalid, the sample must be retested with a new device.

Invalid

REPORTING OF RESULTS

Positive Test	Should be reported as MTb complex (MTbc).
Negative Test	Should be reported as Acid-Fast Bacilli, non-MTbc.
Invalid Test	Do not report results.

LIMITATIONS OF THE PROCEDURE

1. **This test does not rule out the presence of other mycobacterial or mixed bacterial infections.**
2. This test is unable to differentiate between MTb complex (MTbc) organisms.
3. This test should not be used solely for the determination of MTbc infection. The test results are to be used in conjunction with information available from the patient's clinical evaluation and other diagnostic procedures.
4. A negative result does not always rule out the possibility of infection with MTbc. The device is unable to detect MTbc when a mutation arises in the MPT64 gene. The test results are to be used in conjunction with information available from the patient's clinical evaluation and other diagnostic procedures.
5. Some substrains of *M. bovis* BCG among *M. tuberculosis* complex produce no MPT64 antigen and will therefore result in a negative test result with the device.¹⁰

PERFORMANCE CHARACTERISTICS

Cinical Studies

The **BD MGIT** TBc ID test was compared to a molecular identification method (Gen-Probe AccuProbe™ MTbc Identification Test) at one clinical site in a TB high burden country using a total of 247 AFB smear-positive, **BD MGIT** 960 instrument positive samples. AccuProbe identified 121 as MTbc and 126 as NTM. A comparison of these results is shown below in Table 1.

Table 1: Performance Summary of the BD MGIT TBc ID Test compared to AccuProbe with Percent Agreements and Confidence Intervals (CI)

BD MGIT TBc ID	AccuProbe	
	+	–
+	120	9*
–	1	117
(+)= Positive (–)= Negative *Six of the nine positive by TBc ID and negative by AccuProbe were positive by the Hain GenoType™ MTBDR Plus. Positive Percent Agreement (95% CI): 99% (95.5%, 100%) Negative Percent Agreement (95% CI): 93% (86.9%, 96.7%) Overall Percent Agreement (95% CI): 96% (92.7%, 98%)		

Reproducibility

Reproducibility testing was performed at the clinical site using 6 negative samples and 6 positive samples (containing recombinant MPT64 antigen) for a total of 108 replicates. The overall reproducibility for the **BD MGIT** TBc ID test was 100%.

Analytical Studies

Analytical Specificity (Cross Reactivity): 756 samples of 23 non-MTb complex mycobacteria species were tested with no cross reactivity. Fifty-three samples of 10 non-mycobacterial species (common respiratory contaminants) were tested; 3 of the non-mycobacterial samples caused interference at very high levels of **MGIT** tube contamination (>1.0 McFarland).

Internal Studies

The **BD MGIT** TBc ID test was evaluated using mycobacterial strains that were seeded into **BD MGIT** 960 (7 mL) and Manual **MGIT** (4 mL) tubes. A total of 828 samples consisting of 23 MTbc strains (414 devices) and 23 NTM strains (414 devices) were tested. One device produced an invalid result.

Table 2: Performance Summary of the BD MGIT TBc ID Test Evaluated with MTbc and NTM Strains with Percent Agreements

BD MGIT TBc ID	Mycobacteria Strains	
	MTbc	NTM
+	413	0
–	0	414
(+)= Positive (–)= Negative Positive Percent Agreement: 100% (413 / 413) Negative Percent Agreement: 100% (414 / 414) Overall Percent Agreement: 100% (827 / 827)		

AVAILABILITY

Cat. No.	Description
245159	BD MGIT™ TBc Identification Test 25 Test Devices
245111	BD BBL™ MGIT™ 4 mL tube 25 tubes
245113	BD BBL™ MGIT™ 4 mL tube 100 tubes
245122	BD BBL™ MGIT™ 7 mL tube 100 tubes

REFERENCES

- Centers for Disease Control and Prevention. 2008. Morbidity and Mortality Weekly Report. Vol. 57, No.11.
- World Health Organization. 2008. WHO Report 2008, Global Tuberculosis Control, Surveillance, Planning, Financing. Geneva.
- World Health Organization. 2007. WHO Report 2007, Global Tuberculosis Control, Surveillance, Planning, Financing. Geneva.
- Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. U.S. Department of Public Health and Human Services, Atlanta, GA.

- World Health Organization. The use of liquid medium for culture and DST. WHO <http://www.who.int/tb/dots/laboratory/policy/en/index3.html>
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Office Journal L262*, 17/10/2000, p. 0021-0045.
- Hirano, K., A. Aono, M. Takahashi and C. Abe. 2004. Mutations including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB-negative *Mycobacterium tuberculosis* isolates. *J. Clin. Microbiol.* 42:390-392.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долоң лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najnižia dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimitiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / Контроль / 对照

CONTROL +	Positive control / Положительный контроль / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiveve kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívá kontrolé / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂
CONTROL -	Negative control / Отрицательный контроль / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativeve kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativá kontrolé / Negativeve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂
STERILE ED	Method of sterilization: ethylene oxide / Metoda na sterilizacija: этиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringmetode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseerimismetod: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacji: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Sterilizacija edici – этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringmetode: etylenoksid / Metoda sterylizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Method sterilizacijos: этиленоксид / Metodăa sterilizácie: etylenoxid / Metoda sterilizacije: etilen oksid / Steriliseringmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Metod sterilizacji: etilenoksydum / 灭菌方法: 环氧乙烷
STERILE R	Method of sterilization: irradiation / Metoda na sterilizacija: iradijacija / Způsob sterilizace: záření / Steriliseringmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismetod: kiirgus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Sterilizacija edici – сәуле түсіпу / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Method sterilizacijos: apšvitinimas / Metodăa sterilizácie: ožarovenie / Metoda sterilizacije: ozračevanje / Steriliseringmetod: strålning / Sterilizasyon yöntemi: irradyasyon / Metod sterilizacji: opróśnieniam / 灭菌方法: 辐照
	Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologilised riskid / Risques biologiques / Biološki rizik / Βιολογίαλλ veszélyes / Rischio biologico / Биологиялык төуекелдер / 생물학적 위험 / Biologinis pavojus / Biologisk riski / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологічна небезпека / 生物学风险
	Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatus! Lugeka kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristite prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동반된 설명서 참조 / Démesio, žiūrėkite pridėdamus dokumentus / Piesardžība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Výstraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увага: див. супутню документацию / 小心, 请参阅附带文档。
	Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның ықпал етпген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Górná granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sicaklık üst sınırı / Максимальна температура / 温度上限
	Keep dry / Пазете сухо / Skladujte v suchém prostredí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Manterner seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausai / Uzglabāt sausu / Droog houden / Hoides tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Berrepti við volongi / 请保持干燥
	Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ωρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жынау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamani / Час забору / 采集时间
	Peel / Обелете / Otvēfete zde / Abn / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Устіңгі қабатын алып таста / 벗기기 / Plēsti āia / Atlīmēt / Schillen / Trekken av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odrhňnite / Oljuštiti / Dra isår / Ayırma / Відклеїти / 撕下
	Perforation / Перфорация / Perforace / Perforering / Διάτρηση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Tesik тесу / 질취선 / Perforacija / Perforatie / Perforacja / Perforação / Perforare / Перфорация / Perforăcia / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Erep paket buzıylgan bolca, pайдаланба / पैकिङ्गकै 손상लत क्कल सुल सुल कुकु / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / He использовать при повреждении упаковки / Nepoužívejte, ak je obal poškozený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包裝破損, 請勿使用

Keep away from heat / Пазете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Сақын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / He нагрєвать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тепла / 请远离热源

Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Reži / Vágja ki / Tagliare / Keciňz / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciąć / Cortar / Decupați / Отрезать / Odstrihnite / Iseći / Klipp / Kesme / Розрізати / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paėmimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colecta / Data colectării / Дата сбора / Dátum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

AccuProbe is a trademark of Gen-Probe.

GenoType is a trademark of Hain Lifescience GmbH.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.