Part Number:	8840221JAA	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, BD BBL™ Taxo™ Anaerobic Discs	Rev from: 02	Rev to: 03
		Job Number:	745-18

Catalog Number: 231651, 231562

Blank (Sheet) Size: Length: 13.75" Width: 21.5"


Number of Pages: 8 Number of Sheets: 1

Page Size: Length: 13.75" Width: 5.4375" Final Folded Size: 1.25" X 5.4375"

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒

Online / In House Printed: ☐

Web Printed: ☐

See Specification control no. VS-88-4022-1PR for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Shanthi Sagayamary Ramu at 5:46 pm, Oct 08, 2018

Proofing
Approved By: **PROOFING APPROVED BY**
By Sonia F.Thompson at 3:43 pm, Oct 08, 2018

Third Eye By: **THIRD EYE BY**
By Natalie Morio at 1:22 pm, Oct 09, 2018

BD BBL™ Taxo™ Discs for Presumptive Identification of Gram-Negative Anaerobic Bacilli

English: pages 1 – 2 Italiano: pagine 5 – 6
Français: pages 2 – 3 Español: páginas 6 – 8
Deutsch: Seiten 4 – 5

8840221JAA(03)
2018-10

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinize temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

These discs are recommended for use in the presumptive identification of gram-negative anaerobic bacilli based on differences in susceptibility to antimicrobial agents. *This procedure is not for use for therapeutic purposes.*

SUMMARY AND EXPLANATION

Sutter and Finegold reported on the value of special-potency antibiotic discs in differentiating *Bacteroides* and *Fusobacterium* spp. on the basis of antibiotic inhibition patterns.¹ The patterns, along with colonial and microscopic observations and a limited number of biochemical tests, permitted the rapid presumptive identification of anaerobic bacilli. A similar differentiation protocol was adopted by the Centers for Disease Control (CDC).²

Sutter et al. recommend the use of kanamycin 1 mg, vancomycin 5 µg and colistin 10 µg discs with pure cultures streaked on Brucella Blood Agar with Hemin and Vitamin K₁, with a 10 mm zone size breakpoint for reporting susceptibility or resistance.³ The CDC protocol utilizes penicillin 2 units, rifampin 15 µg and kanamycin 1 mg discs on CDC Anaerobe Blood Agar, an enriched medium supplemented with hemin and Vitamin K₁, with larger zone size criteria.² The BD BBL test protocol utilizes these five discs on either medium and the 10 mm zone size breakpoint.

PRINCIPLES OF THE PROCEDURE

These discs are used for the differentiation of gram-negative anaerobic bacilli on the basis of *in vitro* testing of pure cultures against selected antimicrobial agents using the agar diffusion method. Growth of the microorganisms may be inhibited by the antimicrobial agent diffusing out from the disc, resulting in a zone of inhibition around the disc. Differentiation is based on the presence of zones of inhibition of equal to or greater than 10 mm (susceptible) or less than 10 mm (resistant).

PRODUCT DESCRIPTION

These 1/4" discs are made from high quality absorbent paper impregnated with accurately determined amounts of antimicrobial agents. They are furnished in cartridges containing 50 discs each. The last disc in each cartridge is marked "X" and contains the drug as coded. Cartridges are for use in BD BBL single disc dispensers.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

The BD BBL Taxo Kanamycin 1 mg and BD BBL Taxo Rifampin 15 µg discs are not for use in susceptibility testing for therapeutic purposes.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, test plates and other contaminated materials must be sterilized by autoclaving.

Storage Instructions:

- On receipt, store containers of discs at -20 to +8 °C. If the laboratory refrigerator is frequently opened and closed, and a suitable temperature is **not** maintained, place there a supply sufficient only for use within a week.
- Allow containers to come to room temperature before opening. Return unused discs to the refrigerator when application of discs has been completed.
- Use the oldest discs first.
- Discard expired discs. Also discard cartridges from which discs have been frequently removed during a week. Discard discs left out overnight in the laboratory, or test the discs for performance.
- If the discs form incorrect zones with the recommended control organism, the entire procedure should be checked; a faulty zone may be due to the disc, the inoculation, the preparation or depth of the medium or other factors.

The expiration date applies to discs in the intact container, stored as directed.

SPECIMENS

These discs are not for use directly with clinical specimens or other sources containing mixed flora. The organism to be presumptively identified first must be isolated as separate colonies by streaking the specimen onto appropriate culture media.

PROCEDURE

Material Provided: In this package are discs for the differentiation of gram-negative anaerobic bacilli as labeled.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure:

- Preparation of inoculum with test and control cultures
 - Prior to inoculation, loosen the cap of a tube of inoculum broth (e.g., Schaedler Broth with Vitamin K₁ or suitable alternative medium) and heat in a boiling water bath to drive off oxygen, retighten the cap and cool to room temperature, or reduce the inoculum broth overnight in an anaerobic jar.
 - Perform a Gram stain. Use only pure cultures.
 - Select three or four similar colonies and inoculate the inoculum broth.
 - Incubate anaerobically at 35 °C for 18 to 24 h or until heavy turbidity is observed; adjust the turbidity to be comparable to a McFarland No. 1 barium sulfate standard.
- Inoculation
 - Pre-warm to room temperature plates of CDC Anaerobe Blood Agar or Brucella Blood Agar with Hemin and Vitamin K₁. Agar surfaces should be smooth and moist.
 - Within 15 min of adjusting the turbidity of the inoculum, dip a sterile swab into the properly diluted inoculum and rotate it firmly several times against the upper inside wall of the tube to express the excess fluid.
 - Inoculate the entire agar surface three times, turning the plate 60° between streakings to obtain even inoculation.
 - Replace the lid of the plate and hold for 3 to 5 min at room temperature to allow absorption of the moisture from the inoculum.
 - Apply the discs evenly spaced to the inoculated agar surface with sterile forceps or single disc dispenser and tamp with a sterile needle or forceps to make complete contact with the medium surface.
 - Incubate immediately at 35 °C under anaerobic conditions.
 - Examine the plates after 48 h. Measure and record each zone of inhibition to the nearest millimeter: less than 10 mm as resistant; equal to or greater than 10 mm as susceptible.

User Quality Control: Control tests using known cultures should be included each time a susceptibility test is performed. The following organisms are recommended:

Antimicrobial Agent and Code					
Organism	Kanamycin K	Rifampin RA	Penicillin P-2	Vancomycin Va-5	Colistin CL-10
<i>Bacteroides fragilis</i> ATCC® 25285	R	S	R	R	R
<i>Fusobacterium mortiferum</i> ATCC 9817	S	R	S	R	S

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Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Compare reactivity patterns with those in the table.

Inhibition of <i>Bacteroides</i> and <i>Fusobacterium</i> Species with Antibiotic Discs†					
Antimicrobial Agent and Code					
Organism	Kanamycin K	Rifampin RA	Penicillin P-2	Vancomycin Va-5	Colistin CL-10
<i>B. fragilis</i> group*	R	S	R	R	R
<i>P. melaninogenica</i>	R ^S	S	S	R ^S	S ^r
<i>P. buccalis</i> and <i>P. veroralis</i>	R	S	S	R ^S	S ^r
<i>B. ureolyticus</i>	S	S	S	R	S
<i>F. mortiferum</i>	S	R	S ^r	R	S
<i>F. varium</i>	S	R	S ^r	R	S
<i>F. necrophorum</i>	S	S	S	R	S
<i>F. nucleatum</i>	S	S	S	R	S

S = Susceptible (zones ≥ 10 mm); S^r = majority of strains susceptible.
R = Resistant (zones < 10 mm); R^S = majority of strains resistant.
† Data were obtained on CDC Anaerobe Blood Agar and Brucella Blood Agar with Hemin and Vitamin K₁.
* The *B. fragilis* group includes *B. fragilis*, *B. distasonis*, *B. ovatus*, *B. thetaiotaomicron*, *B. vulgatus* and *B. uniformis*.

LIMITATIONS OF THE PROCEDURE

The references should be consulted for the complete schema of morphological and biochemical characteristics that, along with these antibiotic inhibition patterns, enable a presumptive identification to be made.¹⁻⁴ Following presumptive identification, definitive identification procedures should be performed to establish the validity of the presumptive identification. Appropriate references should be consulted for additional information.²⁻⁸

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BD BBL Taxo Anaerobe Differentiation Discs are tested to verify specific product characteristics. Samples are assayed for the potency of each of the antimicrobial agents contained in the kit. Additionally, plates of CDC Anaerobe 5% Sheep Blood Agar and Brucella Blood Agar with Hemin and Vitamin K₁ are swab-inoculated with cultures of *Bacteroides fragilis* (ATCC 25285) and *Fusobacterium mortiferum* (ATCC 9817). Samples are placed on the inoculated plates and incubated at 35 ± 2 °C under anaerobic conditions. After 48 h incubation, the plates are then read for zones of inhibition around each disc. *B. fragilis* is resistant (< 10 mm) to kanamycin, penicillin, vancomycin and colistin and susceptible (≥ 10 mm) to rifampin. *F. mortiferum* is resistant (< 10 mm) to rifampin and vancomycin and susceptible (≥ 10 mm) to kanamycin, penicillin and colistin.

AVAILABILITY

Cat. No.	Description
231651	BD BBL™ Taxo™ Anaerobe Differentiation Discs Set, 50 tests Set components – one cartridge of 50 discs of each of the following antimicrobial agents: BD BBL™ Taxo™ Kanamycin, 1 mg BD BBL™ Taxo™ Rifampin, 15 µg BD BBL™ Sensi-Disc™ Penicillin, 2 units BD BBL™ Sensi-Disc™ Vancomycin, 5 µg BD BBL™ Sensi-Disc™ Colistin, 10 µg
231562	BD BBL™ Taxo™ Kanamycin, 1 mg Discs, Single cartridge of 50 discs

REFERENCES

1. Sutter, V.L., and S.M. Finegold. 1971. Antibiotic disc susceptibility tests for rapid presumptive identification of gram-negative anaerobic bacilli. Appl. Microbiol. 21:13-20.
2. Dowell, V.R., Jr., and G.L. Lombard. 1977. Presumptive identification of anaerobic nonsporeforming gram-negative bacilli. Center for Disease Control, Atlanta.
3. Sutter, V.L., D.M. Citron, M.A.C. Edelstein, and S.M. Finegold. 1985. Wadsworth anaerobic bacteriology manual, 4th ed. Star Publishing Co., Belmont, Ca.
4. Finegold, S.M., and M.A.C. Edelstein. 1985. Gram-negative, nonsporeforming anaerobic bacilli, p. 450-460. In E.H. Lennette, A. Balows, W.J. Hausser, Jr., and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
5. Holdeman, L.V., E.P. Cato, and W.E.C. Moore (ed.). 1977. Anaerobic laboratory manual, 4th ed. Virginia Polytechnic Institute and State University, Blacksburg, Va.
6. Dowell, V.R., and T.M. Hawkins. 1979. Laboratory methods in anaerobic bacteriology, CDC manual. HEW publication (CDC) 79-8272. Center for Disease Control, Atlanta.
7. Kreig, N.R., and J.G. Holt (ed.). 1984. Bergey's manual of systematic bacteriology, vol. 1. Williams & Wilkins, Baltimore.
8. Summamen, P., E.J. Baron, D.M. Citron, C. Strong, H.M. Wexler, S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Ca.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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ATCC is a trademark of the American Type Culture Collection.
© 2018 BD. The BD Logo and all other trademarks are property of Becton, Dickinson and Company.