

Using BD® PhaSeal™ Optima System

for preparation and syringe pump administration
with small volume syringes (<10 mL)

NOTE:

- Follow aseptic technique and local guidelines for safe handling of hazardous drugs.

1 Pharmacy

Prepare the hazardous drug dose for a small volume syringe by drawing up through the BD® PhaSeal™ Optima Injector attached to an extension set.



2 Pharmacy

Clamp the extension set and **disengage** the PhaSeal™ Optima Injector from the vial assembly.



3 Nursing

Add the BD® PhaSeal™ Optima Connector to syringe infusion set. Use a syringe unit to prime the infusion set with solution and then connect to patient.



4 Nursing

Connect the PhaSeal™ Optima Injector with the medication dose to the primed syringe infusion set and PhaSeal™ Optima Connector. **Load** the small volume syringe into the syringe pump and administer per facility protocol.

